
Wasserdestillierapparate aus Glas

Glass Water Stills

2202 - 2304

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**Operating
Instructions**

2202 – 2304 **Wasserdestillierapparate aus Glas** **Glass Water Stills**

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Fully automatic **GFL Glass Water Stills** produce ultra-pure, bacteria and pyrogen free distillate with a very low conductivity. The distillate is in conformity with DAB regulations and the regulations of many international pharmacopeia.

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Before assembly, please check whether contents of package are in good order and complete.

Should you note any damages or have any reasons for complaint, please contact your supplier or directly:

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1. Before installation

The information given in the present manual must **by all means** be carefully read and observed. Only then can a perfect functioning of the Shaking Water Bath be guaranteed.

**Vital information within the manual are emphasized in bold letters.
Safety precautions are additionally marked with the following symbols**

Warning of hot liquids and steam

Warning of hot surfaces

Warning of dangerous electrical voltage

A free of charge guarantee repair cannot be granted for defects due to improper installation or handling.

In order to protect the glass parts of the Water Stills during transport, the interior of the still is filled with packing material, all free openings are covered with tape. These transport protections have to be removed before installation.

2. Location of the Water Still

The unit is suitable for both bench and wall mounting. It is only suitable for **indoor** use.

If the unit is to be placed on a bench, place on solid, even and level surfaces only.

If the unit is to be mounted on a wall, check the carrying capacity of the wall in connection with the total weight of the unit (appliance plus water filling, see technical data).

There are two keyholes for wall mounting at the back of the water still. Place two fastening screws in the wall in the distance of the two keyholes and nest the still on the keyholes.

The unit is **not** suitable for use in explosion endangered surroundings, eg during anaesthesia with inflammable gas or steam types!

3. Voltage

Main switch of the water still must be in **OFF** position.

Mains voltage and voltage stated on the name plate at the left-hand side of the unit must be identical.

If this applies, unit can be connected.

The water still must only be installed connected to a properly installed power connection with neuter (PE), according to the local regulations.

It must be secured that power can be cut all-pole (switch or socket).

Please also refer to item 6 of these operating instructions "Examples for power connections".

4. Water connections

All water connections of the water still are at the right-hand side of the unit and are marked with inscriptions.

4.1 Distillate outlet

Connect a heat-resistant laboratory hose to the spout marked "distillate outlet" and lead the hose to an external distillate storage tank.

Caution! Distilled water leaves the water still with a temperature of more than 90°C.
Danger of scaldings!

4.2 Tap water inlet

The tap water inlet supplies the still with water through a solenoid valve. Use a ½ " (inch, inner diameter 12.7 mm) pressure hose with blockable water connection (stop cock) to connect the hose connection with the water mains.

By all means secure both hose connections with hose clips!

4.3 Inlet separate water supply (only if accessory separate water supply was ordered)

The water inlet for the separate water supply supplies the boiler of the water still with softened or desalinated water through a solenoid valve. Use a ½ " pressure hose with blockable water connection to connect the hose connection with the water mains for pretreated water. After installation of the separate water supply, cooling water will be supplied via the tap water inlet.

By all means secure both hose connections with hose clips!

4.4 Cooling water outlet

Depending on the model, cooling water will flow out of the unit through one or two hose connections. The hoses for cooling water have to be led to a drain on a lower level. The hoses must have a slope on the complete length of the hoses, the cooling water must flow out without back draughts. The hose for cooling water outlet must also be ½ ", it should have a maximum length of approx 1 metre and must be heat resistant.

**Caution! Cooling water leaves
the still with a temperature of up to 70°C.
Danger of scaldings!**

4.5 Waste water outlet


When the function "Cleaning" (water exchange) is activated, water from the boiler will flow out through the connection waste water outlet. Connect a ½ " hose with a maximum length of approx 1 metre. The hose for waste water outlet **must by all means** be led to a drain on a lower level. It must have a slope on the complete length and must not be obstructed. It must be secured that the hose can draw air at the end of the cleaning procedure in order to cut off the water contact between two electrodes in the water still.

4.6 Drainage of the bi-stage boiler (only applicable for models 2302 and 2304)

Open the drain cock to withdraw mono distillate from the bi-stage boiler. During normal operation, the drain cock must be closed. Only then will the bi-stage boiler be filled and started.

**Caution! Distilled water leaves
the still with a temperature
of more than 90°C.
Danger of scaldings!**

5. Installation

After assembling all power and water connections, the cut-off valves of the mains supply (Pos 4.2 and 4.3) must be opened and the main switch must be switched on to start operation.

All cut-off valves of the mains supply should be closed when the still is not in operation.

6. Functional description

6.1 Glass water stills –standard version

After switching the unit on via the main switch, not only the green pilot lamp in the main switch will glow, but also the pilot lamp for "low water". The built-in solenoid valve is opened and releases the water supply to the boiler via the cooling coil of the condenser. The boiler is filled with water until the heating element is below the surface. The water height in the boiler is controlled by a float switch that switches the heating element on. At the same time, the pilot lamp "low water" will switch off.

After a short time, the water in the boiler will start to boil, steam will rise into the glass condenser and condenses on the cooling coil. The first of the distilled water can then be with-drawn from the distillate outlet tap of the **Mono Water Stills** (Models 2202, 2204 and 2208).

In **Bi Water Stills** (models 2302 and 2304) the distillate produced in the mono stage is used to filled the boiler of the bi stage until here the heating element is below water level, too. In the second distillation stage, too, there is a float switch that controls the water level and causes the heating elements to be switched on. When the water in the second stage boils and the steam condenses on the cooling coil, bi distillate can be withdrawn from the distillate outlet tap of **Bi Water Stills** (models 2302 and 2304).

6.2 Glass water stills with the accessory "separate water supply" (article no 2902)

The accessory "**separate water supply**" serves to feed the boiler with softened or desalinated water and the cooling coil with phosphatized or normal tap water. The quantity of inflowing industrial water is controlled via a float switch and a solenoid valve; the amount of inflowing water is limited to the amount of water that evaporates.

Owing to the fact that the industrial water is in this case not preheated in the condensers, the efficiency of the still is reduced by approx 10-15%.

6.3 Accessory "level control switch" (article no 2908)

A float switch is mounted to a plastic tube and has to be placed into a storage vessel (to be provided by the user). When the storage vessel is filled, the level control switch in connection with the electronic control of the water still switches off power and water supplies.

6.4 Function "sterilising"

When pressing the key "sterilising", the raw water supply (cooling water) is cut off for a five minute sterilising period. The steam produced by the heating element is used to sterilise the pure water area (condenser) and is allowed to escape through a vent in the lid of the water still.

Caution!
Heavy vaporization!
Danger of scaldings!

After the sterilising period the water still restarts automatically.

6.5 Function "cleaning"

Depending on the contamination level of the inflowing water and the increasing contamination of the water in the boiler, a fully automatic rinsing process should be started by pressing key **"cleaning"**. This function is also triggered when boiling contaminated water frothes in the boiler and comes into contact with the electrode in the boiler of the mono stage.

During the rinsing process, the boiler is filled with water to its upper rim and then drained completely. The contaminated water flows off through the hose connection marked **"contaminated water drain"** (please also see item 4.5 of these operating instructions). The boiler is automatically refilled with water until the heating element is below water level. The distillation process then restarts.

The purpose of the function "cleaning" is to exchange contaminated water in the mono stage boiler against clean water.

Caution! If the **separate water supply** is installed, the function "cleaning" will only work partially as the low conductivity of the industrial water (e.g. desalinated water) will prevent the water exchange being switched off. In these cases, the water still must be switched off (main switch) after the mono stage boiler is filled with water. The still must not be switched on again before the boiler has been completely drained through the contaminated water drain.

7. Servicing and Maintenance

Caution!
**Before servicing and maintenance always let
the glass water still cool down!**
Danger of scaldings!

Descaling During the rinsing process contaminated water is only exchanged against clean water. Scalings in the boiler are not removed.

Descaling the boiler.

First press key "cleaning" to start the rinsing process. When, at the end of this process, the boiler is drained, switch off the water still on the main switch.

Remove the front screen (lift by the grip holes, tilt to the front and pull down) and pour approx 100 ml commercial descaling agent (e.g. "tin-be" by Messrs Dr Otto Hartmann, 71665 Vaihingen/Enz, Germany) into the front glass funnel which is on the right-hand side next to the boiler.

After that switch on the water still and wait until the boiler is filled with water and the heating element is switched on. When the solution has reached a temperature of approx 50°C to 70°C, switch off the water still. The solution must be allowed to react for approx 1 to 2 hours with the scalings. After this time, switch on the water still again and press key "cleaning" to start the rinsing process. The solution is rinsed and exchanged against clean water, and the water still will restart normal operation.

Descaling and rinsing processes can be repeated at any time, if required.

After descaling, the first few litres of distillate should not be used as it might contain traces of the evaporated descaling agent.

As for bi water stills (models 2302 and 2304) the contents of the back boiler (second stage) should also be discarded because in this case, too, the quality of the distillate might be reduced by traces of descaling agent.

2202 – 2304 **Wasserdestillierapparate aus Glas** **Glass Water Stills**

GFL Glass Water Stills are produced with first class materials and are made to withstand even rough service conditions. Nevertheless, the units should only be subjected to rough conditions within sensible limits.

**For technical support concerning GFL Glas Water Stills,
please do not hesitate to contact our technical staff.**

Before returning a defective unit, please contact your dealer or us.

Please do not return the unit before consulting us.

If we agree to the unit being returned, arrange for careful packing and send the unit to

**GFL Gesellschaft für Labortechnik mbH
Schulze-Delitzsch-Strasse 4
D-30938 Burgwedel
Federal Republic of Germany**

Caution!
Before opening and cleaning the water still
always cut off power!
(Pull the plug or switch off main switch)
Danger of electrical shocks!

Please make sure that no fluids come into contact with cable connections or the electrical parts of the inside unit!

The pulver-coated surfaces and the front screen of the unit can be cleaned with mild detergents, if necessary.

Servicings, repairs or modifications must be carried out according to the commonly recognised Technical Rules and Regulations by competent electricians only.

Only original spare parts must be used. Always demand a detailed confirmation of the carried out tasks by the person in charge (company, date, signature).

2202 – 2304 Wasserdestillierapparate aus Glas Glass Water Stills

8. Technical Data

Exterior (W x D x H)	Models 2202/2204 Models 2208/2302/2304	650 x 200 x 390 mm 650 x 365 x 390 mm
Distillate quality	Models 2202/2204/2208 Models 2302/2304	In conformity with DAB, bacteria and pyrogen free, low gas content. approx 2.2 μS / cm at 20 °C approx 1.6 μS / cm at 20 °C
Distillation capacity	Model 2202 Model 2204 Model 2208 Model 2302 Model 2304	2 l / h mono distillate 4 l / h mono distillate 8 l / h mono distillate 2 l / h bi distillate 2 l / h bi distillate
Cooling water required	Model 2202 Models 2204 / 2302 Models 2208 / 2304	approx 48 l/h approx 96 l/h approx 144 l/h >3,5 bar / 8 bar
Water pressure min / max		
Electrical connection / Mains connection	Typ 2202	230V +/-10%, 50....60Hz, 1,5kW shock proof plug
	Typ 2204	230V +/-10%, 50....60Hz, 3,0kW shock proof plug
	Typ 2208	220V / 3 / PE, +/-10% or 400V / 3 / N / PE, +/-10% 50....60Hz, 6,0kW Mains connection cable without plug
	Typ 2302	230V +/-10%, 50....60Hz, 2,9kW shock proof plug
	Typ 2304	220V / 3 / PE, +/-10% or 400V / 3 / N / PE, +/-10% 50....60Hz, 5,8kW Mains connection cable without plug Caution! Any deviation of the mains power, even within the above stated limits, has effects on the quantity of distilled water produced.
Protection type / protection class		I / IP20
Surrounding conditions		Only inside buildings (not in explosion endangered surroundings) up to 2000 m MSL +10 °C to +40 °C max 80% rel humidity, max 31°C, decreasing to 40% rel. humidity at 40°C
Height Ambient temperature Humidity		
Weight (net, with water filling)	Models 2202/2204 Models 2208/2302/2304	16 kg / 20 kg 24 kg / 32 kg

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







Examples for connection to the mains supply.

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Colour decoding of the individual leads of the mains connection cables for models 2208 and 2304

Farbkennung - Colour decoding	Stromnetz-Mains supply 220V / 3 / PE 50/60Hz	Stromnetz - Mains supply 400V / 3 / N / PE 50/60Hz
ge/gr – gelb/grün - yellow/green	PE	PE
bl – blau - blue	L1	N
sw – schwarz - black		L1
sw – schwarz - black	L2	L2
br – braun - brown	L3	L3

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Glass Water Stills can be supplied in different versions for connection to different mains supplies.


Models 2202, 2204 and 2302 **for connection to 230V** (see nameplate) can be connected **to all mains supplies with 220V or 230V**.

Models 2208 and 2304 **for connection to 220/3** (see nameplate) can **only** be connected to mains supplies with **220V / 3 / PE**.


Models 2208 and 2304 **for connection to 400/3** (see nameplate) can **only** be connected to mains supplies with **400V / 3 / N / PE**.

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


When ordering spare parts, please **always state type and serial number** of the unit in question.

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Dechlorite Filter eliminates chlorine particles in tap water added by the local waterworks. Complete with connection for pressure hose ½ inch, with first filling.

Spare filling for dechlorite filter.

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
Phosphate cartridge prevents scale formation in the condenser by phosphatising of tap water. Complete with connection for pressure hose ½ inch, with first filling.

Spare filling for phosphate cartridge.

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The following modifications concern pages 21 – 22 of the operating instructions

(The **new** markings of the leads are set in brackets in the circuit diagrams)

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**Caution,
modification of the
lead markings
of the mains cables.**

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Effective as of January 2003, the standards (HD 308 S2 – VDE 0293-308) to mark the leads of low tension mains connection cables has been amended.

During a transition period until 01 April 2006, however, cables and leads with the **former** markings may be used.

The individual leads of **GFL Water Stills'** mains connection cables may be marked either according to the **former** or to the **new** standard. We will do our utmost to use only mains connection cables with the new markings within the shortest possible period.

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bioblock.belgium@thermofisher.com
bioblock.ventas@thermofisher.com
bioblock.vente@thermofisher.com
fisher.vendite@thermofisher.com
nl.info@thermofisher.com
info.ch@thermofisher.com

www.bioblock.be
www.es.fishersci.com
www.bioblock.com
www.it.fishersci.com
www.fishersci.nl
www.ch.fishersci.com

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