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Safety Instructions for Borosilicate 3.3 Glassware

Cleaning

Before the first use of glass-filters, they should be cleaned with diluted, warm hydrochloric acid, followed by several rinses with distilled water to remove dust and loose particles. Glassfilters should be cleaned immediately after their use. The surfaces can be wiped with a soft squeegee or brush. Surface rinsing with a wash bottle or back-flushing under the water tap is often sufficient.

Pressure and Vacuum

Even small scratches or internal abrasions of the glass surface can greatly reduce the vacuum- or pressure resistance. Filter pores may clog during filtration and cause a pressure drop. In any case 1 kp/cm² may not be exceeded. Never subject glassware to sudden pressure changes.

Temperature Resistance

The maximum allowable service temperature is approx. 500 degrees Centigrade. To avoid permanent stress, the heating and cooling rate should not exceed 5 C/min. Wet glass-filters should slowly be heated up to 80 degrees C and dried for one hour, before temperature is increased. Permanent stress can reduce the mechanical and thermal resistance. Never subject glassware to sudden temperature changes.

Chemical Resistance

High concentrations of hydrofluoric acid or hot alkali solutions attack the glass with increasing temperature and can cause a deterioration in the filter and a destruction of the pore size.

Safety Advice

When working with glassware always wear protective glasses and protective gloves to avoid injuries. Use a safety screen, -hood or similar protection when working with glassware subject to pressure or vacuum. Be prepared, that reagents may leak from a broken vessel. Please heed these guidelines in combination with the respective state-specific regulations for the use of glassware in the laboratory.

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