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**Wasser-
Destillationsgerät**

Operating Instructions

**Water Distilling
Apparatus**

Mode d'emploi

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Manual de instrucciones

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Operating Instructions Page 7 11

Important notes: Before initial operation of the Water Distilling Apparatus please read and observe carefully the operating instructions. For safety reasons the Water Distilling Apparatus may only be used for the purposes described in these present operating instructions.

Please also observe the operating instructions for the units to be connected.

All specifications in this instruction manual are guidance values which are valid at the time of printing. However, for technical or commercial reasons or in the necessity to comply with the statutory stipulations of various countries, SCHOTT may perform additions to the Water Distilling Apparatus without changing the described properties.

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Operating Instructions**Water Distilling
Apparatus**

TABLE OF CONTENTS	PAGE
Warning and safety instructions	8
Figures	9, 10
Putting into operation	11
Operation	11
Cleaning	11

Important instructions: Read these operating instructions carefully before initial operation of the device. For security reasons, the device must only be used for purposes described in these instructions. This product is subject to a constant adjustment to the state of art. In spite of the utmost accuracy it could be possible that these operating instructions do not describe the device's characteristics to its full extent. In case of doubt, please consult our technical department.

Warnings and safety instructions

To enable trouble-free operation, it is necessary to observe the following directions:

- a) Water Distilling Apparatus whose current carrying capacity is **greater than 16 amps**, are supplied without wall plugs. Connection to the mains supply is to be carried out by an authorized electrician, and the power consumption indicated on the identification plate should be observed.
- b) Before **switching on** Water Distilling Apparatus, fill boiling vessel with water (see para 3 of these Operating Instructions).
- c) Press **reset button** of overheating protection only after Water Distilling Apparatus has **cooled down** (see para 6 of these Operating Instructions).

For safety and functional reasons, the Water Distilling Apparatus may be opened only by authorized persons; e. g. work at the electric appliance must only be carried out by trained specialists. With unauthorized operation of the device as well as careless or deliberate damaging, the warranty is no longer valid.

The Water Distilling Apparatus meets the requirements of Protective Class I. It has been built and tested in accordance with DIN VDE 0700, Part 1/04.88 (HD 251 S3 and Changes 1-3) and DIN VDE 0700, Part 6/11.89 (HD 275 S1, Change 1) in connection with the draft of DIN VDE 0700, Part 6A8/04.89, Protective Measures for Electric Heating Appliances. It has left the factory in perfect technical and safety-related condition. In order to preserve this condition and ensure operation without danger, the user must observe the notices and warnings which are contained in these operating instructions.

Before switching on, it is to be ensured that the operating voltage of the Water Distilling Apparatus and the network voltage are in agreement. The network plug connector may only be introduced into a socket with a protective contact. The protective effect may not be nullified through the use of an extension cable without a protective conductor. Any interruption of the protective conductor within or external to the Water Distilling Apparatus or the loosening of the protective conductor connection can lead to a situation in which the Water Distilling Apparatus represents a danger. Intentional interruption is not allowed.

The Water Distilling Apparatus may **never** be submerged into water.

All specifications in this instruction manual are guidance values which are valid at the time of printing. However, SCHOTT-GERÄTE may make changes at any time for technical or commercial reasons or in the case of necessity to comply with the statutory stipulations of various countries.


Wasser-Destillationsgerät
Water Distilling Apparatus
Appareil à distiller l'eau
Aparato para producir agua destilada

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
7 Netzschalter
 Master switch
 Commutateur principal
 Interruptor de corriente

Adjustment of DZ 8103 Control Component

The marking lines drawn show the correct adjustment of the DZ 8103 control component, required for water level control in the water distillation apparatus.

The water outlet lower edge is required to be located at a level 1 cm higher than the heater coil upper edge. This ensures the heater coil is completely immersed in water.

In the event of the water level falling to below the heater coil upper edge, the apparatus is disconnected.

Putting into operation

1. Connect Water Distilling Apparatus as shown in picture.
2. Adjust regulating part DZ 8108 according fig. 2
3. Remove venting cap (6). Check by hand for tight fit all other plastic screw caps of Water Distilling Apparatus and, if necessary, re-tighten.
4. Open tap water valve (water faucet). As soon as water flows from water outlet (4) (after approx. 2 minutes) turn on heating (master switch (7)).
5. Distillation will start after approx. 5 minutes.

Operation

6. The Water Distilling Apparatus will have attained optimum working conditions when temperature of outflowing water (4) is of approx. 60 °C. Distillate running off is then free from CO₂ and loss of energy is small.
7. The Water Distilling Apparatus is provided with overheating protection. In case of water shortage in boiling vessel, heating turns off automatically. After restoring normal working conditions, distiller will again be ready for use after pressing **hard** on reset button (5).
8. To obtain uniformly good water quality, it is necessary to allow the salt-enriched water in the boiling vessel to run off daily. To this effect, secure again venting cap (6), after having turned off the heating, open water faucet a little more for a few seconds and then shut off completely. Boiling vessel content will siphon off. Cleaning with a detergent solution is necessary after 1 to 4 days of operation according to the salt content of the tap water.

Cleaning

9. Drain Water Distilling Apparatus as under 7 above.
10. Remove venting cap (6) and insert funnel.
11. Slowly fill in 0,5 l detergent solution (1 l detergent solution Type No. DZ 8050, 5 l detergent solution Type No. DZ 8053). Avoid heavy foaming. If necessary, interrupt filling operation.
12. Re-place in position venting cap (6) and tighten again.
13. Close off drain hose (4) (with pinch cock, for instance).
14. Refill boiling vessel with tap water up to approx. 1 cm above heating coil. Heat up to boiling point temporarily. Avoid heavy foaming.
15. Wait until water incrustation (boiler scale) has dissolved up to filling level.
16. Fill boiling vessel completely with tap water, so that the last boiler scale residue may also dissolve.
17. Open drain hose (4). Boiling vessel content will siphon off.
18. Close off drain hose (4). Fill again boiling vessel completely with tap water and allow to siphon off. We recommend to repeat the washing process five times.
19. Do not use distillate after putting the Water Distilling Apparatus into operation, until the conductivity or the pH-value of the distillate have attained the value required. Without measuring, we recommend not to use the distillate of the first hour of operation.

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SCHOTT-GERÄTE GmbH