

# Fast Analysis of Decitabine Using a Solid Core HILIC Column

Prakash Chander and Tushar N. Mehta, Centre of Excellence for Asia Pacific Laboratory  
Thermo Fisher Scientific, Ahmedabad, India

## Key Words

Decitabine,  $\alpha$ -anomer impurity, Accucore HILIC column

## Abstract

This application note demonstrates the use of the Thermo Scientific™ Accucore™ HILIC column for the determination of decitabine and its  $\alpha$ -anomer impurity by HPLC-UV.

## Introduction

Decitabine (Figure 1) is an anti-cancer (antineoplastic or cytotoxic) chemotherapy drug. Decitabine is classified as an antimetabolite and a hypomethylation agent. This application note demonstrates a simple and rapid method for the determination of decitabine and its  $\alpha$ -anomer impurity using an Accucore HILIC column.

Accucore HPLC columns use Core Enhanced Technology™ to facilitate fast and highly efficient separations. The 2.6  $\mu\text{m}$  diameter particles are not totally porous, but instead have a solid core and a porous outer layer. The optimized phase bonding creates a series of high coverage, robust phases. This coverage results in a significant reduction in secondary interactions and thus yields highly efficient peaks with very low tailing. The tightly controlled 2.6  $\mu\text{m}$  diameter of Accucore particles results in much lower backpressures than typically seen with sub-2  $\mu\text{m}$  materials.





Figure 1: Decitabine



## Experimental Details

Consumables	Part Number
Fisher Scientific™ Optima™ LC/MS grade acetonitrile	A/0626/17
Ammonium acetate (HPLC grade)	
Water, from a water purification system	
Decitabine system suitability solution	
All solutions were provided by the customer.	

Vials and Closures		Part Number
Thermo Scientific borosilicate glass vials (2 mL, 12 mm x 32 mm) with 8 mm black screw cap fitted with a silicone/PTFE seal )		60180-600
Separation Conditions		Part Number
Instrumentation:	Thermo Scientific™ Dionex™ UltiMate™ 3000	
Column:	Accucore HILIC 2.6 μm 150 mm x 3.0 mm	17526-153030
Mobile phase:	20 mM ammonium acetate in water / acetonitrile (5:95 v/v)	
Flow rate:	0.5 mL/min	
Column temperature:	40 °C	
Autosampler temperature:	10 °C	
UV detector wavelength:	244 nm	
Injection details:	5 μL partial loop	
Run time:	10 minutes	
Back pressure:	Approximately 270 bar	
Data Processing		
Data were acquired and processed using Thermo Scientific™ Chromeleon™ 7 software		

## Results

The analysis was performed on an Accucore HILIC 2.6 μm, 150 mm x 3.0 mm column. As shown in Figure 2, decitabine and α-anomer impurity were analyzed in less than 10 minutes. The results are shown in Table 1.




Figure 2: Chromatogram of decitabine and its α-anomer impurity

Compound	Retention time (min)	Retention time (% CV)	Area (% CV)
<b>Decitabine</b>	5.30	0.14	1.72
<b>α-anomer</b>	6.22	0.13	2.44

Table 1: Retention time for decitabine and α-anomer impurity

## Conclusion

The Accucore HILIC column provides enhanced retention of polar and hydrophilic analytes. Decitabine and its  $\alpha$ -anomer are easily separated using an Accucore HILIC column, which provides excellent resolution ( $R_s > 5$ ) for these compounds. This demonstrates that the Accucore HILIC column is an excellent choice of column for the rapid analysis of decitabine and related substances.

## References

Hepsiba, G.; Teja, B.B.; Kumar, K.A.; Reddy, Y.R.; Panigrahi, U.K.; Reddy, T.K. *International Journal of PharmTech Research* 2011, 3(1), 237-243.

Appleton, K.; Mackay, H.J.; Judson, I.; Plumb, J.A.; McCormick, C.; Strathdee, G.; Lee, C.; Barrett, S.; Reade, S.; Jadayel, D.; Tang, A.; Bellenger, K.; Mackay, L.; Setanoians, A.; Schatzlein, A.; Twelves, C.; Kaye, S.B.; Brown, R. Phase I and pharmacodynamic trial of the DNA methyltransferase inhibitor decitabine and carboplatin in solid tumors. *J Clin Oncol.* 2007, 25(29), 4603-4609.

For Research Use Only.

[thermoscientific.com/accucore](http://thermoscientific.com/accucore)

© 2013 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. This information is presented as an example of the capabilities of Thermo Fisher Scientific Inc. products. It is not intended to encourage use of these products in any manners that might infringe the intellectual property rights of others. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

**USA and Canada** +1 800 332 3331  
**France** +33 (0)1 60 92 48 34  
**Germany** +49 (0) 2423 9431 20 or 21  
**United Kingdom** +44 (0)1928 534110  
**Japan** +81 3 5826 1615

**China** +86 21 68654588 +86 10 84193588  
+86 20 83145199 800 810 5118  
**India** +91 22 6742 9494 +91 27 1766 2352  
**Australia** 1 300 735 292 (free call domestic)  
**New Zealand** 0800 933 966 (free call domestic)  
**All Other Enquiries** +44 (0) 1928 534 050

**Technical Support**  
**North America** +1 800 332 3331  
**Outside North America** +44 (0) 1928 534 440

**Thermo**  
SCIENTIFIC

Part of Thermo Fisher Scientific