

Fast and easy separation of 23 drugs of abuse using the Accucore Biphenyl Column

Authors: CCS Center of Excellence (CoE)
Application Scientists, Thermo Fisher Scientific

Keywords: Accucore, Biphenyl, drugs of abuse, urine

Introduction

There is an urgent need for robust analytical methods for performing research on the measurement of drugs of abuse. This short protocol describes a method for the separation of 23 drugs of abuse including high, stable resolution of isobaric opioids from human urine by UHPLC-MS/MS using the Thermo Scientific™ Accucore™ Biphenyl HPLC Column. The new Accucore Biphenyl column shows a well-tuned balance of efficiency and selectivity and is a powerful and robust tool for the determination of drugs of abuse.



Important notes

- This method was optimized for urine samples that are diluted before analysis. The same chromatographic conditions may be used for more complex samples and lower limits of detection with additional sample preparation.¹
- Low recovery for basic active compounds, like buprenorphine, may be due to interaction with vials. We recommend using Thermo Scientific™ Chromacol™ GOLD grade vials for the best performance.
- Reference the entire application.²

Materials required

- UHPLC-MS instrument such as Thermo Scientific™ Vanquish™ Horizon UHPLC System with Thermo Scientific™ TSQ Quantiva™ Triple Quadrupole Mass Spectrometer
- Thermo Scientific™ Accucore™ Biphenyl HPLC columns, 50 x 2.1 mm, 2.6 μm (P/N 17826-052130)
- Water, UHPLC-MS grade
- Methanol, UHPLC-MS grade
- Formic acid, analytical grade

Protocol

1. Dilute 100 μL urine sample with 900 μL of 0.1% formic acid in water.




Figure 1. Separation of 23 drugs of abuse under 4 minutes on the Accucore Biphenyl column

Related products

| Description | Part number |
|--|--------------|
| Accucore Biphenyl HPLC Column, 50 x 2.1 mm, 2.6 μm | 17826-052130 |
| Chromacol SureStop GOLD grade vial, clear 2 mL, screw thread | 2-SVWGK |

References

1. CS Center of Excellence (CoE) Application Scientists, Thermo Fisher Scientific, *Drugs of abuse in human plasma using SOLAμ solid phase extraction*, Application Brief 22007
2. Kean Woodmansey, Jon Bardsley and Stacy Tremintin, *Fast and easy separation of 23 drugs of abuse including high, stable resolution of isobaric opioids from human urine by UHPLC-MS/MS*, Technical Note 21883

Current versions of product instructions are available at thermofisher.com/chromexpert

Find out more at thermofisher.com/biphenyl and thermofisher.com/appslab

For Research Use Only. Not for use in diagnostic procedures. © 2020 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific. This information is presented as an example of the capabilities of Thermo Fisher Scientific products. It is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others. Specifications, terms and pricing are subject to change. Not all products are available in all locations. Please consult your local sales representative for details. **AB22004-EN 0120C**

2. Analyze samples using the following HPLC conditions:

| | | | |
|--------------------|---|----|-----|
| Mobile phase A | 0.1% formic acid in water | | |
| Mobile phase B | 0.1% formic acid in methanol | | |
| Flow rate | 0.75 mL/min | | |
| Gradient program | Time (min) | %A | %B |
| | 0.0 | 95 | 5 |
| | 0.15 | 95 | 5 |
| | 4.0 | 0 | 100 |
| | 4.0 | 0 | 100 |
| | 4.6 | 0 | 100 |
| | 4.6 | 95 | 5 |
| 5.5 | 95 | 5 | |
| Column temperature | 45 °C, with active pre-heating and forced air | | |
| Injection volume | 2 μL | | |
| Detection | HESI-MRM | | |

ThermoFisher
SCIENTIFIC