

Dual pulse labeling of S-phase population using click chemistry to measure changes in cell proliferation

Jolene A. Bradford and Scott T. Clarke

Molecular Probes® - Invitrogen Detection Technologies • 29851 Willow Creek Rd • Eugene, Oregon 97402 • USA

Introduction

The measurement of cellular proliferation is fundamental to the assessment of cell health, genotoxicity, and in the evaluation of drug efficacy. Changes in cell proliferation can give insights into mechanisms of cell growth, development, and cell cycle kinetics. A common method for detection of cell proliferation is using the incorporation of a thymidine analog during DNA synthesis. Incorporation of two analogs at different time points can further define cell cycle kinetics. Traditionally the dual pulse method has been done by combining 5-bromo-2'-deoxyuridine (BrdU) immunocytochemistry and 3H-thymidine radiography, or by combining BrdU with iododeoxyuridine (IdU) or chlorodeoxyuridine (CldU) using multiple BrdU antibodies⁴ that cross-react with either IdU or CldU for detection.

With the introduction of the new thymidine analog, 5-ethynyl-2'-deoxyuridine (Edu) with detection via click chemistry^{1,2,3}, dual pulse labeling is simplified. Using sequential pulses of the thymidine analogs Edu and BrdU, a straightforward and reliable method of distinguishing BrdU- and Edu-labeled cells by flow cytometry is introduced. This dual pulse method uses a mouse monoclonal BrdU antibody for the detection of the incorporated BrdU, which shows no cross-reactivity with Edu, clone MoBU-1. This is combined with click-chemistry detection of the incorporated Edu, which is bio-orthogonal and does not react with the incorporated BrdU. This method has broad applications and demonstrates the value of click chemistry combined with other antibody based methods.

Figure 1 – Click Chemistry Based Detection of Edu

The click reaction^{5,6} is a copper(I)-catalyzed variant of the Huisgen [3+2] cycloaddition between an azide and a terminal alkyne: a covalent bond results. In this application, the incorporated Edu contains a terminal alkyne which reacts with a fluorescently labeled azide dye, to covalently label the double stranded DNA.



Click chemistry is employed in this application because the reaction between an alkyne and an azide is bio-orthogonal; it occurs selectively, without the interference of any other functional groups present within complex biological systems. Alkyne and azide compounds are stable, inert functional groups which typically are not present in biological systems. The click reaction will detect the incorporated Edu, but does not cross-react with the incorporated BrdU.

References

- Salic, A., et al. Proc Natl Acad Sci U S A. 2008 Feb 19;105(7):2415-20.
- Buck, S., et al. BioTechniques. June 2008, 44:927-929.
- Diermeier-Daucher, S., et al. Cytometry Part A. EPUB 20 Feb 2009
- Conboy, M., et al. PLoS Biology. May 2007 5(5):1120-1126.
- Rostovtsov, V.V., et al. Angew. Chem. Int. Ed. 2001, 41, 2596.
- Tornøe, C.W., et al. J. Org. Chem. 2001, 67, 3057.
- Current Protocols in Cytometry: vol 1; Robinson, J.P., ed: John Wiley & Sons Inc. 2007.

Figure 2 – Labeling with single and dual pulses demonstrated with Jurkat T-cell leukemia cells




Figure 3 – Dual pulse labeling

Jurkat cells were treated with 10 μ M Edu for one hour and then 10 μ M BrdU is added to the media for one hour, without removal of the Edu or washing of the cells. The cells were then fixed in ETOH, and an acid denaturation method was used before labeling with anti-BrdU-Alexa Fluor® 488 dye, Click-IT®-Alexa Fluor® 647 azide, and FxCycle™ Violet stain for DNA content.





Figure 4 – TF-1 erythroleukemia cells with longer pulses



All the cells were labeled with anti-BrdU-Alexa Fluor® 488 dye, Click-IT® Alexa Fluor® 647 azide, and FxCycle™ Violet stain for DNA content. Plots show dual negative cells green, dual positive cells blue, Edu positive & BrdU negative cells red, BrdU positive and Edu negative cells yellow.

Figure 5 – THP-1 monocytic leukemia cells

Figure 5: THP-1 cells were treated with 10 μ M Edu for one hour, followed by 10 μ M BrdU for varying times of 15 minutes to 2 hours. All the cells were labeled with anti-BrdU-Alexa Fluor® 488 dye, Click-IT® Alexa Fluor® 647 azide, and FxCycle™ Violet stain for DNA content. Dual parameter plots show Edu positive cells orange, BrdU positive cells pink.

