Non-cytotoxic Near-IR DNA Stain for Cell Cycle Analysis in Living Cells

Jolene A. Bradford¹, Chris Langsdorf¹, Gayle Buller¹, Tina Luke² and Derek Davies² ¹Flow Cytometry Systems, Life Technologies, Eugene, Oregon 97402 USA ²Cancer Research UK, London, United Kingdom WC2A3PX

Wavelength (nm

Background: Cell cycle describes the progression of a cell through a cycle of division, a process resulting in cell growth and separation into two daughter cells. Live cell studies of cellular DNA content and cell cycle distribution are useful to detect variations of growth patterns due to a variety of physical, chemical, or biological means, to monitor apoptosis, and to study tumor behavior and suppressor gene mechanisms. In a given population, cells will be distributed among three major phases of cell cycle: G₀/G₁ phase (2N), S phase (DNA synthesis with variable amount of DNA), and G₂/M phase (4N). These applications require dyes that bind to DNA in a stoichiometric manner. Recent advances have produced several live cell dves useful for DNA content analysis using violet (405 nm), blue (488 nm) and green (532 nm) excitation sources with emission in visible wavelengths. We present a new dye with near-infrared emission for DNA content analysis in living cells, Vybrant® DyeCycle™ Ruby stain. This dye can be excited by the commonly available 488 nm and 633 nm excitation sources and has an emission >670 nm.

Methods: Jurkat T cells were labeled with 5µM Vybrant® DyeCycle™ Ruby stain, 5µM Hoechst 33342 or 5µM DRAQ5® (Biostatus) to compare cell health post staining with unstained control cells. Cell vitality and viability were measured after 72 hours using Calcein AM and ethidium homodimer-1. HL-60 cells were labeled with 5µM DyeCycle™ Ruby stain and sorted based on G, DNA content. Post-sort cell count and viability staining were performed to examine any cytotoxic effects. Finally, GFP-expressing cells were labeled with Vybrant® DyeCycle™ Ruby to observe compatibility.

Results: Vybrant® DyeCycle™ Ruby stain demonstrated equivalent cell viability and vitality to the gold standard Hoechst 33342, while DRAQ5™ demonstrated a significant cytotoxic effect. After sorting, no evidence of cytotoxicity was observed with Vybrant® DyeCycle™ Ruby labeled cells; good cell recovery and cell growth were observed. Vybrant® DyeCycle™ Ruby labeling proved compatible with GFPexpressing cells.

Conclusions: Vybrant® DyeCycle™ Ruby stain is a useful stain for DNA content analysis in living cells, using near-infrared emission after either 488 nm or 633 nm excitation. It demonstrates significantly less cytotoxicity than DRAQ5™ dye, allowing it to be used to determine the proliferation rate of living cells and offers the possibility of cell sorting based on DNA content. Additionally Vybrant® DyeCycle™ Ruby stain is compatible with GFP fluorescence.

Figure 1: Spectral Characteristics

The fluorescence excitation and emission spectra the Vybrant® DyeCycle[™] Ruby stain are shown. These spectra were obtained from samples of the dye bound to DNA. The Vybrant® DyeCycle[™] Ruby stain/DNA complex has fluorescence excitation and emission maxima of 638/686nm respectively. This is compatible with excitation lasers from 488 - 640 nm.

Figure 2: DNA content analysis with various configurations

Plots A and B: Histogram of live Jurkat cells stained with Vybrant® DyeCycle[™] Ruby stain showing DNA content distribution, G₂/G₁ and G₂/M phase histogram peaks are separated by the S-phase distribution. **Panel A** shows the distribution of this population of cells when 488 nm excitation was used with a 695/40 bandpass filter. **Panel B** shows the same population when 633 nm excitation was used with a 695/40 bandpass filter.



Plots C-G: Live HL-60 cells labelled with 5uM DyeCycle[™] Ruby where dead cells are excluded using DAPI, run on a BD[™] LSRII looking at fluorescence in several detectors: Panel C uses 488 nm excitation with 660/20 bandpass emission, **Panel D** uses 488 nm excitation with 695/50 bandpass emission, **Panel E** uses 635 nm excitation with 660/20 bandpass emission, **Panel F** uses 635 nm excitation with 735/45 bandpass emission, and **Panel G** uses 635 nm excitation with 736/45 bandpass emission, and **Panel G** uses 635 nm excitation with 736/45 bandpass emission, Panel Panel



Figure 3: Viability and Vitality of cells 72 hours after labeling

Unlabeled control Jurkat T cells (A) were compared to cells labeled with either SµM hoechst 3324 (B), SµM BRAG5TH (C), or SµM Vybrant® DyeCycleTH Ruby stain (D) and allowed to grow for 72 hours in RPMI media + 10% FBS at 37 °C with 5% CO₂. After 72 hours, cell health was measured using Calcein AM as a measure of cell vitality and ethicium homodimer-1 as a measure of cell vibality. Samples were acquired using an BDTH LSRII flow cytometer using 488 nm excitation with 530/30 bandpass for Calcein fluorescence and 610/20 bandpass for ethicidum homodimer-1 fluorescence. Dual parameter pits are displayed on the left, and a brightfield image of each cell culture is displayed on the right. Cells labeled with the gold standard Hoechst 33342 while cells labeled with DyeCyteTH Ruby stain have viability and vitality comparable to cells labeled with the gold standard





DyeCycle™ Ruby Live cells = 71.2 % Dead cells = 26.4 %

References:

 Current Protocols in Cytometry, 7.0.1–7.27.7 (2004); 2. Practical Flow Cytometry, 4th Ed., Shapiro HM, Ed. (2003); 3. Methods Mol Biol 281, 301 (2004); 4. Cytometry A 58, 21 (2004).

Figure 4: Post sort recovery after Vybrant® DyeCycle™ Ruby stain labeling

Both Jurkat and HL-60 cells were labeled using 5 X10^e total cells, suspended in 10mL media and added 10µI DyeCycle[™] Ruby stain, incubated for 15 minutes at 37°. C. Cells were pelleted and resuspended in 1mL media, adding 5µI DAPI at 200µg/ml. Live cells were gated on DAPI-negative cells, followed by a singlet gate based on area and height before sorting the G/G, population and collecting 1.6 million cells of each. Cells were then cultured and evaluated for cell count (left graph) and viability (right graph). No evidence of cytotoxicity was seen and good recovery and growth post sort was observed.



Figure 5: Compatibility with GFP

DT40 chicken B cells, either wild-type or stably transfected with a GFP reporter construct, were labelled with DyeCycleTM Ruby stain. Data was acquired on a BDTM FACSCaliburTM instrument using 488mm excitation and >670 nm emission. The staining is compatible with GFP.



