



**Thermo Scientific
Barnstead Nanopure®
Ultrapure Water Systems**

Pure water. Pure research. Total confidence.

Thermo Scientific Barnstead Nanopure Ultrapure Water Systems

For more than 125 years, Barnstead pure water systems have been a trusted resource for science and industry. Now the most comprehensive portfolio of water systems is sold as Thermo Scientific and enhanced by our worldwide capabilities delivering:

Expertise: our unique Thermo Scientific H₂O SELECT™ testing program will analyze your feed water to match the best water system to your application

Value: key features come standard and low cost-of-ownership means reduced equipment and consumables costs

Quality: water quality validated by independent labs to meet your most stringent application requirements



Barnstead Nanopure Systems

Critical applications in today's laboratories require exceptional performance.

Our Nanopure solutions deliver the expertise, quality and value you need – offering four models with advanced purification technologies to suit your specific research needs.



Nanopure Analytical

Our analytical unit produces ultrapure water suitable for general laboratory applications – such as reagent preparation, high-purity rinsing, as well as for use in AA, ICP and ICP/MS.

- 18.2 MΩ-cm
- 5-10 ppb TOC



Nanopure Analytical UV

Dual wavelength UV oxidation is incorporated into our analytical UV system to reduce organic carbon to the absolute lowest levels. Typical applications include ion chromatography, gas chromatography, HPLC, LC/MS, GC/MS, and TOC analysis.

- 18.2 MΩ-cm
- 1-5 ppb TOC



Nanopure Biological UF

A hollow-fiber ultrafilter (UF) provides water for use in applications requiring pyrogen-free water, including cell and tissue culture, monoclonal antibody production and electrophoresis.

- 18.2 MΩ-cm
- 5-10 ppb TOC
- < 0.001 Eu/ml pyrogens



Nanopure Life Science UV/UF

The hollow-fiber ultrafilter and UV oxidation technologies are combined into one system for applications requiring the lowest levels of organic carbon, pyrogens and nucleases, including PCR, 2-D electrophoresis and cell culture.

- 18.2 MΩ-cm
- 1-5 ppb TOC
- < 0.001 Eu/ml pyrogens

Pure water.

Digital Display

High-visibility display provides an accurate reading of resistivity, conductivity, temperature, and TOC.

Removable Display and Controls

Removable display can be mounted up to 10 feet (3 meters) away from the unit for added flexibility.

Automatic Cleaning Cycle

Simply introduce the cleaner into the system, press a button and walk away; our Nanopure will take care of the rest.

Whisper-Quiet Pump

Provides optimum flow with virtually silent operation.



Choose Your Method of Dispensing

Intermittent, constant or volumetric dispensing of water, based on your needs.

Hollow-Fiber Final Filter

The 0.2 µm final filter requires less frequent replacements than comparable units.

Thermo Scientific Durapak™ Cartridge Pack

High-capacity for long life. Easy to install or change cartridges.

Flexibility in Mounting Location

Locate your system on the bench or wall. Wall mounting bracket included as standard equipment.

Built-in Pressure Regulating Valve

Allows feed water pressures up to 100 psig (6.9 bar) with no additional equipment.

Accudispense Volumetric Dispenser

Patented volumetric technology remotely dispenses a programmed volume of water (250 ml - 60 L) from up to eight feet (2.5 meters) away. Complete recirculation up to the tip of the dispenser eliminates contamination of your water. Easily mounts anywhere in your lab.



Remote Dispenser

Conveniently dispenses ultrapure water from up to eight feet (2.5 meters) away. Unique design virtually eliminates any dead volume, and coiled tube allows for easy use while taking up a minimal amount of space.



Pure research. Total confidence.

Complete Control

- **Continuous recirculation** maintains the highest water purity during periods of inactivity.
- **Built-in pressure regulating valve** ensures optimum system pressures and flows.
- **Automatic system notification** indicates when water quality falls below the programmable set point. Simply program the value you prefer, and the alert will appear on the display when the purity drops below the set value.
- **One-touch automatic dispensing**, volumetric or timed – your choice. Our Nanopure system allows automatic dispensing of volumes between 0.25 and 60 liters.

Intuitive Display

- **Push-button controls** for resistivity, conductivity, TOC, and temperature. Readings are compensated to 25°C. In addition, users can comply with USP specifications and display uncompensated resistivity.
- **Multilingual display** provides information in English, Spanish, German, and French. Japanese language available on select units.

Easy to Maintain

- **Diagnostic menu** evaluates important system components and alerts if a problem exists.
- **Durapack cartridge pack** simplifies cartridge installation and replacement.
- **NIST traceable calibration module** avoids costly downtime and service with on-site performance verification of resistivity meter.



**Durapack
Cartridge Pack**

**Remote Display
and Controls**

Unmatched Expertise and Quality.

FREE H₂O SELECT Water Analysis

Our unique H₂O SELECT testing program provides expert analysis of your feed water. Based on the results, we will identify the best water system for your application. You will receive a comprehensive report including equipment and consumables recommendations.

3 EASY STEPS:

- 1** Request a FREE H₂O SELECT* kit from your sales representative or visit www.thermo.com/select
- 2** Fill the H₂O SELECT bottle with your feed water sample and complete the brief form
- 3** Send your water sample and completed form back to our water experts for a FREE in-depth analysis

Our H₂O SELECT analysis program is free and with no obligation.



Your custom recommendation includes:


- **Summary report** of your feedwater test results
- **System recommendation** based on your feedwater quality, laboratory applications, daily water demand, and budget
- **Estimated life expectancy** of consumables to help you project system operation costs

*H₂O SELECT program may not be available in all countries. Please check with your local sales representative.

Exceptional Value.

Maximizing value is critical for today's laboratories. Our Nanopure systems deliver significant savings – both in capital investment and ongoing operating costs.


Capital Equipment Comparison



Up to 20% savings on initial purchase price with key accessories included as standard equipment

- > **Remote display controls** provide flexibility in choosing a mounting location for the system or control panel
- > **Pressure-regulating valve** allows feed-water pressures up to 100 psig with no additional equipment
- > **Ready-to-mount wall bracket** enables easy bench- or wall-mounting without the purchase of additional parts

Annual Operating Cost Comparison



Up to 50% savings on annual operating costs

- > **High-capacity cartridge pack** has 20% more capacity than alternatives, delivering fewer change-outs and lower costs
- > **UV lamp and ultrafilter replacements** cost 50% less per year than alternative water systems

* Savings based on internal performance data, when compared with alternative systems

Thermo Scientific Barnstead Nanopure Technology Advantage

Intelligent technology to maximize water quality and optimize system operation.

Durapak Cartridge Pack

Our high-capacity cartridge packs are easily installed with one connection. The pack uses two types of activated carbon to absorb both volatile and non-volatile organic compounds and chlorine, and the highest purity semiconductor-grade ion-exchange resins for removing the last traces of ionic contaminants.

Systems fed with deionized water (DI) have their own pack – designed to remove problematic organic and colloidal compounds associated with central DI feed.

Final Filtration

Final filtration is performed with an absolute 0.2-micron gamma irradiated filter. The naturally hydrophilic hollow fibers provide high surface area for longer filter life and fast rinse up to desired purity levels. The fibers are made of cellulose acetate and are encased in a clear polycarbonate housing.


Simple Sanitization

All Nanopure systems are easily sanitized. A liquid cleaning solution is injected by syringe and the software manages the rest of the sanitization process.

Other Materials of Construction

All wetted parts beyond the purification cartridge pack are constructed of low-extractable, corrosion-resistant materials, including the tubing and fittings, which are made of fluoropolymers. The cell well and cartridge pack are constructed of virgin polypropylene, and the cell is constructed of corrosion-resistant titanium.

*Nanopure Analytical
system flow diagram*



UltraFiltration

For applications requiring undetectable levels of pyrogens, such as cell and tissue culture, where ultrafiltration is essential.

Ultrafiltration removes the last traces of pyrogens (bacterial endotoxins) and nucleases. The encapsulated ultrafilter uses unique, polysulfone hollow fibers to remove particulates, nucleases and pyrogens. The filter is located downstream of the purification media and UV chamber for maximum benefit.


Fine hollow UF fibers are folded and secured on the product side of the capsule. Water entering the capsule flows from the outside of the fibers into their hollow cores. The channels merge to become the purified product stream. The system periodically sends water from the outside of the fibers to the drain – thereby removing filtered contaminants and extending filter life.

Fast Rinse-up

The filters are shipped dry without chemical preservatives. This allows fast rinse-up to high-purity water.

Low Extractables

There is almost no organic or inorganic addition to the water as demonstrated by low-TOC and high-resistivity measurements.



UF Hollow Fiber

Water flows from the outside of the fiber into the fiber lumen.

UV Oxidation

For ultra-low TOC (Total Organic Carbon) concentrations and germicidal action. The addition of UV oxidation is ideal for analytical applications, which require the absolute lowest levels of organic carbon.

UV Chamber

A UV oxidation chamber is placed between the third and fourth (final) cartridges. The chamber is constructed of highly polished 316 stainless steel and a UV-transparent quartz sleeve. Within the sleeve resides the mercury vapor lamp.

Germicidal Action


The mercury vapor lamp creates a high fraction of light at 254 nm, which keeps bacterial levels in the system extremely low.

Organic Carbon Oxidation

The lamp also creates a high fraction of 185 nm light which, in conjunction with the 254 nm radiation, produces hydroxyl free radicals ($\cdot\text{OH}$). The radicals quickly oxidize residual organic carbon to CO_2 and water.

Organic Removal Cartridge

A specially formulated media removes oxidation by-products (carbon dioxide and organic intermediates) to produce water virtually free of TOC. A final fraction of semiconductor-grade mixed-bed resin is the last step in producing the highest purity water.



The UV Oxidation Process

Thermo Scientific Barnstead Nanopure Technology Advantage

Total Organic Carbon (TOC) Analyzer

By monitoring TOC and resistivity, you can be sure that the product water exceeds the most stringent requirements of analytical or biological applications.

To ensure accurate measurement of product-water quality before it is dispensed from the system, the TOC analyzer is positioned after all purification technologies, next to the resistivity cell.

Specifications

Volume per measurement: 10–15 ml per reading dispensed to drain

Cycle interval: approximately 3.5-minute cycle between readings

Control: TOC monitoring can be turned on or off via system keypad


Range: 1–250 ppb (@ ≥ 10 M Ω -cm resistivity)

Resolution: 1 ppb

Accuracy: ± 1 ppb or 15% of reading, whichever is greater


Water temperature range: 4–40°C system (40–104°F)

Ambient temperature range: 10–40°C



Analyzer Principles of Operation

Water from our Nanopure passes through a filter (1) and enters the UV reactor (2). The resistivity of the water is measured initially. The high-intensity (185 and 254 nm) UV lamp oxidizes organic compounds present in the water, producing CO₂. The solenoid valve (3) opens and moves the water across the conductivity cell (4), where resistivity is measured again. The difference in resistivity resulting from the increased CO₂ concentration is used to calculate the TOC value, which is then displayed.



*Nanopure Life Science (UV/UF)
with TOC analyzer*

Specifications & Ordering Information

	Nanopure Analytical	Nanopure UV	Nanopure UF	Nanopure Life Science (UV/UF)
Model #				
Standard Unit	D11901	D11911	D11921	D11931
TOC Unit	D11941	D11951	D11961	D11971
Cartridge Pack:				
Distilled/Reverse Osmosis Feed	D50282	D50280	D50282	D50280
Deionized Feed	D50283	D50281	D50283	D50281
0.2µm Final Filter:	included	included	included	included
UV Lamp:	n/a	included	n/a	included
Ultrafilter:	n/a	n/a	included	included
Specifications:				
Resistivity @25°C	18.2 MΩ-cm	18.2 MΩ-cm	18.2 MΩ-cm	18.2 MΩ-cm
TOC (typical)	5-10 ppb	1-5 ppb	5-10 ppb	1-5 ppb
Pyrogen	NA	NA	<0.001 EU/ml	<0.001 EU/ml
Bacteria	<1 cfu/ml	<1 cfu/ml	<1 cfu/ml	<1 cfu/ml
Particles (0.2 micron)	< 1/ml	< 1/ml	< 1/ml	< 1/ml
Silica	<1 ppb	<1 ppb	<1 ppb	<1 ppb
UV Absorbance @254 nm	<0.001 AU	<0.001 AU	<0.001 AU	<0.001 AU
Flow Rate	1.8 lpm	1.8 lpm	1.8 lpm	1.8 lpm

Feed-Water Requirements:

Required Pretreatment	Reverse Osmosis, Deionization, Distillation
Total Organic Carbon	< 1 ppm
Turbidity	< 1.0 NTU
Pressure	100 psig (Gravity feed to 6.9 bar)
Temperature	40–104°F (4–40°C)
Total Dissolved Solids	< 70 ppm as CaCO ₃
Silica	< 1 ppm
Silt Density Index	< 5%

Automatic Volumetric Dispense:

Range	0.25–60 L
Accuracy	± 5% (1 L or greater)
Repeatability	± 3%

Automatic Timed Dispense:

1–40 minutes

Dimensions:

13.5" W x 19.5" H x 17.0" D (34.3 cm x 49.5 cm x 43.2 cm)

Mounting:

Bench or wall (brackets included)

Clearances:

Sides	> 9" (22.9 cm) for servicing
Above	> 3" (7.6 cm) for removal of the top cover
Front	> 4.75" (12.1 cm) for opening front door

Electrical:

90–240 VAC, 100 VA, 47–63 Hz, 1 phase

Tubing Connections:

Inlet	0.375" (11.43 cm) O.D. tubing (supplied)
Atmospheric Drain	0.25" (7.62 cm) O.D. tubing (supplied)

ACCESSORIES



Accudispense



Remote Dispenser



Printer

Model No.	Description
D13661	Accudispense – Remotely dispenses programmed volumes of water (0.25 and 60 liters) from up to eight feet (2.5 meters) away
D11981	Remote Dispenser – Easily attaches to any Nanopure, allowing you to deliver water up to eight feet (2.5 meters) away from the unit
AY1137X1	Printer, 120 or 240V – Date/Time/Resistivity and Temperature (TOC on TOC models) printed every 12 minutes
CMX25	Cleaning Syringe
E896X5	NIST Calibration module - Automatic calibration verification of the purity meter

Our complete line of water purification technologies includes solutions for your most critical and everyday application needs, from tap water pretreatment, deionization, distillation, filtration, and reverse osmosis, as well as storage reservoirs and other accessories.

Exceptional Pretreatment with Thermo Scientific Barnstead TII

Ideal for any daily application needs for primary grade water, our Barnstead TII is fed directly with tap water. It is the most intelligent and flexible system for producing ASTM Type II water. It combines pretreatment, reverse osmosis, deionization, and UV oxidation in a single system that is simple to operate and easy to maintain.



© 2008 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

North America: USA/Canada + 1 800-553-0039

Europe: Austria +43 1 801 40 0, Belgium +32 2 482 30 30, Finland +358 9 329 100, France +33 2 2803 2000, Germany national toll free 08001-536 376, Germany international +49 6184 90 6940, Italy +39 02 95059 1, Netherlands +31 76 571 4440, Russia/CIS +7 (812) 703 42 15, Spain/Portugal +34 93 223 09 18, Switzerland +41 44 454 12 12, UK/Ireland +44 870 609 9203

Asia: China +86 21 6865 4588 or +86 10 5850 3588, India +91 22 5542 9494,

Japan +81 45 453 9220, Other Asian countries +852 2885 4613

Countries not listed: +49 6184 90 6940 or +33 2 2803 2000

www.thermo.com/purewater

Thermo
SCIENTIFIC