Dual-Light[®] Chemiluminescent Reporter Gene Assay for Firefly Luciferase and β-Galactosidase

- Detect firefly luciferase and β-galactosidase in a single sample
- · Rapid combined assay format
- Achieve greater precision in normalization of transfection efficiency



Figure 1. Detection of Firefly Luciferase and β-Galactosidase with Dual-Light® System

Application

- Reporter quantitation and transfection normalization for transiently transfected cell lines and primary cells
- Multiplex analysis of two experimental promoters
- In combination with the pMIR- REPORT[™] miRNA Expression Reporter Vector System to evaluate the interaction between miRNAs and their target sites

Dual-Light[®] Chemiluminescent Reporter Gene Assay system provides rapid and sensitive detection of luciferase and β-galactosidase in a single sample extract. Two reporter genes-an experimental reporter (typically luciferase) and a constitutively expressed reporter used for transfection control (typically β-galactosidase)—are often used to accurately quantitate activity from experimental reporter constructs. The Dual-Light system combines both assays into one convenient simple assay using only one cell extract sample, which gives greater precision in normalization of transfection efficiency than can be achieved with two separate assays.

The assays are performed sequentially. First the luciferase is quantitated with an enhanced luciferase reaction, then β -galactosidase reporter enzyme activity is measured with Galacton-Plus[®] chemiluminescent β -galactosidase substrate. The light signal is measured in a luminometer. The entire assay is completed in less than one hour.

Product Configuration

Dual-Light [®] System Standard Size	Dual-Light [®] System Large Size	Dual-Light [®] System Screening Size
Г1003	T1005	T1004
Capacity: 200 combined assays	Capacity: 600 combined assays	Capacity: 4,000 combined assays
Contents:	Contents:	Contents:
0.2 mL Galacton-Plus® substrate	 0.6 mL Galacton-Plus[®] substrate 	 4 mL Galacton-Plus[®] substrate
5 mL Buffer A*	• 3 x 5 mL Buffer A*	• 20 x 5 mL Buffer A*
22 mL Buffer B*	• 3 x 22 mL Buffer B*	• 20 x 22 mL Buffer B*
25 mL Light Emission Accelerator-II	 75 mL Light Emission Accelerator-II 	 500 mL Light Emission Accelerator-II
70 mL Lysis Solution	 210 mL Lysis Solution 	 1.4 L Lysis Solution
		*Supplied lyophilized

70 mL Lysis Solution	• 210 mL Lysis Solution	 1.4 L Lysis Solution *Supplied lyophilized
Buffers A and B are supplied lyophilized		
Features	Benefits	
Combined assay format	Convenience	
	Precision	
	Speed	
High sensitivity	Detect low levels of reporter enzyme	
	produced by weak promoters or cells with low	
	transfection efficiency	
Wide dynamic range:	Eliminates need for multiple sample dilutions	
Femtogram to nanogram		

Femtogram to nanogram	
	Enables efficient use of sample and assay reagents
Chemiluminescent assay format	High sensitivity—detect femtograms of enzyme

ORDERING INFORMATION

Package Size	P/N	Number of Combined
		Microplate Assays
Standard Size	T1003	400
Large Size	T1005	1200
Screening Size	T1004	10,000

For Research Use Only. Not for use in diagnostic procedures.

© Copyright 2007, Applied Biosystems. All rights reserved. AB (Design), Applied Biosystems, Dual-Light, Galacton-Plus, and Tropix are registered trademarks and pMIR-REPORT is a trademark of Applera Corporation or its subsidiaries in the U.S. and/or certain other countries. All other trademarks are the sole property of their respective owners. Printed in the USA, 09/2007 Publication 120PB19-01



Headquarters

850 Lincoln Centre Drive | Foster City, CA 94404 USA Phone 650.638.5800 | Toll Free 800.345.5224 www.appliedbiosystems.com

International Sales

For our office locations please call the division headquarters or refer to our Web site at www.appliedbiosystems.com/about/offices.cfm