



S ensitive array-based comparative genomic hybridization (CGH) starts with a robust labeling system. The BioPrime® Array CGH Genomic Labeling System is a fast, efficient labeling kit that allows you to generate high yields of labeled genomic samples with strong signal intensities, improving the sensitivity of your CGH experiments.

## Genomic Profiling with the BioPrime<sup>®</sup> Array CGH Genomic Labeling System

The recent emergence of array-based CGH enables researchers to perform genome-wide studies of chromosomal aberrations that lead to a variety of diseases. By differentially labeling normal and tumor DNA samples prior to hybridization to a BAC or cDNA array (Figure 1), you can detect variations in gene copy number that underlie diseases such as cancer (1,2). The ability to detect variations in gene copy number depends largely upon the sensitivity of the labeling system used. The BioPrime<sup>®</sup> Array CGH Genomic Labeling System is optimized for use in array CGH. Using this system, you can expect:

- High yields of labeled genomic DNA
- Strong signal intensities
- Reproducible results

With robust yields and strong signal intensities, you'll obtain the sensitivity you need to perform accurate genomic profiling experiments. And with an optimized kit, you'll make sure you obtain consistent performance.



#### Figure 1 - Array CGH sample preparation overview



# **Brighter signals**

The exo-Klenow polymerase in the BioPrime<sup>®</sup> Array CGH Genomic Labeling System incorporates fluorescently modified nucleotides more effectively than standard Klenow, enabling you to obtain stronger signal intensities (Figure 2).

Figure 2 - Increased signal-to-noise ratios using the BioPrime® Array CGH Genomic Labeling System



With strong signal intensities, you'll obtain the sensitivity you need to detect variations in gene copy number reliably. Figures 3a and 3b show the resulting ratios when genomic DNA from a tumor cell line with four copies of the X chromosome and genomic DNA from a normal female tissue are labeled using the BioPrime<sup>®</sup> Array CGH Genomic Labeling System and co-hybridized to a BAC array. The

Comparison of signal to noise ratios of hybridizations to BAC arrays using 1 µg of genomic DNA labeled using current methods with standard Klenow compared to the BioPrime<sup>®</sup> Array CGH Genomic Labeling System.

> fluorescent ratio of the 4X and 2X samples for chromosome 1 is approximately 1.0 (expected), indicating equivalent copy numbers of chromosome 1 in tumor and normal female genomic DNA (Figure 3a). The median ratio of the samples for the X chromosome is 1.5, indicating an increase in copy number between the tumor cell line and normal genomic DNA (Figure 3b).





Determination of differential gene copy number with a tumor cell line DNA on BAC arrays. Analysis of the fluorescent ratios of normal female genomic DNA and cell line with four copies of the X chromosome on BAC arrays using the BioPrime<sup>®</sup> Array CGH Genomic Labeling System.

### Improved yields

Current labeling methods for array CGH use a high concentration of the large fragment (Klenow) of DNA polymerase I, which lacks 5'-3' exonuclease activity (3,4). This enzyme incorporates fluorescent nucleotides while replicating the genomic DNA sample, providing a simple and efficient means to label genomic targets prior to hybridization. The BioPrime<sup>®</sup> Array CGH Genomic Labeling System uses a high concentration of exo- Klenow, lacking both 5'-3' and 3'-5' exonuclease activity. Using exo- Klenow enables you to obtain higher yields of labeled sample than possible with standard Klenow (Figure 4), increasing the sensitivity of your genomic profiling experiments.

# Figure 4 - Higher yields of labeled sample using BioPrime<sup>®</sup> Array CGH Genomic Labeling System



Average total DNA recovered after labeling 2  $\mu$ g of digested genomic DNA using current protocol (standard Klenow) and the BioPrime\* Array CGH Genomic Labeling System. Duplicate reactions were incubated at 37° C for 2 hours.

Figure 5 - Overview of BioPrime® Array CGH Genomic

Labeling System protocol

# Simple protocol

The BioPrime<sup>®</sup> Array CGH Genomic Labeling System is based on a simple random priming protocol and takes less than three hours to complete, requiring minimal "hands on" time. Briefly, sheared or digested genomic DNA is denatured and random octamers annealed to the template strands. The exo-Klenow fragment then extends from the random priming sites, incorporating fluorescent nucleotides during polymerization. The labeled samples are purified on spin columns to remove free nucleotides. The fluorescently labeled material is then heat denatured briefly before hybridization to a microarray (Figure 5).

### Sheared or digested genomic DNA. Denature double-stranded genomic DNA and anneal random primers. Exo-Klenow fragment extends from random priming sites, incorporating fluorescent nucleotides. Remove free nucleotides and other contaminants using spin columns. Denature and cohybridize labeled genomic targets. along with blocking agents, to microarray. $\mathbf{1}$ Fluorescent nucleotide Exo-Klenow polymeras random primer www.invitrogen.com

# **Complete labeling solution**

The BioPrime<sup>®</sup> Array CGH Genomic Labeling System contains all the major reagents (except fluorescent nucleotides) necessary for fast and efficient labeling of

genomic DNA samples (Table 1). All reagents are functionally tested so you'll obtain quality performance with every labeling reaction.

Exo-Klenow Fragment (40 U/µl)	Purification Buffers (A and B)
2.5X Random Primer Solution	Purification Columns
Control DNA (Salmon Sperm)	Stop Buffer
10X dCTP Nucleotide Mix	Distilled Water
10X dUTP Nucleotide Mix	Collection Tubes

Table 1	- BioPrime®	Array CGH	Genomic	Labeling	System	components
---------	-------------	-----------	---------	----------	--------	------------

### Sensitivity made simple

Genomic target preparation is an integral part of obtaining accurate results during genomic profiling experiments. Using the BioPrime<sup>®</sup> Array CGH Genomic

#### Description

BioPrime<sup>®</sup> Array CGH Genomic Labeling System Human Cot-1 DNA<sup>®</sup> Mouse Cot-1 DNA<sup>®</sup> Yeast tRNA

#### **References:**

- 1. Pinkel D., et. al. (1998) Nature Genetics 20: 207-211.
- 2. Pollack J.R., et. al. (1999) Nature Genetics 23: 41-46.
- 3. http://cmgm.stanford.edu/pbrown/protocols/4\_genomic.html (Stanford)
- 4. http://www.sanger.ac.uk/HGP/methods/cytogenetics/Labelling.shtml (Sanger Institute)

successful array CGH studies. Call and order today.

Labeling System, you'll get the sensitivity you need for

Quantity	Cat. no.		
30 rxns	18095-011		
500 µg	15279-011		
500 µg	18440-016		
50 mg	15401-011		





These products may be covered by one or more Limited Use Label Licenses (See the Invitrogen catalog or our web site, www.invitrogen.com). By the use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses.

For research use only. Not intended for any animal or human therapeutic or diagnostic use. Printed in the U.S.A. ©2003 Invitrogen Corporation. All rights reserved. Reproduction forbidden without permission.

Corporate headquarters:

1600 Faraday Avenue • Carlsbad, CA 92008 USA • Tel: 760 603 7200 • Fax: 760 602 6500 • Toll Free Tel: 800 955 6288 • E-mail: tech\_service@invitrogen.com www.invitrogen.com European headquarters:

Invitrogen Ltd • Inchinnan Business Park • 3 Fountain Drive • Paisley PA4 9RF, UK • Tel: +44 (0) 141 814 6100 • Fax: +44 (0) 141 814 6260 • E-mail: eurotech@invitrogen.com