

## eBioscience™ Fura-2 AM Dye

Catalog Number: 65-0858

For Research Use Only. Not for use in diagnostic procedures.

### Product Information

**Contents:** Fura-2 AM Dye  
**Catalog Number:** 65-0858



**Formulation:** Lyophilized yellow solid. Purity ≥95% as determined by HPLC.

**Temperature Limitation:** Store at -20°C. Protect from light and moisture.

**Batch Code:** Refer to vial

**Use By:** Refer to vial



### Description

Fura-2 is a calcium imaging dye that binds to free  $\text{Ca}^{2+}$ . Fura-2 AM is the cell-permeable acetoxymethyl (AM) ester form of Fura-2. Fura-2 is the preferred dye for ratiometric imaging microscopy with digital image analysis, especially when the alteration of excitation wavelengths is more practical than the detection of multiple emission wavelengths. Upon binding  $\text{Ca}^{2+}$ , the excitation spectrum of Fura-2 shifts to shorter wavelengths between 300 and 400 nm, while the peak emission remains steady around 510 nm. The  $K_d$  of Fura-2 is highly dependent on pH, temperature, ionic strength and viscosity of the cytosol, thus great care should be taken when non-standard conditions are used.

Molecular weight: 1001.86

Peak excitation: variable depending on the concentration of free  $\text{Ca}^{2+}$ , between 300 and 400 nm

Peak emission: 510 nm

Fura-2 AM should be reconstituted in high-quality, freshly opened DMSO. Once reconstituted, it should be protected from light and stored at -20°C. Avoid freeze-thawing.

### Applications Reported

Fura-2 AM has been reported for use in fluorescence microscopy. It is recommended that the concentration used be carefully determined by each investigator for optimal performance in the assay of interest.

### References

Barreto-Chang OL, Dolmetsch RE. Calcium imaging of cortical neurons using Fura-2 AM. J Vis Exp. 2009 Jan 19;(23).

Kang JJ, Toma I, Sipos A, Peti-Peterdi J. From in vitro to in vivo: imaging from the single cell to the whole organism. Curr Protoc Cytom. 2008 Apr;Chapter 12:Unit 12.12.

Malgareoli A, Milani D, Meldolesi J, Pozzan T. Fura-2 measurement of cytosolic free  $\text{Ca}^{2+}$  in monolayers and suspensions of various types of animal cells. J Cell Biol. 1987 Nov;105(5):2145-55.

G Grynkiewicz, M Poenie, RY Tsien. A new generation of  $\text{Ca}^{2+}$  indicators with greatly improved fluorescence properties. J. Biol. Chem. Mar 25, 1985; 260: 3440-3450.

Not for further distribution without written consent.

Copyright © 2016 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.

Tel: 888.999.1371 or 858.642.2058 • Fax: 858.642.2046 • thermofisher.com/ebioscience •  
info@ebioscience.com