






Human IFN gamma ELISPOT

Catalog Number 88-7386

Pub. No. MAN0017483 Rev. A.0 [30]

 **WARNING!** Read the Safety Data Sheets (SDSs) and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves. Safety Data Sheets (SDSs) are available from thermofisher.com/support.

Product information

Symbol	Contents	Human IFN gamma ELISPOT
	Catalog number	88-7386
	Temperature limitation	Store at 2–8°C
	Batch code	Refer to vial
	Use by	Refer to vial
	Caution	Contains sodium azide

Description

This Human IFN gamma ELISPOT reagent set contains the necessary reagents for performing enzyme-linked immunosorbent spot (ELISPOT) assays for high-resolution frequency analysis of IFN gamma-secreting cells. This ELISPOT reagent set is pre-titrated for optimal spot development.

Time requirements

- 1 overnight incubation
- 1–2 days cell activation
- 3–5 hours washing, antibody incubations, and color development

Components

- **Capture Antibody:** Pre-titrated, functional grade (low endotoxin) purified antibody
- **Detection Antibody:** Pre-titrated, biotin-conjugated antibody
- **ELISA/ELISPOT Coating Buffer**
ELISA/ELISPOT Coating Buffer Powder or 10X PBS ELISPOT Coating Buffer
- **5X ELISA/ELISPOT Diluent**
- **Enzyme**
Pre-titrated Avidin-HRP Concentrate
- **Certificate of Analysis:** Lot-specific instructions for dilution of antibodies and enzyme

Other materials needed

- **Reagents:**
 - 96-Well PVDF Membrane ELISPOT Plates (Milipore, Cat. No. MAIPS4510)
 - AEC (3-amino-9-ethyl carbazole) Substrate (Sigma, Cat. No. A-5754)
 - Distilled water (dH₂O)
 - ELISPOT Wash Buffer: ELISA/ELISPOT Wash Buffer Powder (Cat. No. 00-0400) or 1X PBS with 0.05% Tween™ 20
 - Complete RPMI-1640

– 1X PBS

• Instruments:

- Pipettes and pipettors
- Refrigerator
- Incubator
- Laminar flow hood
- Plate washer: Wash bottle or automated wash machine
- ELISPOT plate reader or dissecting microscope for visual inspection

Reagent preparation

• ELISA/ELISPOT Coating Buffer (1X)

Make a 1:10 dilution of 10X PBS ELISPOT Coating Buffer in dH₂O or reconstitute ELISA/ELISPOT Coating Buffer Powder to 1L in dH₂O and sterilize by filtering using a 0.22-μm filter.

• Complete RPMI-1640

Complete RPMI-1640 with 10% Fetal Bovine Serum and 1% Penicillin/Streptomycin/L-Glutamine.

• 5X ELISA/ELISPOT Diluent

Dilute 5X ELISA/ELISPOT Diluent 1:5 in dH₂O.

• ELISA/ELISPOT Wash Buffer

Reconstituted ELISA/ELISPOT Wash Buffer Powder (Cat. No. 00-0400) or 1X PBS with 0.05% Tween™ 20 (0.5 mL Tween™ 20 in 1 L PBS).

• 0.1M Acetate Solution (pH 5.0)

- Combine 148 mL of 0.2 M acetic acid (11.55 mL glacial acetic acid per 1 L of dH₂O) with 352 mL of 0.2 M sodium acetate (27.2 g sodium acetate per 1 L of dH₂O).
- QS to 1 L with dH₂O, then adjust pH to 5.0.

• AEC (3-amino-9-ethyl carbazole) Substrate Solution

- Prepare AEC Stock Solution by dissolving 100 mg of AEC in 10 mL of N, N Dimethylformamide (DMF; Cat. No. 20672).
- Add 333 μL of AEC Stock Solution to 10 mL of 0.1M Acetate Solution (pH 5.0), then filter through a 0.45-μm filter.
- Add 5 μL of 30% H₂O₂ just before use. Mix and use immediately.

Experimental procedures

Aseptic steps

Note: Use sterile buffers and aseptic techniques and perform all steps in a laminar flow hood.

- Dilute Functional Grade purified capture antibody in sterile ELISA/ELISPOT Coating Buffer, as noted on the Certificate of Analysis.
- Coat ELISPOT plate with 100 μL/well of capture antibody solution, then incubate at 2–8°C overnight.
- Decant or aspirate coating antibody from plate.
- Wash plates 2 times with 200 μL/well of sterile ELISA/ELISPOT Coating Buffer, then decant.
- Block plate with 200 μL/well of complete RPMI-1640 at room temperature for 1 hour, then decant or aspirate plate.
- Aliquot mitogen, antigen, or controls diluted in complete RPMI-1640 to appropriate wells at 100 μL/well.

7. Aliquot cells at desired densities (e.g., 1×10^5 – 2×10^6 cells/mL) at 100 μ L/well and incubate at 37°C, in a 5% CO₂ humidified incubator for 24–48 hours.

Note: Optimal kinetics and cell densities vary with target cytokine, treatment, and cell type and must be empirically determined. Cells can be diluted in a sterile tissue culture plate starting at 2×10^6 cells/well in triplicate wells with a series of 1:3 or 1:4 serial dilutions down the plate, and then transferred to the ELISPOT plate.

Non-aseptic steps

1. Decant cells and medium from plates.
2. Wash plate 3 times with ELISA/ELISPOT Wash Buffer.
3. Dilute biotinylated detection antibody in Assay Diluent according to instructions on the Certificate of Analysis.
4. Add 100 μ L/well of detection antibody solution, then incubate at room temperature for 2 hr (or at 2–8°C overnight).
5. Decant antibody solution, then wash 4 times with ELISA/ELISPOT Wash Buffer allowing wells to soak for 1 minute for each wash.
6. Dilute Avidin-HRP reagent in Assay Diluent according to instructions on the Certificate of Analysis.
7. Add 100 μ L/well of Avidin-HRP solution, then incubate at room temperature for 45 minutes.
8. Decant Avidin-HRP solution.
9. Wash plate 3 times with ELISA/ELISPOT Wash Buffer.
10. Wash 2 times with 1X PBS without Tween™ 20.
11. Add 100 μ L/well of freshly-prepared AEC Substrate Solution, then develop at room temperature for 10–60 minutes monitoring development of spots.

12. Stop the substrate reaction by washing wells 3 times with 200 μ L/well of dH₂O.
13. Air-dry the plate, then count spots using a dissecting microscope or automated ELISPOT plate reader.
Store plates in the dark prior to reading.

Customer and technical support

Visit thermofisher.com/support for the latest in services and support, including:

- Worldwide contact telephone numbers
- Product support, including:
 - Product FAQs
 - Software, patches, and updates
 - Training for many applications and instruments
- Order and web support
- Product documentation, including:
 - User guides, manuals, and protocols
 - Certificates of Analysis
 - Safety Data Sheets (SDSs; also known as MSDSs)

Note: For SDSs for reagents and chemicals from other manufacturers, contact the manufacturer.

Limited product warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale found on Life Technologies' website at www.thermofisher.com/us/en/home/global/terms-and-conditions.html. If you have any questions, please contact Life Technologies at www.thermofisher.com/support.



Manufacturer: Bender MedSystems GmbH | Campus Vienna Biocenter 2 | 1030 Vienna, Austria

The information in this guide is subject to change without notice.

DISCLAIMER: TO THE EXTENT ALLOWED BY LAW, LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE, OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.

Important Licensing Information: This product may be covered by one or more Limited Use Label Licenses. By use of this product, you accept the terms and conditions of all applicable Limited Use Label Licenses.

©2017 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.