Thermo

PRODUCT INFORMATION Thermo Scientific EpiJET DNA Methylation Analysis Kit (Taql/HpyF30I)

#K1451	200 rxns
Lot	Expiry Date _
5'тс с	a3'
3' a g c	T 5'

Store at -20°C

www.thermoscientific.com/onebio

Kit components

Component	#K1451 200 rxns
Epi HpyF30I	200 µL
Epi Taql	200 µL
10X EpiBuffer 2	1.2 mL
Control DNA Unmethylated (pUC19/Smal) (0.5 μg/μL)	20 µL
Control DNA CpG Methylated (mpUC19/Smal) (0.5 µg/µL)	20 µL

CERTIFICATE OF ANALYSIS

The kit is tested in Epi HpyF30I and Epi Taql digestion of methylated and unmethylated pUC19/Smal DNA mixed with human blood genomic DNA. Unmethylated plasmid DNA is completely cleaved in one hour with both enzymes in presence of genomic DNA. Methylated plasmid DNA is not cleaved by Epi HpyF30I, but cleaved with Epi Taql in presence of genomic DNA.

```
Quality authorized by: Jurgita Zilinskiene
```

Rev.4

Description

5-methylcytosine is a prominent epigenetic DNA modification which plays an important role in up- and down-regulation of gene expression. In mammalian cells 5-methylcytosines are mainly found in 5'-CpG-3' dinucleotides and are enriched in CpG islands, the GCrich regions most often located in gene promoters. Mspl/Hpall enzyme pair with 5'-CCGG-3' recognition sequence has multiple targets in such regions and is routinely used for cytosine methylation analysis. Outside the CpG islands the CCGG motives are much less abundant and limit the usability of the aforementioned enzyme pair. The Thermo Scientific[™] EpiJET[™] DNA Methylation Analysis Kit (Taql/HpyF30I) uses another pair of restriction enzymes to analyze the DNA methylation status in gene bodies and other low-GC content genomic regions. Epi HpyF30I and Epi Tagl are isoschizomers with TCGA specificity. Methylation of cytosine in the 5'-TCGA-3' tetranucleotide blocks the cleavage with Epi HpyF30I but does not affect the cleavage with Epi Tagl.

The Epi HpyF30I and Epi TaqI enzymes are specially formulated for the epigenetic studies to complete genomic DNA digestion in 1 hour.

Application

CpG methylation analysis at the 5'-TCGA-3' loci.

Important Notes

- To minimize the possibility of pipetting errors during the DNA digestion, prepare enough reaction master mix (including DNA sample, but not including Epi HpyF30I and Epi TaqI) for three reactions plus 10% extra. Mix thoroughly and dispense the mixture into three sterile tubes (19 μ L each). Add 1 μ L of nuclease-free water into tube 1, add 1 μ L of Epi HpyF30I into tube 2, add 1 μ L of Epi TaqI into tube 3 and follow protocol recommendations.
- For control reactions prepare two distinct reaction master mixes (one with Control DNA Unmethylated (pUC19/Smal) and the other with Control DNA CpG Methylated (mpUC19/Smal)). Mix thoroughly and dispense each mixture into three sterile tubes (19 μ L each). Add 1 μ L of nuclease-free water into tubes 1 and 4, add 1 μ L of Epi HpyF30l into tubes 2 and 5, add 1 μ L of Epi Taql into tubes 3 and 6 and follow protocol recommendations.
- DNA sample quality is important for efficient digestion by restriction endonucleases. We recommend use of spin column based kits such as the Thermo Scientific GeneJET[™] Genomic DNA Purification Kit (#K0721/2) for genomic DNA purification.

• Reactions can be scaled up, incrementally increasing reaction volume by 10 μ L and enzyme amount by 1 μ L for every microgram of DNA used. For example, to digest 2 μ g of genomic DNA use 2 μ L of Epi HpyF30l or Epi Taql in 30 μ L reaction volume. If digesting less than 1 μ g, or if DNA concentration is high enough (\geq 125 ng/ μ L) to accommodate 1 μ g into total volume of 8 μ L, reaction volume can be scaled down to 10 μ L.

Protocol

• Read the Important Notes section before starting.

1. Mix the following components in sterile

microcentrifuge tubes:

Component	Undigested DNA	Digested with Epi HpyF30I	Digested with Epi Taql
	1	2	3
Sample DNA	up to 1 µg*	up to 1 µg*	up to 1 µg*
10X Epi Buffer 2	2 µL	2 µL	2 µL
Epi HpyF30I	-	1 µL	-
Epi Taql	-	-	1 µL
Water, nuclease-free (#R0581)	to 20 µL	to 20 µL	to 20 µL
Total volume	20 µL	20 µL	20 µL

2. Incubate all samples at 60°C for 1 hour.

- 3. Terminate the reactions by incubation at 90°C for 10 min.
- To prepare protein free DNA, extract with phenol/chloroform.

*For methylation analysis by qPCR use 0.04-0.4 μ g of sample DNA.

Control Reaction

- To assess DNA sample quality and potential inhibitory effects, perform control reaction using provided unmethylated and methylated control DNA.
- Read the Important Notes section before starting.

1. Mix the following components in sterile microcentrifuge tubes:

Component	Undigested DNA	Digested with Epi HpyF30I	Digested with Epi Taql
	1	2	3
Sample DNA	0.5 µg	0.5 µg	0.5 µg
Control DNA Unmethylated (pUC19/Smal)	0.5 µg	0.5 µg	0.5 µg
10X Epi Buffer 2	2 µL	2 µL	2 µL
Epi HpyF30I	-	1μL	-
Epi Taql	-	-	1 µL
Water, nuclease- free (#R0581)	to 20 µL	to 20 µL	to 20 µL
Total volume	20 µL	20 µL	20 µL

Component	Undiges ted DNA	Digested with Epi HpyF30I	Digested with Epi Taql
	4	5	6
Sample DNA	0.5 µg	0.5 µg	0.5 µg
Control DNA CpG Methylated (mpUC19/Smal)	0.5 µg	0.5 µg	0.5 µg
10X Epi Buffer 2	2 µL	2 µL	2 µL
Epi HpyF30I	-	1 µL	-
Epi Taql	-	-	1 µL
Water, nuclease- free (#R0581)	to 20 µL	to 20 µL	to 20 µL
Total volume	20 µL	20 µL	20 µL

2. Incubate all samples at 60°C for 1 hour.

3. Terminate the reactions by incubation at 90°C for 10 min.

4. Analyze the reaction products by DNA electrophoresis on an agarose gel.

Typical results are presented in Fig 1. Jurkat genomic DNA (uppermost band) is only partially digested by Epi HpyF30I and fully fragmented by Epi Taql that is not sensitive to CpG methylation. Unmethylated plasmid DNA is digested by both enzymes, while methylated plasmid DNA is susceptible only to Epi Taql digestion. Any inhibition of control DNA cleavage by restriction endonucleases indicates the presence of inhibitors. In such a case repurification of sample DNA is recommended



pUC 19 DNA/Smal mpUC 19 DNA/Smal

Figure 1. Digestion of genomic and control DNA by Epi Taql and Epi HpyF30I

Lanes 1 and 4 – undigested DNA,

Lanes 2 and 5 – DNA digested with Epi HpyF30I, Lanes 3 and 6 – DNA digested with Epi TaqI, Lanes 1-3 contain Control DNA Unmethylated (pUC19/Smal), Lanes 4-6 contain Control DNA CpG Methylated (mpUC19/Smal). Genomic DNA in all lanes is Jurkat oDNA.

Guidelines for gPCR

If CpG methylation at a particular TCGA locus is analyzed by qPCR, follow the recommendations provided below.

- For qPCR set primers that flank the TCGA site of interest.
- Use 1-2 μL (2-40 ng) of template prepared according to the Protocol and follow the manufacturer's recommendations for qPCR reaction set-up and cycling conditions.
- During qPCR setup, it is important to avoid DNA cross-contamination. We recommend using a dedicated set of pipettes for qPCR to minimize contamination.
- The accuracy of qPCR is highly dependent on accurate pipetting and thorough mixing of individual reaction components. Take extra care to avoid pipetting errors during qPCR set up and when preparing templates for qPCR. Use of 2-3 technical replicates is highly recommended.
- For calculation of 5-mC percentage use the formula provided in the Calculation of 5-hmC Percentage section. To determine the PCR efficiency value prepare a standard calibration curve by diluting the "Undigested DNA" sample.

Recommendations for qPCR

- The parameters below are recommended for qPCR using Thermo Scientific Maxima[™] SYBR Green/ROX qPCR Master Mix (#K0221/2/3).
- To minimize the possibility of pipetting errors, prepare a reaction master mix by adding the following components (excluding template DNA) for each 20 µL reaction to a tube at room temperature:

10 µL
0.3 µM
0.3 µM
1-2 µL (2-40
ng)
to 20 µL
20 µL

- Mix the reaction master mix thoroughly and dispense the desired volume into PCR tubes or plates.
- Add template DNA (1-2 µL (2-40 ng)) prepared according to the Protocol to the individual PCR tubes or wells containing master mix. Centrifuge briefly if needed.

Note. To minimize inaccuracies associated with pipetting smaller volumes we recommend to dilute template DNA 3-fold and use $3-6 \mu L$.

• Program the thermal cycler according to the recommendations below:

Step	Temperature, °C	Time	Number of Cycles
Initial denaturation	95	10 min	1
Denaturation	95	15 sec	
Annealing	60*	30 sec	40
Extension	72	30 sec	
* use the optimal appealing temperature for your primer pair			

* use the optimal annealing temperature for your primer pair

 Use two-step cycling protocol, if the optimal annealing temperature for your primer pair is 60°C:

Step	Temperature, ℃	Time	Number of Cycles
Initial denaturation	95	10 min	1
Denaturation	95	15 sec	
Annealing/ Extension	60	60 sec	40

Note. Melting curve analysis may be performed to verify the specificity of the PCR product. Primer-dimers may occur during PCR if the primer design is not optimal. The dimers are distinguished from the specific product through their lower melting point.

Interpretation of qPCR Results

For interpretation of qPCR results see the picture provided below.







Cq1 - threshold cycle of "Undigested DNA" sample (continuous curve)

Cq2 - threshold cycle of "Digested with Epi HpyF30I" DNA sample (dashed curve) Cq3 - threshold cycle of "Digested with Epi Tagl" DNA

sample (dotted curve)

* Cq3 – Cq1 < 4.5 means the incomplete digestion of sample DNA by Epi Taql. Repurify sample DNA or extend the digestion time. ** \approx stands for Δ Cq \leq 0.5

Calculation of 5-mC Percentage:

• If Cq3 – Cq1 ≥ 4.5, % of 5-mC is calculated using the formula below:

% of 5-mC =100 / (1 + E)Cq2-Cq1

Where:

Cq1 is the threshold cycle of "Undigested DNA" sample Cq2 is the threshold cycle of "Digested with Epi HpyF30I" sample Cq3 is the threshold cycle of "Digested with Epi TaqI" sample E is the PCR efficiency value (%)

- % of 5-mC can not be accurately estimated when:
- Cq3 Cq1 < 4.5.

This result demonstrates the incomplete digestion of sample DNA by Epi Taql. Repurify sample DNA or extend the digestion time. For high quality genomic DNA purification use commercially available gDNA purification kits, such as GeneJET Genomic DNA Purification Kit (#K0721/2).

• Cq1 \approx Cq2 \approx Cq3.

This result demonstrates the absence of Epi HpyF30I/Epi Taql recognition sequence.

Patent pending

PRODUCT USE LIMITATION

This product is developed, designed and sold exclusively for research purposes and in vitro use only. The product was not tested for use in diagnostics or for drug development, nor is it suitable for administration to humans or animals. Please refer to <u>www.thermoscientific.com/onebio</u> for Material Safety Data Sheet of the product.

O 2013 Thermo Fisher Scientific, Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.