

# InsectSelect<sup>™</sup> BSD System

For Stable Expression of Heterologous Proteins in Lepidopteran Insect Cell Lines using pIB/V5-His

Catalog nos. K820-01, K825-01, V8020-01

**Version I** 07 October 2008 25-0330

## **Table of Contents**

Kit Contents and Storage	iv
Accessory Products	vi
Introduction	1
Overview	1
Methods	3
Culturing Insect Cells	3
Cloning into pIB/V5-His	4
Transforming E. coli	6
Transient Expression in Insect Cells	7
Selecting Stable Cell Lines	12
Scale-Up and Purification	16
Troubleshooting	18
Troubleshooting	18
Appendix	19
Recipes	19
Blasticidin S	21
Map of pIB/V5-His	22
Features of pIB/V5-His	23
Map of pIB/V5-His/CAT	24
OpIE2 Promoter	25
OpIE1 Promoter	26
Technical Support	27
Purchaser Notification	28
References	30

## **Kit Contents and Storage**

#### **Types of Kits**

This manual is supplied with the products listed below.

Product	Catalog no.
InsectSelect <sup>™</sup> BSD System with Sf9 Cells	K820-01
InsectSelect <sup>™</sup> BSD System with High Five <sup>™</sup> Cells	K825-01
pIB/V5-His Vector Kit	V8020-01

#### Shipping/Storage

See the table below for shipping and storage information.

Product	Shipping	Storage
InsectSelect <sup>™</sup> BSD System	Dry ice	Vectors, primers, blasticidin: -20°C
with Sf9 Cells		Cellfectin® Reagent: <b>4°C</b>
InsectSelect <sup>™</sup> BSD System		Cells: Liquid nitrogen
with High Five™ Cells		Medium: 4°C, protected from light
pIB/V5-His Vector Kit	Wet ice	-20°C

#### **Vectors and** Primer

The vectors and primers are supplied with all of the kits listed above. Store at -20°C.

Product	Amount Supplied	Composition
pIB/V5-His	20 μg	40 μl of 0.5 μg/μl pIB/V5-His in 10 mM Tris–HCl, 1 mM EDTA, pH 8.0
pIB/V5-His/CAT	20 μg	40 μl of 0.5 μg/μl pIB/V5-His/CAT control vector in 10 mM Tris–HCl, 1 mM EDTA, pH 8.0
OpIE2 Forward Sequencing Primer	2 μg	Lyophilized in TE, pH 8.0
OpIE2 Reverse Sequencing Primer	2 μg	Lyophilized in TE, pH 8.0

**Primer Sequences** The sequence of each primer is provided below:

Primer	Sequence	pMoles Supplied
OpIE2 Forward	5′-CGCAACGATCTGGTAAACAC-3′	329
OpIE2 Reverse	5′-GACAATACAAACTAAGATTTAGTCAG-3′	250

### Kit Contents and Storage, Continued

#### **Blasticidin**

Blasticidin is supplied with the InsectSelect<sup>™</sup> BSD System kits only.

Blasticidin is shipped as 50 mg (powder). Upon receipt, **store Blasticidin at room temperature**.

Blasticidin is available to order separately (see page vi).

#### Cellfectin<sup>®</sup> Reagent

Cellfectin® Reagent is supplied with the InsectSelect™ BSD System kits only.

Cellfectin® Reagent is shipped as 1 ml (composition: 1 mg/ml lipid in membrane-filtered water).

Upon receipt, store Cellfectin® Reagent at 4°C.

Cellfectin® Reagent is available to order separately (see page vi).

#### **Cells and Medium**

**Supplied with the InsectSelect**<sup>™</sup> **BSD System kits only.** Additional cells and other cell lines are available to order separately (see page vii).

Store the cells in liquid nitrogen.

Store the medium at 4°C, protected from light.

Different cells and media are included, depending on which InsectSelect  $^{\text{\tiny{TM}}}$  BSD System kit you ordered. Refer to the table below. For guidelines and instructions to culture Sf9 and High Five  $^{\text{\tiny{TM}}}$  cells, refer to the Insect Cell Lines manual included with each kit.

Kit	Cells	Medium
InsectSelect™ BSD System with Sf9 Cells	Sf9	<ul> <li>Sf-900 II SFM (1X) or</li> <li>Grace's Insect Cell Culture Medium, Unsupplemented (contains L-glutamine)</li> </ul>
	High Five <sup>™</sup>	Express Five® Serum-Free Medium

#### **Manuals**

The following manuals are supplied with each kit.

Kit	Manual
InsectSelect <sup>™</sup> BSD System with Sf9 Cells	InsectSelect <sup>™</sup> BSD System manual
	Insect Cell Lines manual
InsectSelect <sup>™</sup> BSD System with High	InsectSelect <sup>™</sup> BSD System manual
Five <sup>™</sup> Cells	Insect Cell Lines manual
pIB/V5-His Vector Kit	InsectSelect <sup>™</sup> BSD System manual

## **Accessory Products**

#### Introduction

The products listed in this section are intended for use with the InsectSelect  $^{\text{TM}}$  BSD System and the pIB/V5-His Vector Kit. For more information, refer to our web site at <u>www.invitrogen.com</u> or contact **Technical Support** (page 27).

## Additional Products

The following products are available separately from Invitrogen.

Product	Quantity	Catalog no.
Sf9 Cells, frozen	1 ml vial, $1 \times 10^7$ cells/ml	B825-01
Sf21 Cells, frozen	1 ml vial, $1 \times 10^7$ cells/ml	B821-01
High Five™ Cells, frozen	1 ml vial, $3 \times 10^6$ cells/ml	B855-02
TOP10 Electrocomp <sup>™</sup> Cells	5 × 80 μl	C664-55
	$10 \times 80 \mu l$	C664–11
	$30 \times 80 \mu l$	C664-24
One Shot™ TOP10 (chemically competent cells)	21 × 50 μl	C4040-03
Electrocomp™ TOP10	5 × 80 μl	C664–55
Grace's Insect Cell Culture Medium, Unsupplemented	500 ml	11595-030
Sf-900 II SFM	1 liter	10902-088
Express Five® SFM	1 liter	10486-025
Cellfectin® Reagent	1 ml	10362-010
Blasticidin S	50 mg	R210-01
S.N.A.P. Miniprep Kit	100 reactions	K1900-01
PureLink™ HiPure Plasmid Miniprep Kit	100 preps	K2100-03
Positope <sup>™</sup> Control Protein	5 μg	R900-50

### **Accessory Products, Continued**

#### Additional InsectSelect<sup>™</sup> Kits

Several other kits that allow you to clone and stably express your gene of interest using the InsectSelect<sup>™</sup> technology are available from Invitrogen. These kits include InsectSelect<sup>™</sup> vectors with different antibiotic resistance genes and those that allow secreted protein expression. In addition, the pIZT/V5-His Vector Kit enables expression of a gene of interest and a cycle 3-GFP/Zeocin<sup>™</sup> fusion gene. This allows both visual monitoring of transfection efficiency and generation of a stable cell line. For more information about the various InsectSelect<sup>™</sup> vector kits available from Invitrogen, visit our web site at <a href="www.invitrogen.com">www.invitrogen.com</a> or contact **Technical Support** (page 27). See the table below for ordering information.

Product	Catalog no.
pIB/V5-His TOPO® TA Expression Kit	K890-01
pIZ/V5-His Vector Kit	V8000-01
InsectSelect <sup>™</sup> System with Sf9 Cells	K800-01
InsectSelect <sup>™</sup> System with High Five <sup>™</sup> Cells	K805-01
pIZT/V5-His Vector Kit	V8010-01
pMIB/V5-His Vector Kit	V8030-01

#### Polyacrylamide Gel Electrophoresis

To facilitate separation of your recombinant protein by polyacrylamide gel electrophoresis, a wide range of pre-cast NuPAGE® and Novex® Tris-Glycine polyacrylamide gels are available from Invitrogen. In addition, Invitrogen also carries a large selection of molecular weight protein standards and staining kits. For more information about the appropriate gels, standards, and stains to use to visualize your recombinant protein, refer to our web site at <a href="https://www.invitrogen.com">www.invitrogen.com</a> or contact **Technical Support** (page 27).

### **Accessory Products, Continued**

#### Detecting Recombinant Proteins

Expression of your recombinant fusion protein can be detected using an antibody to the appropriate epitope. The table below describes the antibodies available for detection of C-terminal fusion proteins expressed using pIB/V5-His. Horseradish peroxidase (HRP) or alkaline phosphatase (AP)-conjugated antibodies allow one-step detection using colorimetric or chemiluminescent detection methods.

The amount of antibody supplied is sufficient for 25 Westerns.

Product	Epitope	Catalog no.
Anti-V5 Antibody	Detects 14 amino acid epitope	R960-25
Anti-V5-HRP Antibody	derived from the P and V	R961–25
Anti-V5-AP Antibody	proteins of the paramyxovirus, SV5 (Southern <i>et al.</i> , 1991)	R962-25
	GKPIPNPLLGLDST	
Anti-His (C-term) Antibody	Detects the C-terminal	R930–25
Anti-His(C-term)-HRP Antibody	polyhistidine (6xHis) tag (requires the free carboxyl group	R931–25
Anti-His(C-term)-AP Antibody	for detection (Lindner <i>et al.,</i> 1997) HHHHHH-COOH	R932–25

#### Purifying Recombinant Protein

The metal binding domain encoded by the polyhistidine tag allows simple, easy purification of your recombinant protein by Immobilized Metal Affinity Chromatography (IMAC) using ProBond™ Resin available from Invitrogen. The ProBond™ Purification System and Ni-NTA Purification System are also available from Invitrogen. See the table below for ordering information.

Product	Quantity	Catalog no.
ProBond™ Purification System	6 purifications	K850-01
ProBond <sup>™</sup> Purification Kit with Anti-V5-HRP Antibody	1 kit	K854-01
ProBond <sup>™</sup> Purification Kit with Anti-His(C-term)-HRP Antibody	1 kit	K853-01
ProBond <sup>™</sup> Nickel-Binding Resin	50 ml	R801-01
	150 ml	R801-15
Purification Columns	50	R640-50
(10 ml polypropylene columns)		
Ni-NTA Purification System	6 purifications	K950-01
Ni-NTA Purification System with Anti-V5-HRP Antibody	1 kit	K954-01
Ni-NTA Purification System with Anti-His(C-term)-HRP Antibody	1 kit	K953-01
Ni-NTA Agarose	10 ml	R901-01
	25 ml	R901–15

#### Introduction

#### **Overview**

#### Introduction

The InsectSelect™ System allows you to express your protein of interest in insect cell lines either transiently or stably. The system utilizes a single expression vector, pIB/V5-His to express your gene of interest. This 3.5 kb vector has the following features:

- *OpIE2* promoter for constitutive expression of the gene of interest (Theilmann & Stewart, 1992)
- *OpIE1* promoter for expression of the blasticidin resistance gene (see next bullet) (Theilmann & Stewart, 1991)
- Blasticidin resistance gene for selection of stable cell lines (Takeuchi *et al.*, 1958; Yamaguchi *et al.*, 1965)
- EM7 promoter for expression of ampicillin (or blasticidin) resistance in *E. coli*
- Ampicillin resistance gene for selection of transformants in *E. coli*
- Optional C-terminal peptide containing the V5 epitope and 6xHis tag for detection and purification of your protein of interest

For more information and a map of the vector, see pages 22–23.

#### Description of System

The gene of interest is cloned into pIB/V5-His and transfected into Sf9 or High Five™ cells using lipid-mediated transfection. After transfection, cells can be assayed for expression of the gene of interest. Once you have confirmed that your gene expresses, you can select for a stable polyclonal population or stable clonal cell lines using blasticidin as a selection agent. Stable cell lines can be used to express the protein of interest in either adherent culture or suspension culture.

## Description of Promoters

Baculovirus immediate-early promoters utilize the host cell transcription machinery and do not require viral factors for activation. Both the *OpIE2* and *OpIE1* promoters are from the baculovirus *Orgyia pseudotsugata* multicapsid nuclear polyhedrosis virus (*OpMNPV*). The virus' natural host is the Douglas fir tussock moth; however, the promoters allow protein expression in *Lymantria dispar* (LD652Y), *Spodoptera frugiperda* cells (Sf9) (Hegedus *et al.*, 1998; Pfeifer *et al.*, 1997), Sf21 (Invitrogen), *Trichoplusia ni* (High Five™) (Invitrogen), *Drosophila* (Kc1, S2) (Hegedus *et al.*, 1998; Pfeifer *et al.*, 1997), and mosquito cell lines (unpublished data). The *OpIE2* promoter has been shown to be about 5- to 10-fold stronger than the *OpIE1* promoter (Pfeifer *et al.*, 1997). Both promoters have been sequenced and analyzed. For more detailed information on the *OpIE2* and *OpIE1* promoters, see page 25 and page 26, respectively.

#### **Expression Levels**

The *OpIE2* promoter provides relatively high levels of constitutive expression, although not all proteins will express as high as might be expected from baculovirus late promoters such as polyhedrin or very late promoters such as p10 (Jarvis *et al.*, 1996). However, some researchers have found that the InsectSelect System expresses some proteins better than baculovirus systems. To date, reported expression levels range from 1–2  $\mu$ g/ml (human IL-6; Invitrogen) to 8–10  $\mu$ g/ml (human melanotransferrin) (Hegedus *et al.*, 1999).

### Overview, Continued

#### Blasticidin Resistance

Blasticidin S HCl is a nucleoside antibiotic isolated from *Streptomyces griseochromo-genes* which inhibits protein synthesis in both prokaryotic and eukaryotic cells (Takeuchi *et al.*, 1958; Yamaguchi *et al.*, 1965). Resistance is conferred by expression of either one of two blasticidin S deaminase genes: *BSD* from *Aspergillus terreus* (Kimura *et al.*, 1994) or *bsr* from *Bacillus cereus* (Izumi *et al.*, 1991). These deaminases convert blasticidin S to a non-toxic deaminohydroxy (Kimura & Yamaguchi, 1996; Yamaguchi *et al.*, 1975).

## Experimental Outline

The table below describes the general steps needed to clone and express your gene of interest using the InsectSelect $^{\text{\tiny{M}}}$  BSD kit of your choice. For more details, refer to the manual and the pages indicated.

Step	Action	Source
1	Establish culture of Sf9 or High Five <sup>™</sup> cells <b>Note:</b> Other cell lines ( <i>e.g.</i> Sf21) may be used.	Refer to the Insect Cell Lines manual or use your own protocols.
2	Develop a cloning strategy to ligate your gene of interest into pIB/V5-His.	Pages 4–6, this manual
3	Transform your ligation reactions into a <i>rec</i> A, <i>end</i> A <i>E. coli</i> strain ( <i>e.g.</i> TOP10). Select on LB plates containing 50–100 µg/ml ampicillin.	Page 6, this manual
4	Isolate plasmid DNA and sequence your recombinant expression vector to confirm that your protein is in frame with the C-terminal peptide.	Page 6, this manual
5	Transfect Sf9 or High Five <sup>™</sup> cells.	Pages 7–9, this manual
6	Assay for transient expression of your protein.	Pages 10–11, this manual
7	Create stable cell lines expressing the protein of interest by selecting with blasticidin.	Pages 12–15, this manual
8	Scale-up expression for purification.	Page 16, this manual
		Insect Cell Lines manual
9	Purify your recombinant protein by chromatography on metal-chelating resin (e.g. ProBond $^{\text{TM}}$ ).	Pages 16–17, this manual

#### **Methods**

## **Culturing Insect Cells**

#### Introduction

Before you start your cloning experiments, be sure to have cell cultures of either Sf9 or High Five<sup>™</sup> cells growing and have frozen master stocks available. If you purchased one of the InsectSelect<sup>™</sup> BSD System kits, you will receive either Sf9 cells or High Five<sup>™</sup> cells along with the Insect Cell Lines manual. Use this manual as a guide to initiate cell culture. This manual is available for downloading from our web site at <a href="https://www.invitrogen.com">www.invitrogen.com</a> or by contacting **Technical Support** (page 27).

#### Culturing Sf9 and High Five<sup>™</sup> Cells

To culture Sf9 or High Five<sup>™</sup> cells, refer to the Insect Cell Lines manual. This manual covers the following topics:

- Thawing frozen cells
- Maintaining and passaging cells
- Freezing cells
- Using serum-free medium
- Growing cells in suspension
- Scaling up cell culture



For the best recovery and viability, thaw High Five<sup>™</sup> cells into Express Five<sup>®</sup> Serum-Free medium and thaw Sf9 cells into Sf-900 II SFM (1X) (see page 8)

#### Sf21 Cells

You may also use Sf21 cells as a host for pIB/V5-His. Sf21 cells are larger and we have found that they may produce more protein than Sf9 cells. Refer to the Insect Cell Lines manual for more information.

#### Cells for Transfection

You will need log-phase cells with >95% viability to perform a successful transfection. Review pages 7–10 to determine how many cells you will need for transfection.

## Cloning into pIB/V5-His

#### Introduction

This chapter provides information to help you clone your gene of interest into pIB/V5-His. A diagram is provided on page 5 to help you ligate your gene of interest in frame with the C-terminal peptide sequence.

- For information on transformation into *E. coli*, see page 6.
- For information on transfection into Sf9 or High Five<sup>™</sup> cells see pages 7–9.

#### General Molecular Biology Techniques

For help with DNA ligations, *E. coli* transformations, restriction enzyme analysis, DNA sequencing, and DNA biochemistry, refer to *Molecular Cloning: A Laboratory Manual* (Sambrook *et al.*, 1989) or *Current Protocols in Molecular Biology* (Ausubel *et al.*, 1994).

# Propagation and Maintenance of pIB/V5-His

 $20~\mu g$  of pIB/V5-His vector is supplied in suspension at a concentration of  $40~\mu l$  of  $0.5~\mu g/\mu l$  pIB/V5-His in 10~mM Tris–HCl, 1~mM EDTA, pH 8.0. If you wish to propagate and maintain the pIB/V5-His, we recommend using 10~ng of the vector to transform a recA, endA E. coli strain like TOP10, DH5 , or equivalent using your method of choice.

Select transformants on LB plates containing 50–100 μg/ml ampicillin (see page 6).

## Translation Initiation

Your insert should contain a Kozak consensus sequence with an ATG initiation codon for proper initiation of translation (Kozak, 1987; Kozak, 1990; Kozak, 1991). An example of a Kozak consensus sequence is provided below. Other sequences are possible, but the G or A at position –3 and the G at position +4 (shown in bold) illustrates the most commonly occurring sequence with strong consensus. Replacing one of the two bases at these positions provides moderate consensus, while having neither results in weak consensus. The ATG initiation codon is shown underlined.

#### (G/A)NNATGG

## Fusion to the C-terminal Peptide

If you wish to include the C-terminal peptide for detection with either the V5 or His(C-term) antibodies or purification using the 6xHis tag, you must clone your gene in frame with the peptide. Be sure that your gene does not contain a stop codon upstream of the C-terminal peptide.

If you do not wish to include the C-terminal peptide, include the native stop codon for your gene of interest.

#### Secretion of Recombinant Protein


If your protein of interest is normally secreted, try expressing the protein using the native secretion signal. To date, all mammalian secretion signals tested have functioned properly in insect cells. We have successfully expressed human interleukin-6 (IL6) using the native secretion signal to levels of  $1-2 \,\mu\text{g/ml}$ .

In addition, we recommend that you create a construct to express your protein intracellularly in the event that your protein is not secreted.

### Cloning into pIB/V5-His, Continued

### MCS of pIB/V5-His

The TATA box, start of transcription, and the polyadenylation signal are marked as described in Theilmann and Stewart, 1992. Restriction sites are labeled to indicate the actual cleavage site. Potential stop codons are shown underlined. The multiple cloning site has been confirmed by sequencing and functional testing. The complete sequence of pIB/V5-His is available for downloading from our web site at <a href="https://www.invitrogen.com">www.invitrogen.com</a> or from Technical Support (page 27). For a map and a description of the features of pIB/V5-His, refer to pages 22–23.



## Transforming E. coli

#### Introduction

Once you have completed your ligation reactions, you are ready to transform into *E. coli*. Many strains and transformation protocols are suitable. General recommendations are provided below.

#### E. coli Host

Many *E. coli* strains are suitable for transformation of pIB/V5-His including TOP10 or DH5. We recommend that you propagate vectors containing inserts in *E. coli* strains that are recombination deficient (*rec*A) and endonuclease A deficient (*end*A). For your convenience, TOP10 is available as electrocompetent or chemically competent cells from Invitrogen (see page vi).

#### Transformation Method

You may use your method of choice to transform *E. coli*. To select transformants, use LB plates containing 50–100 µg/ml ampicillin.

**Note**: It is possible to select *E. coli* transformants on  $100-150\,\mu g/ml$  blasticidin. You will need to prepare Low Salt LB (5 g NaCl per liter) and adjust the pH to below 7 before autoclaving. After autoclaving and cooling, add blasticidin. Plates are stable for 1–2 weeks. Since blasticidin is sensitive to salt and pH, we recommend that you use ampicillin for selection of *E. coli* transformants.



We recommend that you sequence your construct to confirm that your gene is fused in frame with the V5 epitope and the polyhistidine tag. Use the OpIE2 Forward and Reverse sequencing primers included in your kit or a primer to your gene of interest to sequence your insert.

**Note**: Resuspend each primer in 20  $\mu$ l sterile water to prepare a 0.1  $\mu$ g/ $\mu$ l stock solution.

#### Long-Term Storage

Once you have confirmed that you have the correct clone, prepare a glycerol stock for long-term storage. It is also a good idea to keep a stock of plasmid DNA at -20°C.

#### To prepare a glycerol stock:

- 1. Grow the *E. coli* strain containing the plasmid overnight.
- 2. Combine 0.85 ml of the overnight culture with 0.15 ml of sterile glycerol.
- Vortex and transfer to a labeled cryovial.
- 4. Freeze the tube in liquid nitrogen or dry ice/ethanol bath and store at -80°C.

## **Transient Expression in Insect Cells**

#### Introduction

Once you have cloned your gene of interest into pIB/V5-His, you are ready to transfect your construct into Sf9 or High Five<sup> $^{\text{IM}}$ </sup> cells using lipid-mediated transfection and test for expression of your protein.

## Plasmid Preparation

Plasmid DNA for transfection into insect cells must be very clean and free from phenol and sodium chloride. Contaminants will kill the cells, and salt will interfere with lipid complexing, decreasing transfection efficiency. We recommend isolating plasmid DNA using the S.N.A.P.™ Miniprep Kit or the PureLink™ HiPure Plasmid Miniprep Kit (page vi). The purified plasmid can be used directly for transfection of insect cells.

#### Method of Transfection

We recommend lipid-mediated transfection with Cellfectin® Reagent. Note that other lipids may be substituted, although transfection conditions may have to be optimized.

#### **Expected Transfection Efficiency using Cellfectin® Reagent:**

- 40-60% for Sf9 cells
- 40–60% for High Five<sup>™</sup> cells

**Note**: Other transfection methods (*e.g.* calcium phosphate and electroporation (Mann & King, 1989) have also been tested with High Five<sup>TM</sup> cells.

## Control of Plasmid Quality

To test the quality of a plasmid DNA preparation, include a mock transfection using DNA only in all transfection experiments. At about 24–48 hours post–transfection, compare the DNA only mock transfection with cells transfected with plasmid. If the plasmid preparation contains contaminants, then the cells will appear unhealthy and start to lyse.

#### **Materials Needed**

You will need the following for each transfection experiment:

- 1–10 µg of highly purified plasmid DNA ( $\sim$ 1 µg/µl in TE buffer)
- Either log-phase Sf9 cells (1.6–2.5 × 10<sup>6</sup> cells/ml, >95% viability) or log-phase High Five<sup>™</sup> cells (1.8–2.3 × 10<sup>6</sup> cells/ml, >95% viability), growing in serumfree medium

**Note:** You may transfect Sf9 cells in Grace's Medium without supplements or FBS. The proteins in the FBS and supplements will interfere with the liposomes, causing the transfection efficiency to decrease.

- Serum-free medium (see next page)
- 60 mm tissue-culture dishes
- 1.5 ml sterile microcentrifuge tubes
- Rocking platform only (NOT orbital)
- 27°C incubator
- Inverted Microscope
- Paper towels and air-tight bags or containers
- 5 mM EDTA, pH 8.0

#### Serum-Free Media

Several serum-free media are available from Invitrogen for use in transfection experiments with pIB/V5-His. Express Five® SFM (page vi) is recommended for use with High Five™ cells while Sf-900 II SFM (1X) (page vi) is optimized for use with Sf9 cells and Sf21 cells. Other serum-free media may be used, although you may have to optimize conditions for transfection and selection. You will need to adapt the cells to serum-free medium before transfection (see the Insect Cell Lines manual for a protocol).

**Note:** You may transfect Sf9 cells in Grace's Medium without supplements or FBS. The proteins in the FBS and supplements will interfere with the liposomes, causing the transfection efficiency to decrease.

#### **Preparing Cells**

For each transfection, use log-phase cells with greater than 95% viability. We recommend that you set up enough plates to perform a time course for expression of your gene of interest. Test for expression 2, 3, and 4 days post-transfection. You will need at least one 60 mm plate for each time point.

- 1. For Sf9 cells or High Five<sup>™</sup> cells, seed 1 × 10<sup>6</sup> cells in appropriate serum-free medium in a 60 mm dish.
  - **Note:** You may transfect Sf9 cells in Grace's Medium without supplements or FBS. The proteins in the FBS and supplements will interfere with the liposomes, causing the transfection efficiency to decrease.
- 2. Rock gently from side to side for 2–3 minutes to evenly distribute the cells. Cells should be 50–60% confluent.
- 3. Incubate the cells for at least 15 minutes without rocking to allow the cells to fully attach to the bottom of the dish to form a monolayer of cells.
- 4. Verify that the cells have attached by inspecting them under an inverted microscope.

## Positive and Negative Controls

We recommend that you include the following controls:

- pIB/V5-His/CAT vector as a positive control for transfection and expression
- Lipid only as a negative control
- DNA only to check for DNA contamination



- If you use another lipid besides Cellfectin® Reagent, review the protocol on the next page and consult the manufacturer's instructions to adapt the protocol for your use. You may have to empirically determine the optimal conditions for transfection.
- <u>Do not linearize</u> the plasmid prior to transfection. Linearizing the plasmid appears to decrease protein expression.

## Transfection Procedure

Plasmid DNA and Cellfectin® Reagent are mixed together in the appropriate medium (see **Serum-Free Media**, previous page) and incubated with freshly seeded insect cells. The amount of cells, liposomes, and plasmid DNA has been optimized for 60 mm culture plates. It is important that you optimize transfection conditions if you use plates or flasks other than 60 mm plates.

- 1. To prepare each transfection mixture, add the following reagents to a 1.5-ml microcentrifuge tube:
  - a. 1 ml Appropriate serum-free medium (see previous page)
  - b. 1–10 µl pIB/V5-His plasmid or construct (~1 µg/µl in TE, pH 8.0)
  - c. 6 µl Cellfectin® Reagent (mix well before use and always add last)
- 2. Gently mix the transfection mixture for 10 seconds.
- 3. Incubate the transfection mixture at room temperature for 15 minutes. While the transfection mixture is incubating, proceed to Step 4.
- 4. Carefully remove the medium from the cells without disrupting the monolayer. If the medium contained serum, wash the cells by carefully adding 2 ml of fresh serum-free medium (see previous page for appropriate media to use).
  - **Note:** This will remove trace amounts of serum that will decrease the efficiency of liposome transfection).
- 5. Again, carefully remove all of the medium from the monolayer and add the entire transfection mix drop-wise into the 60 mm dish. **Repeat** for all transfections
  - **Note:** Distribute the drops evenly over the monolayer. This method reduces the chances of disturbing the monolayer).
- 6. Incubate the dishes at room temperature for 4 hours on a side-to-side, rocking platform. Adjust speed to ~2 side to side motions per minute. **Note**: If you do not have a rocker, manually rock the dishes periodically.
- 7. After incubation, add 1–2 ml of appropriate serum-free medium (see previous page) to each 60 mm dish.
- 8. Place the dishes in a sealed plastic bag with moist paper towels to prevent evaporation and incubate at 27°C.
  - **Note**: It is not necessary to remove the transfection solution as Cellfectin® Reagent is not toxic to the cells. If you are using a different lipid and observe loss of viability, then remove the transfection solution after 4 hours, rinse twice with medium, and replace with 1–2 ml of fresh serum-free medium.
- 9. Harvest the cells 2\*, 3, and 4 days post–transfection and assay for expression of your gene (see next page). There's no need to add fresh medium if the cells are sealed in an airtight plastic bag with moist paper towels.
  - \*Note: To create stable cell lines, proceed to **Stable Transfection Procedure** (page 12) at 48-hours post-transfection. Otherwise, continue to harvest cells at days 3 and 4 post-transfection.

#### Detecting Recombinant Proteins

To detect expression of your recombinant fusion protein by Western blot analysis, you may use the Anti-V5 antibodies or the Anti-His(C-term) antibodies available from Invitrogen (see page viii for ordering information) or an antibody to your protein of interest. In addition, the Positope™ Control Protein (page vi) is available from Invitrogen for use as a positive control for detection of fusion proteins containing a V5 epitope or a 6xHis tag. WesternBreeze® Chromogenic Kits and WesternBreeze® Chemiluminescent Kits are available from Invitrogen to facilitate detection of antibodies by colorimetric or chemiluminescent methods. For more information, refer to our Web site (www.invitrogen.com) or contact **Technical Support** (see page 27).

## Testing for Expression

Use the cells from one 60 mm plate for each expression experiment.

**Before starting** prepare Cell Lysis Buffer and SDS-PAGE sample buffer. Recipes are provided on page 19 for your convenience, but other recipes are suitable. If you are using pre-cast polyacrylamide gels (see below), refer to the manufacturer's instructions to prepare the appropriate sample buffer.

- 1. Prepare an SDS-PAGE gel that will resolve your expected recombinant protein.
- 2. Remove the medium from the cells. If your protein is predicted to be secreted, be sure to save and assay both the medium and the cell pellet.
- 3. Add 100 µl Cell Lysis Buffer to the plate and slough (or scrape) the cells into a microcentrifuge tube. Vortex the cells to ensure they are completely lysed.
- Centrifuge at maximum speed for 1-2 minutes to pellet nuclei and cell membranes. Transfer the supernatant to a new tube.
   Note: If you are expressing a membrane protein, it may be located in the pellet. Be sure to assay the pellet (see below).
- 5. Assay the lysate for protein concentration. You may use the Bradford, Lowry, or BCA assays.
- 6. To assay your samples, mix them with SDS-PAGE sample buffer as follows:
  - Lysate: 30 μl lysate with 10 μl **4X SDS-PAGE** sample buffer.
  - Pellet: Resuspend pellet in 100 µl 1X SDS-PAGE sample buffer.
  - Medium: 30 μl medium with 10 μl **4X SDS-PAGE** sample buffer.

**Note**: Because of the volume of medium, it is difficult to normalize the amount loaded on an SDS-PAGE gel. If you are concerned about normalization, concentrate the medium.

- 7. Boil the samples for 5 minutes. Centrifuge briefly.
- 8. Load approximately 3–30 µg protein per lane. For the cell pellet sample, load the same volume as the lysate. Amount to load depends on the amount of your protein produced.
- 9. Electrophorese your samples, blot, and probe with a suitable antibody (see above).
- 10. Visualize proteins using your desired method.



The C-terminal tag containing the V5 epitope and 6xHis tag will increase the size of your protein by ~3 kDa. Note that any additional amino acids between your protein and the tags are not included in this molecular weight calculation.

#### Polyacrylamide Gel Electrophoresis

To facilitate separation of your recombinant protein by polyacrylamide gel electrophoresis, a wide range of pre-cast NuPAGE® and Novex® Tris-Glycine polyacrylamide gels are available from Invitrogen. In addition, Invitrogen also carries a large selection of molecular weight protein standards and staining kits. For more information about the appropriate gels, standards, and stains to use to visualize your recombinant protein, refer to our web site at <a href="https://www.invitrogen.com">www.invitrogen.com</a> or contact **Technical Support** (page 27).

#### **Assay for CAT**

If you use pIB/V5-His/CAT as a positive control vector, you may assay for CAT expression using your method of choice. Commercial kits to assay for CAT protein are available. There is also a novel, rapid radioactive assay (Neumann *et al.*, 1987).

CAT can be detected by Western blot using antibodies against the C-terminal fusion tag (page viii or an antibody against CAT. The CAT/V5-His protein fusion migrates around 34 kDa on an SDS-PAGE gel.

## **Selecting Stable Cell Lines**

#### Introduction

Once you have demonstrated that your protein is expressed in Sf9 or High Five $^{\text{TM}}$  cells, you may wish to create stable expression cell lines for long-term storage and large-scale production of the desired protein.

## Nature of Stable Cell Lines

Note that stable cell lines are created by multiple copy integration of the vector. Amplification as in the case with calcium phosphate transfection and hygromycin resistance in *Drosophila* is generally not observed.

#### **Before Starting**

Review the information on blasticidin S on page 21. Prepare a stock solution of blasticidin S as described in the protocol provided.

#### Effect of Blasticidin on Sensitive and Resistant Cells

Cytopathic effects should be visible within 3–5 days depending on the concentration of blasticidin in the medium. Sensitive cells will enlarge and become filled with vesicles. The outer membrane will show signs of blebbing, and cells will eventually detach from the plate.

Blasticidin-resistant cells should continue to divide at regular intervals to form distinct colonies. There should not be any distinct morphological changes between blasticidin-resistant cells compared to cells not under selection with blasticidin.

#### Suggested Blasticidin Concentrations

In general, concentrations around 10  $\mu g/ml$  will kill Sf9 (in complete Sf-900 II SFM medium) and concentrations around 20  $\mu g/ml$  will kill High Five<sup>TM</sup> cells (in Express Five® SFM) within one week, although a few cells may remain that exclude trypan blue. To obtain faster and more thorough killing, we recommend using 50–80  $\mu g/ml$  blasticidin. Once blasticidin-resistant clones have been obtained, cells may be maintained in lower concentrations of blasticidin (*e.g.* 10–20  $\mu g/ml$ ). If you are using other media or have trouble selecting cells using the concentrations above, we recommend that you perform a kill curve (see below).

#### Blasticidin Selection Guidelines

If you wish to test your cell line for sensitivity to blasticidin, perform a kill curve as described below. Assays can be done in 24-well tissue culture plates.

- 1. Prepare Sf-900 II SFM medium (1X) or the serum-free medium of choice, supplemented with concentrations ranging from 0–100  $\mu$ g/ml blasticidin. Generally, concentrations that effectively kill lepidopteran insect cells within a week are in the 50–80  $\mu$ g/ml range.
  - Note: While  $10-20~\mu g/ml$  blasticidin will kill cells within a week, higher concentrations will result in faster and more thorough killing. In addition, using higher concentrations of blasticidin may result in enrichment of clones containing multiple integrations of your gene of interest.
- 2. Test varying concentrations of blasticidin on the cell line to determine the concentration that kills your cells within a week (kill curve).
- 3. Use the concentration of drug that kills your cells within a week.

### Selecting Stable Cell Lines, Continued



Stable Transfection Procedure <u>**Do not linearize**</u> the plasmid prior to transfection. Linearizing the plasmid appears to decrease protein expression.

For stable transfections, follow the steps below. Include a mock transfection and a positive control (pIB/V5-His/CAT).

- 1. Complete the **Transfection Procedure** (page 9), through Step 9 until 2-days post-transfection.
- 2. At 2-days, post-transfection, remove the transfection solution and add fresh medium (without blasticidin).
- 3. Split cells 1:5 (20% confluent) and let cells attach overnight before adding selective medium.
- 4. Remove medium and replace with medium containing blasticidin at the appropriate concentration. Incubate cells at 27°C.
- 5. Replace selective medium every 3 to 4 days until you observe foci (colonies) forming. At this point you may use cloning cylinders or dilution to isolate clonal cell lines (next page) or you can let resistant cells grow out to confluence for a polyclonal cell line (2 to 3 weeks).
- 6. To isolate a polyclonal cell line, let the resistant cells grow to confluence and split the cells 1:5 and test for expression.
  - **Important**: Always use medium **without** blasticidin when splitting cells. Let the cells attach before adding selective medium.
- 7. Expand resistant cells into flasks to prepare frozen stocks. **Important:** Always use medium containing blasticidin when maintaining stable lepidopteran cell lines. You may lower the concentration of blasticidin to  $10~\mu g/ml$  for maintenance.

### Selecting Stable Cell Lines, Continued

#### Isolation of Clonal Cell Lines Using Cloning Cylinders

If you elect to select clonal cell lines, try to isolate as many foci (colonies) as possible for expression testing. As in mammalian cell culture, the location of integration may affect expression of your gene.

**Tip**: Perform selections in small plates or wells. When you remove the medium, you must work quickly to prevent the cells from drying out. Using smaller plates or wells limits the number of colonies you can choose at a time. To select more colonies, increase the number of plates or wells, not the size.

#### To select colonies:

- 1. Examine the closed plate under a microscope and mark the location of each colony on the top of the plate. Transfer the markings to the bottom of the plate. Be sure to include orientation marks.
  - **Note**: Each colony will contain 50 to 200 cells. Sf9 cells tend to spread more than High Five™ cells.
- 2. Move the culture dish to the sterile cabinet and remove the lid.
- 3. Apply a thin layer of sterile silicon grease to the bottom of the cloning cylinder (Scienceware, Catalog no. 378747–00 or Belco, Catalog no. 2090–00608), using a sterile cotton-tipped wooden applicator. The layer should be thick enough to retard the flow of liquid from the cylinder, without obscuring the opening on the inside.
  - **Tip**: Cloning cylinders and silicon grease can be sterilized together by placing a small amount of grease in a glass petri dish and placing the cloning cylinders upright in the grease. After autoclaving, the grease will have spread out in a thin layer to coat the bottom of the cylinders.
- 4. Aspirate the culture medium and place the cylinder firmly and directly over the marked area. Use a microscope if it is available to help you direct placement of the cylinder.
- 5. Use 20–100 µl of medium (without blasticidin) to slough the cells. Try to hold the pipette tip away from the sides of the cloning cylinder to avoid the grease (this will take a little practice).
- Remove the cells and medium and transfer to a microtiter plate and let the cells attach. Remove medium and replace with selective medium for culturing. Expand the cell line and test for expression of your gene of interest.

**Important**: Always use medium **without** blasticidin when splitting cells. Let the cells attach before adding selective medium.

### Selecting Stable Cell Lines, Continued

# Isolation of Clonal Cell Lines Using a Dilution Method

You may also select clonal cell lines using a quick dilution method. The objective of this method is to dilute the cells so that under selective pressure only one stable viable cell per well is achieved.

**Note:** The higher your transfection efficiency, the more you should dilute out your cells. The protocol below works well with cells transfected at 5–10% efficiency.

- 1. Forty-eight hours after transfection, dilute the cells to  $1 \times 10^4$  cells/ml in medium **without** blasticidin.
  - **Note**: Other dilutions of the culture should also be used as transfection efficiency will determine how many transformed cells there will be per well.
- 2. Add  $100 \,\mu$ l of the cell solution from Step 1 to 32 wells of a 96-well microtiter plate (8 rows by 4 columns).
- 3. Dilute the remaining cells 1:1 with medium **without** blasticidin and add  $100 \, \mu l$  of this solution to the next group of 32 wells  $(8 \times 4)$ .
- 4. Once again, dilute the remaining cells 1:1 with medium **without** blasticidin and add  $100 \mu l$  of this solution to the last group of 32 wells.
  - **Note**: Although the cells can be diluted to low numbers, cell density is critical for viability. If the density drops below a certain level, the cells will not grow.
- 5. Let the cells attach overnight, then remove the medium and replace with medium **containing** blasticidin.
  - **Note**: Removing and replacing medium may be tedious. If you slough the cells gently, it is possible to dilute the cells directly into selective medium.
- 6. Wrap the plate and incubate at 27°C for 1 week. It is not necessary to change the medium or place in a humid environment.
- 7. Check the plate after a week and mark the wells that have only one colony.
- 8. Continue to incubate the plate until the colony fills most of the well.
- 9. Harvest the cells and transfer to a 24-well plate with 0.5 ml of fresh medium containing blasticidin.
- 10. Continue to expand the clone to 12- and 6-well plates, and finally to a T-25 flask.

#### Assay for Expression

Assay each of your cell lines for yield of the desired protein and select the one with the highest yield for scale-up and purification of recombinant protein. **If your protein is secreted, remember to assay the cell pellet as well as the medium.** You may wish to compare the yield of protein in the cells and medium.

#### Yield of Expressed Protein

In general, the level of secreted protein is comparable to that obtained with viral expression systems in insect cells. We have obtained stable cell lines that express and secrete human interleukin-6 to levels of 1–2  $\mu$ g/ml. Human melanotransferrin has been expressed to levels of 8–10  $\mu$ g/ml (Hegedus *et al.*, 1999).



Remember to prepare master stocks and working stocks of your stable cell lines prior to scale-up and purification. Refer to the Insect Cell Lines manual for information on freezing your cells and scaling up for purification.

## **Scale-Up and Purification**

#### Introduction

Once you have obtained stable cell lines expressing the protein of interest and prepared frozen stocks of your cell lines, you are ready to purify your protein. General information for protein purification is provided below. Eventually, you may expand your stable cell line into larger flasks, spinners, shake flasks, or bioreactors to obtain the desired yield of protein. If your protein is secreted, you may culture cells in serum-free medium to simplify purification.



As you expand your stable cell line, you can maintain the concentration of blasticidin at  $10 \mu g/ml$ .

## Serum-Free Medium

If your protein is secreted, use serum-free medium to facilitate expression and purification (see page 8).

#### Adapting Cells to Different Medium

Cells can be switched from complete medium to serum-free medium during passage. Refer to the Insect Cell Lines manual for more information on how to adapt cells to different medium.



If you plan to use a metal-chelating resin such as ProBond™ to purify your secreted protein from serum-free medium, note that adding serum-free medium directly to the column will strip the nickel ions from the resin. See the information below in Purification of 6xHis-tagged Proteins from Medium for a general recommendation to address this issue.

## Purifying Proteins from Medium

Many protocols are suitable for purifying proteins from the medium. The choice of protocol depends on the nature of the protein being purified. Note that the culture volume needed to purify sufficient quantities of protein is dependent on the expression level of your protein and the method of detection. To purify 6xHis-tagged proteins from the medium, see below.

#### Purification of 6xHis-tagged Proteins from Medium

To purify 6xHis-tagged recombinant proteins from the culture medium, we recommend that you perform dialysis or ion exchange chromatography prior to affinity chromatography on metal-chelating resins.

#### Dialysis allows:

- Removal of media components that strip Ni<sup>+2</sup> from metal-chelating resins **Ion exchange chromatography allows:**
- Removal of media components that strip Ni<sup>+2</sup> from metal-chelating resins
- Concentration of your sample for easier manipulation in subsequent purification steps

Conditions for successful ion exchange chromatography will vary depending on the protein. For more information, refer to *Current Protocols in Protein Science* (Coligan *et al.*, 1998), *Current Protocols in Molecular Biology*, Unit 10 (Ausubel *et al.*, 1994) or the *Guide to Protein Purification* (Deutscher, 1990).

### Scale-Up and Purification, Continued

#### Metal-chelating Resin

You may use the ProBond<sup>™</sup> Purification System or Ni-NTA Purification System (page viii), or a similar product to purify your 6xHis-tagged protein. Both purification systems contain a metal-chelating resin specifically designed to purify 6xHis-tagged proteins. Before starting, be sure to consult the ProBond<sup>™</sup> Purification System manual or Ni-NTA Purification System manual to familiarize yourself with the buffers and the binding and elution conditions. If you are using another resin, consult the manufacturer's instructions.



Many insect cell proteins are naturally rich in histidines, with some containing stretches of six histidines. When using the ProBond™ Purification System or other similar products to purify 6xHis-tagged proteins, these histidine-rich proteins may co-purify with your protein of interest. The contamination can be significant if your protein is expressed at low levels. We recommend that you add 5 mM imidazole to the binding buffer prior to addition of the protein mixture to the column. Addition of imidazole may help to reduce background contamination by preventing proteins with low specificity from binding to the metal-chelating resin.

#### Purification of Intracellularly Expressed Proteins

If you are expressing your 6xHis-tagged protein intracellularly, you may lyse the cells and add the lysate directly to the ProBond<sup>TM</sup> column. You will need  $5 \times 10^6$  to  $1 \times 10^7$  cells for purification of your protein on a 2 ml ProBond<sup>TM</sup> column (see ProBond<sup>TM</sup> Purification System manual).

- 1. Seed  $2 \times 10^6$  cells in two or three 25 cm<sup>2</sup> flasks.
- 2. Grow the cells in selective medium until they reach confluence ( $4 \times 10^6$  cells).
- 3. Wash cells once with PBS.
- 4. Harvest the cells by sloughing.
- 5. Transfer the cells to a sterile centrifuge tube.
- 6. Centrifuge the cells at  $1000 \times g$  for 5 minutes. You may lyse the cells immediately or freeze in liquid nitrogen and store at  $-80^{\circ}$ C until needed.

#### Scale-Up

To scale up insect cell culture, refer to the Insect Cell Lines manual.

## **Troubleshooting**

Problem	Cause	Solution	
Cells Growing Too Slowly (Or not at all)	Refer to the Insect Cell Lines manual	For troubleshooting guidelines regarding cell culture, refer to the Insect Cell Lines manual included with the InsectSelect™ BSD kits.	
Low Transfection Efficiency	Impure DNA	Transfected cells will appear unhealthy when compared to the negative control (DNA only). Use clean, pure DNA isolated by resin based DNA isolation kits ( <i>e.g.</i> S.N.A.P. Midiprep Kit).	
	Poor Cell Viability	Be sure to test cells for viability and make sure you use log-phase cells. Refer to the Insect Cell Lines manual to troubleshoot cell culture.	
	Method of Transfection	Optimize transfection.	
Low or No Protein Expression	Gene not cloned in frame with the C-terminal sequence	If it is not in frame with the C-terminal peptide sequence, expression will not be detected using the antibody to the V5 epitope or the C-terminal histidine tag.	
	No Kozak sequence for proper initiation of translation	Translation will be inefficient and the protein will not be expressed at its optimal level.	
	Optimize expression	If you've tried a time course to optimize expression, try switching cell lines. Proteins may express better in a different cell line.	
	Proteins are degraded	Include protease inhibitors in the Cell Lysis buffer to prevent degradation of recombinant protein.	
	Poor secretion	Check the cell pellet as well as the medium when analyzing secreted expression. Protein may be trapped in the cell and not secreted. To improve secretion, try a different cell line $(e.g. \text{ High Five}^{\text{\tiny{TM}}})$ .	

## **Appendix**

## **Recipes**

#### LB (Luria-Bertani) Medium and Plates

#### **Composition:**

1.0% Tryptone 0.5% Yeast Extract 1.0% NaCl pH 7.0

- 1. For 1 liter, dissolve 10 g tryptone, 5 g yeast extract, and 10 g NaCl in 950 ml deionized water.
- 2. Adjust the pH of the solution to 7.0 with NaOH and bring the volume up to 1 liter.
- 3. Autoclave on liquid cycle for 20 minutes. Allow solution to cool to  $\sim$ 55°C and add antibiotic if needed.
- 4. Store at room temperature or at 4°C.

#### LB agar plates

- 1. Prepare LB medium as above, but add 15 g/L agar before autoclaving.
- 2. Autoclave on liquid cycle for 20 minutes.
- 3. After autoclaving, cool to ~55°C, add antibiotic and pour into 10 cm plates.
- 4. Let harden, then invert and store at 4°C, in the dark.

## Trypan Blue Exclusion Assay

- 1. Prepare a 0.4% stock solution of trypan blue in phosphate buffered saline, pH 7.4.
- 2. Mix 0.1 ml of trypan blue solution with 1 ml of cells and examine under a microscope at low magnification.
- 3. Dead cells will take up trypan blue while live cells will exclude it. Count live cells versus dead cells. Cell viability should be at least 95–99% for healthy log-phase cultures.

## Recipes, Continued

#### **Cell Lysis Buffer**

50 mM Tris, pH 7.8

150 mM NaCl

1% Nonidet P-40

1. This solution can be prepared from the following common stock solutions. For 100 ml, combine

1 M Tris base 5 ml 5 M NaCl 3 ml Nonidet P-40 1 ml

- 2. Bring the volume up to 90 ml with deionized water and adjust the pH to 7.8 with HCl.
- 3. Bring the volume up to 100 ml. Store at room temperature.

Add 1 mM PMSF, 1  $\,$  M leupeptin, and 0.1  $\mu M$  aprotinin to prevent proteolysis, before use.

## **4X SDS-PAGE** Sample Buffer

Combine the following reagents to obtain 10 ml:

0.5 M Tris-HCl, pH 6.8	5 ml
Glycerol (100%)	4 ml
β-mercaptoethanol	0.8 ml
Bromophenol Blue	0.04 g
SDS	0.8 g

Aliquot and freeze at –20°C until needed.

### **Blasticidin S**

#### Molecular Weight, Formula, and Structure

Merck Index: 12: 1350

MW: 458.9

Formula:  $C_{17}H_{26}N_8O_5$ -HCl

NH2

#### Handling Blasticidin

Always wear gloves, mask, goggles, and protective clothing (*e.g.* a laboratory coat) when handling blasticidin. Weigh out blasticidin and prepare solutions in a hood.

To inactivate blasticidin for disposal, add sodium bicarbonate.


# Preparing and Storing Stock Solutions

- Blasticidin S is soluble in water and acetic acid. Water is generally used to prepare stock solutions of 5 to 10 mg/ml.
- Dissolve blasticidin S in sterile water and filter-sterilize the solution.
- Blasticidin S is unstable in solutions with a pH greater than 8.0. Be sure the pH of the solution is 7.0.
- Aliquot in small volumes (see below) and freeze at –20°C for long-term storage or store at 4°C for short term storage.
- Aqueous stock solutions are stable for 1–2 weeks at  $4^{\circ}$ C and 6–8 weeks at  $-20^{\circ}$ C.
- Do not subject stock solutions to freeze/thaw cycles (do not store in a frost-free freezer).
- Upon thawing, use what you need and store at 4°C. Discard after 1–2 weeks.

### Map of pIB/V5-His

#### Map of pIB/V5-His

The figure below summarizes the features of the pIB/V5-His vector (3521 bp). For a more detailed explanation of each feature, see the next page. The complete sequence of pIB/V5-His is available from our web site at <a href="https://www.invitrogen.com">www.invitrogen.com</a> or from Technical Support (page 27).



OpIE2 promoter: bases 1-548

OpIE2 Forward priming site: bases 511-530

Multiple cloning site: bases 557-652

V5 epitope: bases 659-700 6xHis tag: bases 710-727

OpIE2 Reverse priming site: bases 737-762 OpIE2 polyadenylation sequence: bases 745-874

pUC origin: bases 943-1616 (c)

OpIE1 promoter: bases 1665-1956

EM7 promoter: bases 1971-2029

Blasticidin resistance gene (bsd): bases 2048-2447

Ampicillin resistance gene (bla): 2566-3426

## Features of pIB/V5-His

## Features of pIB/V5-His

The features of pIB/V5-His (3521 bp) are described below. All features have been functionally tested. The multiple cloning site has been tested by restriction analysis.

Features	Function		
OpIE2 promoter	Provides constitutive expression of the gene of interest in lepidopteran insect cells (Theilmann & Stewart, 1992).		
OpIE2 Forward priming site	Allows sequencing of the insert from the 5´ end.		
Multiple cloning site (13 unique sites)	Allows insertion of the gene of interest for expression.		
V5 epitope (Gly-Lys-Pro-Ile-Pro-Asn-Pro-Leu-Leu-Gly-Leu-Asp-Ser-Thr)	Allows detection of your recombinant protein with the Anti-V5 Antibodies (Southern <i>et al.,</i> 1991).		
6xHis tag	Allows purification of your recombinant protein on metal-chelating resin such as $ProBond^{M}$ .		
	In addition, the C-terminal 6xHis tag is the epitope for Anti-His(C-term) Antibodies (Lindner <i>et al.</i> , 1997).		
OpIE2 Reverse priming site	Allows sequencing of the insert from the 3´ end.		
OpIE2 polyadenylation sequence	Efficient transcription termination and polyadenylation of mRNA (Theilmann & Stewart, 1992).		
pUC origin	Replication, maintenance, and high copy number in <i>E. coli</i> .		
OpIE1 promoter	Provides constitutive expression of the blasticidin resistance gene in lepidopteran insect cells (Theilmann & Stewart, 1991).		
EM7 promoter	Allows efficient expression of the blasticidin and ampicillin resistance genes in <i>E. coli</i> .		
Blasticidin resistance gene (bsd)	Allows generation of stable insect cell lines (Kimura <i>et al.</i> , 1994).		
Ampicillin resistance gene (bla)	Selection of transformants in <i>E. coli</i> .		
	<b>Note</b> : The native promoter has been removed. Transcription is assumed to start from the EM7 promoter.		


## Map of pIB/V5-His/CAT

#### **Description**

pIB/V5-His/CAT is a 4257 bp control vector expressing chloramphenicol acetyltransferase (CAT). The CAT gene was amplified using PCR and TOPO® Cloned into pIB/V5-His-TOPO®. CAT is expressed as a fusion to the V5 epitope and 6xHis tag. The molecular weight of the protein is 34 kDa.

#### Map of pIB/V5-His/CAT

The figure below summarizes the features of the pIB/V5-His/CAT vector. The complete nucleotide sequence for pIB/V5-His/CAT is available for downloading from our web site at <a href="https://www.invitrogen.com">www.invitrogen.com</a> or by contacting Technical Support (see page 27).



*OpIE1* promoter: bases 2401-2692 EM7 promoter: bases 2707-2765

pUC origin: bases 1679-2352 (c)

Blasticidin resistance gene (*bsd*): bases 2784-3183 Ampicillin resistance gene (*bla*): bases 3302-4162

OpIE2 polyadenylation sequence: bases 1481-1610

### **OpIE2** Promoter

#### Description

The *OpIE2* promoter has been analyzed by deletion analysis using a CAT reporter in both *Lymantria dispar* (LD652Y) and *Spodoptera frugiperda* (Sf9) cells. Expression in Sf9 cells was much higher than in LD652Y cells. Deletion analysis revealed that sequence up to -275 base pairs from the start of transcription is necessary for maximal expression (Theilmann & Stewart, 1992). Additional sequence beyond -275 may broaden the host range expression of this plasmid to other insect cell lines (Tom Pfeifer, personal communication).

In addition, an 18 bp element appears to be required for expression. This 18 bp element is repeated almost completely in three different locations and partially at six other locations. These are marked in the figure below. Elimination of the three major 18 bp elements reduces expression to basal levels (Theilmann & Stewart, 1992). The function of these elements is not known.

Primer extension experiments revealed that transcription initiates equally from either the C or the A indicated. These two transcriptional start sites are adjacent to a CAGT sequence motif that has been shown to be conserved in a number of early genes (Blissard & Rohrmann, 1989).

	6611) 8611	(2110001101 01		2,0,,.		
1	GGATCATGAT	GATAAACAAT	GTATGGTGCT	AATGTTGCTT	CAACAACAAT	TCTGTTGAAC
61	TGTGTTTTCA	TGTTTGCCAA	CAAGCACCTT	TATACTCGGT	GGCCTCCCCA	CCACCAACTT
121	TTTTGCACTG	CAAAAAAACA	CGCTTTTGCA	CGCGGGCCCA	TACATAGTAC	AAACTCTACG
181	TTTCGTAGAC	TATTTTACAT	AAATAGTCTA	CACCGTTGTA	TACGCTCCAA	ATACACTACC
241	ACACATTGAA	CCTTTTTGCA	GTGCAAAAAA	GTACGTGTCG	GCAGTCACGT	AGGCCGGCCT
301	TATCGGGTCG		CGTACGAATC		ACCGGACGAG	TGTTGTCTTA
361	TCGTGACAGG	ACGCCAGCTT	CCTGTGTTGC	TAACCGCAGC	CGGACGCAAC	TCCTTATCGG
421	AACAGGACGC	GCCTCCATAT	CAGCCGCGCG	TTATCTCATG	CGCGTGACCG	GACACGAGGC
481	GCCCGTCCCG	CTTATCGCGC	TATA CTATAAATAC	AGCCCGCAAC	GATCTGGTAA	Start of Transcription  ACACAGTTGA
541	ACAGCATCTG	TTCGAATTTA				

### **OpIE1** Promoter

#### **Description**

The *OpIE1* promoter has been analyzed by deletion analysis using a CAT reporter in both *Lymantria dispar* (LD652Y) and *Spodoptera frugiperda* (Sf9) cells. Deletion analysis revealed that sequence between -186 and -106 is important for maximum transcription in Sf9 cells (Theilmann & Stewart, 1991).

This region contains a canonical CCAAT site (underlined) (Johnson & McKnight, 1989) and an element (R4) that is homologous to the proposed binding site of the *Drosophila* transcription factor Adf-1 (England *et al.*, 1990). Three other Adf-1-like elements are found at three other distal locations. These elements are referred to as R1, R2, R3, and R4. R3 and R4 are marked in the figure below. R1 and R2 are not present in pIB/V5-His but do not appear to be important for expression in Sf9 cells. The function of these elements has not been determined.

Primer extension experiments revealed that transcription initiates from the A in the CAGT sequence. This CAGT sequence motif has been shown to be conserved in a number of early genes (Blissard & Rohrmann, 1989).

			R3			
1661	TTGGTCATGC	GAAACACGCA	CGGCGCGCGC	ACGCAGCTTA	GCACAAACGC	GTCGTTGCAC
1721	GCGCCCACCG	CTAACCGCAG	GCCAATCGGT	CGGCCGGCCT	CATATCCGCT	CACCAGCCGC
	R4					
1781	GTCCTATCGG	GCGCGGCTTC	CGCGCCCATT	TTGAATAAAT	AAACGATAAC	GCCGTTGGTG
						TATA
1841	GCGTGAGGCA	TGTAAAAGGT	TACATCATTA	TCTTGTTCGC	CATCCGGTTG	GTATAAATAG
			Start of tran	scription		
1901	ACGTTCATGT	TGGTTTTTGT	TT <u>CAGT</u> TGCA	AGTTGGCTGC	GGCGCGCGCA	GCACCTTTGC
1961	CGGGATCTGC	CGGGCTGCAG	G CACGTGTTGA	A CAATTAATC	A TCGGCATAG	ST

### **Technical Support**

#### **Web Resources**



Visit the Invitrogen web site at www.invitrogen.com for:

- Technical resources, including manuals, vector maps and sequences, application notes, MSDSs, FAQs, formulations, citations, handbooks, etc.
- Complete technical support contact information
- Access to the Invitrogen Online Catalog
- Additional product information and special offers

#### **Contact Us**

For more information or technical assistance, call, write, fax, or email. Additional international offices are listed on our web site (<a href="www.invitrogen.com">www.invitrogen.com</a>).

Corporate Headquarters:	Japanese Headquarters:	European Headquarters:	
Invitrogen Corporation	Invitrogen Japan	Invitrogen Ltd	
5791 Van Allen Way	LOOP-X Bldg. 6F	Inchinnan Business Park	
Carlsbad, CA 92008 USA	3-9-15, Kaigan	3 Fountain Drive	
Tel: 1 760 603 7200	Minato-ku, Tokyo 108-0022	Paisley PA4 9RF, UK	
Tel (Toll Free): 1 800 955 6288	Tel: 81 3 5730 6509	Tel: +44 (0) 141 814 6100	
Fax: 1 760 602 6500	Fax: 81 3 5730 6519	Tech Fax: +44 (0) 141 814 6117	
E-mail: tech support@invitrogen.com	E-mail: jpinfo@invitrogen.com	E-mail: <u>eurotech@invitrogen.com</u>	

#### **MSDS**

MSDSs (Material Safety Data Sheets) are available on our web site at <a href="https://www.invitrogen.com/msds">www.invitrogen.com/msds</a>.

#### Certificate of Analysis

Product qualification is described in the Certificate of Analysis (CofA), available on our website by product lot number at <a href="https://www.invitrogen.com/cofa">www.invitrogen.com/cofa</a>.

#### **Limited Warranty**

Invitrogen is committed to providing our customers with high-quality goods and services. Our goal is to ensure that every customer is 100% satisfied with our products and our service. If you should have any questions or concerns about an Invitrogen product or service, contact our Technical Support Representatives.

Invitrogen warrants that all of its products will perform according to specifications stated on the certificate of analysis. The company will replace, free of charge, any product that does not meet those specifications. This warranty limits Invitrogen Corporation's liability only to the cost of the product. No warranty is granted for products beyond their listed expiration date. No warranty is applicable unless all product components are stored in accordance with instructions. Invitrogen reserves the right to select the method(s) used to analyze a product unless Invitrogen agrees to a specified method in writing prior to acceptance of the order.

Invitrogen makes every effort to ensure the accuracy of its publications, but realizes that the occasional typographical or other error is inevitable. Therefore, Invitrogen makes no warranty of any kind regarding the contents of any publications or documentation. If you discover an error in any of our publications, please report it to our Technical Support Representatives.

Invitrogen assumes no responsibility or liability for any special, incidental, indirect or consequential loss or damage whatsoever. The above limited warranty is sole and exclusive. No other warranty is made, whether expressed or implied, including any warranty of merchantability or fitness for a particular purpose.

#### **Purchaser Notification**

Limited Use Label License No. 68: InsectSelect<sup>™</sup> Technology The InsectSelect<sup>™</sup> System (the "Expression Kit") was developed into an expression system by scientists at the University of British Columbia (UBC) for high-level expression of recombinant proteins. Components of the InsectSelect<sup>™</sup> System are covered by one or more U.S. patents NUMBERS or patent applications and corresponding foreign patents or patent applications owned and/or licensed by UBC and others.

Invitrogen Corporation ("Invitrogen") has an exclusive license to sell the Expression Kit to scientists for academic research or one year commercial evaluation only, under the terms described below. Use of the Expression Kit for any Commercial Purpose (as defined below) other than evaluation requires the user to obtain a commercial license as detailed below. Before using the Expression Kit, please read the terms and conditions set forth below. Your use of the Expression Kit shall constitute acknowledgment and acceptance of these terms and conditions. If you do not wish to use the Expression Kit pursuant to these terms and conditions, please contact Invitrogen's Technical Supports to return the unused and unopened Expression Kit for a full refund. Otherwise, please complete the Product User Registration Card and return it to Invitrogen.

Invitrogen grants you a non-exclusive license to use the enclosed Expression Kit for academic research or for commercial evaluation purposes only. The Expression Kit is being transferred to you in furtherance of, and reliance on, such license. You may not use the Expression Kit, or the materials contained therein, for any Commercial Purpose without a license for such purpose from Research Corporation Technologies (RCT).

If you are a commercial entity, your right to use the Expression Kit expires after one year. Any commercial entity that wishes to use the Expression Kit beyond this one-year period, must obtain a commercial license from RCT. Commercial entities will be contacted by RCT during this one-year period regarding their desire to obtain a commercial license.

You may terminate your use of the Expression Kit at any time by destroying all InsectSelect™ expression products in your control. Your right to use the Expression Kit will also terminate automatically if you fail to comply with the terms and conditions set forth herein. You shall, upon such termination of your rights, destroy all Expression Kits in your control, and notify Invitrogen of such in writing.

Commercial Purpose include: any use of Expression Products in a Commercial Product; any use of Expression Products in the manufacture of a Commercial Product; any sale of Expression Products; any use (other than evaluation) of Expression products or the Expression Kit to facilitate or advance research or development of a Commercial Product; and any use (other than evaluation) of Expression Products or the Expression Kit to facilitate or advance any research or development program the results of which will be applied to the development of Commercial Products. "Expression Products" means products expressed with the Expression Kit, or with the use of any vectors or host strains in the Expression Kit. "Commercial Product" means any product intended for commercial use.

Access to the Expression Kit must be limited solely to those officers, employees and students of your entity who need access to perform the aforementioned research or evaluation. Each such officer, employee and student must be informed of these terms and conditions and agree, in writing, to be bound by same. You may not distribute the Expression Kit or the vectors or host strains contained in it to others. You may not transfer modified, altered, or original material from the Expression Kit to a third party without written notification to, and written approval from Invitrogen. You may not assign, sublicense, rent, lease or otherwise transfer any of the rights or obligations set forth herein, except as expressly permitted by Invitrogen and RCT.

Inquiries for commercial use should be directed to: Research Corporation Technologies, 101 North Wilmot Road, Suite 600, Tucson, AZ 85711-3335, Tel: 1-520-748-4400 Fax: 1-520-748-0025.

### Purchaser Notification, Continued

Limited Use Label License No. 22: Vectors and Clones Encoding Histidine Hexamer This product is licensed under U.S. Patent Nos. 5,284,933 and 5,310,663 and foreign equivalents from Hoffmann-LaRoche, Inc., Nutley, NJ and/or Hoffmann-LaRoche Ltd., Basel, Switzerland and is provided only for use in research. Information about licenses for commercial use is available from QIAGEN GmbH, Max-Volmer-Str. 4, D-40724 Hilden, Germany.

Limited Use Label License No. 36: Cellfectin<sup>®</sup> Transfection Reagent This product is the subject of one or more of U.S. Patent Nos. 5,674,908, 5,834,439, 6,110,916, and foreign equivalents owned by Invitrogen. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. The buyer may transfer information or materials made through the use of this product to a scientific collaborator, provided that such transfer is not for any Commercial Purpose, and that such collaborator agrees in writing (a) to not transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for research and not for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. Invitrogen Corporation will not assert a claim against the buyer of infringement of the above patents based upon the manufacture, use or sale of a therapeutic, clinical diagnostic, vaccine or prophylactic product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. If the purchaser is not willing to accept the limitations of this limited use statement, Invitrogen is willing to accept return of the product with a full refund. For information on purchasing a license to this product for purposes other than research, contact Licensing Department, Invitrogen Corporation, 5791 Van Allen Way, Carlsbad, California 92008. Phone (760) 603-7200. Fax (760) 602-6500.

Limited Use Label License No. 51: Blasticidin and the Blasticidin Selection Marker Blasticidin and the blasticidin resistance gene (*bsd*) are the subject of U.S. Patent No. 5,527,701sold under patent license for research purposes only. Inquiries for commercial use should be directed to: Kaken Pharmaceutical Company, Ltd., Bunkyo Green Court, Center Office Building, 19-20 Fl, 28-8 Honkomagome 2-chome, Bunkyo-ku, Tokyo 113-8650, Japan, Tel: 81 3-5977-5008; Fax: 81 3-5977-5008.

Limited Use Label License No. 66 High Five<sup>™</sup> Cells The High Five<sup>™</sup> cell line is patented by the Boyce Thompson Institute for Plant Research, Ithaca, New York, and is covered under U.S. patent no. 5,300,435. The High Five<sup>™</sup> cell line is sold for research purposes only. Commercial use requires a license from Boyce Thompson. For more information, please contact: Joyce L. Frank, Tel: 607-254-1220; Fax: 607-254-1242; E-mail: jf51@cornell.edu.

#### References

- Ausubel, F. M., Brent, R., Kingston, R. E., Moore, D. D., Seidman, J. G., Smith, J. A., and Struhl, K. (1994) *Current Protocols in Molecular Biology*, Greene Publishing Associates and Wiley-Interscience, New York
- Blissard, G. W., and Rohrmann, G. F. (1989) Location, Sequence, Transcriptional Mapping, and Temporal Expression of the gp64 Envelope Glycoprotein Gene of the *Orgyia pseudotsugata* Multicapsid Nuclear Polyhedrosis Virus. Virology 170, 537-555
- Coligan, J. E., Dunn, B. M., Ploegh, H. L., Speicher, D. W., and Wingfield, P. T. (1998) *Current Protocols in Protein Science*. Current Protocols (Chanda, V. B., Ed.), John Wiley and Sons, Inc., New York
- Deutscher, M. P. (ed) (1990) *Guide to Protein Purification* Vol. 182. Methods in Enzymology. Edited by Abelson, J. N., and Simon, M. I., Academic Press, San Diego, CA.
- England, B. P., Heberlien, U., and Tjian, R. (1990) Purified *Drosophila* Transcription Factor, ADH Distal Factor-1 (Adf-1), Binds to Sites in Several *Drosophila* Promoters and Activates Transcription. J. Biol. Chem. 265, 5086-5094
- Hegedus, D. D., Pfeifer, T. A., Hendry, J., Theilmann, D. A., and Grigliatti, T. A. (1998) A Series of Broad Host Range Shuttle Vectors for Constitutive and Inducible Expression of Heterologous Proteins in Insect Cell Lines. Gene 207, 241-249
- Hegedus, D. D., Pfeifer, T. A., Theilmann, D. A., Kennard, M. L., Gabathuler, R., Jefferies, W. A., and Grigliatti, T. A. (1999) Differences in the Expression and Localization of Human Melanotransferrin in Lepidopteran and Dipteran Insect Cell Lines. Protein Expression and Purification *15*, 296-307
- Izumi, M., Miyazawa, H., Kamakura, T., Yamaguchi, I., Endo, T., and Hanaoka, F. (1991) Blasticidin S-Resistance Gene (*bsr*): A Novel Selectable Marker for Mammalian Cells. Exp. Cell Res. 197, 229-233
- Jarvis, D. L., Weinkauf, C., and Guarino, L. A. (1996) Immediate-Early Baculovirus Vectors for Foreign Gene Expression in Transformed or Infected Insect Cells. Protein Expression and Purification 8, 191-203
- Johnson, P. F., and McKnight, S. L. (1989) Eukaryotic Transcriptional Regulatory Proteins. Ann. Rev. Biochem. 58, 799-839
- Kimura, M., Takatsuki, A., and Yamaguchi, I. (1994) Blasticidin S Deaminase Gene from *Aspergillus terreus* (*BSD*): A New Drug Resistance Gene for Transfection of Mammalian Cells. Biochim. Biophys. ACTA 1219, 653-659
- Kimura, M., and Yamaguchi, I. (1996) Recent Development in the Use of Blasticidin S, a Microbial Fungicide, as a Useful Reagent in Molecular Biology. Pesticide Biochem. Physiol. *56*
- Kozak, M. (1987) An Analysis of 5´-Noncoding Sequences from 699 Vertebrate Messenger RNAs. Nucleic Acids Res. 15, 8125-8148
- Kozak, M. (1990) Downstream Secondary Structure Facilitates Recognition of Initiator Codons by Eukaryotic Ribosomes. Proc. Natl. Acad. Sci. USA 87, 8301-8305
- Kozak, M. (1991) An Analysis of Vertebrate mRNA Sequences: Intimations of Translational Control. J. Cell Biology 115, 887-903
- Lindner, P., Bauer, K., Krebber, A., Nieba, L., Kremmer, E., Krebber, C., Honegger, A., Klinger, B., Mocikat, R., and Pluckthun, A. (1997) Specific Detection of His-tagged Proteins With Recombinant Anti-His Tag scFv-Phosphatase or scFv-Phage Fusions. BioTechniques 22, 140-149
- Mann, S. G., and King, L. A. (1989) Efficient Transfection of Insect Cells with Baculovirus DNA Using Electroporation. J. Gen. Virol. 70, 3501-3505
- Neumann, J. R., Morency, C. A., and Russian, K. O. (1987) A Novel Rapid Assay for Chloramphenicol Acetyltransferase Gene Expression. BioTechniques *5*, 444-447

## References, Continued

- Pfeifer, T. A., Hegedus, D. D., Grigliatti, T. A., and Theilmann, D. A. (1997) Baculovirus Immediate-Early Promoter-Mediated Expression of the Zeocin™ Resistance Gene for Use as a Dominant Selectable Marker in Dipteran and Lepidopteran Insect Cell Lines. Gene 188, 183-190
- Sambrook, J., Fritsch, E. F., and Maniatis, T. (1989) *Molecular Cloning: A Laboratory Manual*, Second Ed., Cold Spring Harbor Laboratory Press, Plainview, New York
- Southern, J. A., Young, D. F., Heaney, F., Baumgartner, W., and Randall, R. E. (1991) Identification of an Epitope on the P and V Proteins of Simian Virus 5 That Distinguishes Between Two Isolates with Different Biological Characteristics. J. Gen. Virol. 72, 1551-1557
- Takeuchi, S., Hirayama, K., Ueda, K., Sakai, H., and Yonehara, H. (1958) Blasticidin S, A New Antibiotic. The Journal of Antibiotics, Series A 11, 1-5
- Theilmann, D. A., and Stewart, S. (1991) Identification and Characterization of the IE-1 Gene of *Orgyia* pseudotsugata Multicapsid Nuclear Polyhedrosis Virus. Virology 180, 492-508
- Theilmann, D. A., and Stewart, S. (1992) Molecular Analysis of the trans-Activating IE-2 Gene of *Orgyia pseudotsugata* Multicapsid Nuclear Polyhedrosis Virus. Virology *187*, 84-96
- Yamaguchi, H., Yamamoto, C., and Tanaka, N. (1965) Inhibition of Protein Synthesis by Blasticidin S. I. Studies with Cell-free Systems from Bacterial and Mammalian Cells. J. Biochem (Tokyo) 57, 667-677
- Yamaguchi, I., Shibata, H., Seto, H., and Misato, T. (1975) Isolation and Purification of Blasticidin S Deaminase from *Aspergillus terreus*. J. Antibiotics 28, 7-14

©1999–2008 Invitrogen Corporation. All rights reserved.

For research use only. Not intended for any animal or human therapeutic or diagnostic use.

## Notes

## Notes

## Notes



#### **Corporate Headquarters**

5791 Van Allen Way Carlsbad, CA 92008 T: 1 760 603 7200

F: 1 760 602 6500

E: tech\_support@invitrogen.com

For country-specific contact information, visit our web site at **www.invitrogen.com**