# Image-iT<sup>®</sup> Lipid Peroxidation Kit \*for live cell analysis\*

Catalog number C10445

Table 1 Contents and storage

Material	Amount	Concentration	Storage	Stability
Image-iT <sup>®</sup> Lipid Peroxidation Sensor (Component A)	5 × 25 μL	10 mM solution in DMSO	<ul> <li>≤-20°C</li> <li>Protect from light</li> </ul>	When stored as directed, the product is stable for 6 months from the date of receipt.
Cumene hydroperoxide (Component B)	100 μL	5.4 M	Room temperature	

**Excitation/Emission:** In reduced state, the excitation and emission maxima of the Image-iT<sup>®</sup> Lipid Peroxidation Sensor is 581/591 nm; after oxidation, the probe shifts the excitation and emission to 488/510 nm.

## Introduction

Lipid peroxidation generally refers to the oxidative degradation of cellular lipids by reactive oxygen species. Peroxidation of unsaturated lipids affects cell membrane properties<sup>1</sup>, signal transduction pathways<sup>2</sup>, apoptosis, and the deterioration of foods and other biological compounds<sup>3</sup>. Lipid hydroperoxides have been reported to accumulate in oxidatively stressed individuals, including HIV-infected patients<sup>4</sup>. Lipid peroxidation may play a role in aging, as well as pathological processes such as drug-induced phototoxicity and atherosclerosis<sup>4</sup>, and it is often the cause of free radical–mediated damage in cells.

The Image-iT<sup>®</sup> Lipid Peroxidation Kit is based on BODIPY<sup>®</sup> 581/591 C11 reagent and is a sensitive fluorescent reporter for lipid peroxidation. Upon oxidation in live cells, fluorescence shifts from red to green of the phenylbutadiene segment of the fluorophore<sup>5</sup>, providing a ratiometric indication of lipid peroxidation by traditional and high content microscopy as well as flow cytometry. This oxidation-dependent emission shift enables ratiometric fluorescence imaging of lipid peroxidation in live cells<sup>6</sup>. Other common applications of BODIPY<sup>®</sup> 581/591 C11 reagent include fluorometric assays of antioxidant efficacy in plasma<sup>7</sup> and in lipid vesicles<sup>8</sup>. The oxidation and nitro-oxidation products of this BODIPY<sup>®</sup> fatty acid have been characterized by mass spectrometry<sup>9</sup>.

- Optimized kit for ratiometric detection of lipid peroxidation in live cells.
- Upon oxidation, the reagent shifts fluorescence emission peak from ~590 nm (red) to ~510 nm (green).
- Cumene hydroperoxide is included as a positive control to induce lipid peroxidation.
- The reagent is provided as a stable, ready-to-use DMSO solution with a simple protocol compatible with standard workflows in fluorescence microscopy.
- The lipid peroxidation sensor is provided as five single-use vials and contains enough reagent for five 96-well plates or 100 coverslips.

Figure 1 Fluorescence emission spectra of the Image-iT<sup>®</sup> Lipid Peroxidation Sensor, before and after lipid peroxidation.







**Figure 3** Human osteosarcoma cells (U-2 OS) were plated on 35-mm glass bottom dishes (MatTek) and stained with 10  $\mu$ M Lipid Peroxidation Sensor for 30 minutes in complete growth medium at 37°C. For tocopherol treatment, the cells were pre-treated with 150  $\mu$ M  $\alpha$ -tocopherol for 30 minutes. The cells were then treated with vehicle (DMSO), 100  $\mu$ M menadione, 200  $\mu$ M tert-butyl hydroperoxide (TBHP), or 100  $\mu$ M cumene hydroperoxide (CH) for 2 hours at 37°C. Cells were stained with Hoechst 33342 during last 30 minutes of compound incubation. The cells were then washed 3X with phosphate buffered saline (PBS) and then imaged on a Zeiss Axiovert inverted microscope using a 40X objective using filters for Hoechst, FITC, and Texas Red® channels. The signal was then quantifyl ipid peroxidation in cells. In control cells, most of the signal is in red channel and the ratio of 590/510 is high. When the cells are treated with ROS producing reagents (i.e., menadione, TBHP, and CH), the ratios are lower. Tocopherol pre-treatment decreases the lipid peroxidation in cells as seen from increased ratios in tocopherol pre-treated samples.



### **Before Starting**

Materials required but not provided	<ul> <li>Cells and culture medium</li> <li>Phosphate buffered saline (PBS, pH 7.2–7.6)</li> </ul>
Caution	DMSO is known to facilitate the entry of organic molecules into tissues. Handle reagents containing DMSO using equipment and practices appropriate for the hazards posed by such materials. Dispose of the reagents in compliance with all pertaining local regulations. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. Always wear suitable laboratory protective clothing and gloves when handling this reagent. Cumene hydroperoxide (Component B) is toxic. Use appropriate precautions when using this compound.
Prepare stock solution of cumene hydroperoxide	To prepare a 100 mM stock solution of cumene hydroperoxide, add 1 $\mu L$ of cumene hydroperoxide (Component B, 5.4 M) to 54 $\mu L$ of 100% ethanol.

	The following protocol was developed with BPAE, HepG2, U-2OS cells with an optimized Image-iT <sup><math>\circ</math></sup> Lipid Peroxidation Sensor concentration of 10 $\mu$ M, but the assay can be adapted to any cell type. Growth media, cell density, cell type variations may influence labeling. In initial experiments, we recommend testing a concentration range of Image-iT <sup><math>\circ</math></sup> Lipid Peroxidation Sensor to determine the optimal dye concentration for your cell type.
Labeling and detection	<b>1.</b> Plate cells at a desired density and incubate them overnight at 37°C.
	<b>2.</b> Treat the cells with the compound of interest and incubate for the recommended time.
	3. Add Image-iT <sup><math>\circ</math></sup> Lipid Peroxidation Sensor (Component A) at a final concentration of 10 $\mu$ M to the cells. Then incubate for another 30 minutes at 37°C.
	4. Remove media and wash cells three times with PBS.
	<b>5.</b> Read the fluorescence at to separate wavelengths; one at excitation/emission of 581/591 nm (Texas Red <sup>®</sup> filter set) for the reduced dye, and the other at excitation/emission of 488/510 nm (traditional FITC filter set) for the oxidized dye (see <b>Imaging and Analysis</b> , below).
	<b>6.</b> The ratio of the emission fluorescence intensities at 590 nm to 510 nm gives the read-out for lipid peroxidation in cells.
Inducing oxidative stress with cumene hydroperoxide	1. Plate cells at a desired density and incubate them overnight at 37°C.
	2. Add cumene hydroperoxide (Component B) to the cells at a final concentration of 100 $\mu M$ and incubate at 37°C for 2 hours.
	3. Stain with Image-iT <sup>*</sup> Lipid Peroxidation Sensor (Component A) by adding the reagent at a final concentration of 10 $\mu$ M for the last 30 minutes of cumene hydroperoxide incubation. Add Hoechst 33342 dye at this point if performing high content imaging.
	<b>4.</b> Remove media and wash cells three times with PBS, and image cells within 2 hours of staining.
	<b>Note:</b> We recommend that you keep the cells in Live Cell Imaging Solution (Cat. no A14291DJ) to keep them healthier.
Imaging and analysis	Perform fluorescence imaging using traditional Texas Red <sup>®</sup> (590 nm) and FITC (510 nm) emission filters. Determine lipid peroxidation by quantitating the fluorescence intensities and calculating the ratio of intensity in Texas Red <sup>®</sup> channel to the intensity in FITC channel.
	If using high content imaging, use an additional channel (NucBlue <sup>™</sup> /Hoechst 33342) to identify cells. Automatic calculation of ratios can be set up. Refer to the appropriate instrument user guide to set up the ratio calculation.

Biochim Biophys Acta 1487, 61 (2000);
 Circ Res 84, 229 (1999);
 Plant Physiol 125, 1591 (2001);
 J Biol Chem 269, 798 (1994);
 Methods Enzymol 319, 603 (2000);
 FEBS Lett 453, 278 (1999);
 Anal Biochem 354, 290 (2006);
 Biochemistry 45, 8135 (2006);
 Free Radic Biol Med 36, 1635 (2004).

## Product List Current prices may be obtained at www.invitrogen.com or from our Customer Service Department.

<b>Catalog no.</b> C10445	<b>Product Name</b> Image-iT® Lipid Peroxidation Kit *for live cell analysis*	<b>nit Size</b> 1 kit
Related Pro		
A14291DJ		500 mL
C10422	CellROX <sup>®</sup> Deep Red Reagent *for oxidative stress detection*	× 50 uL
C10423	CellEvent® Caspase-3/7 Green Detection Reagent *2 mM solution in DMSO*	100 µL
C10443		× 50 μL
C10444	CellROX <sup>®</sup> Green Reagent *for oxidative stress detection*	× 50 μL
C10446	Click-iT® Lipid Peroxidation Imaging Kit - Alexa Fluor® 488	1 kit
C10447	Click-iT® LAA (Linoleamide alkyne) *for lipid peroxidation detection*	× 20 μL
C10448	CellROX® Reagent Variety Pack *for oxidative stress detection*	1 kit
C6827	CM-H <sub>2</sub> DCFDA (5-(and-6)-chloromethyl- 2',7'-dichlorodihydrofluorescein diacetate, acetyl ester) *mixed isomers*	
	*special packaging*	× 50 µg
D399	H <sub>2</sub> DCFDA (2',7'-dichlorodihydrofluorescein diacetate (2',7'-dichlorofluorescin diacetate)	100 mg
D11347	dihydroethidium (hydroethidine) *special packaging*10	×1 mg
D23844	DAF-FM diacetate (4-amino-5-methylamino- 2',7'-difluorofluorescein diacetate) *special packaging*	< 50 mg
110291	Image-iT® DEAD Green™ viability stain	25 µL
136007	Image-iT <sup>®</sup> LIVE Green Reactive Oxygen Species Detection Kit	1 kit
M36008	MitoSOX® Red mitochondrial superoxide indicator *for live-cell imaging*	× 50 µg
T10096	ThiolTracker <sup>™</sup> Violet (Glutathione Detection Reagent) *for 5 microplates*	each
R37603	BackDrop™ Background Suppressor *for live cells*	1 kit
R37605	NucBlue™ Live Cell Stain *Hoechst 33342 special formulation*	1 kit

## **Purchaser Notification**

#### **Corporate Headquarters**

5791 Van Allen Way Carlsbad, CA 92008 USA Phone: +1 760 603 7200 Fax: +1 760 602 6500 Email: techsupport@lifetech.com

#### **European Headquarters**

Inchinnan Business Park 3 Fountain Drive Paisley PA4 9RF UK Phone: +44 141 814 6100 Toll-Free Phone: 0800 269 210 Toll-Free Tech: 0800 838 380 Fax: +44 141 814 6260 Tech Fax: +44 141 814 6117 Email: euroinfo@invitrogen.com Email Tech: eurotech@invitrogen.com

#### Japanese Headquarters

LOOP-X Bldg. 6F 3-9-15, Kaigan Minato-ku, Tokyo 108-0022 Japan Phone: +81 3 5730 6509 Fax: +81 3 5730 6519 Email: jpinfo@invitrogen.com

Additional international offices are listed at www.lifetechnologies.com

These high-quality reagents and materials must be used by, or directly under the supervision of, a technically qualified individual experienced in handling potentially hazardous chemicals. Read the Safety Data Sheet provided for each product; other regulatory considerations may apply.

#### **Obtaining Support**

For the latest services and support information for all locations, go to www.lifetechnologies.com.

At the website, you can:

- · Access worldwide telephone and fax numbers to contact Technical Support and Sales facilities
- · Search through frequently asked questions (FAQs)
- Submit a question directly to Technical Support (techsupport@lifetech.com)
- Search for user documents, SDSs, vector maps and sequences, application notes, formulations, handbooks, certificates of analysis, citations, and other product support documents
- Obtain information about customer training
- Download software updates and patches

#### SDS

Safety Data Sheets (SDSs) are available at www.lifetechnologies.com/sds.

#### **Certificate of Analysis**

The Certificate of Analysis provides detailed quality control and product qualification information for each product. Certificates of Analysis are available on our website. Go to **www.lifetechnologies.com/support** and search for the Certificate of Analysis by product lot number, which is printed on the product packaging (tube, pouch, or box).

#### For research use only. Not for use in diagnostic procedures.

LIFE TECHNOLOGIES CORPORATION AND/OR ITS AFFILIATE(S) DISCLAIM ALL WARRANTIES WITH RESPECT TO THIS DOCUMENT, EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO THOSE OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, OR NON-INFRINGEMENT. TO THE EXTENT ALLOWED BY LAW, IN NO EVENT SHALL LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) BE LIABLE, WHETHER IN CONTRACT, TORT, WARRANTY, OR UNDER ANY STATUTE OR ON ANY OTHER BASIS FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING BUT NOT LIMITED TO THE USE THEREOF.

LIFE TECHNOLOGIES CORPORATION AND/OR ITS AFFILIATES MAKE NO REPRESENTATION WHATSOEVER THAT SERVICES PROVIDED BY LIFE TECHNOLOGIES AND/OR ITS AFFILIATES SATISFY OR WILL SATISFY ANY REQUIREMENTS OF ANY GOVERNMENTAL BODY OR OTHER ORGANIZATION, INCLUDING BUT NOT LIMITED TO, ANY REQUIREMENT OF THE UNITED STATES FOOD AND DRUG ADMINISTRATION OR THE INTERNATIONAL ORGANIZATION FOR STANDARDIZATION. CUSTOMER AGREES THAT IT IS CUSTOMER'S RESPONSIBILITY TO ENSURE THAT SUCH SERVICES ARE ADEQUATE TO MEET ITS REGULATION/CERTIFICATION REQUIREMENTS AND THAT ALL REQUIREMENTS OF ANY GOVERNMENTAL BODY OR OTHER ORGANIZATION, INCLUDING, BUT NOT LIMITED TO, ANY REQUIREMENT OF THE UNITED STATES FOOD AND DRUG ADMINISTRATION OR THE INTERNATIONAL ORGANIZATION FOR STANDARDIZATION ARE THE RESPONSIBILITY OF THE CUSTOMER.

#### Limited Warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale found on Life Technologies' website at www.lifetechnologies.com/termsandconditions. If you have any questions, please contact Life Technologies at www.lifetechnologies.com/support.

#### Limited Use Label License: Research Use Only

The purchase of this product conveys to the purchaser the limited, non-transferable right to use the purchased amount of the product only to perform internal research for the sole benefit of the purchaser. No right to resell this product or any of its components is conveyed expressly, by implication, or by estoppel. This product is for internal research purposes only and is not for use in commercial services of any kind, including, without limitation, reporting the results of purchaser's activities for a fee or other form of consideration. For information on obtaining additional rights, please contact **outlicensing@lifetech.com** or Out Licensing, Life Technologies Corporation, 5791 Van Allen Way, Carlsbad, California 92008.

The trademarks mentioned herein are the property of Life Technologies Corporation or their respective owners. ArravScan<sup>®</sup> is a registered trademark of Cellnomics. Inc.

©2012 Life Technologies Corporation. All rights reserved.

