One Shot[®] MAX Efficiency[®] DH10B-T1^R Competent Cells

Catalog no. 12331-013

Version C 6 April 2004 25-0424



www.invitrogen.com tech_service@invitrogen.com

Table of Contents

Kit Contents and General Information	2
Storage	2
Contents	2
Genotype	2
Quality Control Procedures	2
Overview	3
Description	3
Applications	3
General Handling	3
Transforming Chemically Competent Cells	4
Materials Supplied by the User	4
Before Starting	4
Procedure	4
Transforming Control Plasmid DNA	5
Calculation	5
Technical Service	6
World Wide Web	6
Contact Us	6
MSDS Requests	6
Limited Warranty	7
References	8

Kit Contents and General Information

Storage Upon receipt, store the cells at -70°C or -80°C. Do not store in liquid nitrogen. Contents Each kit (20 reactions, Catalog no. 12331-013) contains the following: Component Composition Amount Chemically Competent DH10B-T1^R $21 \times 50 \ \mu l$ --SOC Medium 2% Tryptone 6 ml 0.5% Yeast Extract 10 mM NaCl 2.5 mM KCl 10 mM MgCl₂ 10 mM MgSO₄ 20 mM glucose pUC19 10 pg/µl in 5 mM Tris-HCl, 50 µl 0.5 mM EDTA, pH 8 Genotype F-mcrA Δ (mrr-hsdRMS-mcrBC) ϕ 80dlacZ Δ M15 lacX74 recA1 endA1 araD139 Δ (ara*leu*)7697 galU galK rpsL nupG λ - tonA (confers resistance to phage T1) Untransformed cells are tested for the following: **Quality Control** Procedures Leucine auxotrophy (Leu⁻) • Lack of endogenous β -galactosidase activity (Lac⁻) Inability to grow on galactose (Gal⁻) . Inhibited growth on nitrofurantoin (recA) Antibiotic sensitivity (ampicillin, kanamycin, and tetracycline) Antibiotic resistance (streptomycin) (rpsL) Resistance to phage T5, a standard test for resistance to phage T1 (tonA) Competent cells are transformed with either supercoiled pUC19 plasmid, methylated pBR322, or non-methylated pBR322 (50 pg each). Transformation efficiencies should be[.] $>1 \times 10^9$ cfu/µg DNA for pUC19 Relative transformation efficiency of methylated to non-methylated pBR322 should be >20%.

Overview

Description	DH10B-T1 ^R Competent <i>E. coli</i> is resistant to the lytic bacteriophages T1 and T5 (Killmann <i>et al.</i> , 1996). Since bacteriophage T1 transmits easily by aerosolization, it is a hazard in high throughput laboratories and genomic centers. This strain contains a mutation in the <i>ton</i> A locus, conferring resistance to both T1 and T5 bacteriophages.			
Applications	DH10B-T1 ^R is an <i>Escherichia coli</i> K-12 strain suitable for:			
	 Cloning DNA that contains methylcytosine and methyladenine (i.e. genomic DNA) (Blumenthal, 1989; Raleigh <i>et al.</i>, 1988; Upcroft and Healey, 1987; Woodcock <i>et al.</i>, 1989) 			
	Construction of gene banks			
	Generation of cDNA libraries using plasmid-derived vectors			
	• Blue/white screening of transformants on selective plates containing Bluo-gal or X-gal			
	• Plasmid rescue (Grant et al., 1990)			
	Transformation of large plasmids			
	High quality plasmid preparation			
	Lastly, it is resistant to the effects of ligase and ligation buffers.			
Important Important	One Shot [®] DH10B-T1 ^R <i>E. coli</i> does not require IPTG to induce expression from the <i>lac</i> promoter.			
	If blue/white screening is required to select for transformants spread 40 μ l of 40 mg/ml X-Gal in dimethylformamide on top of the agar. Let the X-Gal diffuse into the agar for approximately 1 hour.			
General Handling	Be extremely gentle when working with competent cells. Competent cells are highly sensitive to changes in temperature or mechanical lysis caused by pipetting. Transformation should be started immediately following the thawing of the cells on ice. Mix by swirling or tapping the tube gently, not by pipetting or vortexing.			

Transforming Chemically Competent Cells

Materials Supplied	You will need the following items for transformation:				
by the User	•	37°C shaking and non-shaking incubator			
	٠	10 cm diameter LB agar plates with appropriate antibiotic			
	•	Ice bucket with ice			
	•	42°C water bath			
Before Starting	•	Equilibrate a water bath to 42°C			
	٠	Warm the vial of SOC medium to room temperature			
	٠	Spread X-Gal onto LB agar plates with antibiotic, if desired			
	•	Warm the plates in a 37°C incubator for 30 minutes			
	•	Obtain a test tube rack that will hold all transformation tubes so that they can all be put into a 42°C water bath at once.			
Procedure	The instructions provided below are for general use. Plasmid DNA should be free of phenol, ethanol, protein, and detergents for maximum transformation efficiency.				
	1.	Briefly centrifuge the ligation reaction and place on ice.			
	2.	Thaw, on ice, one 50 μl vial of One Shot^ ${\ensuremath{\mathbb R}}$ cells for each ligation/transformation.			
	3.	Pipet 1 to 5 μ l of each ligation reaction directly into the competent cells and mix by tapping gently. Do not mix by pipetting up and down. Store the remaining ligation reaction at -20°C.			
	4.	Incubate the vial on ice for 30 minutes.			
	5.	Incubate for exactly 30 seconds in the 42°C water bath. Do not mix or shake.			
	6.	Remove vial from the 42°C bath and place on ice.			
	7.	Add 250 μ l of pre-warmed SOC medium to each vial. (SOC is a rich medium; sterile technique must be practiced to avoid contamination.)			
	8.	Place the vial in a microcentrifuge rack on its side and secure with tape to avoid loss of the vial. Shake the vial at 37°C for exactly 1 hour at 225 rpm in a shaking incubator.			
	9.	Spread 20 μ l to 200 μ l from each transformation vial on separate, labeled LB agar plates. We recommend that you plate two different volumes. Note : You may have to dilute cells 1:10 to obtain well-spaced colonies.			
	10.	Store the remaining transformation reaction at +4°C and plate out the next day, if desired.			
	11.	Invert the plates and incubate at 37°C overnight.			
	12.	Select colonies and analyze by plasmid isolation, PCR, or sequencing.			

Transforming Control Plasmid DNA

	If you do not obtain the expected number of colonies, we recommend that you test the efficiency of the competent cells. Transform DH10B-T1 ^R with the supercoiled pUC19 plasmid supplied with the kit as described below.			
	1. Prepare LB agar plates containing 100 μg/ml ampicillin.			
	2. Transform 5 μl (50 pg) pUC19 into 50 μl of competent cells.			
	3. Follow Steps 4-8 on page 4.			
	4. Dilute the transformation reaction 1:100 and plate 30 μ l on selective plates.			
	5. Incubate overnight at 37°C and count colonies. Calculate transformation efficiency using the formula below.			
Calculation	Calculate the transformation efficiency as transformants per 1 μ g of plasmid DNA. For chemically competent cells, use the formula below to calculate transformation efficiency:			
# of colonies 50 pg plasmid DNA	$- \times \frac{10^{6} \text{ pg}}{\mu \text{g}} \times \frac{300 \mu \text{l total transformation volume}}{30 \mu \text{l plated}} \times 100 = \frac{\# \text{ transformants}}{\mu \text{g plasmid DNA}}$			
	Expected transformation efficiency: $\geq 1 \ge 1 \ge 10^9$ cfu/µg supercoiled plasmid			
Note	Transformation efficiencies for cDNA and ligation of inserts to vectors will be lower than for a supercoiled control plasmid such as pUC19. For cDNA, transformation efficiencies may be 10- to 100-fold lower. For ligation of inserts to vectors, transformation efficiencies may be 10-fold lower.			

Technical Service

World Wide Web



Visit the Invitrogen Web Resource using your World Wide Web browser. At the site, you can:

- Get the scoop on our hot new products and special product offers
- View and download vector maps and sequences
- Download manuals in Adobe[®] Acrobat[®] (PDF) format
- Explore our catalog with full color graphics
- Obtain citations for Invitrogen products
- Request catalog and product literature

Once connected to the Internet, launch your Web browser (Internet Explorer 5.0 or newer or Netscape 4.0 or newer), then enter the following location (or URL):

http://www.invitrogen.com

...and the program will connect directly. Click on underlined text or outlined graphics to explore. Don't forget to put a bookmark at our site for easy reference!

Contact Us For more information or technical assistance, call, write, fax, or email. Additional international offices are listed on our Web page (www.invitrogen.com).

Corporate Headquarters:	Japanese Headquarters:	European Headquarters:
Invitrogen Corporation	Invitrogen Japan K.K.	Invitrogen Ltd
1600 Faraday Avenue	Nihonbashi Hama-Cho Park Bldg. 4F	Inchinnan Business Park
Carlsbad, CA 92008 USA	2-35-4, Hama-Cho, Nihonbashi	3 Fountain Drive
Tel: 1 760 603 7200	Tel: 81 3 3663 7972	Paisley PA4 9RF, UK
Tel (Toll Free): 1 800 955 6288	Fax: 81 3 3663 8242	Tel: +44 (0) 141 814 6100
Fax: 1 760 602 6500	E-mail: jpinfo@invitrogen.com	Tech Fax: +44 (0) 141 814 6117
E-mail:		E-mail: eurotech@invitrogen.com
tech_service@invitrogen.com		<u> </u>

MSDS Requests

To request an MSDS, visit our Web site at www.invitrogen.com. On the home page, go to 'Technical Resources', select 'MSDS', and follow instructions on the page.

Continued on next page

Technical Service, Continued

Limited Warranty Invitrogen is committed to providing our customers with high-quality goods and services. Our goal is to ensure that every customer is 100% satisfied with our products and our service. If you should have any questions or concerns about an Invitrogen product or service, please contact our Technical Service Representatives. Invitrogen warrants that all of its products will perform according to the specifications stated on the certificate of analysis. The company will replace, free of charge, any product that does not meet those specifications. This warranty limits Invitrogen Corporation's liability only to the cost of the product. No warranty is granted for products beyond their listed expiration date. No warranty is applicable unless all product components are stored in accordance with instructions. Invitrogen reserves the right to select the method(s) used to analyze a product unless Invitrogen agrees to a specified method in writing prior to acceptance of the order. Invitrogen makes every effort to ensure the accuracy of its publications, but realizes that the occasional typographical or other error is inevitable. Therefore Invitrogen makes no warranty of any kind regarding the contents of any publications or documentation. If you discover an error in any of our publications, please report it to our Technical Service Representatives. Invitrogen assumes no responsibility or liability for any special, incidental, indirect or consequential loss or damage whatsoever. The above limited warranty is sole and exclusive. No other warranty is made, whether expressed or implied, including any warranty of merchantability or fitness for a particular purpose.

References

Blumenthal, R. M. (1989). Cloning and Restriction of Methylated DNA in Escherichia coli. Focus 11, 41-46.

- Grant, S. G., Jessee, J., Bloom, F. R., and Hanahan, D. (1990). Differential Plasmid Rescue from Transgenic Mouse DNAs into Escherichia coli Methylation-restriction Mutants. Proc. Natl. Acad. Sci. USA 87, 4645-4649.
- Killmann, H., Benz, R., and Braun, V. (1996). Properties of the FhuA channel in the Escherichia coli outer membrane after deletion of FhuA portions within and outside the predicted gating loop. J Bacteriol 178, 6913-20.
- Raleigh, E. A., Murray, N. E., Revel, H., Blumenthal, R. M., Westaway, D., Reith, A. D., Rigby, P. W., Elhai, J., and Hanahan, D. (1988). McrA and McrB Restriction Phenotypes of Some *E. coli* Strains and Implications for Gene Cloning. Nuc. Acids Res. 16, 1563-1575.
- Upcroft, P., and Healey, A. (1987). Rapid and Efficient Method for Cloning of Blunt-Ended DNA Fragments. Gene *51*, 69-75.
- Woodcock, D. M., Crowther, P. J., Doherty, J., Jefferson, S., DeCruz, E., Noyer-Weidner, M., Smith, S. S., Michael, M. Z., and Graham, M. W. (1989). Quantitative Evaluation of *Escherichia coli* Host Strains for Tolerance to Cytosine Methylation in Plasmid and Phage Recombinants. Nuc. Acids Res. 17, 3469-3478.

©2001-2004 Invitrogen Corporation. All rights reserved.

For research use only. Not intended for any animal or human therapeutic or diagnostic use.