

pCEP4

Catalog no. V044-50

Rev. date: 28 December 2010 Manual part no. 28-0038 MAN0000615

User Manual

Table of Contents

| Kit Contents and Storage. | iv |
|------------------------------|----|
| Introduction | 1 |
| Product Overview . | 1 |
| Methods . | 2 |
| Cloning into pCEP4 | 2 |
| Transfection. | 5 |
| Creating Stable Cell Lines . | 6 |
| Appendix | 8 |
| pCEP4 Vector . | 8 |
| pCEP4/CAT. | |
| - Accessory Products . | |
| Technical Support | |
| References | |

Kit Contents and Storage

Shipping andpCEP4 vectors are shipped on wet ice. Upon receipt, store vectors at -20°C.Storage

Kit Contents All vectors are supplied as detailed below. **Store the vectors at –20°C.**

| Vector | Composition | Amount |
|-----------|--|--------|
| pCEP4 | 40 μL of 0.5 μg/μL vector in 10 mM Tris- HCl, 1 mM EDTA, pH 8.0 | 20 µg |
| pCEP4/CAT | 20 μL of 0.5 μg/μL vector in 10 mM Tris- HCl, 1 mM EDTA, pH 8.0 | 10 µg |

Introduction

Product Overview

| pCEP4 | pCEP4 is an episomal mammalian expression vector that uses the cytomegalovirus (CMV) immediate early enhancer/promoter for high level transcription of recombinant genes inserted into the multiple cloning site. The Epstein-Barr Virus replication origin (oriP) and nuclear antigen (encoded by the EBNA-1 gene) is carried by this plasmid to permit extrachromosomal replication in human, primate, and canine cells. pCEP4 also carries the hygromycin B resistance gene for stable selection in transfected cells. pCEP4/CAT is provided as a positive control for the relative level of expression of recombinant proteins in a cell line of interest. It expresses the chloramphenicol acetyl transferase (CAT) protein from the CMV enhancer/promoter. Like its | |
|-------------------------|---|--|
| | parent vector pCEP4, pCEP4/CAT contains the hygromycin B resistance gene for selection. | |
| Experimental Outline | Use the following outline to clone and express your gene of interest in pCEP4.1. Consult the multiple cloning site described on page 3 to design a strategy to clone your gene into pCEP4. | |
| | Ligate your insert into the appropriate vector and transform into <i>E. coli</i>. Select transformants on LB plates containing 50 to 100 μg/mL ampicillin. | |
| | 3. Analyze your transformants for the presence of insert by restriction digestion. | |
| | 4. Select a transformant with the correct restriction pattern and use sequencing to confirm that your gene is cloned in the proper orientation. | |
| | 5. Transfect your construct into the mammalian cell line of interest using your own method of choice. Generate a stable cell line, if desired. | |
| | 6. Test for expression of your recombinant gene by western blot analysis or functional assay. | |
| | | |

Methods

Cloning into pCEP4

| General Molecular Biology Techniques | For help with DNA ligations, <i>E. coli</i> transformations, restriction enzyme analysis, purification of single-stranded DNA, DNA sequencing, and DNA biochemistry, refer to <i>Molecular Cloning: A Laboratory Manual</i> (Sambrook <i>et al.</i> , 1989) or <i>Current Protocols in Molecular Biology</i> (Ausubel <i>et al.</i> , 1994). | | | |
|--|--|--|--|--|
| <i>E. coli</i> Strain | Many <i>E. coli</i> strains are suitable for the propagation of this vector. We recommend that you propagate vectors containing inserts in <i>E. coli</i> strains that are recombination deficient (<i>rec</i> A) and endonuclease A-deficient (<i>end</i> A). | | | |
| | For your convenience, TOP10 is available as chemically competent or electrocompetent cells from Invitrogen (see page 11). | | | |
| Maintaining pCEP4 | To propagate and maintain pCEP4, use a small amount of the supplied $0.5 \ \mu g/\mu L$ stock solution in TE, pH 8.0 to transform a <i>recA</i> , <i>endA E</i> . <i>coli</i> strain like TOP10 or equivalent. Select transformants on LB plates containing 50 to 100 $\mu g/mL$ ampicillin. Be sure to prepare a glycerol stock of your plasmid-containing <i>E</i> . <i>coli</i> strain for long-term storage (see page 4). | | | |
| Important Important | Your insert should contain a Kozak consensus sequence with an ATG initiation codon for proper initiation of translation (Kozak, 1987; Kozak 1990). An example of a Kozak consensus sequence is provided below. Other sequences are possible, but the G or A at position -3 and the G at position $+4$ (shown in bold) illustrates the most commonly occurring sequence with strong consensus. Replacing one of the two bases at these positions provides moderate consensus, while having neither results in weak consensus. The ATG initiation codon is shown underlined. $(G/A)NNATGG$ | | | |
| | Your insert should also contain a stop codon for proper termination of your gene. | | | |
| | | | | |

Cloning into pCEP4, Continued

| | ple Cloning of pCEP4 | indicate the cle sequencing and downloading | eavage site. The r d functional testi f rom our websit | ite for pCEP4. Re nultiple cloning s ng. The sequenc e (<u>www.invitrog</u> ap and a descript | site has been con e of pCEP4 is av <u>en.com</u>) or from | firmed by ailable for Technical |
|-----|-----------------------------------|--|---|---|--|---------------------------------------|
| 1 | GTTGACATTG | ATTATTGACT | AGTTATTAAT | AGTAATCAAT | TACGGGGTCA | TTAGTTCATA |
| | | enh | ancer region (5' end) |) | | |
| 61 | GCCCATATAT | GGAGTTCCGC | GTTACATAAC | TTACGGTAAA | TGGCCCGCCT | GGCTGACCGC |
| 121 | CCAACGACCC | CCGCCCATTG | ACGTCAATAA | TGACGTATGT | TCCCATAGTA | ACGCCAATAG |
| 181 | GGACTTTCCA | TTGACGTCAA | TGGGTGGAGT | ATTTACGGTA | AACTGCCCAC | TTGGCAGTAC |
| 241 | ATCAAGTGTA | TCATATGCCA | AGTCCGCCCC | CTATTGACGT | CAATGACGGT | AAATGGCCCG |
| 301 | CCTGGCATTA | TGCCCAGTAC | ATGACCTTAC | GGGACTTTCC | TACTTGGCAG | TACATCTACG |
| 361 | | CGCTATTACC | ATGGTGATGC | GGTTTTGGCA | GTACACCAAT | GGGCGTGGAT |
| 421 | AGCGGTTTGA enhancer region (3' | CTCACGGGGA end) | TTTCCAAGTC | TCCACCCCAT | TGACGTCAAT | GGGAGTTTGT |
| 481 | TTTGGCACCA | AAATCAACGG | GACTTTCCAA | AATGTCGTAA | TAACCCCGCC pCEP Forward | |
| 541 | I AAATGGGCGG | TAGGCGTGTA | CGGTGGGAGG | I I I TCTATATAAG | CAGAGCTCGT | TTAGTGAACC |
| | | ptional start | Kpn Pvu Nh | | Nhe I* Not I | Xho I |
| 601 | J GTCAGATCTC | TAGAAGCTGG | · · · · | CTAGCAAGCT | TGCTAGCGGC | CGCTCGAGGC |
| 661 | | BamHI GGATCCAGAC | ATGATAAGAT | | EBV Reverse primer | АССАСААСТА |

Cloning into pCEP4, Continued

| <i>E. coli</i> Transformation | Transform your ligation mixtures into a competent <i>recA</i> , <i>endA E</i> . <i>coli</i> strain (<i>e.g.</i> TOP10) and select transformants on LB plates containing 50 to 100 μ g/mL ampicillin. Select 10–20 clones and analyze for the presence and orientation of your insert. | | | |
|---|--|--|--|--|
| | Transformation Method : You may use any method of your choice for transformation. Chemical transformation is the most convenient for most researchers. Electroporation is the most efficient and the method of choice for large plasmids. | | | |
| Sequencing Your Construct | Several primers are available separately that you may use to sequence your construct. These are marked in the multiple cloning site diagram on page 3. For ordering information, see page 11. Alternatively, you may design your own primer for sequencing. | | | |
| Preparing a Glycerol Stock | Once you have identified the correct clone, purify the colony and make a glycerol stock for long-term storage. You should keep a DNA stock of your plasmid at -20°C. | | | |
| | Streak the original colony out on an LB plate containing 50 μg/mL ampicillin. Incubate the plate at 37°C overnight. | | | |
| | 2. Isolate a single colony and inoculate into $1-2$ mL of LB containing 50 μ g/mL ampicillin. | | | |
| | 3. Grow the culture to mid-log phase ($OD_{600} = 0.5-0.7$). | | | |
| Mix 0.85 mL of culture with 0.15 mL of sterile glycerol and transferred cryovial. | | | | |
| | 5. Store at -80° C. | | | |
| | | | | |

Transfection

| Introduction | Once you have verified that your gene is cloned in the correct orientation and contains an initiation ATG and a stop codon, you are ready to transfect your cell line of choice. We recommend that you include the positive control vector and a mock transfection (negative control) to evaluate your results. | | |
|----------------------------|--|--|--|
| Plasmid Preparation | Plasmid DNA for transfection into eukaryotic cells must be clean and free from phenol and sodium chloride. Contaminants will kill the cells, and salt will interfere with lipid complexing, decreasing transfection efficiency. We recommend isolating plasmid DNA using the PureLink [™] HQ Mini Plasmid Purification Kit or PureLink [™] HiPure Miniprep or Midiprep Kits available from Invitrogen (see page 11). | | |
| Methods of Transfection | For established cell lines (<i>e.g.</i> HeLa), consult original references or the supplier of your cell line for the optimal method of transfection. We recommend that you follow exactly the protocol for your cell line. Pay particular attention to medium requirements, when to pass the cells, and at what dilution to split the cells. Further information is provided in <i>Current Protocols in Molecular Biology</i> (Ausubel <i>et al.</i> , 1994). | | |
| | Methods for transfection include calcium phosphate (Chen and Okayama, 1987; Wigler <i>et al.</i> , 1977), lipid-mediated (Felgner <i>et al.</i> , 1989; Felgner and Ringold, 1989) and electroporation (Chu <i>et al.</i> , 1987; Shigekawa and Dower, 1988). Invitrogen offers the Lipofectamine [™] 2000 Reagent for mammalian cell transfection (see page 11). | | |
| Positive Control | pCEP4/CAT is provided as a positive control vector for mammalian transfection and expression (see page 10) and may be used to optimize transfection conditions for your cell line. The gene encoding chloramphenicol acetyl transferase (CAT) is expressed in mammalian cells under the control of the CMV promoter. A successful transfection will result in CAT expression that can be easily assayed (see below). | | |
| Assay for CAT Protein | You may assay for CAT expression using your method of choice. Invitrogen offers a CAT assay kit for detection of the protein (page 11). | | |

Creating Stable Cell Lines

| Introduction | pCEP4 contains the hygromycin resistance gene for selection of stable cell lines using hygromycin B. We recommend that you test the sensitivity of your mammalian host cell to hygromycin B, as natural resistance varies among cell lines. General information and guidelines are provided in this section for your convenience. | | | | |
|--|---|--|--|--|--|
| Hygromycin B Activity | Hygromycin B (527.5 MW) is an aminocyclitol that inhibits protein synthesis by disrupting translocation and promoting mistranslation. Hygromycin B-phospho-transferase detoxifies hygromycin-B by phosphorylation. | | | | |
| CAUTION | • Hygromycin is light sensitive. Store the liquid stock solution at +4°C protected from exposure to light. | | | | |
| | • Hygromycin is toxic. Do not ingest solutions containing the drug. | | | | |
| | Wear gloves, a laboratory coat, and safety glasses or goggles when handling hygromycin and hygromycin-containing solutions. | | | | |
| Determining Antibiotic Sensitivity | To successfully generate a stable cell line expressing your gene of interest from pCEP4, you need to determine the minimum concentration of hygromycin B required to kill your untransfected host cell line. Typically, concentrations ranging from 10 to $400 \mu g/mL$ hygromycin are sufficient to kill most untransfected mammalian cell lines. We recommend that you test a range of concentrations (see protocol below) to ensure that you determine the minimum concentration necessary for your host cell line. | | | | |
| | 1. Plate or split a confluent plate so the cells will be approximately 25% confluent. Prepare a set of 7 plates. Allow cells to adhere overnight. | | | | |
| | The next day, substitute culture medium with medium containing varying concentrations of hygromycin (0, 10, 25, 50, 100, 200, 400 μg/mL hygromycin). | | | | |
| | 3. Replenish the selective media every 3–4 days, and observe the percentage of surviving cells. | | | | |
| | Count the number of viable cells at regular intervals to determine the appropriate concentration of hygromycin that prevents growth within 2–3 weeks after addition of hygromycin. | | | | |
| | Note: Cells will divide once or twice in the presence of lethal doses of hygromycin, so the effects of the drug may take several days to become apparent. Complete inhibition of cell growth can take 2–3 weeks of growth in selective medium. | | | | |

Creating Stable Cell Lines, Continued


| Selecting Stable Integrants | Once you have determined the appropriate hygromycin concentration to use for selection in your host cell line, you can generate a stable cell line expressing your gene of interest. | | |
|-------------------------------------|---|--|--|
| | 1. Transfect your mammalian host cell line with your pCEP4 construct using the desired protocol. Remember to include a plate of untransfected cells as a negative control and the pCEP4/CAT plasmid as a positive control. | | |
| | 2. 24 hours after transfection, wash the cells and add fresh medium to the cells. | | |
| | 48 hours after transfection, split the cells into fresh medium containing hygromycin at the pre-determined concentration required for your cell line. Split the cells such that they are no more than 25% confluent. | | |
| | 4. Feed the cells with selective medium every 3–4 days until hygromycin- resistant foci can be identified. | | |
| | 5. Pick and expand colonies in 96- or 48-well plates. | | |
| Maintaining Stable Transfectants | p CEP4 is an episomally-maintained plasmid (Reisman and Sugden, 1986; Yates <i>al.</i> , 1985). Transfected cells may lose the pCEP4 plasmid if they are not maintained under selection or are continuously cultured for long periods of time (over six months). To prevent loss of pCEP4 from transfected cells, we recommend that you follow these guidelines when working with the cells: | | |
| | • Always use early-passage cells. Grow and freeze multiple vials of transfected cells to ensure that you have an adequate supply of early-passage cells. | | |
| | • Always maintain cells in medium containing 50 µg/mL hygromycin. | | |
| | • Do not maintain cells in continuous culture for longer than 6 months. | | |

Appendix

pCEP4 Vector

Map of pCEP4

The figure below summarizes the features of the pCEP4 vector. **The sequence for** pCEP4 is available for downloading from our website (<u>www.invitrogen.com</u>) or from Technical Support (see page 12).



pCEP4 Vector, Continued

Features of pCEP4 (10,186 bp) contains the following elements. All features have been functionally tested.

| Feature | Benefit |
|--|--|
| Human cytomegalovirus (CMV) immediate-early promoter/enhancer | Allows efficient, high-level expression of your recombinant protein (Andersson et al., 1989; Boshart et al., 1985; Nelson et al., 1987) |
| Multiple cloning site | Allows insertion of your gene and facilitates cloning |
| SV40 polyadenylation signal | Efficient transcription termination and polyadenylation of mRNA |
| EBV origin of replication (oriP) and nuclear antigen (EBNA-1) | High-copy episomal replication in primate and canine cell lines (Reisman and Sugden, 1986; Yates et al., 1985) |
| Ampicillin resistance gene | Selection of vector in <i>E. coli</i> |
| (β-lactamase) | |
| pUC origin | High-copy number replication and growth in <i>E. coli</i> |
| Herpes Simplex Virus thymidine kinase (TK) promoter | Allows efficient, high-level expression of the hygromycin resistance gene (McKnight, 1980) |
| Hygromycin resistance gene | Selection of stable transfectants in mammalian cells (Gritz and Davies, 1983; Palmer et al., 1987) |
| Herpes Simplex Virus thymidine kinase (TK) promoter polyadenylation signal | Efficient transcription termination and polyadenylation of mRNA |


pCEP4/CAT

10891 nucleotides

Map of

pCEP4/CAT

The figure below summarizes the features of the pCEP4/CAT vector. The nucleotide sequence for pCEP4/CAT is available for downloading from our website (<u>www.invitrogen.com</u>) or by contacting Technical Support (page 12).



CMV promoter: bases 1-588 Chloramphenicol acetyl transferase (CAT) gene: bases 675-1334 SV40 polyadenylation signal: bases 1390-1631 OriP: bases 2049-4024 EBNA-1 gene (complementary strand): bases 4325-6250 Ampicillin resistance gene: bases 6876-7736 pUC origin: bases 7745-8520 TK promoter: bases 8888-9050 Hygromycin resistance gene: bases 9114-10124 TK polyadenylation signal: bases 10136-10407

Accessory Products

Introduction

The products listed below are designed for use with pBudCE4.1. For details, visit www.invitrogen.com or contact **Technical Support** (page 12).

| Item | Quantity | Catalog no. |
|--|------------------------|-------------|
| One Shot® TOP10 Chemically Competent Cells | $21 \times 50 \ \mu L$ | C4040-03 |
| One Shot [®] TOP10 Electrocomp [™] Cells | $21 \times 50 \ \mu L$ | C4040-52 |
| PureLink™ HiPure Plasmid Miniprep Kit | 100 preps | K2100-03 |
| PureLink™ HiPure Plasmid Midiprep Kit | 25 preps | K2100-04 |
| PureLink [™] HQ Mini Plasmid Purification Kit | 100 preps | K2100-01 |
| Lipofectamine [™] 2000 Reagent | 1.5 mL | 11668-019 |
| Hygromycin B | 20 mL | 10687-010 |
| Fast Cat [®] Chloramphenicol Acetyltransferase Assay Kit | 1 kit | F2900 |

Primers

For your convenience, Invitrogen offers a custom primer synthesis service. Visit www.invitrogen.com for more details.

Technical Support

Web Resources



- Visit the Invitrogen website at <u>www.invitrogen.com</u> for:
- Technical resources, including manuals, vector maps and sequences, application notes, MSDSs, FAQs, formulations, citations, handbooks, etc.
- Complete technical support contact information
- Access to the Invitrogen Online Catalog
- Additional product information and special offers

Contact Us

For more information or technical assistance, call, write, fax, or email. Additional international offices are listed on our website (<u>www.invitrogen.com</u>).

| Corporate Headquarte 5791 Van Allen Way Carlsbad, CA 92008 US Tel: 1 760 603 7200 Tel (Toll Free): 1 800 95 Fax: 1 760 602 6500 E-mail: <u>tech_support@</u> | 55 6288 | Japanese Headquarters: LOOP-X Bldg. 6F 3-9-15, Kaigan Minato-ku, Tokyo 108-0022 Tel: 81 3 5730 6509 Fax: 81 3 5730 6519 E-mail: jpinfo@invitrogen.com | European Headquarters: Inchinnan Business Park 3 Fountain Drive Paisley PA4 9RF, UK Tel: +44 (0) 141 814 6100 Tech Fax: +44 (0) 141 814 6117 E-mail: <u>eurotech@invitrogen.com</u> |
|--|---|--|---|
| MSDS | Material Safety Data Sheets (MSDSs) are available on our website at www.invitrogen.com/msds. | | |
| Certificate of Analysis | The Certificate of Analysis provides detailed quality control and product qualification information for each product. Certificates of Analysis are available on our website. Go to www.invitrogen.com/support and search for the Certificate of Analysis by product lot number, which is printed on the box. | | |
| with high-quality satisfied with our about an Invitrog All Invitrogen pro- certificate of analy meet those specifi <u>product.</u> No warr warranty is applied instructions. The of product unless the the order. Invitrogen makes occasional typogr warranty of any k discover an error Representatives. Life Technologie incidental, indire warranty is sole a | | of Life Technologies Corporation) is c goods and services. Our goal is to ens products and our service. If you shoul en product or service, contact our Tech oducts are warranted to perform accor- ysis. The Company will replace, free of ications. <u>This warranty limits the Comp</u> anty is granted for products beyond the cable unless all product components an Company reserves the right to select the e Company agrees to a specified methor every effort to ensure the accuracy of raphical or other error is inevitable. The cind regarding the contents of any pub in any of our publications, please repo s Corporation shall have no responsib- ter or consequential loss or damage ward exclusive. No other warranty is marranty of merchantability or fitness for the section of the section of the section of the section of the | ure that every customer is 100% Id have any questions or concerns unical Support Representatives. ding to specifications stated on the charge, any product that does not pany's liability to only the price of the neir listed expiration date. No re stored in accordance with ne method(s) used to analyze a od in writing prior to acceptance of its publications, but realizes that the erefore the Company makes no lications or documentation. If you rt it to our Technical Support bility or liability for any special, hatsoever. The above limited ade, whether expressed or implied, |

References

- Andersson, S., Davis, D. L., Dahlbäck, H., Jörnvall, H., and Russell, D. W. (1989). Cloning, Structure, and Expression of the Mitochondrial Cytochrome P-450 Sterol 26-Hydroxylase, a Bile Acid Biosynthetic Enzyme. J. Biol. Chem. 264, 8222-8229.
- Ausubel, F. M., Brent, R., Kingston, R. E., Moore, D. D., Seidman, J. G., Smith, J. A., and Struhl, K. (1994). Current Protocols in Molecular Biology (New York: Greene Publishing Associates and Wiley-Interscience).
- Boshart, M., Weber, F., Jahn, G., Dorsch-Häsler, K., Fleckenstein, B., and Schaffner, W. (1985). A Very Strong Enhancer is Located Upstream of an Immediate Early Gene of Human Cytomegalovirus. Cell *41*, 521-530.
- Chen, C., and Okayama, H. (1987). High-Efficiency Transformation of Mammalian Cells by Plasmid DNA. Molec. Cell. Biol. 7, 2745-2752.
- Chu, G., Hayakawa, H., and Berg, P. (1987). Electroporation for the Efficient Transfection of Mammalian Cells with DNA. Nucleic Acids Res. *15*, 1311-1326.
- Felgner, P. L., Holm, M., and Chan, H. (1989). Cationic Liposome Mediated Transfection. Proc. West. Pharmacol. Soc. 32, 115-121.
- Felgner, P. L. a., and Ringold, G. M. (1989). Cationic Liposome-Mediated Transfection. Nature 337, 387-388.
- Gritz, L., and Davies, J. (1983). Plasmid-Encoded Hygromycin-B Resistance: The Sequence of Hygromycin-B-Phosphotransferase Gene and its Expression in *E. coli* and *S. Cerevisiae*. Gene 25, 179-188.
- Kozak, M. (1987). An Analysis of 5´-Noncoding Sequences from 699 Vertebrate Messenger RNAs. Nucleic Acids Res. 15, 8125-8148.
- Kozak, M. (1991). An Analysis of Vertebrate mRNA Sequences: Intimations of Translational Control. J. Cell Biology 115, 887-903.
- Kozak, M. (1990). Downstream Secondary Structure Facilitates Recognition of Initiator Codons by Eukaryotic Ribosomes. Proc. Natl. Acad. Sci. USA *87*, 8301-8305.
- McKnight, S. L. (1980). The Nucleotide Sequence and Transcript Map of the Herpes Simplex Virus Thymidine Kinase Gene. Nucleic Acids Res. *8*, 5949-5964.
- Nelson, J. A., Reynolds-Kohler, C., and Smith, B. A. (1987). Negative and Positive Regulation by a Short Segment in the 5´-Flanking Region of the Human Cytomegalovirus Major Immediate-Early Gene. Molec. Cell. Biol. 7, 4125-4129.
- Palmer, T. D., Hock, R. A., Osborne, W. R. A., and Miller, A. D. (1987). Efficient Retrovirus-Mediated Transfer and Expression of a Human Adenosine Deaminase Gene in Diploid Skin Fibroblasts from an Adenosine-Deficient Human. Proc. Natl. Acad. Sci. U.S.A. 84, 1055-1059.
- Reisman, D., and Sugden, B. (1986). *trans* Activation of an Epstein-Barr Viral Transcriptional Enhancer by the Epstein-Barr Viral Nuclear Antigen 1. Mol. Cell. Biol. *6*, 3838-3846.
- Sambrook, J., Fritsch, E. F., and Maniatis, T. (1989). Molecular Cloning: A Laboratory Manual, Second Edition (Plainview, New York: Cold Spring Harbor Laboratory Press).
- Shigekawa, K., and Dower, W. J. (1988). Electroporation of Eukaryotes and Prokaryotes: A General Approach to the Introduction of Macromolecules into Cells. BioTechniques *6*, 742-751.
- Wigler, M., Silverstein, S., Lee, L.-S., Pellicer, A., Cheng, Y.-C., and Axel, R. (1977). Transfer of Purified Herpes Virus Thymidine Kinase Gene to Cultured Mouse Cells. Cell *11*, 223-232.

Yates, J. L., Warren, N., and Sugden, B. (1985). Stable Replication of Plasmids Derived from Epstein-Barr Virus in Various Mammalian Cells. Nature *313*, 812-815.

©2009, 2010 Life Technologies Corporation. All rights reserved.

For research use only. Not intended for any animal or human therapeutic or diagnostic use.

invitrogen

Corporate Headquarters 5791 Van Allen Way Carlsbad, CA 92008 T: 1 760 603 7200 F: 1 760 602 6500 E: tech_support@invitrogen.com

For country-specific contact information, visit our web site at www.invitrogen.com