

## SuperScript® One-Cycle cDNA Kit for Use with Affymetrix® One-Cycle Assays

Cat. no. A10752-030

Size: 30 reactions

Store at  $-20^{\circ}\text{C}$

### Kit Components

The following components are provided with the kit. Store all reagents at  $-20^{\circ}\text{C}$ .

Component	Amount
T7-Oligo(dT), 50 $\mu\text{M}$	60 $\mu\text{l}$
5X 1st-Strand Reaction Mix	120 $\mu\text{l}$
DTT, 0.1 M	60 $\mu\text{l}$
dNTP Mix, 10 mM (10 mM each dATP, dCTP, dGTP, dTTP)	120 $\mu\text{l}$
SuperScript® II RT (200 U/ $\mu\text{l}$ )	60 $\mu\text{l}$
5X 2nd-Strand Reaction Mix	900 $\mu\text{l}$
<i>E. coli</i> DNA Ligase (10 U/ $\mu\text{l}$ )	30 $\mu\text{l}$
<i>E. coli</i> DNA Polymerase I (10 U/ $\mu\text{l}$ )	120 $\mu\text{l}$
RNase H (2 U/ $\mu\text{l}$ )	30 $\mu\text{l}$
T4 DNA Polymerase (5 U/ $\mu\text{l}$ )	60 $\mu\text{l}$
RNase-free water	3.1 ml
EDTA, 0.5 M	300 $\mu\text{l}$


### Description

The SuperScript® One-Cycle cDNA Kit contains all of the reagents necessary to make double-stranded cDNA from purified total RNA or mRNA for use with the GeneChip® One-Cycle Target Labeling System from Affymetrix. This kit is designed to fit into the standard one-cycle labeling protocol for GeneChip® arrays from Affymetrix, as described in the GeneChip® Expression Analysis Technical Manual.

It is important to use high-quality RNA isolated from tissue or cells for your cDNA preparation. High-quality RNA can be obtained by isolating total RNA with the PureLink™ Micro-to-Midi Total RNA Purification System or TRIzol® Reagent. See page 4 for ordering information.

### Workflow Overview

Following isolation of total RNA or mRNA from tissue or cells, use this kit to generate double-stranded cDNA as shown in the workflow below. Then proceed directly to purification using the GeneChip® Sample Cleanup Module from Affymetrix.



## Before Starting

- Briefly centrifuge all reagents in the kit before use.
- For best results, use a thermocycler for all incubation steps.
- RNA can be quantified using UV absorbance at 260 nm or a quantitation kit such as the Quant-iT™ RNA Assay Kit. RNA quality can be assessed by agarose gel electrophoresis or on an Agilent 2100 Bioanalyzer.

## First-Strand Synthesis

The 20- $\mu$ l reaction described in this protocol is designed to convert 1–15  $\mu$ g of total RNA or 0.2–2  $\mu$ g of mRNA into first-strand cDNA. The amount of SuperScript® II RT added to the reaction will depend on the amount of starting RNA.

- **Total RNA:** We recommend using 200 units of SuperScript® II RT for <8  $\mu$ g total RNA or 400 units for 8–15  $\mu$ g total RNA.
  - **mRNA:** We recommend using 200 units of SuperScript® II RT per  $\mu$ g of mRNA.
  - **Poly-A RNA controls:** The GeneChip® Eukaryotic Poly-A RNA Control Kit provides poly-A RNA controls specifically designed for the Affymetrix® GeneChip® One-Cycle Target Labeling system. If you are using these controls, dilute them as described in the GeneChip® Expression Analysis Technical Manual, then add 2  $\mu$ l of the diluted controls per reaction into the first-strand synthesis reaction as shown in Step 1 below.
1. For each reaction, add T7 oligo(dT) primer, RNA, RNase-free water, and optional controls according to the table below to a RNase-free 1.5-ml microcentrifuge tube:

Component	Volume if using <8 $\mu$ g Total RNA or <1 $\mu$ g mRNA	Volume if using 8–15 $\mu$ g Total RNA or 1–2 $\mu$ g mRNA
Sample RNA ( $\mu$ l)	variable	variable
<i>Optional:</i> Diluted GeneChip® Eukaryotic Poly-A RNA Controls	2 $\mu$ l	2 $\mu$ l
T7-oligo(dT), 50 $\mu$ M	2 $\mu$ l	2 $\mu$ l
RNase-free water	to 12 $\mu$ l	to 11 $\mu$ l

2. Flick the tube gently to mix, then centrifuge briefly to collect the contents. Heat the mixture to 70°C for 10 minutes, and then chill on ice.
3. In a separate tube, prepare a master mix of the following reagents. Multiply the number of reactions by the volume per reaction listed in the table below; be sure to prepare extra master mix to allow for pipetting variations.

Component	Volume per reaction
5X 1st-Strand Reaction Mix	4 $\mu$ l
DTT, 0.1 M	2 $\mu$ l
dNTP Mix, 10 mM	1 $\mu$ l
Total volume	7 $\mu$ l

4. Mix the master mix gently by flicking the tube and centrifuge briefly to collect the contents. Add 7  $\mu$ l of the first-strand master mix to each reaction tube.
5. Vortex gently and then centrifuge briefly to collect the tube contents. Immediately incubate each tube at 42°C for 2 minutes.
6. Add SuperScript® II RT to each tube as shown in the following table, to bring the total volume to 20  $\mu$ l:

Component	Volume if using <8 $\mu$ g Total RNA or <1 $\mu$ g mRNA	Volume if using 8–15 $\mu$ g Total RNA or 1–2 $\mu$ g mRNA
SuperScript® II RT	1 $\mu$ l	2 $\mu$ l

7. Mix gently, and incubate at 42°C for 1 hour.
8. Place the tube on ice to terminate the reaction, and proceed to **Second-Strand Synthesis**.

## Second-Strand Synthesis

1. In a separate tube, prepare a master mix of the following reagents. Multiply the number of reactions by the volume per reaction listed in the table below; be sure to prepare extra master mix to allow for pipetting variations.

Component	Volume per reaction
RNase-free water	91 $\mu$ l
5X 2nd-Strand Reaction Mix	30 $\mu$ l
dNTP Mix, 10 mM	3 $\mu$ l
<i>E. coli</i> DNA Ligase (10 U/ $\mu$ l)	1 $\mu$ l
<i>E. coli</i> DNA Polymerase I (10 U/ $\mu$ l)	4 $\mu$ l
RNase H (2 U/ $\mu$ l)	1 $\mu$ l
Final volume	130 $\mu$ l

2. Mix the master mix gently by flicking the tube and centrifuge briefly to collect the contents. On ice, add 130  $\mu$ l of the master mix to each 20- $\mu$ l first-strand reaction tube.
3. Vortex gently to mix, and incubate for 2 hours at 16°C. Do not allow the temperature to rise above 16°C.
3. Add 2  $\mu$ l (10 units) of T4 DNA Polymerase, and continue to incubate at 16°C for 5 minutes.
4. Place the tube on ice, and add 10  $\mu$ l of EDTA, 0.5 M, to stop the reaction.
5. Proceed to cleanup of the double-stranded cDNA product using the GeneChip® Sample Cleanup Module from Affymetrix, followed by *in vitro* transcription and target labeling as described in the GeneChip® Expression Analysis Technical Manual. If you are not immediately proceeding to cleanup, store the double-stranded cDNA at -20°C.

## Product Qualification

The Certificate of Analysis (CofA) provides detailed quality control information for this product. You can search for the Certificate of Analysis on our product support page at [www.invitrogen.com/support](http://www.invitrogen.com/support). Enter the product lot number in the search field and select the **Certificates of Analysis** search. Note that the lot number is printed on each box.

## Additional Products

Related products are available separately from Invitrogen. Ordering information is provided below. For more information, visit our website at [www.invitrogen.com](http://www.invitrogen.com) or contact Technical Support.

Product	Amount	Catalog no.
Quant-iT™ RNA Assay Kit	1 kit	Q-33140
PureLink™ Micro-to-Midi™ Total RNA Purification System	50 rxns	12183-018
TRIzol® Reagent	100 ml	15596-026
	200 ml	15596-018

### Limited Use Label License No. 5: Invitrogen Technology

The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. The buyer may transfer information or materials made through the use of this product to a scientific collaborator, provided that such transfer is not for any Commercial Purpose, and that such collaborator agrees in writing (a) not to transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for research and not for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. For products that are subject to multiple limited use label licenses, the most restrictive terms apply. Invitrogen Corporation will not assert a claim against the buyer of infringement of patents owned or controlled by Invitrogen Corporation which cover this product based upon the manufacture, use or sale of a therapeutic, clinical diagnostic, vaccine or prophylactic product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. If the purchaser is not willing to accept the limitations of this limited use statement, Invitrogen is willing to accept return of the product with a full refund. For information on purchasing a license to this product for purposes other than research, contact Licensing Department, Invitrogen Corporation, 1600 Faraday Avenue, Carlsbad, California 92008. Phone (760) 603-7200. Fax (760) 602-6500. Email: [outlicensing@invitrogen.com](mailto:outlicensing@invitrogen.com).

### Trademarks of Other Companies

Affymetrix® and GeneChip® are registered trademarks of Affymetrix, Inc.

The uses herein should not imply a sponsorship or endorsement by, or affiliation with, Affymetrix, Inc.

TRIzol® is a registered trademark of Molecular Research Center, Inc.

©2008 Invitrogen Corporation. All rights reserved.

For research use only. Not intended for any animal or human therapeutic or diagnostic use.