DEDICATED TO MICROBIOLOGY

Oxoid Limited, Wade Road, Basingstoke, Hants RG24 8PW, UK

Made in the United Kingdom

EN

Oxoid™ Antimicrobial Susceptibility Disc Dispenser

REF ST6090, ST6100, ST8090, ST8100

INTENDED USE

The Thermo Scientific™ Oxoid™ Antimicrobial Susceptibility Disc Dispenser is a 6 or 8 place height adjustable disc dispenser (with a 90 mm or a 100 mm skirt) used as an accessory to the Oxoid Antimicrobial Susceptibility Test (AST) disc range. The Disc Dispenser applies discs equidistance onto agar plates with a simple, one-handed operation. The integrity of the discs when cartridges are opened is ensured by a low moisture environment provided by a desiccant located inside the base of the device.

The device is for professional use only, is not automated, and nor is it a companion diagnostic.

MATERIALS PROVIDED

Product name	Product Code
Thermo Scientific™ Oxoid™ (6 Places, 90 mm skirt) Antimicrobial Susceptibility Disc Dispenser	ST6090
Thermo Scientific™ Oxoid™ (6 Places, 100 mm skirt) Antimicrobial Susceptibility Disc Dispenser	ST6100
Thermo Scientific™ Oxoid™ (8 Places, 90 mm skirt) Antimicrobial Susceptibility Disc Dispenser	ST8090
Thermo Scientific™ Oxoid™ (8 Places, 100 mm skirt) Antimicrobial Susceptibility Disc Dispenser	ST8100

MATERIALS REQUIRED BUT NOT SUPPLIED

1. 6mm Antimicrobial susceptibility test discs
2. Cleaning solution - diluted (2%) non-toxic quaternary ammonium salt disinfectant or diluted (1%) sodium hypochlorite
3. Swabs
4. Petri dishes

DISC DISPENSER USER GUIDE


The OXOID Disc Dispenser is a hand-operated device for placing Antimicrobial Susceptibility Test Discs on the surface of agar in a Petri dish. The device is available in 4 formats, combining either 6-place or 8-place presentation with an option of base ring suitable for either 90mm or 100mm Petri dish use.


The dispenser will accept OXOID Antimicrobial Susceptibility Test Disc cartridges only.

FEATURES

Easily adjustable base ring. No tools required.


Plunger-lock system indicates that one or more cartridges are empty and also prevents contamination of tamp pins.

Sealable storage container complete with desiccant.

Labels on dispenser and container for laboratory's own identification markings.


Precision positioning of discs equally spaced on 54mm circle.

Discs are firmly placed on to agar surface to allow optimum diffusion.

STORAGE

Store dispenser in its container with desiccant provided to protect from moisture. When not in use, store dispenser in its container at 2°C - 8°C. On removal from refrigerator, allow to attain room temperature before use.

The desiccant is active whilst orange in colour. When desiccant is green, discard and replace with fresh desiccant (ST9100).


CLEANING

Before cleaning, always remove all cartridges.

Only use diluted (2%) non-toxic quaternary ammonium salt disinfectant or diluted (1%) sodium hypochlorite solution. Use a swab to clean underside of dispenser and ejection holes.

Never use phenolic compounds, concentrated sodium hypochlorite or concentrated detergent solutions. These may degrade the surface of the dispenser.

OPERATING

Adjust base ring to suit agar depth.

Load cartridges: ensure that they locate into place with a positive 'click'.

Place dispenser over Petri dish.

Depress plunger in full and smooth continuous downward movement and release. Do not exert undue pressure on the plunger.


Avoid excessive pressure if dispenser appears to lock – check for empty cartridges or blockage (see troubleshooting guide).

DECONTAMINATION

Under normal circumstances no component of the dispenser should become contaminated. However, should this occur, decontaminate as follows:

Before commencing decontamination, always remove *all* cartridges.


Partially immerse dispenser to top of base ring, as shown, in a diluted (2%) non-toxic quaternary ammonium salt disinfectant overnight.

Rinse with distilled water.

Fully dry dispenser at ambient temperature before use.

TROUBLESHOOTING GUIDE

If the dispenser fails to dispense, please follow the next 7 steps in the order presented:

1. Do not depress plunger again.

2. Check for empty cartridges (replace as necessary).


3. Check that the cartridge is fitted correctly.

4. Remove all cartridges.

5. Check for discs trapped in dispenser. Avoid use of metal implements.


6. Check for faulty discs remaining in cartridges (replace cartridge as necessary).

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

Symbol Legend

Symbol	Definition
	Catalogue number
	In Vitro Diagnostic Medical Device
	Batch Code
	Temperature limitation (storage temp.)
	Consult instructions for use or consult electronic instructions for use
	Manufacturer
	Authorized representative in the European Community/ European Union
	European Conformity Assessment
	UK Conformity Assessment
	Unique device identifier
	Made in the United Kingdom
	Importer - To indicate the entity importing the medical device into the locale. Applicable to the European Union

© 2025 Thermo Fisher Scientific Inc. All rights reserved.

All other trademarks are the property of Thermo Fisher Scientific and its subsidiaries. This information is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others. Please contact your local sales representative for all enquiries.

For technical assistance please contact your local distributor.

Revision Information

Version	Date of issue and modifications introduced
3.0	11 July 2025 Typographical error corrected

DEDICATED TO MICROBIOLOGY

Oxoid Limited, Wade Road, Basingstoke, Hants RG24 8PW, UK

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