Sf-900[™] Medium 1.3X

Description

Sf-900[™] Medium 1.3X is a complete serum-free medium (SFM) developed for the convenient and reproducible formulation of 1% agarose overlays used for plaque assays of baculovirus in *Spodoptera frugiperda* (Sf9, Sf21) cells. When mixed as directed with agarose, such as 4% Agarose Gel, this medium concentrate allows convenient formulation of an overlay containing Sf-900[™] Serum-Free Insect Cell Culture Medium of proper osmolarity and 1X nutrient complement. Sf-900[™] SFM 1.3X is supplemented with L-glutamine.

Product	Catalog no.	Amount	Storage	Shelf life*
Sf-900™ Medium 1.3X, liquid	10967-032 100 mL		2°C to 8°C; Protect from light	9 months

* Shelf life duration is determined from Date of Manufacture.

Product use

For Research Use Only. Not for use in diagnostic procedures.

Safety information

Read the Safety Data Sheets (SDSs) and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves.

Use

When mixed 3:1 with 4% Agarose Gel according to directions, a 1X Sf-900TM SFM, 1% LMP (low melting point) agarose solution is produced. Either the nature or concentration of the agarose component may be easily varied by substituting a 4X concentrate of the desired agarose in water in the formulation of the overlay.

Prepare medium

- Sf-900[™] Medium 1.3X is a complete liquid medium. Supplementation with L-glutamine is not required.
- Antibiotics are not recommended; however 2.5–5 mL/L of Penicillin-Streptomycin may be used when required.

Prepare agarose overlay

The following plaque assay protocol is for preparing 60 mL (sufficient for 5 × 6-well plates, 2 mL per well) of a typical 1% agarose overlay using Sf-900[™] Medium 1.3X and 4% Agarose Gel. Using strict aseptic technique:

- 1. Pre-warm Sf-900[™] Medium 1.3X and a sterile 100-mL bottle at 37°C.
- 2. Dispense 45 mL of pre-warmed Sf-900[™] Medium 1.3X into the sterile 100-mL bottle and return to a 37°C water bath.
- 3. Melt 4% Agarose Gel in a 70°C water bath (approximately 10 minutes) and move to the 37°C water bath.
- 4. Working quickly to prevent the agarose from gelling, dispense 15 mL of the liquefied 4% Agarose Gel into the warm Sf-900[™] Medium 1.3X, and mix by gentle swirling. Return to 37°C water bath.
- 5. After the appropriate viral inoculation incubation period, completely remove inoculum fluid from the insect cultures and gently replace with 37°C agarose overlay.
- 6. Allow overlay to harden (5–20 minutes), then place plates in a humidified environment and incubate at 28°C until plaques develop (4–7 days).

Note: The visual, immunologic, chemical and microscopic characteristics of the plaque formation process, and the methods for calculation of the inoculum titer and purification of individual plaques vary with the nature of the strain of baculovirus in use. An excellent review is presented by V. Luckow (see **References**).

Related products

Product	Catalog no.
4% Agarose Gel	18300
Sf9 Cells Adapted in Sf-900™ III SFM	12659
Sf21 Cells Adapted in Sf-900™ III SFM	12682
Sf9 Cells Adapted in Sf-900™ II SFM	11496
Sf21 Cells Adapted in Sf-900™ II SFM	11497
Penicillin-Streptomycin, liquid	15070
Sf-900™ III SFM, (1X)	12658
Sf-900™ II SFM, (1X)	10902
Grace's Insect Medium, Supplemented (1X)	11605
Gibco® Bottle, 100 mL	10339
BaculoDirect™ N-Term Expression Kit	12562-054
BaculoDirect [™] N-Term Transfection Kit	12562-062
BaculoDirect™ C-Term Expression Kit	12562-013
BaculoDirect [™] C-Term Transfection Kit	12562-039
Bac-N-Blue™ Transfection Kit	K855-01
Bac-to-Bac® Baculovirus Expression System	10359
Bac-to-Bac® Vector Kit	10360

Explanation of symbols and warnings The symbols present on the product label are explained below:

X		LOT				REF
Temperature Limitation	Manufacturer	Batch code		Use By:		Catalog number
\triangle	i		×		STERILE A	
Caution, consult accompanying documents	Consult instructions for use		Keep away from light		Sterilized using aseptic processing techniques	

Limited product warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale found on Life Technologies' website at **www.lifetechnologies.com/termsandconditions**. If you have any questions, please contact Life Technologies at **www.lifetechnologies.com/support**.

Important licensing information

This product may be covered by one or more Limited Use Label Licenses. By use of this product, you accept the terms and conditions of all applicable Limited Use Label Licenses.

References

Luckow, Verne A. in *Recombinant DNA Technology and Applications* (ed. Ales Prokop, Rakesh K. Bajpai, Chester S.Ho; New York, McGraw Hill), 4:97-153, "Cloning and Expression of Heterologous Genes in Insect Cells with Baculovirus Vectors." ISBN: 0-07-029075-X (1991).

For additional technical information such as Safety Data Sheets (SDS), Certificates of Analysis, visit www.lifetechnologies.com/support For further assistance, email **techsupport@lifetech.com**

©2014 Thermo Fisher Scientific Inc. All rights reserved.

All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.

DISCLAIMER - LIFE TECHNOLOGIES CORPORATION AND/OR ITS AFFILIATE(S) DISCLAIM ALL WARRANTIES WITH RESPECT TO THIS DOCUMENT, EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO THOSE OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, OR NON-INFRINGEMENT. TO THE EXTENT ALLOWED BY LAW, IN NO EVENT SHALL LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) BE LIABLE, WHETHER IN CONTRACT, TORT, WARRANTY, OR UNDER ANY STATUTE OR ON ANY OTHER BASIS FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING BUT NOT LIMITED TO THE USE THEREOF.