# remel

## BactiDrop™ Indole (Kovacs')

EN

REF R21522

50

### INTENDED USE

BactiDrop Indole (Kovacs') is a reagent for use in qualitative procedures for the detection of indole production by microorganisms grown in indole broth. The device is used in a diagnostic workflow to aid clinicians in the treatment options for patients suspected of having bacterial infections.

The device is not automated, is for professional use only and is not a companion diagnostic.

### SUMMARY AND EXPLANATION

Previous tests for the detection of indole production required a xylene extraction procedure to detect small amounts of indole produced from the breakdown of tryptophan. Kovacs' utilized amyl alcohol in place of ethyl alcohol to bypass the extraction procedure involved in the Ehrlich-Boehme reagent.<sup>1</sup> Gadebusch and Gabriel reported that substitution of isoamyl alcohol also resulted in greater stability of the Kovacs' reagent.<sup>2</sup>

### PRINCIPLE

Intracellular enzymes, collectively called "tryptophanase", mediate the production of indole by hydrolytic activity against the amino acid, tryptophan. Indole combines with dimethylaminobenzaldehyde to form a red compound. The reaction occurs by a condensation process formed by an acid splitting of the protein.<sup>3</sup>

### REAGENTS (CLASSICAL FORMULA)\*

p-Dimethylaminobenzaldehyde (CAS 100-10-7) 50.0 g  
Hydrochloric Acid, Conc (CAS 7647-01-0) ..... 250.0 ml  
Amyl Alcohol (CAS 71-41-0) ..... 750.0 ml

\*Adjusted as required to meet performance standards.

### PRECAUTIONS

**DANGER!** Causes burns by all exposure routes. Flammable liquid and vapor. Breathing vapors may cause drowsiness and dizziness. May be absorbed through intact skin. Aspiration hazard if swallowed. Can enter lungs and cause damage. May cause central nervous system depression. Repeated or prolonged exposure may cause erosion of exposed teeth. Corrosive to metal.

This product is for *in vitro* diagnostic use and should be used by properly trained individuals. Precautions should be taken against the dangers of microbiological hazards by properly sterilizing specimens, containers, and media after use. Directions should be read and followed carefully. Refer to Safety Data Sheet on company website for detailed information on chemical reagents.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the

competent authority of the Member State in which the user and/or the patient is established.

In the event of malfunction, do not use the device

### STORAGE

This product is ready for use and no further preparation is necessary. Store product in its original container at 20-25°C until used. Do not freeze or overheat. Protect from light.

### PRODUCT DETERIORATION

This product should not be used if (1) there is evidence of dehydration, (2) the color has changed, (3) the expiration date has passed, or (4) there are other signs of deterioration. The expiration date applies to the product in its intact container when stored as directed.

### SPECIMEN COLLECTION, STORAGE, TRANSPORT

Specimens should be collected and handled following recommended guidelines.<sup>4,5</sup>

### MATERIALS PROVIDED

50 BactiDrop Indole (Kovacs') – glass ampule enclosed in a protective plastic sheath with a dropper dispenser tip (0.75/ampule)  
1 Ampule crusher (reusable)

### MATERIALS REQUIRED BUT NOT SUPPLIED

(1) Loop sterilization device, (2) Inoculating loop, swabs, collection containers, (3) Incubators, alternative environmental systems, (4) Supplemental media, (5) Quality control organisms, (6) Indole broth (R061172), (7) Pipettes.

### PROCEDURE

Place dropper in the assembled, reusable ampule crusher provided. Hold the dropper/crusher in an upright position and lightly tap the bottom to dislodge any bubbles that may have formed. Grasp the middle of the dropper/crusher with the thumb and forefinger, and with the tip pointing away, press gently to crush the ampule. Invert dropper and squeeze slightly to dispense in a dropwise fashion. Discard remaining portion of partially used ampule at end of workday.

### Test Procedure:

1. Dispense 5 drops of BactiDrop Indole (Kovacs') down the side of a tube of indole medium that has been inoculated with the test organism and incubated at 35°C for 24 hours. Gently shake the tube.
2. Observe for a pink color development in the form of a ring interfaced between the broth and the reagent.

### INTERPRETATION

Positive Test - A pink-red ring between the broth and reagent

Negative Test - No color development

### QUALITY CONTROL

All lot numbers of BactiDrop Indole (Kovacs') have been tested using the following quality control organisms and have been found to be acceptable. Testing with a

positive and negative control should be performed in accordance with established laboratory quality control procedures. If aberrant quality control results are noted, patient results should not be reported.

<b>CONTROL</b>	<b>RESULTS</b>
<i>Escherichia coli</i> ATCC® 25922	Positive
<i>Salmonella enterica</i> serovar Typhimurium ATCC® 14028	Negative




**LIMITATIONS**

1. Test only colonies cultured on non-glucose containing media, as glucose inhibits the production of indole.<sup>3</sup>
2. Organisms from MacConkey Agar and EMB Agar cannot be tested; indicators in these media may cause a positive reaction.<sup>3</sup>
3. Certain strains of *Proteus vulgaris*, *Providencia* and *Aeromonas* will give a false negative reaction.<sup>6</sup>
4. Verify that media used contains adequate tryptophan by testing known positive and negative control organisms.
5. The enzyme that degrades tryptophan is diffusible in agar. Weakly false positive reactions may occur if a mixed culture of indole-positive and indole-negative organisms is used. Test only pure cultures.<sup>3</sup>

**BIBLIOGRAPHY**

1. Kovacs, N.Z. 1928. Immunforsch. 55:311.
2. Gadebusch, H.H. and S. Gabriel. 1956. Am. J. Clin. Pathol. 26:1373-1375.
3. MacFaddin, J.F. 2000. Biochemical Tests for Identification of Medical Bacteria. 3<sup>rd</sup> ed. Lippincott Williams & Wilkins, Philadelphia, PA.
4. Versalovic, J., K.C. Carroll, G. Funke, J.H. Jorgensen, M.L. Landry, and D.W. Warnock. 2011. Manual of Clinical Microbiology. 10<sup>th</sup> ed. ASM Press, Washington, D.C.
5. Forbes, B.A., D.F. Sahm, and A.S. Weissfeld. 2007. Bailey and Scott's Diagnostic Microbiology. 12<sup>th</sup> ed. Mosby Elsevier, St. Louis, MO.
6. Blazevic, D.J. and G.M. Ederer. 1975. Principles of Biochemical Tests in Diagnostic Microbiology. John Wiley & Sons, New York, NY.

**Symbol Legend**

	Catalog Number
	In Vitro Diagnostic Medical Device
	Batch Code (Lot Number)
	Consult Instructions for Use (IFU)
	Temperature Limitation (Storage Temp.)
	Use By (Expiration Date)
	Contains sufficient for <n> tests
	Do not use if packaging damaged
	Authorized European Representative
	Manufacturer
	UK Conformity Assessment
	European Conformity Assessment

Remel Inc.  
12076 Santa Fe Trail Drive  
Lenexa, KS 66215, USA  
[www.thermofisher.com/microbiology](http://www.thermofisher.com/microbiology)  
(800) 255-6730  
International: (913) 888-0939

For technical information contact your local distributor.  
©2022 Thermo Fisher Scientific Inc. All rights reserved  
BactiDrop is a trademark of Thermo Fisher Scientific and its subsidiaries

ATCC and ATCC catalogue marks are a trademark of American Type Culture Collection.

CAS (Chemical Abstracts Service Registry No.)

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Remel Inc., 12076 Santa Fe Trail Drive, Lenexa, KS 66215, USA
©2022 Thermo Fisher Scientific Inc. Med ensamrätt.
Alla övriga varumärken tillhör Thermo Fisher Scientific Inc. och dess dotterbolag.