remel

BactiDrop™ Lactophenol Aniline Blue EN

REF R21526

INTENDED USE

BactiDrop Lactophenol Aniline Blue is a stain used for the microscopic detection of fungal elements. The device is used in a diagnostic workflow to aid clinicians in the treatment options for patients suspected of having fungal infections.

The device is not automated, is for professional use only and is not a companion diagnostic.

SUMMARY AND EXPLANATION

Lactophenol was first described by Strasburger-Koernicke as a satisfactory mounting fluid.¹ Linder recommended lactophenol for the examination and preservation of fungi since no plasmolysis occurred with this mounting fluid, and organisms were killed at the same time.² Lactophenol Aniline Blue is recommended by Koneman et al. for direct mounting and staining of fungi.³

PRINCIPLE

The phenol contained in this mounting fluid acts as a gross cytoplasmic poison, precipitating cellular proteins and inactivating enzyme systems. Glycerol permits the examination of the wet mount 18-24 hours after preparation by preventing drying. Aniline blue is a dye that stains the hyaline fungal structures, thereby making them more distinct.³

REAGENTS (CLASSICAL FORMULA)*

Aniline Blue (CAS 28631-66-5).....	0.5	g
Glycerol (CAS 56-81-5)	400.0	ml
Lactic Acid (CAS 50-21-5)	200.0	ml
Phenol (CAS 108-95-2)	200.0	ml
Demineralized Water (CAS 7732-18-5).....	200.0	ml

*Adjusted as required to meet performance standards.

PRECAUTIONS

DANGER! May be fatal if inhaled, absorbed through skin or swallowed. Causes digestive and respiratory tract burns. Causes eye and skin burns. Combustible liquid and vapor. May cause central nervous system depression. May cause liver and kidney damage. Air sensitive. Light sensitive. Hygroscopic (absorbs moisture through the air).

This product is for *in vitro* diagnostic use and should be used by properly trained individuals. Precautions should be taken against the dangers of microbiological hazards by properly sterilizing specimens, containers, and media after use. Directions should be read and followed carefully.

Refer to Safety Data Sheet on company website for detailed information on chemical reagents.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

In the event of malfunction, do not use the device

STORAGE

This product is ready for use and no further preparation is necessary. Store product in its original container at 20-25°C until used. Do not freeze or overheat.

PRODUCT DETERIORATION

This product should not be used if (1) there is evidence of dehydration, (2) the color has changed, (3) the expiration date has passed, or (4) there are other signs of deterioration. The expiration date applies to the product in its intact container when stored as directed.

SPECIMEN COLLECTION, STORAGE, TRANSPORT

Specimens should be collected and handled following recommended guidelines.^{4,5}

MATERIALS PROVIDED

50 BactiDrop Lactophenol Aniline Blue – glass ampule enclosed in a protective plastic sheath with a dropper dispenser tip (0.75/ampule)
1 Ampule crusher (reusable)

MATERIALS REQUIRED BUT NOT SUPPLIED

(1) Loop sterilization device, (2) Inoculating loop, swabs collection containers, (3) Incubators, alternative environmental systems, (4) Supplemental media, (5) Quality control organisms, (6) Glass slides, coverslips, (7) Teasing needle, (8) Microscope, (9) Clear tape.

PROCEDURE

Place dropper in the assembled, reusable ampule crusher provided. Hold the dropper/crusher in an upright position and lightly tap the bottom to dislodge any bubbles that may have formed. Grasp the middle of the dropper/crusher with the thumb and forefinger, and with the tip pointing away, press gently to crush the ampule. Invert dropper and squeeze slightly to dispense in a dropwise fashion. Discard remaining portion of partially used ampule at end of workday.

Test Procedure:

Dispense stain as needed for procedures involving examination of mycological elements. Refer to appropriate reference material or laboratory procedures.

QUALITY CONTROL

All lot numbers of BactiDrop Lactophenol Aniline Blue have been tested using the following quality control organisms and have been found to be acceptable. Testing of control organisms should be performed in accordance with established laboratory quality control procedures. If aberrant quality control results are noted, patient results should not be reported.

CONTROL

Trichophyton mentagrophytes
ATCC® 9533

RESULTS

Hyphae visible,
good morphology

LIMITATIONS

1. Microscopic examination of fungi utilizing lactophenol aniline blue is helpful in recognition and presumptive identification. Consult appropriate references on techniques for definitive identification.³
2. The drawback of the tease mount method is the difficulty in preserving continuity between the spores, fruiting structures, and hyphae after such rigorous treatment. This may be critical when microscopic definition is required for differential identification. In these instances, the use of a slide culture is recommended.³
3. Molds may not sporulate or may do so poorly on media containing antimicrobial agents.⁶
4. Cotton swabs should not be used in slide preparations as the cotton strands may resemble hyphae.⁶

BIBLIOGRAPHY

1. Strasburger-Koernicke. 1913. Botanisches Praktikum. 5:767.
2. Linder, D.H. 1929. Science. 70:430.
3. Koneman, E.W., S. D. Allen, W.M. Janda, P.C. Schreckenberger, and W.C. Winn, Jr. 1997. Color Atlas and Textbook of Diagnostic Microbiology. 5th ed. Lippincott Williams & Wilkins, Philadelphia, PA.
4. Garcia, L.S. 2010. Clinical Microbiology Procedures Handbook. 3rd ed. ASM Press, Washington, D.C.
5. Versalovic, J., K.C. Carroll, G. Funke, J.H. Jorgensen, M.L. Landry, and D.W. Warnock. 2011. Manual of Clinical Microbiology. 10th ed. ASM Press, Washington, D.C.
6. Larone, D.H. 2011. Medically Important Fungi, A Guide to Identification. 5th ed. ASM Press, Washington, D.C.

Symbol Legend

	Catalog Number
	In Vitro Diagnostic Medical Device
	Batch Code (Lot Number)
	Consult Instructions for Use (IFU)

	Temperature Limitation (Storage Temp.)
	Use By (Expiration Date)
	Contains sufficient for <n> tests
	Do not use if packaging damaged
	Authorized European Representative
	Manufacturer
	UK Conformity Assessment
	European Conformity Assessment

Remel Inc
12076 Santa Fe Trail Drive
Lenexa, KS 66215, USA
www.thermofisher.com
(800) 255-6730
International: (913) 888-0939

For technical information contact your local distributor.

©2022 Thermo Fisher Scientific Inc. All rights reserved

BactiDrop is a trademark of Thermo Fisher Scientific and its subsidiaries

ATCC and ATCC catalogue are a trademark of American Type Culture Collection.

CAS (Chemical Abstracts Service Registry No.)

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Remel Inc., 12076 Santa Fe Trail Drive, Lenexa, KS 66215, USA
©2022 Thermo Fisher Scientific Inc. Med ensamrätt.
Alla övriga varumärken tillhör Thermo Fisher Scientific Inc. och dess dotterbolag.