remel

Pathotec™ CO

(Cytochrome Oxidase)

EN

REF R38191.....  100

1. INTENDED USE

The Pathotec CO (cytochrome oxidase), is a reagent impregnated paper strip for use in qualitative procedures for the detection of the enzyme cytochrome oxidase activity in bacteria grown on plated media. The device is used in a diagnostic workflow to aid clinicians in the treatment options for patients suspected of having bacterial infections.

The device is not automated, is for professional use only and is not a companion diagnostic.

2. SUMMARY AND EXPLANATION

The enzyme, cytochrome oxidase (CO), is produced in large amounts by some genera and species of bacteria. However, this enzyme is not produced by members of the family Enterobacteriaceae. Furthermore, typical genera and species of this family of microorganisms are not only cytochrome oxidase negative but also nitrate reductase positive. When Pathotec CO is used with a nitrate reductase test, the results of these two tests generally indicate whether an isolate belongs to the family Enterobacteriaceae (exceptions are *Stenotrophomonas maltophilia*¹ and some members of the genus *Pasteurella*).

Pathotec CO has been found to compare well with conventional procedures for detecting cytochrome oxidase activity in microorganisms and offers the following advantages: (1) Accuracy – complete agreement between Pathotec CO and conventional methods^{2,3,4,5}; (2) Speed – following primary isolation, Pathotec CO can be read 30 seconds after inoculation; (3) Flexibility – blends conveniently with current procedures and requires only a small storage area compared to conventional media; (4) Reliability – reagents which duplicate conventional procedures are applied to paper, dried, and packaged in desiccant-capped, amber vials to ensure stability.

3. PRINCIPLE

The enzyme cytochrome oxidase couples with α -naphthol and N,N-dimethyl-*p*-phenylenediamine in the presence of oxygen to form a bright blue product, indophenol blue, and water.

4. REAGENTS

Pathotec CO (Cytochrome Oxidase) Strip (25 strips/vial)

Reactive ingredients: N,N-dimethyl-*p*-phenylenediamine
 α -naphthol

5. PRECAUTIONS

This product is for *in vitro* diagnostic use and should be used by properly trained individuals. Precautions should be taken against the dangers of microbiological hazards by properly sterilizing specimens, containers, and media after use. Directions should be read and followed closely.

Please refer to the Safety Data Sheet (SDS) on company website and product labelling for information on potentially hazardous components.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

In the event of malfunction do not use device.

6. STORAGE

This product is ready for use and no further preparation is necessary. The product should be stored in its original container at 2-8°C until used. Allow product to come to room temperature before use. Protect product from light.

7. PRODUCT DETERIORATION

The product should not be used if (1) the color has changed, (2) the expiration date has passed, or (3) there are other signs of deterioration. Protect strips from moisture, heat, and light. Promptly replace the cap and store at 2-8°C.

8. SPECIMEN COLLECTION, STORAGE, AND TRANSPORTATION

Specimens should be collected and handled following recommended procedures.^{6,7}

9. MATERIALS REQUIRED BUT NOT SUPPLIED

(1) Loop sterilization device, (2) Inoculating loop, swabs, collection containers, (3) Incubators, alternative environmental systems, (4) Supplemental media, (5) Quality control organisms, (6) Forceps, (7) Platinum loop or needle, or wooden applicator stick.

10. PROCEDURE

Clinical specimens should be inoculated onto primary isolation media as soon as possible after receipt. Specimens must be carefully streaked to obtain well-defined, isolated colonies.^{1,6}

NOTE: Colony sizes of Enterobacteriaceae on primary media vary substantially. One medium-sized colony will contain a sufficient number of cells for the Pathotec CO test.

1. Place Pathotec CO Test Strip on a paper towel or in a petri dish.
2. Select a well-isolated colony from a primary isolation medium.
3. Using a platinum loop or needle, or wooden applicator stick, smear the growth of an 18-24 hour old culture onto the pink-gray reagent band of the Pathotec CO strip.
4. Allow to stand approximately 30 seconds.
5. Observe the color of the inoculated area.

11. INTERPRETATION OF THE TEST

- Positive Test - Blue color development
Negative Test - Absence of blue color

12. QUALITY CONTROL

All lot numbers of Pathotec CO have been tested using the following quality control organisms and have been found to be acceptable. Testing of control organisms should be performed in accordance with established laboratory quality control procedures. If aberrant quality control results are noted, patient results should not be reported.

CONTROL	INCUBATION	RESULTS
<i>Pseudomonas aeruginosa</i> ATCC® 27853	Aerobic, 30 seconds @ Room Temperature	Positive
<i>Alcaligenes faecalis</i> ATCC® 8750	Aerobic, 30 seconds @ Room Temperature	Positive
<i>Escherichia coli</i> ATCC® 25922	Aerobic, 30 seconds @ Room Temperature	Negative

13. PERFORMANCE CHARACTERISTICS

The performance of the Pathotec CO has been confirmed by several investigators. Complete agreement has been reported between Pathotec CO and conventional procedures.^{2,3,4,5}

14. LIMITATIONS

1. Use a fresh culture (18-24 hours old). Older cultures or cultures grown on old or deteriorated media will have decreased levels of enzyme activity and may give false negative test results. Direct use of broth cultures is not recommended.⁸
2. Select colonies for testing which have not been cultured on selective media or media containing glucose, as selective agents and fermentation may inhibit oxidase enzyme activity resulting in possible false-negative results.⁸
3. A platinum inoculating loop or needle, or a wooden applicator stick is recommended for use in oxidase testing. The presence of any trace of iron (nichrome) can catalyze an oxidase reaction, resulting in a false-positive reaction.⁸

15. BIBLIOGRAPHY

1. Edwards, P.R. and W.H. Ewing. 1972. Identification of Enterobacteriaceae. 3rd ed. Burgess Publishing Co., Minneapolis, MN.
2. Blazevic, D.J., P.C. Schreckenberger, and J.M. Matsen. 1973. Appl. Microbiol. 26:886-889.
3. Edberg, S.C., M. Novak, H. Slater, and J.M. Singer. 1975. J. Clin. Microbiol. 2:469-473.
4. Rosner, R. 1973. Appl. Microbiol. 26:890-893.
5. Smith, P.B., D.L. Rhoden, and K.M. Tomfohrde. 1975. J. Clin. Microbiol. 1:359-362.
6. Jorgensen J.H., M.A. Pfaller et al. 2015. Manual of Clinical Microbiology. 11th edition. Washington ASM Press.

7. Tille, P. Bailey and Scott’s Diagnostic Microbiology, 15th Edition.
8. MacFaddin, J.F. 2000. Biochemical Tests for Identification of Medical Bacteria. 3rd ed. Lippincott Williams & Wilkins, Philadelphia, PA.

16. PACKAGING

REF R38191, Pathotec CO 100 Strips/Pk (4 x 25/Vial)

17. SYMBOL LEGEND

	Catalogue Number
	In Vitro Diagnostic Medical Device
	Batch Code (Lot Number)
	Consult Instructions for Use (IFU)
	Temperature Limitations (Storage temp.)
	Contains sufficient for <N> tests
	Do not use if package is damaged
	Do not re-use
	Use By (Expiration Date)
	Keep away from sunlight
	Authorized representative in the European Community
	UK Conformity Assessed
	Manufacturer

Pathotec™ is a registered trademark of Remel Inc.
ATCC® is a registered trademark of American Type Culture Collection.
IFU38191, Revised March, 2022 Printed in UK
Remel Inc., 12076 Santa Fe Trail Drive, Lenexa, KS 66215, USA
www.thermofisher.com/microbiology
Tel: (800) 255-6730 • International: (913) 888-0939
www.oxid.com/IFU
Europe +800 135 79 135 • US 1 855 2360 190
CA 1 855 805 8539 • ROW +31 20 794 7071

remel Pathotec® CO (Cytochrome-Oxidase)

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Reaktive Inhaltsstoffe: N,N-Dimethyl-p-phenyldiamin α -Naphthol

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Pathotec® es una marca registrada de Remel Inc.