Potassium Tellurite 3.5%

EN

REF SR0030J

* This instructions for use (IFU) document are intended to be read in conjunction with the IFU for Hoyle Medium Base (CM0083B).

Intended Use

Potassium Tellurite (SR0030J) is a selective supplement that is intended to be used in the preparation of Hoyle Medium Base (CM0083B) supplemented with Laked Horse Blood (SR0048C) for the isolation and differentiation of *Corynebacterium diphtheriae* types from wound and upper respiratory samples, including throat and nasal swabs.

The device is for professional use only, is not automated, nor is it a companion diagnostic

Summary and Explanation

Corynebacterium species are Gram-positive rods with clubbed end appearance and exist in single or pairs. The medically relevant corynebacteria include *Corynebacterium pseudotuberculosis*, *Corynebacterium ulcerans* and the most toxigenic form *Corynebacterium diphtheriae*, with four biovars: *intermedius*, *gravis*, *mitis* and *belfanti*¹. Diphtheria is an illness caused by the toxigenic strains of *C. diphtheriae* and less often by *C. ulcerans*². The only known reservoir of *C. diphtheriae* is humans³.

The clinical manifestation of diphtheria can be described as respiratory diphtheria which occur as invasive or noninvasive in the respiratory tract or cutaneous diphtheria associated with skin lesions. The mode of transmission of disease can occur either through droplets from symptomatic or asymptomatic patients or from contaminated milk or dairy products. The onset of symptoms starts with a mild fever and sore throat in two to four days. During the infection exotoxin induces the formation of a white glossy necrotic coagulum in the walls of the respiratory tract which then turns into grey to black patches in appearance leading to suffocation and aspiration of membrane which ultimately leads to death in patients of all age groups⁴. However, with a timely administration of antitoxin and antimicrobial therapy, the illness can be treated completely⁵.

Principle of Method

Potassium Tellurite acts as a selective agent and has inhibitory activity against most gram-positive and gram-negative bacteria except *Corynebacterium* species. *C. diphtheriae* reduces potassium tellurite to tellurium and produce grey black-coloured colonies

Typical Formula

Potassium Tellurite 2mL

Materials Provided

SR0030J: 10x 2ml vials. Yield will depend upon formulation.

Materials Required but Not Supplied

- Hoyle Medium Base product code CM0083B
- Laked Horse Blood SR0048C

Storage

- Store product in its original packaging between 2°C and 8°C.
- Keep container tightly closed.
- The product may be used until the expiry date stated on the label.
- Protect from moisture.
- Store away from light.
- Allow reconstituted product to equilibrate to room temperature before use.

Warnings and Precautions

- Each vial is for single use. Do not re-use.
- For in vitro diagnostic use only.
- For professional use only.
- Inspect the product packaging before first use.
- Do not use the product if there is any visible damage to the packaging or vials.
- Do not use the product beyond the stated expiry date.
- Do not use the device if signs of contamination are present.
- Do not use the device if the colour has changed or there are other signs of deterioration.
- It is the responsibility of each laboratory to manage waste produced according to their nature and degree of hazard and to have them treated or disposed of in accordance with any federal, state and local applicable regulations. Directions should be read and followed carefully. This includes the disposal of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products.

Refer to the Safety Data Sheet (SDS) for safe handling and disposal of the product (www.thermofisher.com).

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

Specimen Collection, Handling and Storage

Specimen should be collected and handled following local recommended guidelines, such as the UK Standards for Microbiology Investigations (UK SMI) ID 2 and Q 5.

Procedure

- Suspend 4g of Hoyle Medium Base (CM0083B) in 100ml of distilled water and bring to the boil to dissolve completely.
- Sterilise by autoclaving at 121°C for 15 minutes.
- Cool to 50°C.
- Add 5ml of Laked Horse Blood (SR0048C) and 1ml Potassium Tellurite (SR0030J).
- Mix well and pour into sterile containers

Refer to the IFU for CM0083B. IFUs can be found at www.thermofisher.com

Bibliography


1. Public Health England. 2023 'Identification of Corynebacterium species' ID 2i5
2. Mokrousov, Igor. "Corynebacterium diphtheriae: genome diversity, population structure and genotyping perspectives." Infection, Genetics and Evolution 9, no. 1 (2009): 1-15.
3. Rajamani Sekar, S. K., B. Veeraraghavan, S. Anandan, N. K. Devanga Ragupathi, L. Sangal, and S. Joshi. "Strengthening the laboratory diagnosis of pathogenic Corynebacterium species in the Vaccine era." Letters in applied microbiology 65, no. 5 (2017): 354-365.
4. Zakikhany, Katherina, and Androulla Efstratiou. "Diphtheria in Europe: current problems and new challenges." Future microbiology 7, no. 5 (2012): 595-607.
5. Sharma, Naresh Chand, Androulla Efstratiou, Igor Mokrousov, Ankur Mutreja, Bhabatosh Das, and Thandavarayan Ramamurthy. "Diphtheria (Primer)." Nature Reviews: Disease Primers (2019).

Symbol Legend

Symbol	Definition
	Catalogue number
	In Vitro Diagnostic Medical Device
	Batch code
	Temperature limit
	Use-by date
	Keep away from sunlight
	Do not re-use
	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests
	Do not use if packaging damaged and consult instructions for use
	Manufacturer
	Authorized representative in the European Community/ European Union

	European Conformity Assessment
	UK Conformity Assessment
	Unique device identifier
	Importer - To indicate the entity importing the medical device into the locale. Applicable to the European Union
	Made in the United Kingdom

©2025 Thermo Fisher Scientific Inc. All rights reserved. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

For technical assistance please contact your local distributor.

Revision information

Version	Modifications Introduced
3.0	2025-04-28 Typographical error corrected

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