Brilliant Green Agar (Modified)

EN

REF CM0329B, CM0329R, CM0329T & CM0329K

Intended Use

Brilliant Green Agar Modified (CM0329B, CM0329R, CM0329T and CM0329K) devices are selective media for the isolation of *Salmonella* species (not *Salmonella typhi* or *Salmonella paratyphi A*) from clinical samples such as faecal samples. The devices can also be used for testing food samples. Brilliant Green Agar Modified (CM0329B, CM0329R, CM0329T and CM0329K) devices are used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having *Salmonellae* infections.

The devices are for professional use only, are not automated, nor are they companion diagnostics

Summary and Explanation

Salmonella species are rod shaped motile bacteria with peritrichous flagella with 2 to 5µm length and 0.7 to 1.5µm diameter of Enterobacteriaceae family¹. *Salmonella* species is broadly classified into two major groups that are *Salmonella enterica* and *Salmonella bongori* and the former subdivided into typhoidal and paratyphoid serovar. Almost 70% of the *Salmonella* infections in the UK are caused by the non-typhoidal serovar *Salmonella enteritidis*¹. *Salmonella enterica* often causes gastroenteritis in humans and the mode of transmission can be food like undercooked meat, egg, dairy products, contaminated water, infected people and animals².

The implications of *Salmonella* infections are reported more in immunocompromised patients. In patients with HIV, bacteraemia due to non-*Salmonella* infections is common. Studies report that in immunocompromised HIV patients the *Salmonella* bacteraemia can also spread to other tissues and organs especially lungs³. Consequently, it is important to distinguish *Salmonella* species in the clinical samples as it is the third most common cause of bacterial gastroenteritis in the UK¹.

Principle of Method

Brilliant Green Agar Modified (CM0329B, CM0329R, CM0329T and CM0329K) contains peptones, Lab Lemco and yeast extract to support the growth of *Salmonella* species. The presence of lactose and sucrose within the media act as an energy source, the fermentation of which results in a drop in pH, resulting in the phenol red indicator dye changing from red to yellow.

Typical Formula

	grams per litre
'Lab-Lemco' powder	5.0
Peptone	10.0
Yeast extract	3.0
Di-sodium hydrogen phosphate	1.0
Sodium dihydrogen phosphate	0.6
Lactose	10.0
Sucrose	10.0
Phenol red	0.09
Brilliant green	0.0047
Agar	12.0

Materials Provided

CM0329B: 500g of Brilliant Green Agar Base CM0329R: 2.5Kg of Brilliant Green Agar Base CM0329T: 5Kg of Brilliant Green Agar Base CM0329K: 25Kg of Brilliant Green Agar Base

Materials Required but Not Supplied

- Inoculating loops, swabs, collection containers
- Incubators
- Quality control organisms
- Sterile containers

Storage

- Store product in its original packaging between 10°C and 30°C.
 - Keep container tightly closed.
 - The product may be used until the expiry date stated on the label.
 - Protect from moisture.
 - Store away from light.
 - Allow reconstituted product to equilibrate to room temperature before use.
- Once reconstituted, store media between 2°C and 12°C.

Warnings and Precautions

- Do not inhale. May cause allergy or asthma symptoms or difficulty breathing if inhaled.
- Causes serious eye irritation.
- May cause an allergic skin reaction.

- If on skin wash with plenty of soap and water.
- If in eyes, rinse cautiously with water for several minutes.
- Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists, seek medical advice/attention.
- If inhaled, if breathing is difficult, remove subject to fresh air and keep in a position comfortable for breathing. If experiencing respiratory symptoms, call a POISON CENTER or doctor/physician.
- For in vitro diagnostic use only.
- For professional use only.
- Inspect the product packaging before first use.
- Do not use the product if there is any visible damage to the packaging (pot or cap).
- Do not use the product beyond the stated expiry date.
- Do not use the device if signs of contamination are present.
- It is the responsibility of each laboratory to manage waste produced according to their nature and degree of hazard and to have them treated or disposed of in accordance with any federal, state and local applicable regulations. Directions should be read and followed carefully. This includes the disposal of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products.
- Ensure the lid of the container is kept tightly closed after first opening and between use to minimise moisture ingress, which may result in incorrect product performance.

Refer to the Safety Data Sheet (SDS) for safe handling and disposal of the product (www.thermofisher.com).

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

Specimen Collection, Handling and Storage

Specimen should be collected and handled following local recommended guidelines, such as the UK Standards for Microbiology Investigations (UK SMI) ID 01, ID 24, S 7 and Q 5.

Procedure

Suspend 52g in litre of distilled water. With frequent agitation, bring to the boil to dissolve completely. Cool to 50°C. Mix well and pour into sterile containers. DO NOT AUTOCLAVE. This product is subject to slight colour variations. However, this does not affect the microbiological performance.

Interpretation

Once reconstituted *Salmonella enteritidis* and *Salmonella virchow* will show as 1-4mm red colonies. *Salmonella typhimurium* and *Salmonella nottingham* will show as 0.5- 4mm red colonies. *Pseudomonas aeruginosa* will show as 0.5 – 1.5mm red colonies.

Quality Control

It is the responsibility of the user to perform Quality Control testing taking into account the intended use of the medium, and in accordance with any local applicable regulations (frequency, number of strains, incubation temperature etc.).

The performance of this medium can be verified by testing the following reference strains.

Incubation Conditions: 18h @ 37° ± 2°C aerobic

Positive Controls	
Inoculum level: 10-100 cfu Colony count is ≥ 50% of the control medium count	
<i>Salmonella enteritidis</i> ATCC® 13076™	1-4mm red colonies and medium
<i>Salmonella typhimurium</i> ATCC® 14028™	0.5-4mm red colonies and medium
<i>Salmonella virchow</i> NCTC 5742	1-4mm red colonies and medium
<i>Salmonella nottingham</i> NCTC 7832	0.5-4mm red colonies and medium
<i>Pseudomonas aeruginosa</i> ATCC® 9027™	0.5-1.5mm red colonies and medium
Negative Controls	
Inoculum level: 10 ³ -10 ⁵ cfu	
<i>Escherichia coli</i> ATCC® 11775™	No growth or ppt-1mm yellow/green colonies
<i>Proteus mirabilis</i> ATCC® 29906™	No growth or ppt-1mm pink colonies, no swarming
<i>Proteus mirabilis</i> ATCC® 12453™	No growth or ppt-1mm pink colonies, no swarming

<i>Enterococcus cloacae</i> ATCC® 13047™	No growth or 0.5-1mm yellow/green colonies and medium
Inoculum level: 10 ⁴ -10 ⁶ cfu	
<i>Escherichia coli</i> ATCC® 25922™	No growth or ppt-1mm yellow/green colonies

Limitations

It should be noted that lactose fermenting salmonellae will not produce typical red colonies. Some salmonellae such as *Salmonella typhi* and *Salmonella paratyphi* are inhibited by the selective agents in the medium and will not grow. Not all non-target organisms may be completely inhibited and *Proteus citrobacter* and *Pseudomonas* spp., when able to grow mimic enteric pathogens producing red colonies. Identifications are presumptive and should be confirmed by appropriate serological/biochemical methods.

The medium should be protected from light at all times except during inoculation and after incubation.

Performance Characteristics

Accuracy has been demonstrated through review of the QC data. Correct detection of *Salmonella* strains is confirmed by the inclusion of a well-characterised isolate in the QC processes performed as part of the manufacture of each batch of the devices, which must meet the defined acceptance criteria. The precision of Brilliant Green Agar Modified (CM0329B, CM0329R, CM0329T and CM0329K)

was demonstrated by an overall pass rate of 100% obtained for the product over nearly 2 years of testing (June 2020 – March 2022; 10 batches). This shows that the performance is reproducible.


Brilliant Green Agar Modified (CM0329B, CM0329R, CM0329T and CM0329K) are tested in-house as part of the QC process since the products were launched in 1995. For target organisms, when using 10-100 cfu inoculum of *Salmonella enteritidis* ATCC® 13076™, *Salmonella typhimurium* ATCC® 14028™, *Salmonella virchow* NCTC 5742, *Salmonella Nottingham* NCTC 7832 and *Pseudomonas aeruginosa* ATCC® 9027™ and incubating the device at 37°C for 18 hours, the user can recover organisms with colony size and morphology as listed in this document. For non-target organisms, when using 10³-10⁵ cfu inoculum level of *Escherichia coli* ATCC® 11775™, *Proteus mirabilis* ATCC® 29906™, *Proteus mirabilis* ATCC® 12453™, *Enterococcus cloacae* ATCC® 13047™ or 10⁴-10⁶ cfu inoculum of *Escherichia coli* ATCC® 25922™ and incubating the device at and incubating the device at 37 ± 2°C for 18 hours the user can recover organisms with colony size and morphology as listed in this document.

Bibliography

1. Public Health England. 2021. 'Standards for Microbiology Investigations Identification of *Salmonella* species.' Issue no: 4 | Issue date: 12.03.21. <https://www.gov.uk/government/publications/smi-id-24-identification-of-salmonella-species>
2. Ford, L., Glass, K., Veitch, M., Wardell, R., Polkinghorne, B., Dobbins, T., ... & Kirk, M. D. (2016). Increasing incidence of *Salmonella* in Australia, 2000-2013. PLoS one, 11(10), e0163989.
3. Ridha, A. G., Malbrain, M. L. N. G., Mareels, J., Verbraeken, H., & Zachee, P. (1996). Lung abscess due to nontyphoid salmonella in an immunocompromised host: case report with review of the literature. Acta Clinica Belgica, 51(3), 175-183.

Symbol Legend

Symbol	Definition
	Catalogue number
	In Vitro Diagnostic Medical Device
	Batch code
	Temperature limit
	Use-by date
	Keep away from sunlight
	Do not re-use
	Consult instructions for use or consult electronic instructions for use


	Contains sufficient for <n> tests
	Do not use if packaging damaged and consult instructions for use
	Manufacturer
	Authorized representative in the European Community/ European Union
	European Conformity Assessment
	UK Conformity Assessment
	Unique device identifier
	Importer - To indicate the entity importing the medical device into the locale. Applicable to the European Union
	Made in the United Kingdom

©2022 Thermo Fisher Scientific Inc. All rights reserved. ATCC and ATCC catalogue marks are a trademark of American Type Culture Collection.

NCTC and NCTC catalogue marks are a trademark of National Collection of Type Cultures.

All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

For technical assistance please contact your local distributor.

Revision information

Version	Modifications Introduced
2.0	2023-12-19

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Alla övriga varumärken tillhör Thermo Fisher Scientific Inc. och dess dotterbolag.
Oxoid Limited, Wade Road, Basingstoke,