Campylobacter Growth Supplement

[REF] **SR0232E**

EN

* This instructions for use (IFU) document is intended to be read in conjunction with the IFU for Columbia Blood Agar Base (CM0331B/R/T/K), Blood Agar Base No.2 (CM0271B/R/T/K), Campylobacter Agar Base (CM0689B) and Nutrient Broth No.2 (CM0067B).

Intended Use

Campylobacter Growth Supplement (SR0232E) is a supplement intended to be used in the preparation of Columbia Blood Agar Base (CM0331B/R/T/K), Blood Agar Base No.2 (CM0271B/R/T/K), Campylobacter Agar Base (CM0689B) and Nutrient Broth No.2 (CM0067B) for the enhanced growth and aerotolerance of *Campylobacter* species from faecal samples.

Campylobacter Growth Supplement (SR0232E) is used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having bacterial infections. The device is for professional use only, is not automated and nor is it a companion diagnostic.

Summary and Explanation

Campylobacter species are mobile gram-negative rods, 0.5 to 8µm in length with a curved, spiral or S shape¹. This fastidious group grow under anaerobic or microaerobic conditions (5-10% of O₂ and CO₂), with some strains requiring a hydrogen-enriched atmosphere¹. *Campylobacter* species are thermophilic, with optimum growth temperatures between 37°C and 42°C¹. The group is most commonly associated with diarrheal and sickness infections in humans, with *C. jejuni* and *C. coli* the most commonly isolated species from patients with bacterial gastroenteritis¹. The group has been isolated from a variety of clinical sample types, including blood, faeces, cerebrospinal fluid, the intestinal tract, gall bladder, urine and oral cavities¹. *Campylobacter jejuni* is the most common pathogen isolated from samples, accounting for 90% of campylobacter infections worldwide². It has been estimated that 1 in 10 people are infected by *Campylobacter* spp. each year, with 33 million healthy life years (HLY) lost³. In 2015, *Campylobacter* species officially became a nationally notifiable infection².

Principle of Method

Campylobacter Growth Supplement (SR0232E) incorporates sodium pyruvate to provide an energy source and ferrous sulphate as an iron nutrient source to support the enhanced growth and aerotolerance of *Campylobacter* spp. Sodium metabisulphite and sodium pyruvate have antioxidant properties and are used to improve cell recovery and produce a reducing atmosphere suitable for microaerophiles.

Typical Formula

Sodium pyruvate	0.125g
Sodium metabisulphite	0.125g
Ferrous sulphate	0.125g

Materials Provided

SR0232E: 10x 2 mL vials each for 500 mL of medium

Materials Required but Not Supplied

1. Columbia Blood Agar Base (dehydrated) product code CM0331B/R/T/K
2. Blood Agar Base No.2 (dehydrated) product code CM0271B/R/T/K
3. Campylobacter Agar Base product code CM0689B
4. Nutrient Broth No.2 (dehydrated) product code CM0067B
5. Inoculating loops
6. Swabs
7. Collection containers
8. Incubators
9. Quality control organisms

Storage

- Store product in its original packaging between 2°C and 8°C.
- Keep container tightly closed.
- The product may be used until the expiry date stated on the label.
- Protect from moisture.
- Store away from light.
- Allow reconstituted product to equilibrate to room temperature before use.

Warnings and Precautions

Signal Word: Danger

Hazard Statements

H318 - Causes serious eye damage

Precautionary Statements

P280 - Wear eye protection/ face protection

P305 + P351 + P338 - IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing
P310 - Immediately call a POISON CENTER or doctor/physician

Each vial is for single use. Do not re-use.

For in vitro diagnostic use only.

For professional use only.

Inspect the product packaging before first use.

Do not use the product if there is any visible damage to the packaging or vials.

Do not use the product beyond the stated expiry date.

Do not use the device if signs of contamination are present.

Do not use the device if the colour has changed or there are other signs of deterioration.

It is the responsibility of each laboratory to manage waste produced according to their nature and degree of hazard and to have them treated or disposed of in accordance with any federal, state and local applicable regulations. Directions should be read and followed carefully. This includes the disposal of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products

All waste must be emptied into the HAZARDOUS WASTE container situated in the Microbiology laboratory.

The vial and rubber bung must be rinsed into the HAZARDOUS WASTE container.

The metal lid and rinsed rubber bung should be placed in the waste bin.

The rinsed glass vial should then be placed in the glass bin.

Refer to the Safety Data Sheet (SDS) for safe handling and disposal of the product (www.thermofisher.com).

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

Specimen Collection, Handling and Storage

There are no specimen collection and/or transport materials provided with the accessory device. Specimens should be collected and handled following local recommended guidelines, such as the UK Standards for Microbiology Investigations (UK SMI) ID 23, S7 and Q5.

Procedure

1. Allow supplement to equilibrate to room temperature before use.
2. Aseptically add the vial contents to 500 mL of sterile nutrient medium cooled to 50°C, prepared from the following media.

Base	Supplement
Campylobacter Agar Base – CM0689B	Modified Preston Campylobacter Selective Supplement (SR0204E)
Columbia Blood Agar Base - CM0331B/R/T/K	Blaser-Wang Campylobacter Selective Supplement (SR0098E) OR Butzler Campylobacter Selective Supplement (SR0085E)
Blood Agar Base No.2 - CM0271B/R/T/K	Blaser-Wang Campylobacter Selective Supplement (SR0098E)
Nutrient Broth No.2 – CM0067B	Modified Preston Campylobacter Selective Supplement (SR0204E)

All mediums may be supplemented with 5-7% Lysed Horse Blood (SR0048C), 5-7% Defibrinated Horse Blood (SR0050B/C/D/E/F/R) OR 5-7% Defibrinated Sheep Blood (SR0051B/C/D/E).

3. Mix well and aseptically dispense into sterile containers.

Refer to the IFU for CM0689B, Columbia Blood Agar Base - CM0331B/R/T/K, CM0271B/R/T/K, CM0067B. IFUs can be found at www.thermofisher.com.

Bibliography

1. Public Health England. 2018. 'Identification of Campylobacter species'. UK Standards for Microbiology Investigations ID 23 (3.1). https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/685065/ID_23i3.1.pdf.
2. Centers for Disease Control and Prevention (CDC). 2019. 'Campylobacter (Campylobacteriosis): For Health Professionals'. Available at: <https://www.cdc.gov/campylobacter/technical.html>.
3. World Health Organisation (WHO). 2013. 'The global view of campylobacteriosis: report of an expert consultation'. Available at: <https://www.who.int/publications/i/item/9789241564601>

Symbol Legend

Symbol	Definition
REF	Catalogue number
IVD	In Vitro Diagnostic Medical Device
LOT	Batch code

	Temperature limit
	Use-by date
	Keep away from sunlight
	Do not re-use
	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests
	Do not use if packaging damaged and consult instructions for use
	Manufacturer
	Authorized representative in the European Community/ European Union
	European Conformity Assessment
	UK Conformity Assessment
	Unique device identifier
	Importer - To indicate the entity importing the medical device into the locale. Applicable to the European Union
	Made in the United Kingdom

©2022 Thermo Fisher Scientific Inc. All rights reserved.
All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, UK

For technical assistance please contact your local distributor.

Revision information

Version	Date of issue and modifications introduced
3.0	2024-05-30

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

Litteratur

- Public Health England. 2018. 'Identification of Campylobacter species'. UK Standards for Microbiology Investigations ID 23 (3.1). https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/685065/ID_23i3.1.pdf.
- Centers for Disease Control and Prevention (CDC). 2019. 'Campylobacter (Campylobacteriosis): For Health Professionals'. Available at: <https://www.cdc.gov/campylobacter/technical.html>.
- World Health Organisation (WHO). 2013. 'The global view of campylobacteriosis: report of an expert consultation'. Available at: <https://www.who.int/publications/i/item/9789241564601>

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- Public Health England. 2018. 'Identification of *Campylobacter* species'. UK Standards for Microbiology Investigations ID 23 (3.1). https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/685065/ID_23i3.1.pdf.
- Centers for Disease Control and Prevention (CDC). 2019. 'Campylobacter (Campylobacteriosis): For Health Professionals'. Available at: <https://www.cdc.gov/campylobacter/technical.html>.
- World Health Organisation (WHO). 2013. 'The global view of campylobacteriosis: report of an expert consultation'. Available at: <https://www.who.int/publications/i/item/9789241564601>

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Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, Reino Unido

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©2022 Thermo Fisher Scientific Inc. Med ensamrätt.
Alla övriga varumärken tillhör Thermo Fisher Scientific Inc. och dess dotterbolag.
Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, Storbritannien