## Mycoplasma Broth Base & Mycoplasma Agar Base

EN

**REF** CM0401B, CM0403B & CM0403T

### Intended Use

Mycoplasma Broth Base (CM0403B/T) and Mycoplasma Agar Base (CM0401B) are basic media, which when supplemented, supports the growth of *Mycoplasma* and *Ureaplasma* species from genital and respiratory samples.

The devices are used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having bacterial infections.

The devices are for professional use only, are not automated, nor are they a companion diagnostic.

### Summary and Explanation

*Mycoplasma* and *Ureaplasma* species are a major causative agent of human respiratory tract infections, including both upper and lower respiratory infections in all ages<sup>1</sup>. *M. pneumoniae* is a common species implicated in respiratory infections, with symptoms ranging from malaise to headaches and fever<sup>1</sup>. More severe infections can result in tracheobronchitis, and severe lung infections may require hospital care<sup>2</sup>. *Mycoplasma* species are also aetiological agents of community-acquired pneumonia, where infection is more frequent and widespread in crowded settings- hospitals, nursing homes and schools are common environments for increased transmission rates<sup>1,2</sup>. It is estimated that 2 million cases of *Mycoplasma* infection occur each year in the US, with up to 1 in 5 cases of community acquired pneumonia due to *Mycoplasma* infection alone<sup>2</sup>. More severe infections can result in extra-pulmonary manifestations, including those in central nervous systems, dermatological and cardiac manifestations<sup>1</sup>. Some infections may also result in encephalitis, haemolytic anaemia, renal dysfunction and skin disorders<sup>2</sup>. For those with weakened immune systems, asthma sufferers or those recovering from respiratory illness, these opportunistic pathogens can pose a significant threat and may require hospital treatment on infection<sup>2</sup>.

### Principle of Method

Isolation of *Mycoplasma* and *Ureaplasma* species is achieved through the inclusion of antibacterial agents and nutritional sources to support the growth of target species. Mycoplasma Agar Base (CM0401B) and Mycoplasma Broth Base (CM0403B/T) are comprised of bacterial peptones and Lab-Lemco powder to provide essential nutrition for the enhanced growth of *Mycoplasma* and *Ureaplasma* species. Mineral supplements have been incorporated to improve the growth and colony characteristics of *Mycoplasma*. Sodium chloride salt is required to maintain osmotic equilibrium within the medium to allow for optimum growth of bacteria. Agar is included in Mycoplasma Agar Base (CM0401B) as the gelling agent that allows the product to solidify in the petri dish and incorporates anti-swarming.

### Typical Formula

CM0401B	grams per litre
Bacteriological peptone	10.0
'Lab-Lemco' powder	10.0
Sodium chloride	5.0
Mineral Supplement	0.5
Agar	10.0

CM0403B/T	grams per litre
Bacteriological peptone	10.0
'Lab-Lemco' powder	10.0
Sodium chloride	5.0
Mineral Supplement	0.5

### Materials Provided

CM0401B: 500g of Mycoplasma Agar Base  
 CM0403B: 500g of Mycoplasma Broth Base  
 CM0403T: 5kg of Mycoplasma Broth Base

### Materials Required but Not Supplied

- Inoculating loops, swabs, collection containers
- Incubators
- Quality control organisms
- Petri dish
- Supplements (SR0059C)

### Storage

- Store product in its original packaging between 10°C and 30°C.
- Keep container tightly closed.
- The product may be used until the expiry date stated on the label.
- Protect from moisture.
- Store away from light.

- Allow reconstituted product to equilibrate to room temperature before use. Once reconstituted, store media between 2°C and 10°C.

### Warnings and Precautions

- For *in vitro* diagnostic use only.
- For professional use only.
- Inspect the product packaging before first use.
- Do not use the product if there is any visible damage to the packaging (pot or cap).
- Do not use the product beyond the stated expiry date.
- Do not use the device if signs of contamination are present.
- It is the responsibility of each laboratory to manage waste produced according to their nature and degree of hazard and to have them treated or disposed of in accordance with any federal, state and local applicable regulations. Directions should be read and followed carefully. This includes the disposal of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products.
- Ensure the lid of the container is kept tightly closed after first opening and between use to minimise moisture ingress, which may result in incorrect product performance.

Refer to the Safety Data Sheet (SDS) for safe handling and disposal of the product ([www.thermofisher.com](http://www.thermofisher.com)).

### Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

### Specimen Collection, Handling and Storage

Specimen should be collected and handled following local recommended guidelines, such as the UK Standards for Microbiology Investigations (UK SMI) B5, B9, B19, B28 and B57.

### Procedure

#### CM0401B

Add 35.5g to 1 litre of distilled water. Boil to dissolve the agar and distribute in 80ml volumes. Sterilise by autoclaving at 121°C for 15 minutes. Cool to approximately 50°C and add one vial of Mycoplasma Supplement (SR0059C) reconstituted as directed.

#### CM0403B/T

Add 25.5g to 1 litre of distilled water. Mix well and distribute in 80ml volumes. Sterilize by autoclaving at 121°C for 15 minutes. Cool to 50°C and add the sterile Mycoplasma supplement (SR0059C), reconstituted as directed.

### Interpretation

Once the medium is reconstituted:

Characteristic 'fried egg' colonies indicates *Mycoplasma* spp. or *Acholeplasma laidlawii*.

Examine for colonies at x 40 magnification, using obliquely transmitted light. The colonies are characteristic with the centre of the colony embedded beneath the surface, giving a 'fried egg' appearance.

Growth and colony characteristics shall be comparable to standard.

### Quality Control

It is the responsibility of the user to perform Quality Control testing taking into account the intended use of the medium, and in accordance with any local applicable regulations (frequency, number of strains, incubation temperature etc.).

The performance of this medium can be verified by testing the following reference strains.

**CM0401B:** Incubation Conditions: Up to 21 days @ 37°C under microaerophilic conditions

<b>Positive Controls</b>	
Inoculum level: 10 <sup>4</sup> – 10 <sup>6</sup> cfu	
<i>Mycoplasma hominis</i> ATCC® 23114™	Characteristic 'fried egg' colonies
<i>Mycoplasma pneumoniae</i> ATCC® 15531™	Characteristic 'fried egg' colonies
<i>Mycoplasma bovigenitalium</i> ATCC® 19852™	Characteristic 'fried egg' colonies
<i>Mycoplasma orale</i> type 1 ATCC® 23714™	Characteristic 'fried egg' colonies
<i>Acholeplasma laidlawii</i> ATCC® 23206™	Characteristic 'fried egg' colonies
<b>Negative Controls</b>	
Inoculum level: 10 <sup>4</sup> – 10 <sup>6</sup> cfu	
Incubation conditions: 5 days @ 37°C under microaerophilic conditions	
<i>Streptococcus pyogenes</i> ATCC® 19615™	No Growth
<i>Escherichia coli</i> ATCC® 25922™	Inhibited ≥ 3 Log
<i>Staphylococcus aureus</i> ATCC® 25923™	No Growth

<i>Proteus mirabilis</i> ATCC® 12453™	No Growth
--	-----------

**CM0403B/T:** Incubation Conditions: Up to 21 days @ 37°C under microaerophilic conditions

<b>Positive Controls</b> Inoculum level: 10 <sup>4</sup> – 10 <sup>6</sup> cfu	
<i>Mycoplasma hominis</i> ATCC® 23114™	Characteristic 'fried egg' colonies
<i>Mycoplasma pneumoniae</i> ATCC® 15531™	Characteristic 'fried egg' colonies
<i>Mycoplasma bovigenitalium</i> ATCC® 19852™	Characteristic 'fried egg' colonies
<i>Mycoplasma orale</i> type 1 ATCC® 23714™	Characteristic 'fried egg' colonies
<i>Acholeplasma laidlawii</i> ATCC® 23206™	Characteristic 'fried egg' colonies
<b>Negative Controls</b> Inoculum level: 10 <sup>4</sup> – 10 <sup>6</sup> cfu Incubation conditions: 5 days @ 37°C under microaerophilic conditions	
<i>Escherichia coli</i> ATCC® 25922™	Inhibited ≥ 3 Log

### Limitations

Due to the formulation and inhibitory agents included in Mycoplasma Agar Base (CM0401B), Mycoplasma Broth Base (CM0403B/T) and Mycoplasma Supplement G (SROO059C), some non-target species may grow or may only be partially inhibited. However, the devices are optimised for the growth and isolation of *Mycoplasma* and *Ureaplasma* species. Non-target growth observed, due to its colony size and morphology, would still be distinguishable from target growth.

### Performance Characteristics

Accuracy has been demonstrated through review of the QC data. Correct detection of *Mycoplasma* and *Ureaplasma* strains is confirmed by the inclusion of a well-characterised isolate in the QC processes performed as part of the manufacture of each batch of the devices. The precision of Mycoplasma Agar Base (CM0401B) was demonstrated by an overall pass rate of 100% obtained for the product over 19 months of testing (29.12.2020 – 24.07.2022; 10 batches). The precision of Mycoplasma Broth Base (CM0403B/T) was demonstrated by an overall pass rate of 100% obtained for the product over 17 months of testing (10.02.2021 – 24.07.2022; 10 batches). This shows that the performance is reproducible.


Mycoplasma Agar Base (CM0401B) and Mycoplasma Broth Base (CM0403B/T) are tested in-house as part of the QC process since the product was launched and results are available since 1997. When using 10<sup>4</sup> – 10<sup>6</sup> cfu inoculum of *Mycoplasma hominis*, *Mycoplasma pneumoniae*, *Mycoplasma bovigenitalium*, *Mycoplasma orale* type 1 and *Acholeplasma laidlawii*. And incubating the device at 37°C for 21 days in microaerophilic conditions, the user can recover organisms with the morphology as listed in this document.

### Bibliography

1. European Centre for Disease Prevention and Control. 2012. 'Epidemiological update: Mycoplasma pneumoniae infections – recent increases reported in EU countries.' <https://www.ecdc.europa.eu/en/news-events/epidemiological-update-mycoplasma-pneumoniae-infections-recent-increases-reported-eu>
2. Centers for Disease Control and Prevention. 2020. 'Mycoplasma pneumoniae Infection.' <https://www.cdc.gov/pneumonia/atypical/mycoplasma>

### Symbol Legend

Symbol	Definition
	Catalogue number
	In Vitro Diagnostic Medical Device
	Batch code
	Temperature limit
	Use-by date

	Keep away from sunlight
	Consult instructions for use or consult electronic instructions for use
	Do not use if packaging damaged and consult instructions for use
	Manufacturer
	Authorized representative in the European Community/ European Union
	European Conformity Assessment
	UK Conformity Assessment
	Unique device identifier
	Importer - To indicate the entity importing the medical device into the locale. Applicable to the European Union
	Made in the United Kingdom

©2025 Thermo Fisher Scientific Inc. All rights reserved.  
ATCC and ATCC catalogue marks are a trademark of American Type Culture Collection.  
All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Oxoid Limited, Wade Road, Basingstoke, RG24 8PW, United Kingdom

For technical assistance please contact your local distributor.

### Revision information

Version	Date of modifications introduced
2.0	2025-02-06 QC results amended for E. coli Minor formatting changes

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© 2025 Thermo Fisher Scientific Inc. Med ensamrätt.
Alla övriga varumärken tillhör Thermo Fisher Scientific Inc. och dess dotterbolag.