	https://www.thermofisher.com	
	EU	+800 135 79 135
	US	1 855 236 0910
	CA	1 855 805 8539
	ROW	+31 20 794 7071

Oxoid™ Moxifloxacin Discs 5 µg (MXF5)

REF Code CT1633B

Antimicrobial Susceptibility Test Discs

NOTE: This IFU should be read in conjunction with the AST generic instructions for use supplied with the product and available online.

Moxifloxacin Antimicrobial Susceptibility Test (AST) Discs are 6 mm paper discs which contain specific amounts of the antimicrobial agent moxifloxacin. The discs are labelled on both sides with details of the antimicrobial (MXF) and amount present (µg): MXF5 (5µg).

Discs are supplied in cartridges containing 50 discs. There are 5 cartridges per pack. Each cartridge is individually sealed together with a desiccant tablet in a foil covered transparent blister pack. Moxifloxacin AST Discs can be dispensed using an Oxoid Disc Dispenser (sold separately). Each individual disc should only be used once.

Intended Use

Moxifloxacin Antimicrobial Susceptibility Test Discs are used in the semi-quantitative agar diffusion test method for in vitro susceptibility testing. Used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having a microbial infection, these discs are intended to determine susceptibility against microorganisms for which moxifloxacin has been shown to be active both clinically and in vitro. To be used with a pure, agar grown culture.

The device is not automated, is for professional use only and is not a companion diagnostic.

The test provides information to categorise organisms as either resistant, intermediate or susceptible to the antimicrobial agent.

Further requirements for specimen collection, handling and storage of samples can be found in local procedures and guidelines. There is no specified testing population for these discs.

Published clinical breakpoints in the current version of the FDA,³ CLSI M100^{1ac} or EUCAST² breakpoint tables must be used to interpret the zone size result.

Species with published breakpoints according to current FDA literature:

Gram-Positive

- *Enterococcus faecalis*
- *Staphylococcus* spp.
- *Streptococcus pneumoniae*

Gram-Negative

- *Haemophilus influenzae* and *Haemophilus parainfluenzae*
- Enterobacteriaceae

Species with published breakpoints according to current CLSI literature:

Gram-Positive

- *Staphylococcus* spp.
- *Streptococcus pneumoniae*

Gram-Negative

- *Haemophilus influenzae* and *Haemophilus parainfluenzae*

Species with published breakpoints according to current EUCAST literature:

Gram-Positive

- *Staphylococcus* spp. (*S. aureus*)
- *Streptococcus* groups A, B, C, and G
- *Streptococcus pneumoniae*
- *Corynebacterium* spp.

Gram-Negative

- Enterobacterales
- *Haemophilus influenzae*
- *Moraxella catarrhalis*

Each disc is single use only. The pack contains sufficient test devices for multiple single-use tests.

Principle of Method

Moxifloxacin AST Discs may be used in the semi-quantitative agar diffusion test method for in vitro susceptibility testing. For full instructions relating to the generation and interpretation of the results according to CLSI^{1bc}/EUCAST^{2ab} methodology please refer to the relevant current standards. Tables showing CLSI^{1bc}/EUCAST² compound/concentrations can be found in their documentation referenced below. Pure cultures of clinical isolates are inoculated onto the test medium and the AST disc placed on the surface. The antibiotic within the disc diffuses through the agar to form a gradient. After incubation, the zones of inhibition around the discs are measured and compared against recognised zone diameter ranges for the specific antimicrobial agent(s)/organism(s) under test.

Metrological Traceability of Calibrator and Control Material Values

Metrological traceability of values assigned to calibrators and control materials intended to establish or verify trueness of a method for Moxifloxacin AST Discs is based on internationally recognized procedures and standards.

For recommended concentrations, the zone limits are in accordance with the current performance standards for antimicrobial disc susceptibility tests as detailed by CLSI^{1bc} and/or EUCAST².

Materials Provided

Moxifloxacin AST Discs consist of 6 mm diameter paper discs impregnated with a specific amount of antimicrobial agent. The discs are marked on both sides to indicate the agent and amount. Moxifloxacin AST Discs are supplied in cartridges of 50 discs. There are 5 cartridges in each pack. Cartridges are individually packed in a foil-sealed blister pack with a desiccant.

See Table 1 below for a description of components associated with the device. For a description of active reagents that influence the result of the device, please refer to table 2.

Table 1. Materials Provided with CT1633B	
Component Description	Material Description
Cartridge with spring, cap, and plunger (x5)	Assembly components and plastic cartridge containing 50x AST discs.
Desiccant tablet (x5)	Pale beige to brown, small lozenge shaped tablets. 1 supplied with each cartridge.
Foil	Foil individually sealing each cartridge with its desiccant.
Susceptibility Test Discs (x250)	Absorbent paper individual discs. 6mm. 50 in each cartridge. 5 cartridges per pack.

Table 2. Description of Moxifloxacin AST Disc Reagents	
Reagent	Description of Function
Moxifloxacin	Slightly yellowish to yellow powder that is soluble in water. Moxifloxacin is a fluoroquinolone antimicrobial agent that inhibits chromosomal DNA partitioning during bacterial cell division. Moxifloxacin achieves this by inhibiting topoisomerase II and IV.

The concentration of antibiotic on the AST disc is analysed for every batch and is controlled using internal and external specifications (e.g. FDA³). The actual concentration is detailed on the Certificate of Analysis.

Moxifloxacin AST Discs can be dispensed using an Oxoid Disc Dispenser, which is not included with the device.

Shelf-life and Storage Conditions

Unopened cartridges of Moxifloxacin AST discs have a shelf-life of 36 months if stored under the recommended conditions. Unopened cartridges must be stored at -20°C to 8°C until required.

Once opened, cartridges should be stored within the dispenser in the container provided (with an unsaturated (orange) desiccant), or other suitable opaque airtight container with a desiccant to protect the discs from moisture. Dispensers should be stored within the container at 2° to 8°C and be allowed to come to room temperature before opening to prevent the formation of condensation. Once opened from their desiccant-containing packaging, discs should be used within 7 days and only if stored as described in this IFU.

Analytical Performance Characteristics

Table 3. Raw data of zone sizes, in accordance with CLSI methodology, taken from the 3 most recent batches of Moxifloxacin (5 µg) MXF5 Antimicrobial Susceptibility Test Discs (CT1633B).

CLSI Methodology										
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance) (%)
							1	2	3	
3320580	<i>Escherichia coli</i>	25922	MHA	28	35	31.5	28	28	28	-3.5 ± 0 (CV=0)
	<i>Staphylococcus aureus</i>	25923	MHA	28	35	31.5	28	28	28	-3.5 ± 0 (CV=0)
	<i>Pseudomonas aeruginosa</i>	27853	MHA	17	25	21	20	20	20	-1 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	49247	HTM	31	39	35	31	31	31	-4 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	49619	MHA + 5%SB	25	31	28	26	26	26	-2 ± 0 (CV=0)
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance) (%)
							1	2	3	
3274025	<i>Escherichia coli</i>	25922	MHA	28	35	31.5	29	29	29	-2.5 ± 0 (CV=0)
	<i>Staphylococcus aureus</i>	25923	MHA	28	35	31.5	28	28	28	-3.5 ± 0 (CV=0)
	<i>Pseudomonas aeruginosa</i>	27853	MHA	17	25	21	19	19	19	-2 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	49247	HTM	31	39	35	34	34	34	-1 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	49619	MHA + 5%SB	25	31	28	27	27	27	-1 ± 0 (CV=0)
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance) (%)
							1	2	3	
3232846	<i>Escherichia coli</i>	25922	MHA	28	35	31.5	29	29	29	-2.5 ± 0 (CV=0)
	<i>Staphylococcus aureus</i>	25923	MHA	28	35	31.5	28	28	28	-3.5 ± 0 (CV=0)

	<i>Pseudomonas aeruginosa</i>	27853	MHA	17	25	21	24	24	25	3.33 ± 0.58 (CV=0.17)
	<i>Haemophilus influenzae</i>	49247	HTM	31	39	35	33	35	35	-0.67 ± 1.15 (CV=1.73)
	<i>Streptococcus pneumoniae</i>	49619	MHA + 5%SB	25	31	28	29	29	29	1 ± 0 (CV=0)

Table 4. Raw data of zone sizes, in accordance with EUCAST methodology, taken from the 3 most recent batches of Moxifloxacin (5 µg) MXF5 Antimicrobial Susceptibility Test Discs (CT1633B).

EUCAST Methodology										
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance) (%)
							1	2	3	
3320580	<i>Escherichia coli</i>	25922	MHA	28	35	31.5	28	28	28	-3.5 ± 0 (CV=0)
	<i>Staphylococcus aureus</i>	25913	MHA	25	31	28	30	30	28	1.33 ± 1.15 (CV=0.87)
	<i>Streptococcus pneumoniae</i>	49619	MHA + 5%SB	24	30	27	25	25	25	-2 ± 0 (CV=0)
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance) (%)
							1	2	3	
3274025	<i>Escherichia coli</i>	25922	MHA	28	35	31.5	29	29	29	-2.5 ± 0 (CV=0)
	<i>Staphylococcus aureus</i>	25913	MHA	25	31	28	28	29	29	0.67 ± 0.58 (CV=0.87)
	<i>Streptococcus pneumoniae</i>	49619	MHA + 5%SB	24	30	27	27	27	27	0 ± 0 (CV=0)
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance) (%)
							1	2	3	
3232846	<i>Escherichia coli</i>	25922	MHA	28	35	31.5	29	29	29	-2.5 ± 0 (CV=0)
	<i>Staphylococcus aureus</i>	25913	MHA	25	31	28	27	28	28	-0.33 ± 0.58 (CV=1.73)

	<i>Streptococcus pneumoniae</i>	49619	MHA + 5%SB	24	30	27	29	29	29	2 ± 0 (CV=0)
--	---------------------------------	-------	------------	----	----	----	----	----	----	--------------

Clinical Performance Characteristics

Idelevich et al., 2016 aimed to test the efficacy of the ADAGIO automated disc diffusion reading system (Bio-Rad) in reading clinical isolates.⁴ Antimicrobial susceptibility testing was performed according to EUCAST guidelines. A variety of Oxoid antimicrobial susceptibility test discs were utilised including Oxoid Moxifloxacin AST Discs. The tested clinical isolates were categorised as susceptible, intermediate, or resistant by two investigators (which were blinded to the results and the automated reading) or the ADAGIO system. The results were classed as a very major error (a false susceptible results by the automated system when the isolates are determined as resistant by the standard method), major error (a false resistant results measured by the automated system determined as susceptible by the standard method), minor error (a false intermediate result measured by the automated system), and categorical agreement. The study tested Oxoid Moxifloxacin AST Discs against 29 *S. agalactiae*, 14 *S. dysgalactiae*, 25 *S. pyogenes*, 28 *S. pneumoniae*, and 7 *M. catarrhalis*. All the isolates were susceptible to moxifloxacin and were in complete categorical agreement between the standard method and the automated method. This study determined that the automated ADAGIO system allowed for the automated reading of Oxoid AST discs in fastidious bacteria with good categorical agreement with the manual reading. The study did not report inconsistent or out-of-specification results for Oxoid Moxifloxacin AST Discs and, therefore, demonstrated that the Oxoid device is an effective method for the in vitro analysis of fastidious bacteria.

A cross-sectional study of 95 healthy new-borns investigated the development of aerobic conjunctival flora of neonates who were not given prophylactic antibiotic eye drops after birth and were born via different methods (spontaneous vaginal delivery or elective caesarean section).⁵ One day after delivery a conjunctival sample was taken and plated, cultures were incubated, and the antibiotic susceptibility profile evaluated by antimicrobial susceptibility discs including Oxoid Moxifloxacin AST Discs. *S. aureus* ATCC 25923 and *E. coli* ATCC 25922 were used as the QC strains for the AST discs. The study found that *S. aureus* growth, coagulase negative staphylococcus, and *Klebsiella* was observed in 70%, 20%, and 10% of culture positive swabs in new-borns from caesarean section respectively. Comparatively, culture positive swabs were observed in 17 out of 75 new-borns delivered by spontaneous vaginal delivery, however, these differences are not statistically significant. None of the isolates obtained from new-borns delivered by caesarean section were resistant to moxifloxacin, however, the isolates obtained from new-borns delivered by spontaneous vaginal delivery demonstrated a 21.05% resistance rate. This study demonstrated that new-borns delivered by caesarean section had a higher culture positivity to *S. aureus*, however, lower resistance rates to antimicrobials. Due to the lack of inconsistent/out-of-specification results, the study showed that Oxoid Moxifloxacin AST Discs were effective in the in vitro analysis of a range of bacterial strains.

Ayaz et al., 2016, aimed to screen compounds from *Polygonum hydropiper* against pathogenic bacterial and fungal strains, and compare these compounds to current antimicrobials.⁶ The susceptibility of *S. aureus*, *P. mirabilis*, *S. typhi*, *P. aeruginosa*, *K. pneumoniae*, *E. faecalis*, and *E. coli* was determined by various Oxoid AST Discs including Moxifloxacin Discs. The study showed that all the strains were susceptible to all the tested antibiotics. This study demonstrated that Oxoid Moxifloxacin AST Discs were effective in the in vitro analysis of a range of bacterial strains due to the lack of inconsistent/out-of-specification results.

Serious Incidents


Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

References

- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Idelevich, E. A., Becker, K., Schmitz, J., Knaack, D., Peters, G., & Köck, R. (2016). Evaluation of an automated system for reading and interpreting disk diffusion antimicrobial susceptibility testing of fastidious bacteria. *PLoS ONE*, 11(7), 8–15. <https://doi.org/10.1371/journal.pone.0159183>.
- Kara, M., Kivanç, S. A., Olcaysü, O. O., Akova Budak, B., Özmen, A. T., Kivanç, M., & Hörmet Öz, H. T. (2018). The newborn conjunctival flora at the post delivery 24 hours. *Journal of Current Ophthalmology*, 30(4), 348–352. <https://doi.org/10.1016/j.joco.2018.03.002>.
- Ayaz, M., Junaid, M., Ullah, F., Sadiq, A., Ovais, M., Ahmad, W., ahmad, S., & Zeb, A. (2016). Chemical profiling, antimicrobial and insecticidal evaluations of *Polygonum hydropiper* L. *BMC Complementary and Alternative Medicine*, 16(1), 1–14. <https://doi.org/10.1186/s12906-016-1491-4>.

© 2022 Thermo Fisher Scientific Inc. All rights reserved. ATCC® is a trademark of ATCC. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. This information is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others.

Glossary of Symbols

Symbol/Label	Meaning
	Manufacturer
	In Vitro Diagnostic Medical Device
	Temperature limit
	Batch Code
	Catalog Number
	Do not re-use
	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests
	Use-by date
	Do not use if package is damaged and Consult instructions for use
	Authorized representative in the European Community/ European Union
	Unique device identifier
	USA: Caution: Federal law restricts this device to sale by or on order of a Physician
	European Conformity Mark
	UK Conformity Mark

Revision Information

Version	Date of issue and modifications introduced
1.0	2022-10-26. New document

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Literatura

- Clinical & Laboratory Standards Institute (CLSI). Proučite najnovije izdanje:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Proučite najnoviju verziju:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Idelevich, E. A., Becker, K., Schmitz, J., Knaack, D., Peters, G., & Köck, R. (2016). Evaluation of an automated system for reading and interpreting disk diffusion antimicrobial susceptibility testing of fastidious bacteria. *PLoS ONE*, 11(7), 8–15. <https://doi.org/10.1371/journal.pone.0159183>.
- Kara, M., Kivanç, S. A., Olcaysü, O. O., Akova Budak, B., Özmen, A. T., Kivanç, M., & Hörmüt Öz, H. T. (2018). The newborn conjunctival flora at the post delivery 24 hours. *Journal of Current Ophthalmology*, 30(4), 348–352. <https://doi.org/10.1016/j.joco.2018.03.002>.
- Ayaz, M., Junaid, M., Ullah, F., Sadiq, A., Ovais, M., Ahmad, W., Ahmad, S., & Zeb, A. (2016). Chemical profiling, antimicrobial and insecticidal evaluations of *Polygonum hydropiper* L. *BMC Complementary and Alternative Medicine*, 16(1), 1–14. <https://doi.org/10.1186/s12906-016-1491-4>.

© 2022 Thermo Fisher Scientific Inc. Sva su prava pridržana. Oznaka ATCC® zaštitni je znak Američke zbirke tipskih kultura (ATCC). Svi ostali zaštitni znakovi vlasništvo su društva Thermo Fisher Scientific Inc. i njegovih društava kćeri. Ove informacije nisu namijenjene poticanju upotrebe ovih proizvoda na bilo koji način koji bi mogao dovesti do povrede prava intelektualnog vlasništva drugih.

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Referencer

- Clinical & Laboratory Standards Institute (CLSI). Se den seneste udgave:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Se den seneste version:
 - EUCAST -agar-diffusionsmetode
 - Læsevejledning, EUCAST -agar-diffusionsmetode til antimikrobiel modtagelighedstestning
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Idelevich, E. A., Becker, K., Schmitz, J., Knaack, D., Peters, G., & Köck, R. (2016). Evaluation of an automated system for reading and interpreting disk diffusion antimicrobial susceptibility testing of fastidious bacteria. *PLoS ONE*, 11(7), 8–15. <https://doi.org/10.1371/journal.pone.0159183>.
- Kara, M., Kıvanç, S. A., Olcaysü, O. O., Akova Budak, B., Özmen, A. T., Kıvanç, M., & Hörmet Öz, H. T. (2018). The newborn conjunctival flora at the post delivery 24 hours. *Journal of Current Ophthalmology*, 30(4), 348–352. <https://doi.org/10.1016/j.joco.2018.03.002>.
- Ayaz, M., Junaid, M., Ullah, F., Sadiq, A., Ovais, M., Ahmad, W., ahmad, S., & Zeb, A. (2016). Chemical profiling, antimicrobial and insecticidal evaluations of *Polygonum hydropiper* L. *BMC Complementary and Alternative Medicine*, 16(1), 1–14. <https://doi.org/10.1186/s12906-016-1491-4>.

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Literatur

- Clinical & Laboratory Standards Institute (CLSI). Konsultieren Sie die neueste Ausgabe:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Konsultieren Sie die neueste Version:
 - EUCAST Disk Diffusion Method
 - Lesehilfe, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Teil 460 (2005).
- Idelevich, E. A., Becker, K., Schmitz, J., Knaack, D., Peters, G., & Köck, R. (2016). Evaluation of an automated system for reading and interpreting disk diffusion antimicrobial susceptibility testing of fastidious bacteria. *PLoS ONE*, 11 (7), 8–15. <https://doi.org/10.1371/journal.pone.0159183>.
- Kara, M., Kivanç, S. A., Olcaysü, O. O., Akova Budak, B., Özmen, A. T., Kivanç, M., & Hörmet Öz, H. T. (2018). The newborn conjunctival flora at the post delivery 24 hours. *Journal of Current Ophthalmology*, 30(4), 348–352. <https://doi.org/10.1016/j.joco.2018.03.002>.
- Ayaz, M., Junaid, M., Ullah, F., Sadiq, A., Ovais, M., Ahmad, W., Ahmad, S., & Zeb, A. (2016). Chemical profiling, antimicrobial and insecticidal evaluations of *Polygonum hydropiper* L. *BMC Complementary and Alternative Medicine*, 16(1), 1–14. <https://doi.org/10.1186/s12906-016-1491-4>.

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Bibliografia

1. Clinical & Laboratory Standards Institute (CLSI). Consultare l'ultima edizione:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. Comitato europeo sui test di suscettibilità antimicrobica (European Committee on Antimicrobial Susceptibility Testing, EUCAST). Consultare l'ultima versione:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Idelevich, E. A., Becker, K., Schmitz, J., Knaack, D., Peters, G. e Köck, R. (2016). Evaluation of an automated system for reading and interpreting disk diffusion antimicrobial susceptibility testing of fastidious bacteria. *PLoS ONE*, 11(7), 8–15. <https://doi.org/10.1371/journal.pone.0159183>.
5. Kara, M., Kivanç, S. A., Olcaysü, O. O., Akova Budak, B., Özmen, A. T., Kivanç, M. e Hörmet Öz, H. T. (2018). The newborn conjunctival flora at the post delivery 24 hours. *Journal of Current Ophthalmology*, 30(4), 348–352. <https://doi.org/10.1016/j.joco.2018.03.002>.
6. Ayaz, M., Junaid, M., Ullah, F., Sadiq, A., Ovais, M., Ahmad, W., ahmad, S. e Zeb, A. (2016). Chemical profiling, antimicrobial and insecticidal evaluations of *Polygonum hydropiper* L. *BMC Complementary and Alternative Medicine*, 16(1), 1–14. <https://doi.org/10.1186/s12906-016-1491-4>.

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Referências

1. Clinical & Laboratory Standards Institute (CLSI, Instituto de padrões clínicos e laboratoriais dos EUA). Consultar a edição mais recente:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST, Comité Europeu para o Teste à Suscetibilidade Antimicrobiana). Consultar a versão mais recente:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA, Agência dos alimentos e dos medicamentos dos EUA). CFR Title 21, Volume 5, Part 460 (2005).
4. Idelevich, E. A., Becker, K., Schmitz, J., Knaack, D., Peters, G., & Köck, R. (2016). Evaluation of an automated system for reading and interpreting disk diffusion antimicrobial susceptibility testing of fastidious bacteria. *PLoS ONE*, 11(7), 8–15. <https://doi.org/10.1371/journal.pone.0159183>.
5. Kara, M., Kivanç, S. A., Olcaysü, O. O., Akova Budak, B., Özmen, A. T., Kivanç, M., & Hörmet Öz, H. T. (2018). The newborn conjunctival flora at the post delivery 24 hours. *Journal of Current Ophthalmology*, 30(4), 348–352. <https://doi.org/10.1016/j.joco.2018.03.002>.
6. Ayaz, M., Junaid, M., Ullah, F., Sadiq, A., Ovais, M., Ahmad, W., Ahmad, S., & Zeb, A. (2016). Chemical profiling, antimicrobial and insecticidal evaluations of *Polygonum hydropiper* L. *BMC Complementary and Alternative Medicine*, 16(1), 1–14. <https://doi.org/10.1186/s12906-016-1491-4>.

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Odkazy

1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Hlava 21, Zväzok 5, Časť 460 (2005).
4. Idelevich, E. A., Becker, K., Schmitz, J., Knaack, D., Peters, G., & Köck, R. (2016). Evaluation of an automated system for reading and interpreting disk diffusion antimicrobial susceptibility testing of fastidious bacteria. *PLoS ONE*, 11(7), 8–15. <https://doi.org/10.1371/journal.pone.0159183>.
5. Kara, M., Kivanç, S. A., Olcaysü, O. O., Akova Budak, B., Özmen, A. T., Kivanç, M., & Hörmet Öz, H. T. (2018). The newborn conjunctival flora at the post delivery 24 hours. *Journal of Current Ophthalmology*, 30(4), 348–352. <https://doi.org/10.1016/j.joco.2018.03.002>.
6. Ayaz, M., Junaid, M., Ullah, F., Sadiq, A., Ovais, M., Ahmad, W., Ahmad, S., & Zeb, A. (2016). Chemical profiling, antimicrobial and insecticidal evaluations of *Polygonum hydropiper* L. *BMC Complementary and Alternative Medicine*, 16(1), 1–14. <https://doi.org/10.1186/s12906-016-1491-4>.

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Referencias

1. Clinical & Laboratory Standards Institute (CLSI). Consultar la última edición:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests.
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing.
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consultar la última versión:
 - a. EUCAST Disk Diffusion Method
 - b. «Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing».
3. Food and Drug Administration (FDA). Título 21 del CFR, volumen 5, parte 460 (2005).
4. Idelevich, E. A., Becker, K., Schmitz, J., Knaack, D., Peters, G. y Köck, R. (2016). «Evaluation of an automated system for reading and interpreting disk diffusion antimicrobial susceptibility testing of fastidious bacteria». *PLoS ONE*, 11(7), 8–15. <https://doi.org/10.1371/journal.pone.0159183>.
5. Kara, M., Kivanç, S. A., Olcaysü, O. O., Akova Budak, B., Özmen, A. T., Kivanç, M. y Hörmet Öz, H. T. (2018). «The newborn conjunctival flora at the post delivery 24 hours». *Journal of Current Ophthalmology*, 30(4), 348–352. <https://doi.org/10.1016/j.joco.2018.03.002>.
6. Ayaz, M., Junaid, M., Ullah, F., Sadiq, A., Ovais, M., Ahmad, W., Ahmad, S. y Zeb, A. (2016). «Chemical profiling, antimicrobial and insecticidal evaluations of Polygonum hydropiper L». *BMC Complementary and Alternative Medicine*, 16(1), 1–14. <https://doi.org/10.1186/s12906-016-1491-4>.

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