www.thermofisher.com

Bismuth Sulphite Agar (Modified)

[REF] CM0201B

Intended Use

Bismuth Sulphite Agar (Modified) (CM0201B) is a selective medium intended to be used for the isolation and preliminary identification of *Salmonella* Typhi and other salmonellae from faecal samples.

Bismuth Sulphite Agar (Modified) (CM0201B) is intended to be used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having salmonellosis.

The device is for professional use only, is not automated and nor is it a companion diagnostic.

Summary and Explanation

Salmonella species are rod shaped motile bacteria of the *Enterobacteriaceae* family with peritrichous flagellae, 2 to 5 µm in length and 0.7 to 1.5 µm in diameter¹. *Salmonella* Typhi causes severe disease and is sometimes fatal². In humans *Salmonella* species can cause serious health problems including gastroenteritis, typhoid fever, paratyphoid fever or even death¹. Human salmonellosis is notifiable as a causative agent and food poisoning is a notifiable disease under the Health Protection (Notification) Regulations 2010². *Salmonella* species are broadly classified into two major groups: *Salmonella enterica* and *Salmonella bongori*, and the former subdivided into serovars including Typhi and Paratyphi. Almost 70% of the salmonella infections in the UK are caused by the non-typhoidal serovar *Salmonella* Enteritidis¹.

Bismuth Sulphite Agar (Modified) (CM0201B) is a modification of the original Wilson and Blair selective medium for the isolation and preliminary identification of *Salmonella* Typhi and other salmonellae from pathological material³.

Principle of Method

The peptone present in the medium serves as the nitrogen source and glucose as a source of energy. Di-sodium phosphate maintains the pH of the medium and agar is the solidifying agent. Freshly precipitated bismuth sulphite acts together with brilliant green as a selective agent by suppressing the growth of coliforms, whilst permitting the growth of salmonellae. Sulphur compounds provide a substrate for hydrogen sulphide production, whilst the metallic salts in the medium stain the colony and surrounding medium black or brown in the presence of hydrogen sulphide.

Typical Formula

	<u>grams per litre</u>
Peptone	5.0
'Lab-Lemco' powder	5.0
Glucose	5.0
Disodium Phosphate	4.0
Ferrous Sulphate	0.3
Bismuth Sulphite Indicator	8.0
Brilliant Green	0.016
Agar	12.7

Materials Provided

500 g of dehydrated Bismuth Sulphite Agar (Modified) (CM0201B).

500 g of dehydrated Bismuth Sulphite Agar (Modified) (CM0201B) yields approximately 25 L after reconstitution.

Materials Required but Not Supplied

- Inoculating loops, swabs, collection containers.
- Incubators.
- Quality control organisms.
- Petri dishes.

Storage

- Store product in its original packaging between 10°C and 30°C.
- Keep container tightly closed.
- The product may be used until the expiry date stated on the label.
- Protect from moisture.
- Store away from light.
- Allow reconstituted product to equilibrate to room temperature before use.

Once reconstituted, store media between 2°C and 8°C for a maximum of two days.

Warnings and Precautions

- Storing the poured plates at 4°C for 3 days causes the medium to change colour to green, making it less selective with small numbers of salmonellae being recovered.
- Do not inhale. May cause allergy or asthma symptoms or difficulty breathing if inhaled.
- Causes serious eye irritation.
- May cause an allergic skin reaction.
- If on skin wash with plenty of soap and water.
- If in eyes, rinse cautiously with water for several minutes.
- Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists, seek medical advice/attention.
- If inhaled, if breathing is difficult, remove subject to fresh air and keep in a position comfortable for breathing. If experiencing respiratory symptoms, call a POISON CENTER or doctor/physician.
- For *in vitro* diagnostic use only.
- For professional use only.
- Inspect the product packaging before first use.
- Do not use the product if there is any visible damage to the packaging (pot or cap).
- Do not use the product beyond the stated expiry date.
- Do not use the product if signs of contamination are present.
- It is the responsibility of each laboratory to manage waste produced according to their nature and degree of hazard and to have them treated or disposed of in accordance with any federal, state and local applicable regulations. Directions should be read and followed carefully. This includes the disposal of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products.
- Ensure the lid of the container is kept tightly closed after first opening and between use to minimise moisture ingress, which may result in incorrect product performance.

Refer to the Safety Data Sheet (SDS) for safe handling and disposal of the product (www.thermofisher.com).

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

Specimen Collection, Handling and Storage

Specimen should be collected and handled following local recommended guidelines, such as the UK Standards for Microbiology Investigations (UK SMI) ID 24, S 7 and Q 5.

Procedure

Suspend 20 g of CM0201B in 500 mL of distilled water in a 1 litre flask. Heat gently with frequent agitation until the medium just begins to boil, then simmer for 30 seconds to dissolve the agar. Cool to 50-55°C, mix well to disperse suspension and pour thick plates (25 mL medium per plate). Allow the medium to solidify with the dish uncovered. Dry the plates before use but take care to avoid over drying. Correctly prepared plates should have a smooth, cream-like opacity with a pale straw colour. There should be no sedimentation of the indicator.

Once set, inoculate by streaking the specimen onto the medium using a standard loop. Incubate plates for 48 hours at 37°C. Plates should be examined for typical colonies at 18-24 hours and again, if necessary, at 48 hours. Visually inspect plates to assess colony growth and colour under good lighting.

All negative plates should be incubated for 48 hours.

Interpretation

Salmonella Typhi

Black 'rabbit-eye' colonies with a black zone and metallic sheen surrounding the colony after 18 hours. Uniformly black after 48 hours incubation.

Other Salmonella species

Variable colony appearance after 18 hours, they may be black, green or clear and mucoid. Uniformly black colonies are seen after 48 hours, often with widespread staining of the medium and a pronounced metallic sheen.

Other organisms, e.g. coliform bacteria, Serratia,

Proteus species

Usually inhibited but occasional strains give dull green or brown colonies with no metallic sheen or staining of the surrounding medium.

Quality Control

It is the responsibility of the user to perform Quality Control testing taking into account the intended use of the medium, and in accordance with any local applicable regulations (frequency, number of strains, incubation temperature etc.).

The performance of this medium can be verified by testing the following reference strains.

Incubation Conditions: 48 hours at 37°C

Positive Controls Inoculum level: 10 – 100 cfu (colony forming units) A satisfactory result is a Colony count $\geq$ 70% of the control medium count.	
<i>Salmonella</i> Typhi ATCC® 19430™	0.5 – 2 mm black 'rabbit eye' colonies.

<i>Salmonella typhimurium</i> ATCC® 19430™	0.25 – 2 mm black colonies with sheen.
<i>Salmonella virchow</i> NCTC 5742	0.25 – 2 mm black colonies with sheen.
<i>Salmonella abony</i> NCTC 6017	0.25 – 2 mm black colonies with sheen.
<i>Salmonella poona</i> NCTC 4840	0.25 – 2 mm black colonies with sheen.
<i>Salmonella</i> Enteritidis ATCC® 13076™	0.25 - 1.5 mm green colonies
Inoculum level: 50 – 200 cfu (colony forming units) A satisfactory result is a Colony count equal to or less than 100% of the control medium count.	
<i>Escherichia coli</i> ATCC® 25922™	0.5 - 1.5 mm green colonies.
<i>Escherichia coli</i> ATCC® 8739™	0.5 - 1.5 mm green colonies.
<i>Klebsiella pneumoniae</i> ATCC® 13883™	0.5 - 3.5 mm green colonies.
<i>Citrobacter freundii</i> ATCC® 8090™	0.5- 1.5 mm dark green colonies.
Negative Controls Inoculum level: 10 ⁴ -10 ⁶ cfu. Negative strains are inhibited. For <i>Pseudomonas aeruginosa</i> ATCC® 9027™, a satisfactory result is represented by a negative diagnostic reaction.	
<i>Staphylococcus aureus</i> ATCC® 6538™	No growth
<i>Enterococcus faecalis</i> ATCC® 29212™	No growth
<i>Pseudomonas aeruginosa</i> ATCC® 9027™	No growth to +++ 0.25 - 1.0 mm green colonies.

Limitations

Atypical colonies may appear if the medium is heavily inoculated with organic matter. Such a situation may be prevented by suspending the sample in sterile saline and using the supernatant for inoculation.

Prepared plates should not be stored for longer than two days at 2-8°C after which the brilliant green oxidises to give a green medium which may be inhibitory to some salmonellae.

Identifications are presumptive and colonies should be confirmed using appropriate methods.

Performance Characteristics

Accuracy has been demonstrated through review of the QC data. Correct detection of *Salmonella* Typhi and other salmonellae is confirmed by the inclusion of well-characterised isolates in the QC processes performed as part of the manufacture of each batch of the device, which must meet the defined acceptance criteria. The precision of Bismuth Sulphite Agar (Modified) (CM0201B) was demonstrated by an overall pass rate of 90% obtained for the product over 5 years of testing (10.04.2017 - 30.03.2022; 30 batches). This shows that the performance is reproducible.

An overall pass-rate of 90% (based on data from the most recent 30 batches) provides sufficient data on accuracy

considering the low risk of the device and it is not reflective of the product that is released for clinical use. With the normal variation in the manufacturing process, a small number of failures may occur, however, there is no risk to the end user as batches will not be released for sale if out of specification.

The device is tested in-house as part of the QC process since the product was launched in 1963, with results available from 1996. For target organisms, when using 10-100 cfu inoculum of *Salmonella* Typhi (ATCC® 14028™), *Salmonella* virchow (NCTC 5742), *Salmonella* abony (NCTC 6017), *Salmonella* poona (NCTC 4840) and *Salmonella* Enteritidis (ATCC® 13076™) and incubating the device at 37°C for 48 hours, the user can recover organisms with colony morphology that meets the defined acceptance criteria. For non-target organisms, when using 20-200 cfu inoculum of *Escherichia coli* (ATCC® 25922™), *Escherichia coli* (ATCC® 8739™), *Klebsiella pneumoniae* (ATCC® 13883™) and *Citrobacter freundii* (ATCC® 8090™) and incubating the device at 37°C for 48 hours, the user can recover organisms with colony morphology that meets the defined acceptance criteria. For negative strains, when using 10⁴ to 10⁶ cfu inoculum of *Staphylococcus aureus* (ATCC® 6538™), *Enterococcus faecalis* (ATCC® 29212™) and *Pseudomonas aeruginosa* (ATCC® 9027™) the user can observe inhibition of organisms, meeting the defined acceptance criteria. For *Pseudomonas aeruginosa* ATCC® 9027™, a satisfactory result is also represented by a negative diagnostic reaction (growth of green colonies), meeting the defined acceptance criteria.

Bibliography

- Public Health England. 2021. 'Identification of Salmonella species.' UK Standards for Microbiology Investigations. Issue no: 4. <https://www.gov.uk/government/publications/smi-id-24-identification-of-salmonella-species>
- Public Health England. 2020. 'Gastroenteritis and diarrhoea'. UK Standards for Microbiology Investigations S 7. Issue 2. <https://www.gov.uk/government/publications/smi-s-7-gastroenteritis-and-diarrhoea>
- Wilson W. J. and Blair E. M. McV (1927) *J. Hyg. Camb.* 26. 374.

Symbol Legend

Symbol	Definition
	Catalogue number
	In Vitro Diagnostic Medical Device
	Batch code
	Temperature limit
	Use-by date
	Keep away from sunlight

	Do not re-use
	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests
	Do not use if packaging damaged and consult instructions for use
	Manufacturer
	Authorized representative in the European Community/ European Union
	European Conformity Assessment
	UK Conformity Assessment
	Unique device identifier
	Made in the United Kingdom

ATCC Licensed
Derivative

©2022 Thermo Fisher Scientific Inc. All rights reserved. ATCC and ATCC catalogue marks are a trademark of American Type Culture Collection. NCTC and NCTC catalogue marks are a trademark of National Collection of Type Cultures. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Oxoid Limited, Wade Road, Basingstoke, RG24 8PW, UK

For technical assistance please contact your local distributor.

Revision information

Version	Date of modifications introduced
2.0	2023-03-30

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