[www.thermofisher.com](http://www.thermofisher.com)

## Staph Strep Selective Supplement

### REF SR0070E

This instructions for use (IFU) document is intended to be read in conjunction with the IFU for Columbia Blood Agar Base (CM0331B).

### Intended Use

Staph/Strep Selective Supplement (SR0070E) is a selective supplement that may be used in the preparation of CNA Staph/Strep Medium for the selective isolation of *Staphylococcus* and *Streptococcus* species from clinical samples (e.g. wounds, throat, genital, nose, groins etc.) or non-clinical samples.

Staph/Strep Selective Supplement (SR0070E) is used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having bacterial infections.

The device is for professional use only, is not intended for self-testing, and is not a companion diagnostic.

### Summary and Explanation

Gram-positive and Gram-negative microorganisms exist in a wide range of ecological niches, including environmental bacteria, commensals, and human pathogens<sup>1</sup>. Prominent Gram-positive human pathogens include *Streptococcus* species, *Staphylococcus* species<sup>2</sup>. The *Streptococcus* genus is comprised of many pathogenic commensals and some of the important *Streptococcus* species are *Streptococcus pneumoniae*, *Streptococcus viridans* and *Streptococcus pyogenes*. and the optimum growth temperature of these species is identified as 37°C<sup>1</sup>.

Other prominent Gram-positive human pathogens include *Staphylococcus aureus*<sup>3</sup>. *S. aureus* is an opportunistic pathogen, entering the body through cuts or abrasions where it can cause a spectrum of disease in both immunocompetent and immunocompromised hosts and can range from superficial skin infections such as cellulitis and abscesses to life-threatening illnesses such as pneumonia, toxic shock syndrome and sepsis. *S. aureus* is one of the leading causes of hospital acquired infections (HAIs), further complicated by the rapid emergence of antibiotic resistance in this pathogen<sup>3</sup>. Gram-negative pathogen that are involved in blood stream related infections (BSI) in patients receiving invasive operations are *Escherichia coli* and is also considered as the second most common bacterium that is isolated from the blood of neonates identified with sepsis<sup>1</sup>. *Neisseria* species is another human pathogen that is nutritionally fastidious and haemolytic in nature and can cause meningococcal disease in humans<sup>4</sup>. Among anaerobic pathogenic bacteria *Clostridium* belongs to Clostridiaceae family and a few strains are known to cause illness in humans<sup>5</sup>.

### Principle of Method

Staph/Strep Selective Supplement (SR0070E) is a selective supplement that may be used in the preparation of CNA

Staph/Strep Medium for the growth of *Staphylococcus* and *Streptococcus* species from clinical samples. Colistin sulphate is a polymyxin antibiotic active against Gram-negative bacteria including *Pseudomonas aeruginosa*. Nalidixic acid is also active against mainly Gram-negative bacteria including *Proteus* spp. Unrestricted growth of Gram-positive cocci (including *Staphylococcus* and *Streptococcus* species) is permitted.

### Typical Formula

|                   | mg/500ml |
|-------------------|----------|
| Nalidixic acid    | 7.5mg    |
| Colistin sulphate | 5.0mg    |

### Materials Provided

SR0070E: 10x freeze dried vials, each for supplementing 500ml of medium.

### Materials Required but Not Supplied

- CM0331B: 500g of dehydrated Columbia Blood Agar Base
- Inoculating loops, swabs, collection containers.
- Petri dishes.

### Storage

- Store product in its original packaging between 2°C and 8°C.
- Keep container tightly closed.
- The product may be used until the expiry date stated on the label.
- Protect from moisture.
- Store away from light.

- Allow reconstituted product to equilibrate to room temperature before use.

## Warnings and Precautions

### Signal Word: Danger

#### Hazard Statements

H302 - Harmful if swallowed

H317 - May cause an allergic skin reaction

H334 - May cause allergy or asthma symptoms or breathing difficulties if inhaled

#### Precautionary Statements

P301 + P312 - IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell

P280 - Wear protective gloves/protective clothing

P302 + P352 - IF ON SKIN: Wash with plenty of soap and water

P333 + P313 - If skin irritation or rash occurs: Get medical advice/attention

P285 - In case of inadequate ventilation wear respiratory protection

P342 + P311 - If experiencing respiratory symptoms: Call a POISON CENTER or doctor/physician

P304 + P340 - IF INHALED: Remove person to fresh air and keep comfortable for breathing

Each vial is for single use. Do not re-use.

For in vitro diagnostic use only.

For professional use only.

Inspect the product packaging before first use.

Do not use the product if there is any visible damage to the packaging or vials.

Do not use the product beyond the stated expiry date.

Do not use the device if signs of contamination are present.

Do not use the device if the colour has changed or there are other signs of deterioration.

It is the responsibility of each laboratory to manage waste produced according to their nature and degree of hazard and to have them treated or disposed of in accordance with any federal, state and local applicable regulations. Directions should be read and followed carefully. This includes the disposal of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products.

Refer to the Safety Data Sheet (SDS) for safe handling and disposal of the product ([www.thermofisher.com](http://www.thermofisher.com)).

## Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

**Specimen Collection, Handling and Storage** There are no specimen collection and/or transport materials provided with this device. Specimens should be collected and handled following local recommended guidelines, such as the UK Standards for Microbiology Investigations (UK SMI) ID 4, ID 7, B 5, B 9, B 11, and B 28.

## Procedure

To one vial add 5ml of 95% Ethanol and mix gently to dissolve. Aseptically add the vial contents to 500ml of Columbia Blood Agar Base (CM0331B) containing 5-7% horse blood (SR0050C) cooled to 50°C. Mix well and pour into sterile containers.

## Bibliography

1. Public Health England. 2021a. 'Identification of Streptococcus species, Enterococcus species and morphologically similar organisms. UK Standards for Microbiology Investigations. ID 4 Issue 4.. <https://www.gov.uk/uk-standards-for-microbiology-investigations-smi-quality-and-consistency-in-clinical-laboratories>
2. Public Health England 2020a. 'Group A Streptococcal Infections: Third Report on Seasonal Activity, 2019/20'. Vol. 14.
3. Public Health England 2020b. 'Identification of Staphylococcus Species, Micrococcus Species and Rothia Species'. UK Standards for Microbiology Investigations.ID 7 (4). <https://www.gov.uk/government/publications/s-mi-id-7-identification-of-staphylococcus-species-micrococcus-species-and-rothia-species>.
4. Public Health England. 2015. 'Identification of Neisseria species. ID 6 <https://www.gov.uk/government/publications/s-mi-id-6-identification-of-neisseria-species>
5. Public Health England. 2016. 'Identification of Clostridium species. ID 8. <https://www.gov.uk/government/publications/s-mi-id-8-identification-of-clostridium-species>


## Symbol Legend

| Symbol     | Definition       |
|------------|------------------|
| <b>REF</b> | Catalogue number |

|  |                                                                                                                  |
|--|------------------------------------------------------------------------------------------------------------------|
|  | In Vitro Diagnostic Medical Device                                                                               |
|  | Batch code                                                                                                       |
|  | Temperature limit                                                                                                |
|  | Use-by date                                                                                                      |
|  | Keep away from sunlight                                                                                          |
|  | Do not re-use                                                                                                    |
|  | Consult instructions for use or consult electronic instructions for use                                          |
|  | Contains sufficient for <n> tests                                                                                |
|  | Do not use if packaging damaged and consult instructions for use                                                 |
|  | Manufacturer                                                                                                     |
|  | Authorized representative in the European Community/<br>European Union                                           |
|  | European Conformity Assessment                                                                                   |
|  | UK Conformity Assessment                                                                                         |
|  | Unique device identifier                                                                                         |
|  | Importer - To indicate the entity importing the medical device into the locale. Applicable to the European Union |
|  | Made in the United Kingdom                                                                                       |

©2022 Thermo Fisher Scientific Inc. All rights reserved. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

For technical assistance please contact your local distributor.

Revision information

| Version | Modifications introduced |
|---------|--------------------------|
| 2.0     | 2024-01-25               |

[www.thermofisher.com](http://www.thermofisher.com)

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## Literaturverzeichnis

1. Public Health England. 2021a. 'Identification of Streptococcus species, Enterococcus species and morphologically similar organisms. UK Standards for Microbiology Investigations. ID 4 Issue 4.. <https://www.gov.uk/uk-standards-for-microbiology-investigations-smi-quality-and-consistency-in-clinical-laboratories>
2. Public Health England 2020a. 'Group A Streptococcal Infections: Third Report on Seasonal Activity, 2019/20'. Vol. 14.
3. Public Health England. 2020b. 'Identification of Staphylococcus Species, Micrococcus Species and Rothia Species'. UK Standards for Microbiology Investigations.ID 7 (4). <https://www.gov.uk/government/publications/s-mi-id-7-identification-of-staphylococcus-species-micrococcus-species-and-rothia-species>.
4. Public Health England. 2015. "Identification of Neisseria species. ID 6 <https://www.gov.uk/government/publications/s-mi-id-6-identification-of-neisseria-species>
5. Public Health England. 2016. "Identification of Clostridium species. ID 8. <https://www.gov.uk/government/publications/s-mi-id-8-identification-of-clostridium-species>

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- Public Health England. 2021a. 'Identification of Streptococcus species, Enterococcus species and morphologically similar organisms. UK Standards for Microbiology Investigations. ID 4 Issue 4. <https://www.gov.uk/uk-standards-for-microbiology-investigations-smi-quality-and-consistency-in-clinical-laboratories>
- Public Health England 2020a. 'Group A Streptococcal Infections: Third Report on Seasonal Activity, 2019/20'. Vol. 14.
- Public Health England 2020b. 'Identification of Staphylococcus Species, Micrococcus Species and Rothia Species'. UK Standards for Microbiology Investigations. ID 7 (4). <https://www.gov.uk/government/publications/s-mi-id-7-identification-of-staphylococcus-species-micrococcus-species-and-rothia-species>
- Public Health England. 2015. 'Identification of Neisseria species. ID 6 <https://www.gov.uk/government/publications/s-mi-id-6-identification-of-neisseria-species>
- Public Health England. 2016. 'Identification of Clostridium species. ID 8. <https://www.gov.uk/government/publications/s-mi-id-8-identification-of-clostridium-species>

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Oxoid Limited, Wade Road, Basingstoke, RG24 8PW, H.B.

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1. Public Health England. 2021a. 'Identification of Streptococcus species, Enterococcus species and morphologically similar organisms. UK Standards for Microbiology Investigations. ID 4 Issue 4. <https://www.gov.uk/uk-standards-for-microbiology-investigations-smi-quality-and-consistency-in-clinical-laboratories>
2. Public Health England 2020a. 'Group A Streptococcal Infections: Third Report on Seasonal Activity, 2019/20'. Vol. 14.
3. Public Health England 2020b. 'Identification of Staphylococcus Species, Micrococcus Species and Rothia Species'. UK Standards for Microbiology Investigations. ID 7 (4). <https://www.gov.uk/government/publications/s-mi-id-7-identification-of-staphylococcus-species-micrococcus-species-and-rothia-species>.
4. Public Health England. 2015. "Identification of Neisseria species. ID 6 <https://www.gov.uk/government/publications/s-mi-id-6-identification-of-neisseria-species>
5. Public Health England. 2016. "Identification of Clostridium species. ID 8. <https://www.gov.uk/government/publications/s-mi-id-8-identification-of-clostridium-species>

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Oxoid Limited, Wade Road, Basingstoke, RG24 8PW, Reino Unido

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### Bibliographie

1. Public Health England. 2021a. 'Identification of Streptococcus species, Enterococcus species and morphologically similar organisms. UK Standards for Microbiology Investigations. ID 4 Issue 4. <https://www.gov.uk/uk-standards-for-microbiology-investigations-smi-quality-and-consistency-in-clinical-laboratories>
2. Public Health England 2020a. 'Group A Streptococcal Infections: Third Report on Seasonal Activity, 2019/20'. Vol. 14.
3. Public Health England 2020b. 'Identification of Staphylococcus Species, Micrococcus Species and Rothia Species'. UK Standards for Microbiology Investigations. ID 7 (4). <https://www.gov.uk/government/publications/s-mi-id-7-identification-of-staphylococcus-species-micrococcus-species-and-rothia-species>.
4. Public Health England. 2015. 'Identification of Neisseria species. ID 6 <https://www.gov.uk/government/publications/s-mi-id-6-identification-of-neisseria-species>
5. Public Health England. 2016. 'Identification of Clostridium species. ID 8. <https://www.gov.uk/government/publications/s-mi-id-8-identification-of-clostridium-species>

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 Oxoid Limited, Wade Road, Basingstoke, RG24 8PW, Royaume-Uni

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### Bibliografija

1. Public Health England. 2021a. 'Identification of Streptococcus species, Enterococcus species and morphologically similar organisms. UK Standards for Microbiology Investigations. ID 4 Issue 4.. <https://www.gov.uk/uk-standards-for-microbiology-investigations-smi-quality-and-consistency-in-clinical-laboratories>
2. Public Health England 2020a. 'Group A Streptococcal Infections: Third Report on Seasonal Activity, 2019/20'. Vol. 14.
3. Public Health England 2020b. 'Identification of Staphylococcus Species, Micrococcus Species and Rothia Species'. UK Standards for Microbiology Investigations. ID 7 (4). <https://www.gov.uk/government/publications/s-mi-id-7-identification-of-staphylococcus-species-micrococcus-species-and-rothia-species>.
4. Public Health England. 2015. 'Identification of Neisseria species. ID 6 <https://www.gov.uk/government/publications/s-mi-id-6-identification-of-neisseria-species>
5. Public Health England. 2016. 'Identification of Clostridium species. ID 8. <https://www.gov.uk/government/publications/s-mi-id-8-identification-of-clostridium-species>

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1. Public Health England. 2021a. 'Identification of Streptococcus species, Enterococcus species and morphologically similar organisms. UK Standards for Microbiology Investigations. ID 4 Issue 4.. <https://www.gov.uk/uk-standards-for-microbiology-investigations-smi-quality-and-consistency-in-clinical-laboratories>
2. Public Health England 2020a. 'Group A Streptococcal Infections: Third Report on Seasonal Activity, 2019/20'. Vol. 14.
3. Public Health England 2020b. 'Identification of Staphylococcus Species, Micrococcus Species and Rothia Species'. UK Standards for Microbiology Investigations.ID 7 (4). <https://www.gov.uk/government/publications/s-mi-id-7-identification-of-staphylococcus-species-micrococcus-species-and-rothia-species>.
4. Public Health England. 2015. "Identification of Neisseria species. ID 6 <https://www.gov.uk/government/publications/s-mi-id-6-identification-of-neisseria-species>
5. Public Health England. 2016. "Identification of Clostridium species. ID 8. <https://www.gov.uk/government/publications/s-mi-id-8-identification-of-clostridium-species>

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### Bibliografia

1. Public Health England. 2021a. 'Identification of Streptococcus species, Enterococcus species and morphologically similar organisms. UK Standards for Microbiology Investigations. ID 4 wydanie 4.'. <https://www.gov.uk/uk-standards-for-microbiology-investigations-smi-quality-and-consistency-in-clinical-laboratories>
2. Public Health England 2020a. 'Group A Streptococcal Infections: Third Report on Seasonal Activity, 2019/20'. Tom 14.
3. Public Health England 2020b. 'Identification of Staphylococcus Species, Micrococcus Species and Rothia Species'. UK Standards for Microbiology Investigations. ID 7 (4). <https://www.gov.uk/government/publications/s-mi-id-7-identification-of-staphylococcus-species-micrococcus-species-and-rothia-species>.
4. Public Health England. 2015. 'Identification of Neisseria species. ID 6' <https://www.gov.uk/government/publications/s-mi-id-6-identification-of-neisseria-species>
5. Public Health England. 2016. 'Identification of Clostridium species. ID 8'. <https://www.gov.uk/government/publications/s-mi-id-8-identification-of-clostridium-species>

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### Bibliografia

1. Public Health England. 2021a. "Identification of Streptococcus species, Enterococcus species and morphologically similar organisms. UK Standards for Microbiology Investigations. ID 4 Issue 4. <https://www.gov.uk/uk-standards-for-microbiology-investigations-smi-quality-and-consistency-in-clinical-laboratories>
2. Public Health England 2020a. 'Group A Streptococcal Infections: Third Report on Seasonal Activity, 2019/20'. Vol 14.
3. Public Health England, 2020b. 'Identification of Staphylococcus Species, Micrococcus Species and Rothia Species'. UK Standards for Microbiology Investigations.ID 7 (4). <https://www.gov.uk/government/publications/s-mi-id-7-identification-of-staphylococcus-species-micrococcus-species-and-rothia-species>.
4. Public Health England. 2015. "Identification of Neisseria species. ID 6 <https://www.gov.uk/government/publications/s-mi-id-6-identification-of-neisseria-species>
5. Public Health England. 2016. "Identification of Clostridium species. ID 8. <https://www.gov.uk/government/publications/s-mi-id-8-identification-of-clostridium-species>

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### Bibliografie

1. Public Health England. 2021a. 'Identification of Streptococcus species, Enterococcus species and morphologically similar organisms. UK Standards for Microbiology Investigations. ID 4 Issue 4. <https://www.gov.uk/uk-standards-for-microbiology-investigations-smi-quality-and-consistency-in-clinical-laboratories>
2. Public Health England 2020a. 'Group A Streptococcal Infections: Third Report on Seasonal Activity, 2019/20'. Vol. 14.
3. Public Health England 2020b. 'Identification of Staphylococcus Species, Micrococcus Species and Rothia Species'. UK Standards for Microbiology Investigations.ID 7 (4). <https://www.gov.uk/government/publications/s-mi-id-7-identification-of-staphylococcus-species-micrococcus-species-and-rothia-species>.
4. Public Health England. 2015. 'Identification of Neisseria species. ID 6 <https://www.gov.uk/government/publications/s-mi-id-6-identification-of-neisseria-species>
5. Public Health England. 2016. 'Identification of Clostridium species. ID 8. <https://www.gov.uk/government/publications/s-mi-id-8-identification-of-clostridium-species>

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Oxoid Limited, Wade Road, Basingstoke, RG24 8PW, Regatul Unit

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## Referenser

1. Public Health England. 2021a. "Identification of Streptococcus species, Enterococcus species and morphologically similar organisms. UK Standards for Microbiology Investigations. ID 4 Issue 4.. <https://www.gov.uk/uk-standards-for-microbiology-investigations-smi-quality-and-consistency-in-clinical-laboratories>
2. Public Health England 2020a. "Group A Streptococcal Infections: Third Report on Seasonal Activity, 2019/20". Vol. 14.
3. Public Health England 2020b. "Identification of Staphylococcus Species, Micrococcus Species and Rothia Species". UK Standards for Microbiology Investigations. ID 7 (4). <https://www.gov.uk/government/publications/smi-id-7-identification-of-staphylococcus-species-micrococcus-species-and-rothia-species>.
4. Public Health England. 2015. "Identification of Neisseria species. ID 6 <https://www.gov.uk/government/publications/smi-id-6-identification-of-neisseria-species>
5. Public Health England. 2016. "Identification of Clostridium species. ID 8. <https://www.gov.uk/government/publications/smi-id-8-identification-of-clostridium-species>

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Oxoid Limited, Wade Road, Basingstoke, RG24 8PW, Storbritannien