	https://www.thermofisher.com	
		
	EU	+800 135 79 135
	US	1 855 236 0910
	CA	1 855 805 8539
	ROW	+31 20 794 7071

Oxoid™ Cefoperazone Discs 75 µg (CFP75)

REF CT0249B

Antimicrobial Susceptibility Test Discs

NOTE: This IFU should be read in conjunction with the AST generic instructions for use supplied with the product and available online.

Cefoperazone Antimicrobial Susceptibility Test (AST) Discs are 6 mm paper discs which contain specific amounts of the antimicrobial agent cefoperazone. The discs are labelled on both sides with details of the antimicrobial (CFP) and amount present (µg): CFP75 (75µg).

Discs are supplied in cartridges containing 50 discs. There are 5 cartridges per pack. Each cartridge is individually sealed together with a desiccant tablet in a foil covered transparent blister pack. Cefoperazone AST Discs can be dispensed using an Oxoid Disc Dispenser (sold separately). Each individual disc should only be used once.

Intended Use

Antimicrobial Susceptibility Test Discs are used in the semi-quantitative agar diffusion test method for in vitro susceptibility testing. Used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having a microbial infection, these discs are intended to determine susceptibility against microorganisms for which cefoperazone has been shown to be active both clinically and in vitro. To be used with a pure, agar grown culture.

The device is not automated, is for professional use only and is not a companion diagnostic.

The test provides information to categorise organisms as either resistant, intermediate or susceptible to the antimicrobial agent.

Further requirements for specimen collection, handling and storage of samples can be found in local procedures and guidelines. There is no specified testing population for these discs.

Published clinical breakpoints in the current version of the FDA,³ CLSI M100,^{1ac} or EUCAST² breakpoint tables must be used to interpret the zone size result.

Species with published breakpoints according to current FDA literature:

Gram-Positive

- N/A

Gram-Negative

- N/A

Species with published breakpoints according to current CLSI M100:

Gram-Positive

- N/A

Gram-Negative

- Enterobacterales

Species with published breakpoints according tables:

Gram-Positive

- N/A

Gram-Negative

- N/A

Each disc is single use only. The pack contains sufficient test devices for multiple single-use tests.

Principle of Method

Cefoperazone AST Discs may be used in the semi-quantitative agar diffusion test method for in vitro susceptibility testing. For full instructions relating to the generation and interpretation of the results according to CLSI^{1bc}/EUCAST^{2ab} methodology, refer to the relevant current standards. Tables showing CLSI^{1bc}/EUCAST² compound/concentrations can be found in their documentation referenced below. Pure cultures of clinical isolates are inoculated onto the test medium and the AST disc placed on the surface. The antibiotic within the disc diffuses through the agar to form a gradient. After incubation, the zones of inhibition around the discs are measured and compared against recognised zone diameter ranges for the specific antimicrobial agent(s)/organism(s) under test.

Metrological Traceability of Calibrator and Control Material Values

Metrological traceability of values assigned to calibrators and control materials intended to establish or verify trueness of a method for Cefoperazone AST Discs is based on internationally recognized procedures and standards.

For recommended concentrations, the zone limits are in accordance with the current performance standards for antimicrobial disc susceptibility tests as detailed by CLSI^{1bc} and/or EUCAST².

Materials Provided

Cefoperazone AST Discs consist of 6 mm diameter paper discs impregnated with a specific amount of antimicrobial agent. The discs are marked on both sides to indicate the agent and amount. Cefoperazone AST Discs are supplied in cartridges of 50 discs. There are 5 cartridges in each pack. Cartridges are individually packed in a foil-sealed blister pack with a desiccant.

See Table 1 below for a description of components associated with the device. For a description of active reagents that influence the result of the device, refer to Table 2.

Table 1. Materials Provided with CT0249B	
Component Description	Material Description
Cartridge with spring, cap, and plunger (x5)	Assembly components and plastic cartridge containing 50x AST discs.
Desiccant tablet (x5)	Pale beige to brown, small lozenge shaped tablets. 1 supplied with each cartridge.
Foil	Foil individually sealing each cartridge with its desiccant.
Susceptibility Test Discs (x250)	Absorbent paper individual discs. 6mm. 50 in each cartridge. 5 cartridges per pack.

Table 2. Description of Cefoperazone AST Disc Reagent	
Reagent	Description of Function
Cefoperazone	White solid that is soluble in water. Cefoperazone is a broad-spectrum beta-lactam antimicrobial agent that causes bacterial cell lysis through cell wall inhibition.

The concentration of antibiotic on the AST disc is analysed for every batch and is controlled using internal and external specifications (e.g. FDA³). The actual concentration is detailed on the Certificate of Analysis.

Cefoperazone AST Discs can be dispensed using an Oxoid Disc Dispenser, which is not included with the device.

Shelf-life and Storage Conditions

Unopened cartridges of Cefoperazone AST discs have a shelf-life of 36 months if stored under the recommended conditions. Unopened cartridges must be stored at -20°C to 8°C until required.

Once opened, cartridges should be stored within the dispenser in the container provided with an unsaturated (orange) desiccant or other suitable opaque airtight container with a desiccant to protect the discs from moisture. Dispensers should be stored within the container at 2° to 8°C and be allowed to come to room temperature before opening to prevent the formation of condensation. Once opened from their desiccant-containing packaging, discs should be used within 7 days and only if stored as described in this IFU.

Analytical Performance Characteristics

Table 3. Raw data of zone sizes, in accordance with CLSI methodology, taken from the 3 most recent batches of Cefoperazone (75 µg) CFP75 Antimicrobial Susceptibility Test Discs (CT0249B).

CLSI Methodology										
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance) (%)
							1	2	3	
3289635	<i>Escherichia coli</i>	25922	MHA	28	34	31	28	28	28	-3 ± 0 (CV=0)
	<i>Pseudomonas aeruginosa</i>	27853	MHA	23	29	26	25	25	25	-1 ± 0 (CV=0)
	<i>Staphylococcus aureus</i>	25923	MHA	24	33	28.5	32	32	31	3.17 ± 0.58 (CV=0.18)
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance) (%)
							1	2	3	
3232779	<i>Escherichia coli</i>	25922	MHA	28	34	31	32	32	24	-1.67 ± 4.62 (CV=-2.77)
	<i>Pseudomonas aeruginosa</i>	27853	MHA	23	29	26	26	26	26	0 ± 0 (CV=0)
	<i>Staphylococcus aureus</i>	25923	MHA	24	33	28.5	27	28	28	-0.83 ± 0.58 (CV=-0.69)
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance) (%)
							1	2	3	
3116895	<i>Escherichia coli</i>	25922	MHA	28	34	31	31	31	31	0 ± 0 (CV=0)
	<i>Pseudomonas aeruginosa</i>	27853	MHA	23	29	26	25	25	25	-1 ± 0 (CV=0)
	<i>Staphylococcus aureus</i>	25923	MHA	24	33	28.5	31	31	31	-2.5 ± 0 (CV=0)

Clinical Performance Characteristics

A 2009 study aimed to determine the emergence of extended spectrum beta-lactamases (ESBLs) in *Escherichia coli* which are therefore multi-drug resistant in various public hospitals in Malaysia.⁴ This is a retrospective study of *E. coli* isolates collected from intensive care units in 5 different public hospitals in 2004 in Malaysia. The isolates were fingerprinted, genotyped and tested for antimicrobial susceptibility by disc diffusion in accordance with CLSI guidelines – including Oxoid Cefoperazone (75 µg) AST Discs. A total of 47 *E. coli* isolates were collected and determined as genetically diverse and heterogeneous (excluding to two isolates from different patients of the same ward - likely due to nosocomial spread – however, these isolates contained different sized plasmids). Cefoperazone (75 µg) AST Discs demonstrated resistance in 30% of these *E. coli* isolates; 76.5% isolates were multi-drug resistance. A total of 40 *E. coli* were ESBL-producers by genotypic detection of ESBL-encoding genes, double-disk synergy testing and phenotypic confirmatory test. Of the ESBL-producing isolates, 87.5% harbored the *bla*_{TEM} gene. The study demonstrated a 30% resistance of *E. coli* clinical isolates out of 40/47 ESBL-producing *E. coli* using the Oxoid Cefoperazone (75 µg) AST Discs, which was further confirmed by double-disk synergy and phenotypic confirmation. This study did not report any inconsistent or out-of-specification results for the Oxoid Cefoperazone (75 µg) Discs. Therefore, this study demonstrated that the Oxoid Cefoperazone (75 µg) CFP75 AST Discs are an effective method for *in vitro* analysis of *E. coli*.

In 2004, Ellabib et al. aimed to determine the most important causative agents of neonatal bacteraemia causing treatment in the intensive care unit over a 5-year period.⁵ Patients were recruited between January 1997 and January 2001. Kirby-Bauer disc diffusion was performed in accordance with the national committee for clinical laboratory standards utilising numerous antimicrobials including Oxoid Cefoperazone (75 µg) AST Discs. Over the study period 1431 neonatal intensive care unit blood cultures were taken from 1092 cases of suspected bacteraemia. Out of these, 801 were found positive (corresponding to 648 neonates). The susceptibility of the common leading Gram-negative pathogens demonstrated 100% susceptibility of *S. marcescens*, *S. liquefaciens*, *K. pneumoniae*, and *E. cloacae*, and 93% susceptibility of *E. agglomerans*. The study did not report any inconsistent or out of specification results for Oxoid Cefoperazone (75 µg) AST Discs and, therefore, demonstrated that the device is an effective method for the *in vitro* analysis of Gram-negative bacteria causing bacteremia in neonates and for the detection of emerging resistance.

Misbah et al. aimed to identify the antimicrobial susceptibility profile of *Acinetobacter* strains collected from the University Malaya Medical Centre, Kuala Lumpur in 1987 and 1996-1998.⁶ Clinical isolates from patients admitted to the University Malaya Medical Centre were collected in 1987 (n=21) and 1996-1998 (n=88). Antimicrobial susceptibility testing was carried out with Oxoid AST discs including Cefoperazone (75 µg) AST Discs and were interpreted in accordance with National Committee for Clinical Laboratory Standards (NCCLS). The MIC of antimicrobials was determined by agar-dilution. All of the 1996-1998 isolates were resistant to Cefoperazone an increase from 95.2% in 1987. DNA fingerprinting revealed the presence of four different lineages of *Acinetobacter* spp. with two of these lineages being dominant (genotypes 1 and 4). Genotype 4 demonstrated an increased resistance profile, whereas genotype 1 was mainly sensitive, suggesting genotype 4 developed *in situ*. The study did not report any inconsistent or out of specification results for the Oxoid Cefoperazone (75 µg) CFP75 AST Discs and, therefore, demonstrated that the device is an effective method for the *in vitro* analysis of *Acinetobacter* spp. in hospitals and for the detection of emerging resistance.

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

References

- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (–)
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Lim, K.-T., Yasin, R., Yeo, C.-C., Puthucheary, S., and Thong, K.-L. (2009). Characterization of Multidrug Resistant ESBL-Producing *Escherichia coli* Isolates from Hospitals in Malaysia. *Journal of Biomedicine and Biotechnology*, 2009. <https://doi.org/10.1155/2009/165637>.
- Ellabib, M. S., Ordonez, A., Ramali, A., Walli, A., Benayad, T., and Shebrlo, H. (2004). Changing Pattern of Neonatal Bacteremia. Microbiology and Antibiotic Resistance. *Saudi Medical Journal*. 25(12), 1951–1956.
- Misbah, S., AbuBakar, S., Hassan, H., Hanifah, Y. A., & Yusof, M. Y. (2004). Antibiotic Susceptibility and REP-PCR Fingerprints of *Acinetobacter* spp. Isolated from a Hospital Ten Years Apart. *Journal of Hospital Infection*. 58(4), 254-261. <https://doi.org/10.1016/j.jhin.2004.07.007>.

© 2023 Thermo Fisher Scientific Inc. All rights reserved. ATCC® is a trademark of ATCC. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. This information is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others.

Glossary of symbols

Symbol/Label	Meaning
	Manufacturer

	In Vitro Diagnostic Medical Device
	Temperature limit
	Batch Code
	Catalog Number
	Do not re-use
	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests
	Use-by date
	Do not use if package is damaged and Consult instructions for use
	Authorized representative in the European Community/ European Union
	Unique device identifier
	USA: Caution: Federal law restricts this device to sale by or on order of a Physician
	European Conformity Mark
	UK Conformity Mark

Revision Information

Version	Date of issue and modifications introduced
1.0	2023-01-13

	https://www.thermofisher.com	
	EC	+800 135 79 135
	US	1 855 236 0910
	CA	1 855 805 8539
	ROW	+31 20 794 7071

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1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (–)
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Lim, K.-T., Yasin, R., Yeo, C.-C., Puthucheary, S., and Thong, K.-L. (2009). Characterization of Multidrug Resistant ESBL-Producing *Escherichia coli* Isolates from Hospitals in Malaysia. *Journal of Biomedicine and Biotechnology*, 2009. <https://doi.org/10.1155/2009/165637>.
5. Ellabib, M. S., Ordóñez, A., Ramali, A., Walli, A., Benayad, T., and Shebrlo, H. (2004). Changing Pattern of Neonatal Bacteremia. Microbiology and Antibiotic Resistance. *Saudi Medical Journal*. 25(12), 1951–1956.
6. Misbah, S., AbuBakar, S., Hassan, H., Hanifah, Y. A., & Yusof, M. Y. (2004). Antibiotic Susceptibility and REP-PCR Fingerprints of *Acinetobacter* spp. Isolated from a Hospital Ten Years Apart. *Journal of Hospital Infection*. 58(4), 254–261. <https://doi.org/10.1016/j.jhin.2004.07.007>.

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Reference

1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (–)
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Lim, K.-T., Yasin, R., Yeo, C.-C., Puthucheary, S., and Thong, K.-L. (2009). Characterization of Multidrug Resistant ESBL-Producing *Escherichia coli* Isolates from Hospitals in Malaysia. *Journal of Biomedicine and Biotechnology*, 2009. <https://doi.org/10.1155/2009/165637>.
5. Ellabib, M. S., Ordonez, A., Ramali, A., Walli, A., Benayad, T., and Shebrlo, H. (2004). Changing Pattern of Neonatal Bacteremia. Microbiology and Antibiotic Resistance. *Saudi Medical Journal*. 25(12), 1951–1956.
6. Misbah, S., AbuBakar, S., Hassan, H., Hanifah, Y. A., & Yusof, M. Y. (2004). Antibiotic Susceptibility and REP-PCR Fingerprints of *Acinetobacter* spp. Isolated from a Hospital Ten Years Apart. *Journal of Hospital Infection*. 58(4), 254–261. <https://doi.org/10.1016/j.jhin.2004.07.007>.

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1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (–)
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Lim, K.-T., Yasin, R., Yeo, C.-C., Puthucheary, S., and Thong, K.-L. (2009). Characterization of Multidrug Resistant ESBL-Producing *Escherichia coli* Isolates from Hospitals in Malaysia. *Journal of Biomedicine and Biotechnology*, 2009. <https://doi.org/10.1155/2009/165637>.
5. Ellabib, M. S., Ordonez, A., Ramali, A., Walli, A., Benayad, T., and Shebrlo, H. (2004). Changing Pattern of Neonatal Bacteremia. Microbiology and Antibiotic Resistance. *Saudi Medical Journal*. 25(12), 1951–1956.
6. Misbah, S., AbuBakar, S., Hassan, H., Hanifah, Y. A., & Yusof, M. Y. (2004). Antibiotic Susceptibility and REP-PCR Fingerprints of *Acinetobacter* spp. Isolated from a Hospital Ten Years Apart. *Journal of Hospital Infection*. 58(4), 254-261. <https://doi.org/10.1016/j.jhin.2004.07.007>.

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Referințe

1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (–)
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Lim, K.-T., Yasin, R., Yeo, C.-C., Puthucherry, S., and Thong, K.-L. (2009). Characterization of Multidrug Resistant ESBL-Producing *Escherichia coli* Isolates from Hospitals in Malaysia. *Journal of Biomedicine and Biotechnology*, 2009. <https://doi.org/10.1155/2009/165637>.
5. Ellabib, M. S., Ordonez, A., Ramali, A., Walli, A., Benayad, T., and Shebrlo, H. (2004). Changing Pattern of Neonatal Bacteremia. Microbiology and Antibiotic Resistance. *Saudi Medical Journal*. 25(12), 1951–1956.
6. Misbah, S., AbuBakar, S., Hassan, H., Hanifah, Y. A., & Yusof, M. Y. (2004). Antibiotic Susceptibility and REP-PCR Fingerprints of *Acinetobacter* spp. Isolated from a Hospital Ten Years Apart. *Journal of Hospital Infection*. 58(4), 254–261. <https://doi.org/10.1016/j.jhin.2004.07.007>.

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Referencie

- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (–)
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Lim, K.-T., Yasin, R., Yeo, C.-C., Puthucheary, S., and Thong, K.-L. (2009). Characterization of Multidrug Resistant ESBL-Producing *Escherichia coli* Isolates from Hospitals in Malaysia. *Journal of Biomedicine and Biotechnology*, 2009. <https://doi.org/10.1155/2009/165637>.
- Ellabib, M. S., Ordóñez, A., Ramali, A., Walli, A., Benayad, T., and Shebrlo, H. (2004). Changing Pattern of Neonatal Bacteremia. Microbiology and Antibiotic Resistance. *Saudi Medical Journal*. 25(12), 1951–1956.

6. Misbah, S., AbuBakar, S., Hassan, H., Hanifah, Y. A., & Yusof, M. Y. (2004). Antibiotic Susceptibility and REP-PCR Fingerprints of *Acinetobacter* spp. Isolated from a Hospital Ten Years Apart. *Journal of Hospital Infection*. 58(4), 254-261. <https://doi.org/10.1016/j.jhin.2004.07.007>.

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Referencias

1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (–)
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Lim, K.-T., Yasin, R., Yeo, C.-C., Puthucheary, S., and Thong, K.-L. (2009). Characterization of Multidrug Resistant ESBL-Producing *Escherichia coli* Isolates from Hospitals in Malaysia. *Journal of Biomedicine and Biotechnology*, 2009. <https://doi.org/10.1155/2009/165637>.
5. Ellabib, M. S., Ordóñez, A., Ramali, A., Walli, A., Benayad, T., and Shebrlo, H. (2004). Changing Pattern of Neonatal Bacteremia. Microbiology and Antibiotic Resistance. *Saudi Medical Journal*. 25(12), 1951–1956.
6. Misbah, S., AbuBakar, S., Hassan, H., Hanifah, Y. A., & Yusof, M. Y. (2004). Antibiotic Susceptibility and REP-PCR Fingerprints of *Acinetobacter* spp. Isolated from a Hospital Ten Years Apart. *Journal of Hospital Infection*. 58(4), 254-261. <https://doi.org/10.1016/j.jhin.2004.07.007>.

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Oxoid Limited, Wade Road, Basingstoke,
Hampshire, RG24 8PW, Reino Unido

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