Key Code TSMX5240H
www.thermofisher.com

Europe +800 135 79 135 US 1 855 236 0910
CA 1 855 805 8539 ROW +31 20 794 7071

Staphytest Plus

REF	DR0850M.....	▽ 100	EN
	DR0850B.....	▽ 500	

1. INTENDED USE

Staphytest Plus™ kit is a qualitative latex slide agglutination test for the differentiation of *Staphylococcus aureus* from other *Staphylococcus* species isolates grown on agar by the detection of clumping factor, Protein A and certain capsular polysaccharides found in methicillin-resistant *Staphylococcus aureus*. Used in a diagnostic workflow to aid clinicians in treatment options for patients suspected of having staphylococcal infection or colonization. The device is not automated, is for professional use only and is not a companion diagnostic.

2. PRINCIPLE OF THE TEST

Traditionally, differentiation between coagulase-positive and coagulase-negative staphylococci has been performed either with the tube coagulase test that detects extracellular staphylocoagulase or the slide coagulase test that detects the clumping factor (bound coagulase) present on the bacterial cell surface. Several other differentiation tests are also available including the passive haemagglutination test (Oxoid Staphylase – DR595) and the DNase test. It has been reported that approximately 97% of human strains of *Staphylococcus aureus* possess both bound coagulase and extracellular staphylocoagulase.

Protein A is found on the cell surface of about 95% of human strains of *S. aureus* and has the ability to bind the Fc portion of immunoglobulin G (IgG)². It has been observed that certain methicillin-resistant strains of *S. aureus* (MRSA) may express undetectable levels of clumping factor and Protein A^{3,4,5}. It has been shown however that these strains all possess capsular polysaccharide⁶. The capsule can mask both Protein A and the clumping factor thereby preventing agglutination. Staphytest Plus uses blue latex particles coated with porcine fibrinogen and rabbit IgG including specific polyclonal antibodies raised against capsular polysaccharides of *S. aureus*^{7,8}. When the reagent is mixed on a card with colonies of *S. aureus*, rapid agglutination occurs through the reaction between (i) fibrinogen and clumping factor, (ii) Fc portion of IgG and Protein A, (iii) specific IgG and capsular polysaccharide. Agglutination may also occur with other species which possess clumping factor or Protein A such as *Staphylococcus hyicus* and *Staphylococcus intermedius*. If neither clumping factor, Protein A nor specific capsular polysaccharides are present, agglutination will not occur and the result will be regarded as negative. The most frequent coagulase and Protein A negative isolates of staphylococci are *Staphylococcus epidermidis*.

3. COMPONENTS OF THE KIT (DR850M)

DR0851M Staphytest Plus Test Reagent (5.6 ml)
A suspension of synthetic blue latex particles coated with both porcine fibrinogen and rabbit IgG together with specific polyclonal antibodies raised against capsular polysaccharide of *S. aureus*. Each bottle contains sufficient reagent for 100 tests.

DR0852M Staphytest Plus Control Reagent (5.6 ml)

A suspension of unsensitised synthetic blue latex particles. Each bottle contains sufficient reagents for 100 tests.

DR0500G Reaction Cards.

There are 35 disposable reaction cards provided in the kit.

Instructions for Use.

4. MATERIALS REQUIRED BUT NOT PROVIDED

Timer

Microbiological Loop

A suitable laboratory disinfectant

Positive Control: *S. aureus* strain such as ATCC 25923

Negative Control: *S. epidermidis* strain such as ATCC 12228.

5. PRECAUTIONS

IVD This product is for *in vitro* diagnostic use only. Do not freeze.

Reagents contain 0.095% sodium azide as a preservative. Sodium azide is toxic and may react with lead or copper plumbing to produce metal azides which are explosive by contact detonation. To prevent azide accumulation in plumbing flush with copious amounts of water immediately after waste disposal.

Specimen materials may contain pathogenic organisms handle with appropriate precautions.


Directions should be read and followed closely.

Please refer to the Safety Data Sheet (SDS) on company website and product labelling for information on potentially hazardous components.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

In the event of malfunction do not use device.

6. STORAGE

This kit must be stored at 2–8°C. Do not freeze. Under these conditions the reagents will retain their activity until the expiry date shown on the kit box.

7. CONTROL PROCEDURES

On each occasion the kit is used the following control procedures must be performed:

7.1. Positive Control: Use a known *S. aureus* strain such as ATCC 25923 (Culti-Loop™ R4607010). Follow the method given in Test Procedure. Ensure that agglutination occurs within 20 seconds.

7.2. Negative Control: Use a known *S. epidermidis* strain such as

ATCC12228 (Culti-Loop™ R4606500).

Follow the method given in Test Procedure. Ensure that the reagent remains smooth and non-agglutinated for the entire 20 seconds of the test.

Do not use the test if reactions with the control organisms are incorrect.

Important Procedure Notes

Do not allow the test latex to become contaminated by letting the dropper tip touch the specimen on the reaction card.

Ensure that the caps are securely fitted after each use to prevent contamination and drying out of the reagent. After use return the kit to the refrigerator ensuring that the bottle is stored in an upright position. Ensure that sufficient growth is removed from the culture plate; an insufficient quantity may give rise to false negative reactions.

8. SPECIMEN COLLECTION AND PREPARATION

For details of specimen collection and treatment a standard text book should be consulted⁹

Gram positive, catalase positive colonies may be tested from any of the following culture media:

Blood Agar

Nutrient Agar

Tryptone Soya Agar

Tryptone Soya Agar with 5% blood

Mannitol Salt Agar

Columbia Blood Agar

Columbia CNA Agar

Mueller-Hinton Agar with 5% blood

Baird-Parker Agar

Brilliance MRSA 2 Agar

Oxacillin Resistance Screening Agar (ORSA).

The use of fresh cultures grown overnight is recommended (18 to 36 hours incubation).

The tendency of colonies to cause autoagglutinating reactions increases with incubation beyond the recommended 36 hour period.

9. STANDARD TEST METHOD

9.1. Bring the Latex Test and Control Reagents to room temperature. Make sure that the Latex Reagents are mixed by vigorous shaking and expel any latex from the dropper pipette for complete mixing.

9.2. Dispense 1 drop of Test Reagent onto one of the circles on the reaction card and 1 drop of Control Reagent onto another circle.

9.3. Using a loop, pick up and smear the equivalent of 5 average-sized suspect Staphylococcal colonies (equivalent to 2–3 mm diameter of growth) onto a circle from a culture media plate and mix this in the Control Latex reagent. Spread to cover the circle. Discard the loop appropriately.

9.4. Using a separate loop proceed in the same way with the Test Reagent.

9.5. Pick up and rock the card for up to 20 seconds and observe for agglutination under normal lighting conditions. Do not use a magnifying glass.

9.6. When the test is completed dispose of the reaction cards safely into disinfectant.

9.7. Test Method for Oxacillin Resistance Screening Agar

9.7.1 As Standard Test Method.

9.7.2 Using a loop, pick up and smear the equivalent of 5 average-sized suspect Staphylococcal colonies (equivalent to 2–3 mm diameter of growth) from a culture media plate onto a circle. Using the loop, spread the colony material into a thin film.

9.7.3 Dispense 1 drop of Control Latex directly onto the thin film and mix IMMEDIATELY.

9.7.4–6. As Standard Test Method.

10. READING AND INTERPRETATION OF RESULTS

Positive Result

A result is positive if agglutination of the blue Test Reagent particles occurs within 20 seconds. This presumptively identifies the strain as a *S. aureus*.

Negative Result

A negative result is obtained if no agglutination occurs and a smooth blue suspension remains after 20 seconds in the test circle. This presumptively identifies the strain as a non *S. aureus*.

Equivocal Result

Slight graininess of the Test Reagent accompanied by no change in the appearance of the Control Reagent should be interpreted as an equivocal result. Strains should be re-tested following subculture onto non-selective media.

Uninterpretable Result

The test is uninterpretable if the Control Reagent shows agglutination. This indicates that the culture causes autoagglutination.

Granular or Stringy Reactions

Occasionally granular or stringy reactions may be seen due to the particulate nature of the test material. When such reactions are seen to occur they should be interpreted using the following criteria:

The result is **positive** when, using the Test Reagent greater clearing of the blue background is observed when compared with the reaction of the Control Reagent. The result is **negative** when there are no differences between the clearing of the blue background using the Test and Control Reagents.

Reactions occurring after 20 seconds should be ignored.

11. LIMITATIONS OF THE PROCEDURE

11.1. The tendency of isolated colonies to cause autoagglutination increases with incubation times beyond the recommended 36 hour period.

11.2. The antibody used in Staphytest Plus has been optimised to avoid potential cross-reactions with shared antigens from coagulase negative staphylococci.

It should be noted that this has been shown to reduce sensitivity to some type 18 MRSA strains¹⁰.

11.3. Some species of Staphylococci other than *S. aureus* notably *S. hyicus*, *S. intermedius*, *S. lugdunensis*, *S. xylosus*, *S. schleiferi* and *S. haemolyticus*^{11,12,13,14}, may give positive results in coagulase tests and/or rapid latex procedures. If necessary the species may be identified by biochemical

test procedures e.g. using a test for PYRase activity (Oxoid O.B.I.S. PYR ID0580M), *S. aureus* and *S. hyicus* will be PYRase negative and all the other strains named above will be positive^{15,16,17}. *S. hyicus* and *S. intermedius* are rarely encountered in the clinical laboratory.

11.4. Staphylococci isolated from urine specimens which give a weak positive¹⁸ result with Staphytest Plus may be *Staphylococcus saprophyticus*. Further identification of such isolates may be conducted using biochemical tests and novobiocin sensitivity (*S. saprophyticus* is resistant to novobiocin).

11.5. Some streptococci and possibly other organisms possessing immunoglobulin or plasma binding factors may react in the latex test and some species such as *Escherichia coli* are able to agglutinate latex particles^{19,20} non-specifically. To overcome these non-specific results a Gram-stain should be performed to ensure only typical Staphylococci are tested.

12. PERFORMANCE CHARACTERISTICS (DATA ON FILE)

The performance characteristics of Oxoid Staphytest Plus have been determined using data from the studies detailed below. It is important to note, however, that *S. aureus* is known to show considerable antigenic variation with respect to different geographical locations.

Clinical study

Oxoid Staphytest Plus was evaluated at a large French hospital. A total of 299 isolates was tested using tube coagulase as the gold standard method. In the following data analysis, results from strains known to be cross reactors^{11,12,13,14} and results from autoagglutinating strains have been omitted (n=283). The relative sensitivity was 96.5% and the relative specificity was 97.6%.

Industrial study

Oxoid Staphytest Plus was evaluated in food laboratories in a multi-centre study in the United Kingdom. A total of 621 samples was evaluated, these were colonies taken from Baird-Parker Agar, which had been inoculated with food or environmental material. The gold standard method was tube coagulase. In the

following data analysis, results from strains known to be cross reactors^{11,12,13,14} and results from autoagglutinating strains have been omitted (n=604). The relative sensitivity was 97.6% and the relative specificity was 97.1%.

13. REFERENCES:

1. Essers, L. and Radebold, K. (1980).“Rapid and Reliable Identification of Staphylococcus aureus by a Latex Agglutination Test”. J. Clin. Microbiol.12: 641–643.

2. Taussig, M. J. (1984). Processes in Pathology and Microbiology. 2nd Edn. 520–530. Blackwell, Oxford.

3. Ruane, P. J., Morgan, M. A., Citron,D. M. and Mulligan, M. E. (1986). “Failure of Rapid Agglutination Methods to Detect Oxacillin-Resistant Staphylococcus aureus”. J. Clin.Microbiol. 24: 490–492.

4. Roberts, J. I. S. and Gaston, M. A.(1987). “Protein A and coagulase expression in epidemic and nonepidemic Staphylococcus aureus”. J. Clin. Pathol. 40: 837–840.

5. Wanger, A. R., Morris, S. L., Ericsson, C., Singh, K. V. and LaRocco, M. T. (1992). “Latex Agglutination-Negative Methicillin-Resistant Staphylococcus aureus Recovered from Neonates: Epidemiologic Features and Comparison of Typing Methods”. J.Clin. Microbiol. 30: 2583–2588.

6. Fournier, J. M., Boutonnier, A. and Bouvet, A. (1989). “Staphylococcus aureus Strains Which Are Not Identified by Rapid Agglutination Methods Are of Capsular Serotype 5”. J. Clin.Microbiol. 27: 1372–1374.

7. Fournier, J. M., Bouvet, A., Boutonnier, A., Audurier, A., Goldstein, F., Pierre, J.,Bure, A., Lebrun, L. and Hochkeppel, H. K. (1987). Predominance of Capsular Polysaccharide Type 5 among Oxacillin-Resistant Staphylococcus aureus”. J. Clin.Microbiol. 25: 1932–1933.

8. Karakawa, W. W., Fournier, J. M., Vann, W. F., Arbeit, R., Schneerson, R. S. and Robbins, J. B. (1985). “Method for the Serological Typing of the Capsular Polysaccharides of Staphylococcus aureus”. J. Clin. Microbiol. 22:445–447.

9. Kloos, W. E. and Jorgensen, J. H. (1988). Staphylococci. pp. 143–153. In Manual of Clinical Microbiology. 4th Edn. (Eds) Lennette, E. H., Balows, A., Hauser, W. J. and Shadomy, H. J.Assoc. Amer. Microbiol. Washington.

10. Data on file at Oxoid Ltd.

11. Jean-Pierre, H., Darbas, H., Jean-Roussenq, A. and Boyer, G. (1989) “Pathogenicity in Two Cases ofStaphylococcus schleiferi, a Recently Described Species”. J. Clin. Microbiol. 27: 2110–2111.

12. Freney, J., Brun, Y, Bes, M., Heugnier, H., Grimont, F., Grimont, P. A. D., Nervie, C. and Fleurette, J.(1988). “Staphylococcus lugdunensis sp. nov. and Staphylococcus schleiferi sp. nov., Two species from Human Clinical Specimens”. Int. J. Sup.Bacteriol. 38: 168–172.

13. Phillips, W. E. and Kloos, W. E. (1981).“Identification of Coagulase-Positive Staphylococcus intermedius and Staphylococcus hyicus. subsp. hyicus Isolates from Veterinary Clinical Specimens”. J. Clin. Microbiol: 14:671–673.

14. van Griethuysen, A., Bes, M., Etienne, J., Zbinden, R. and Kluytmans, J.(2001). “An International Multicenter Evaluation of a new Latex Agglutination Test for Identification of Staphylococcus aureus”. J. Clin. Microbiol. 39: 86–89.

15. Schnitzler, N., Rainer, M., Conrads, G.,Frank, D. and Haase, G. (1988).“Staphylococcus lugdunensis: Report of a case of Peritonitis and an Easy-To-Perform Screening Strategy”. J. Clin.Microbiol. 26: 1939–1949.

16. Ann-Herbert, A., Crowder, C. G.,Hancock, G. A., Jarvis, W. R. and Thornsberry, C. (1998). “Characteristics of Coagulase-Negative-Staphylococci That Help Differentiate These Species of the Family Micrococcaceae”. J. Clin.Microbiol. 36: 812–813.

17. Gregson, D. B., Low, D. E., Skulnick, M.and Simor, A. E. (1988). “Problems with Rapid Agglutination Methods for Identification of Staphylococcus aureus When Staphylococcus saprophyticus Is Being Tested”. J. Clin. Microbiol. 26:1398–1399.

18. Myhre, E. B. and Kuusela, P. (1983).“Binding of Human Fibronectin to Group A, C and G Streptococci”. Infect. Immun. 40: 29–34.

19. Runehagen, A., Schonbeck, C.,Hedneru, Hessel, B. and Kronvall, G.(1981). “Binding of Fibrinogen Degradation Products to S. aureus and to -Hemolytic Streptococci Group A, C and G”. Acta. path. microbiol.Scand., Sect B. 89: 49–55.

14. SYMBOL DEFINITIONS

	Catalogue Number
	In Vitro Diagnostic Medical Device
	Consult Instructions for Use (IFU)
	Temperature Limitations (Storage temp.)
	Contains sufficient for <N> tests
	Do not use if package is damaged
	Batch Code (Lot Number)
	Use By (Expiration Date)
	Importer
	Unique Device Identifier
	Authorized representative in the European Community
	UK Conformity Assessed
	European Conformity Assessment
	Manufacturer

OXOID Limited, Wade Road, Basingstoke,
Hampshire, RG24 8PW, UK
+44 (0) 1256 841144 thermofisher.com

For technical assistance www.thermofisher.com

Version	Date of modifications introduced
X5240H	April 2024 Updated to conform to IVDR

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6. Fournier, J. M., Boutonnier, A. and Bouvet, A. (1989). “Staphylococcus aureus Strains Which Are Not Identified by Rapid Agglutination Methods Are of Capsular Serotype 5”. J. Clin.Microbiol. 27: 1372-1374.

7. Fournier, J. M., Bouvet, A., Boutonnier, A., Audurier, A., Goldstein, F., Pierre, J.,Bure, A., Lebrun, L. and Hochkeppel, H. K. (1987). Predominance of Capsular Polysaccharide Type 5 among Oxacillin-Resistant Staphylococcus aureus”. J. Clin.Microbiol. 25: 1932-1933.

8. Karakawa, W. W., Fournier, J. M., Vann, W. F., Arbeit, R., Schneerson, R. S. and Robbins, J. B. (1985). “Method for the Serological Typing of the Capsular Polysaccharides of Staphylococcus aureus”. J. Clin. Microbiol. 22:445-447.

9. Kloos, W. E. and Jorgensen, J. H. (1988). Staphylococci. pp. 143–153. In Manual of Clinical Microbiology. 4th Edn. (Eds) Lennette, E. H., Balows, A., Hausser, W. J. and Shadomy, H. J.Assoc. Amer. Microbiol. Washington.

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16. Ann-Herbert, A., Crowder, C. G.,Hancock, G. A., Jarvis, W. R. and Thornsberry, C. (1998). “Characteristics of Coagulase-Negative-Staphylococci That Help Differentiate These Species of the Family Micrococcaceae”. J. Clin.Microbiol. 36: 812-813.

17. Gregson, D. B., Low, D. E., Skulnick, M.and Simor, A. E. (1988). “Problems with Rapid Agglutination Methods for Identification of Staphylococcus aureus When Staphylococcus saprophyticus Is Being Tested”. J. Clin. Microbiol. 26:1398-1399.

18. Myhre, E. B. and Kuusela, P. (1983).“Binding of Human Fibronectin to Group A, C and G Streptococci”. Infect. Immun. 40: 29-34.

19. Runehagen, A., Schonbeck, C.,Hedneru, Hessel, B. and Kronvall, G.(1981). “Binding of Fibrinogen Degradation Products to S. aureus and to -Hemolytic Streptococci Group A, C and G”. Acta. path. microbiol.Scand., Sect B. 89: 49-55.

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6. Fournier, J. M., Boutonnier, A. and Bouvet, A. (1989). "Staphylococcus aureus Strains Which Are Not Identified by Rapid Agglutination Methods Are of Capsular Serotype 5". J. Clin. Microbiol. 27: 1372-1374.

7. Fournier, J. M., Bouvet, A., Boutonnier, A., Audurier, A., Goldstein, F., Pierre, J., Bure, A., Lebrun, L. and Hochkeppel, H. K. (1987). Predominance of Capsular Polysaccharide Type 5 among Oxacillin-Resistant Staphylococcus aureus". J. Clin. Microbiol. 25: 1932-1933.

8. Karakawa, W. W., Fournier, J. M., Vann, W. F., Arbeit, R., Schneerson, R. S. and Robbins, J. B. (1985). "Method for the Serological Typing of the Capsular Polysaccharides of Staphylococcus aureus". J. Clin. Microbiol. 22: 445-447.

9. Kloos, W. E. and Jorgensen, J. H. (1988). Staphylococci. pp. 143-153. In Manual of Clinical Microbiology. 4th Edn. (Eds) Lennette, E. H., Balows, A., Hausser, W. J. and Shadomy, H. J. Assoc. Amer. Microbiol. Washington.

10. Data on file at Oxoid Ltd.

11. Jean-Pierre, H., Darbas, H., Jean-Roussenq, A. and Boyer, G. (1989) "Pathogenicity in Two Cases of Staphylococcus schleiferi, a Recently Described Species". J. Clin. Microbiol. 27: 2110-2111.

12. Freney, J., Brun, Y., Bes, M., Heugnier, H., Grimont, F., Grimont, P. A. D., Nerve, C. and Fleurette, J. (1988). "Staphylococcus lugdunensis sp. nov. and Staphylococcus schleiferi sp. nov., Two species from Human Clinical Specimens". Int. J. Sup. Bacteriol. 38: 168-172.

13. Phillips, W. E. and Kloos, W. E. (1981). "Identification of Coagulase-Positive Staphylococcus intermedius and Staphylococcus hyicus. subsp. hyicus Isolates from Veterinary Clinical Specimens". J. Clin. Microbiol. 14: 671-673.

14. van Griethuysen, A., Bes, M., Etienne, J., Zbinden, R. and Kluytmans, J. (2001). "An International Multicenter Evaluation of a new Latex Agglutination Test for Identification of Staphylococcus aureus". J. Clin. Microbiol. 39: 86-89.

15. Schnitzler, N., Rainer, M., Conrads, G., Frank, D. and Haase, G. (1988). "Staphylococcus lugdunensis: Report of a case of Peritonitis and an Easy-To-Perform Screening Strategy". J. Clin. Microbiol. 26: 1939-1949.

16. Ann-Herbert, A., Crowder, C. G., Hancock, G. A., Jarvis, W. R. and Thornsberry, C. (1998). "Characteristics of Coagulase-Negative-Staphylococci That Help Differentiate These Species of the Family Micrococcaceae". J. Clin. Microbiol. 36: 812-813.

17. Gregson, D. B., Low, D. E., Skulnick, M. and Simor, A. E. (1988). "Problems with Rapid Agglutination Methods for Identification of Staphylococcus aureus When Staphylococcus saprophyticus Is Being Tested". J. Clin. Microbiol. 26: 1398-1399.

18. Myhre, E. B. and Kuusela, P. (1983). "Binding of Human Fibronectin to Group A, C and G Streptococci". Infect. Immun. 40: 29-34.

19. Runehegan, A., Schonbeck, C., Hedner, Hessel, B. and Kronvall, G. (1981).

"Binding of Fibrinogen Degradation Products to S. aureus and to -Hemolytic Streptococci Group A, C and G". Acta. path. microbiol. Scand., Sect B. 89: 49-55.

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Oxoid Limited, Wade Road, Hampshire
Basingstoke, RG24 8PW, UK

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6. Fournier, J. M., Boutonnier, A. and Bouvet, A. (1989). "Staphylococcus aureus Strains Which Are Not Identified by Rapid Agglutination Methods Are of Capsular Serotype 5". J. Clin. Microbiol. 27: 1372–1374.

7. Fournier, J. M., Bouvet, A., Boutonnier, A., Audurier, A., Goldstein, F., Pierre, J., Bure, A., Lebrun, L. and Hochkeppel, H. K. (1987). "Predominance of Capsular Polysaccharide Type 5 among Oxacillin-Resistant *Staphylococcus aureus*". J. Clin. Microbiol. 25: 1932–1933.

8. Karakawa, W. W., Fournier, J. M., Vann, W. F., Arbeit, R., Schneerson, R. S. and Robbins, J. B. (1985). "Method for the Serological Typing of the Capsular Polysaccharides of *Staphylococcus aureus*". J. Clin. Microbiol. 22: 445–447.

9. Kloos, W. E. and Jorgensen, J. H. (1988). Staphylococci. pp. 143–153. In Manual of Clinical Microbiology. 4th Edn. (Eds) Lennette, E. H., Balows, A., Hausser, W. J. and Shadomy, H. J. Assoc. Amer. Microbiol. Washington.

10. Podaci dostupni tvrtki Oxoid Ltd.

11. Jean-Pierre, H., Darbas, H., Jean-Roussienq, A. and Boyer, G. (1989) "Pathogenicity in Two Cases of *Staphylococcus schleiferi*, a Recently Described Species". J. Clin. Microbiol. 27: 2110–2111.

12. Freney, J., Brun, Y., Bes, M., Heugnier, H., Grimont, F., Grimont, P. A. D., Nervie, C. and Fleurette, J. (1988). "Staphylococcus lugdunensis sp. nov. and *Staphylococcus schleiferi* sp. nov., Two species from Human Clinical Specimens". Int. J. Sup. Bacteriol. 38: 168–172.

13. Phillips, W. E. and Kloos, W. E. (1981). "Identification of Coagulase-Positive *Staphylococcus intermedius* and *Staphylococcus hyicus* subsp. *hyicus* Isolates from Veterinary Clinical Specimens". J. Clin. Microbiol. 14: 671–673.

14. van Griethuysen, A., Bes, M., Etienne, J., Zbinden, R. and Kluytmans, J. (2001). "An International Multicenter Evaluation of a new Latex Agglutination Test for Identification of *Staphylococcus aureus*". J. Clin. Microbiol. 39: 86–89.

15. Schnitzler, N., Rainer, M., Conrads, G., Frank, D. and Haase, G. (1988). "Staphylococcus lugdunensis: Report of a case of Peritonitis and an Easy-To-Perform Screening Strategy". J. Clin. Microbiol. 26: 1939–1949.

16. Ann-Herbert, A., Crowder, C. G., Hancock, G. A., Jarvis, W. R. and Thornsberry, C. (1998). "Characteristics of Coagulase-Negative-Staphylococci That Help Differentiate These Species of the Family Micrococcaceae". J. Clin. Microbiol. 36: 812–813.

17. Gregson, D. B., Low, D. E., Skulnick, M. and Simor, A. E. (1988). "Problems with Rapid Agglutination Methods for Identification of *Staphylococcus aureus* When *Staphylococcus saprophyticus* Is Being Tested". J. Clin. Microbiol. 26: 1398–1399.

18. Myhre, E. B. and Kuusela, P. (1983). "Binding of Human Fibronectin to Group A, C and G Streptococci". Infect. Immun. 40: 29–34.

19. Runeheggen, A., Schonbeck, C., Hedneru, Hessel, B. and Kronvall, G. (1981). "Binding of Fibrinogen Degradation Products to *S. aureus* and to Hemolytic Streptococci Group A, C and G". Acta. path. microbiol. Scand., Sect. B. 89: 49–55.

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DR0852M Staphytest Plus Control reagents (5,6 ml);

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Oxoid Limited, Wade Road, Hampshire Basingstoke, RG24 8PW, UK