DrySpot Legionella Test

- DR0200M DrySpot Legionella pneumophila 1** ▼60
- DR0210M DrySpot Legionella pneumophila 2-14** ▼60
- DR0220M DrySpot Legionella Species** ▼60

1. INTENDED USE

DR0200M

The DrySpot Legionella pneumophila 1 Test is a qualitative latex agglutination assay for the confirmatory identification of *Legionella pneumophila* serogroup 1 from isolates grown on agar. Used in a diagnostic workflow to aid clinicians in treatment options for patients suspected of having bacterial infections. The device is not automated, is for professional use only and is not a companion diagnostic.

DR0210M

The DrySpot Legionella pneumophila 2-14 Test is a qualitative latex agglutination assay for the confirmatory identification of *Legionella pneumophila* serogroups 2-14 from isolates grown on agar. Used in a diagnostic workflow to aid clinicians in treatment options for patients suspected of having bacterial infections. The device is not automated, is for professional use only and is not a companion diagnostic.

DR0220M

The DrySpot Legionella Species Test is a qualitative latex agglutination assay for the confirmatory identification of *Legionella* species from isolates grown on agar. Used in a diagnostic workflow to aid clinicians in treatment options for patients suspected of having bacterial infections. The device is not automated, is for professional use only and is not a companion diagnostic.

The test kit allows identification of seven non-pneumophila *Legionella* species that have been implicated in human disease.

DR0200M, DR0210M and DR0220M may also be used for the identification of *Legionella* species isolates from environmental sources.

2. PRINCIPLE OF THE TEST

Legionella pneumophila has been shown to be a major cause of both pneumonia and an acute self limiting febrile disease called Pontiac Fever. *L. pneumophila* strains and other *Legionella* species are isolated predominantly from patients with pneumonia and from the environment (mainly water).

L. pneumophila is the most common cause of Legionnaires' disease. At present, 14 different serotypes are known of which serogroup 1 accounts for 90% of cases.

Rare isolations have also been made from sites other than the respiratory tract such as wound abscesses. The major reservoirs of *Legionella* species appear to be fresh water sites, air-conditioning units, cooling towers and water plumbing fixtures.

The DrySpot Legionella Test uses antibody-sensitised blue polystyrene 'latex' particles dried onto cards. The latex will agglutinate in the presence of specific *Legionella* cell wall

antigens to form visible clumps. The tests provide a fast and simple screening procedure for pathogenic *Legionella* species and serotypes.^{1,2}

3. COMPONENTS OF THE KITS

DR0200M	DrySpot Legionella pneumophila serogroup 1 Test
DR0201M	Legionella pneumophila serogroup 1 Test Cards. A suspension of synthetic blue latex particles with specific rabbit antibody reactive with <i>L. pneumophila</i> serogroup 1 antigen and dried onto cards (Test reaction area). A suspension of synthetic blue latex particles sensitised with non-reactive rabbit globulins (Control reaction area). Two pouches each containing 10 cards and a moisture-absorbent sachet. There are 3 Test and 3 Control reaction areas on each test card – 60 tests in total.
DR0205M	Positive Control Strips (10 Sticks – pink spots). Pink-dyed inactivated polyvalent antigenic extract of <i>L. pneumophila</i> 1 strains dried onto strips.
DR0240M	Negative Control Strips (10 Sticks – green spots). Green-dyed inactivated antigenic extract of <i>L. moravica</i> and <i>L. birminghamensis</i> dried onto strips.
DR0230M	DrySpot Legionella Diluent Buffer Mixing Paddles Plastic Pouch Clip for storage of opened pouches Instructions For Use
DR0210M	DrySpot Legionella pneumophila serogroups 2–14 Test
DR0211M	Legionella pneumophila serogroups 2–14 Test Cards. A suspension of synthetic blue latex particles sensitised with specific rabbit antibody reactive with <i>L. pneumophila</i> serogroups 2–14 antigen and dried onto cards (Test reaction area). A suspension of synthetic blue latex particles sensitised with non-reactive rabbit globulins (Control reaction area). Two pouches each containing 10 cards and a moisture-absorbent sachet. There are 3 Test and 3 Control reaction areas on each test card – 60 tests in total.
DR0215M	Positive Control Strips (10 Sticks – pink spots). Pink-dyed inactivated polyvalent antigenic extract of <i>L. pneumophila</i> 2–14 strains dried onto strips.
DR0240M	Negative Control Strips (10 Sticks – green spots). Green-dyed inactivated antigenic extract of <i>L. moravica</i> and <i>L. birminghamensis</i> dried onto strips.
DR0230M	DrySpot Legionella Diluent Buffer Mixing Paddles Plastic Pouch Clip for storage of opened pouches Instructions For Use
DR0220M	DrySpot Legionella Species Test
DR0221M	Legionella Species Test Cards. A suspension of synthetic blue latex particles sensitised with specific rabbit antibody reactive with the following

species and serotypes: *L. longbeachae* 1 and 2, *L. bozemanii* 1 and 2, *L. dumoffii*, *L. gormanii*, *L. jordanis*, *L. micdadei*, *L. anisa* dried onto cards (Test reaction area).

A suspension of synthetic blue latex particles sensitised with non-reactive rabbit globulins (Control reaction area).

Two pouches each containing 10 cards and a moisture-absorbent sachet. There are 3 Test and 3 Control reaction areas on each test card – 60 tests in total.

DR0225M	Positive Control Strips (10 Sticks – pink spots). Pink-dyed inactivated polyvalent antigenic extract of <i>Legionella</i> species dried onto strips.
DR0240M	Negative Control Strips (10 Sticks – green spots). Green-dyed inactivated antigenic extract of <i>L. moravica</i> and <i>L. birminghamensis</i> dried onto strips.
DR0230M	DrySpot Legionella Diluent Buffer Mixing Paddles Plastic Pouch Clip for storage of opened pouches Instructions For Use

4. MATERIALS REQUIRED BUT NOT PROVIDED:

Sterile microbiological loop.

Suitable laboratory disinfectant e.g. sodium hypochlorite solution >1.3% w/v.

5. PRECAUTIONS

IVD

This product is for *in vitro* diagnostic use only.

In the event of malfunction do not use device.

Specimen materials may contain pathogenic organisms, handle with appropriate precautions.

Reagents contain 0.1% sodium azide as a preservative.

Sodium azide may react with lead or copper plumbing to produce metal azides which are explosive by contact detonation. To prevent azide accumulation in plumbing flush with copious amounts of water immediately after waste disposal

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

6. STORAGE AND OPENING

This kit must be stored between 2°C and 25°C. If stored in a cold environment, allow pouches to reach room temperature before opening to prevent condensation of moisture on the cards. The DrySpot reagents will deteriorate and may give false results if they are allowed to absorb moisture.

Open the pouches by cutting with scissors just below the seal.

Once opened, remove the number of cards required for immediate testing (testing within the next 10 minutes) and re-seal the pouch immediately by clamping the open end of the bag between the two halves of the plastic clip provided.

If fewer tests are required cut the reaction cards along the indicated lines and return the unused portions to the pouch. Do not return used portions because they will cause contamination of remaining cards in the pouch.

The Control Sticks are also provided in a moisture-impermeable pouch. Ensure that the same techniques are used to avoid moisture damage.

Under these conditions the reagents will retain their activity until the expiry date shown on the kit box.

7. CONTROL PROCEDURES

The dried Control Strips provided should be used in the following way to check the correct working of the latex reagent each day before routine tests are performed.

Add a drop of DrySpot Legionella Buffer (DR0230M) to the small ring at the base of the test oval reaction area. Remove a Positive Control Strip (DR0205M, DR0215M or DR0225M as appropriate) from the pouch by tearing one from the others on the strip, taking care to avoid touching the flexible end where the dry spots are located. Re- seal the inner bag and pouch. Turn the stick over so that the coloured spots are at the bottom and place the stick on the card with the spots touching the liquid. Push down so that the end bends at the hinge and mix in a circular manner for 10 seconds to rehydrate the dried reagent. Continue to use the stick to mix the liquid into the appropriate DrySpot Legionella Test Reagent until all the reagents are fully rehydrated and homogeneous. Rock the card and look for agglutination. This procedure should be repeated using a Negative Control Strip (DR0240M).

The **positive control** must show agglutination with the dried reagent (of its respective kit) within 1 minute.

The **negative control** must show no agglutination within 1 minute.

Do not use the test if reactions with the control reagents are incorrect.

8. IMPORTANT PROCEDURE NOTES

Do not touch the circles on the reaction cards as this may cause contamination and affect the reaction.

In a high humidity environment do not leave the pouches open for more than 2 minutes. If there is evidence of moisture in the spots do not use.

Do not place the drop of Buffer directly onto the dry latex spots.

Pouch clips should be retained for future use to allow multiple packs to be opened. Although suitable for room temperature storage the kit or pouches must not be stored near heat sources or where exposure to sunlight may cause increased temperatures.

9. SPECIMEN COLLECTION AND PREPARATION

Isolates derived from environmental and clinical samples may be cultured on standard non-selective or selective *Legionella* culture media. Typical isolation schemes are given in references 3 and 4. To confirm the isolate as a *Legionella* it is necessary to show a requirement for L-cysteine (*Legionella* Agar without cysteine CM0655 + SR0175). This confirmation may be performed before or after the latex test.

The following Oxoid media and supplements may be used for the culture of *Legionella* before performing the latex test.

BCYE (CM0655 + SR0110), BPMAα (CM0655 + SR0110 + SR0111), MWY (CM0655 + SR0110 + SR0118), GVPC (CM0655 + SR0110 + SR0152). For further details of these products please consult your local distributor.

Cultures may be tested at any stage of growth providing that the colonies are of sufficient size. Older cultures, however, may produce stringy reactions making interpretation more difficult.

10. TEST METHOD

1. Add 1 drop of DrySpot Legionella Buffer (DR0230M) to the small ring (at the bottom of each oval) in both the Test and Control reaction areas ensuring that the liquid does not mix with the dried latex reagents.
2. Use a sterile loop (or one of the paddles provided) to pick off a suspect colony (at least a 2 mm equivalent area) from a culture media plate and carefully emulsify in the Buffer. Ensure that the resulting suspension is smooth.
3. Using the paddle mix the suspension of organisms into the dry latex spots until the mixture is homogeneous and spread over the reaction area. Discard the loop or paddle safely.
4. Using a separate loop or paddle, proceed in the same way with the Control Latex.
5. Rock the card for up to 60 seconds and look for agglutination under normal lighting conditions. Do not use a magnifying glass.
6. When the test is completed dispose of the reaction cards safely.

11. READING AND INTERPRETATION OF RESULTS

Positive Result
A result is **positive** if agglutination of the polystyrene 'latex' particles from the respective kit occurs within 1 minute. This indicates the presence of *Legionella*.

Negative Result
A **negative** result is obtained if no agglutination occurs and a smooth blue suspension remains after 60 seconds in the test area.

Reactions occurring after 60 seconds should be ignored.

12. UNINTERPRETABLE RESULT

The test is uninterpretable if the Control Reagent shows agglutination. This indicates that the culture causes autoagglutination.

13. GRANULAR OR STRINGY REACTIONS

Occasionally granular or stringy reactions may be seen due to the particulate nature of the test material. When such reactions are seen to occur they should be interpreted using the following criteria:

The result is **positive** when, using the Test Reagent, greater clearing of the blue background occurs compared with the reaction of the Control Reagent. The result is **negative**, when there is no difference between clearing of the blue background using the Test and Control Reagents.

14. LIMITATIONS

The DrySpot Legionella Test provides a presumptive result. Confirm positive results using biochemical tests⁵.

The test is designed to differentiate between different species and serotypes of *Legionella*. Cultures should be confirmed as Gram-negative rods which do not grow on cysteine-deficient media.

A negative DrySpot Legionella Test does not mean that the culture is not a *Legionella* species. It only indicates that the culture is not *L. pneumophila* serogroup 1 (when tested with the *L. pneumophila* serogroup 1 Test), not *L. pneumophila* serogroups 2–14 (when tested with the *L. pneumophila* serogroups 2–14 Test) and not *L. longbeachae* 1 and 2, *L. bozemanii* 1 and 2, *L. dumoffii*, *L. gormanii*, *L. jordanis*, *L. micdadei*, *L. anisa* when tested with the Legionella Species Test⁶. A cross reaction may occur between *L. pneumophila* serogroup 1 and serogroup 9 due to naturally-occurring group antigens. If both the *L. pneumophila* serogroup

1 and 2–14 reagents agglutinate with the isolate then this cross reaction should be suspected.

Cross reactions with the Legionella Species Test Reagent have been reported to occur occasionally with certain serotypes of other Legionellae (e.g. *L. parisiensis*, *L. sainthelensi*, *L. steigerwaltii*, *L. wadsworthii*, *L. santicrucis*, *L. tusconensis*, *L. gratiana*, *L. cincinnatiensis*)⁶.

15. PERFORMANCE CHARACTERISTICS

The reagents of the three DrySpot Legionella Tests have been tested for cross-reactivity against a panel of organisms listed below. No cross-reactivity was observed with any of the organisms.

<i>L. cherri</i>	<i>L. fairfieldensis</i>
<i>L. birminghamensis</i>	<i>L. brunensis</i>
<i>L. rubrilucens</i>	<i>L. spiritensis</i>
<i>L. maceachernii</i>	<i>Pseudomonas fluorescens</i>
<i>L. oakridgensis</i>	<i>Pseudomonas cepacia</i>
<i>L. erythra</i>	<i>Pseudomonas aeruginosa</i>
<i>L. feelei</i>	<i>Aeromonas hydrophila</i>
<i>L. hackeliae</i>	<i>Stentotrophomonas maltophilia</i>
<i>L. israeliensis</i>	<i>Bacillus subtilis</i>
<i>L. jamestowniensis</i>	<i>Citrobacter freundii</i>
<i>L. quinlivanii</i>	<i>Escherichia coli</i>
<i>L. moravica</i>	<i>Serratia marcescens</i>
<i>L. sainthelensi</i>	<i>Streptococcus faecalis</i>
<i>Staphylococcus aureus</i>	<i>Salmonella typhimurium</i>
<i>Vibrio cholera</i>	

The DrySpot Legionella Test was evaluated in a clinical and an environmental laboratory. A total of 207 isolates were tested, covering *Legionella pneumophila* serogroups 1–14 and non-pneumophila Legionellae. Each isolate was confirmed by serology.

The results of the trial are summarised below.⁷

	DrySpot Legionella Test/Serology		
		Sensitivity	Specificity
	Number	%	%
<i>Legionella pneumophila</i> serogroup 1	140/140	100	100
<i>Legionella pneumophila</i> serogroups 2–14	138/140	100	97.8 ¹
Other Legionellae included in the kit	136/141	100	94.4 ²

- 1 The Legionella pneumophila 2–14 reagent reacted with serogroups 15 and 16. Currently there is only limited data available on these serogroups. It has been demonstrated that a common antigen is shared by serogroups 15 and 16 with *Legionella pneumophila* serogroups 3 and 4 respectively.^{8,9}

L. pneumophila serogroup 15 has not currently been isolated from clinical or environmental samples in Europe¹¹ and has only been isolated once in the USA.¹⁰

This most recent serogroup to be designated contains only one strain (Lansing-3 (ATCC® 35251)). A 16th serogroup was proposed from studies of the Jena-1 isolate^{12,13} until further analysis showed that the strain did not form a unique serogroup but was a member of serogroup 4 *L. pneumophila* (monoclonal group Portland 1).^{14,15}

Additional Internal data has been collected which demonstrates that the DrySpot Legionella pneumophila serogroups 2-14 Test is able to detect the *Legionella pneumophila* serogroup 15 as well as detecting *L. pneumophila* serogroups 2-14.

However, as only one serogroup 15 strain is available for testing, it is not thought appropriate to rename the DrySpot

Legionella pneumophila serogroups 2–14 Test based on the results from a single isolate.

The Legionella DrySpot Tests benefit the user by allowing discrimination of samples into three groups: *L. pneumophila* serogroup 1, *L. pneumophila* serogroups 2-15 (with the 2-14 reagent) and other *Legionella* species in a fast and simple screening procedure.

- 2 The false positives given by the Legionella Species Kit are known cross-reactions reported in the literature.

16. REFERENCES

1. Sedgwick, A. K. and Tilton, R. C. (1983). *J. Clin. Microbiol.*, **17**: 365–368.
2. Ciesielski, C. A., Blaser, M. J. and Wang, W. L. (1986). *Infect. Immun.*, **51**: 397–404.
3. Dennis, P. J. L. (1988). Isolation of Legionellae from Environmental Specimens p. 31–34. In Harrison, T. G. and Taylor, A. G. (eds). A Laboratory Manual for Legionellae. John Wiley & Sons Ltd., Chichester, UK.
4. Dournon, E. (1988). Isolation of Legionellae from Clinical Specimen p. 13–30. In Harrison, T. G. and Taylor, A. G. (eds). A Laboratory Manual for Legionella. John Wiley & Sons Ltd., Chichester, UK.
5. Cowan and Steel's Manual for the Identification of Medical Bacteria 3rd ed. Barrow, G. I. and Feltham, R. K. A. (eds) (1993). p. 161–163. University Press, London.
6. Harrison, T. G. and Taylor, A. G. (1988). Identification of Legionellae by Serological Methods. In Harrison, T. G. and Taylor, A. G. (eds). A Laboratory Manual for *Legionella*. John Wiley & Sons Ltd., Chichester, UK.
7. Data on file, Oxoid Ltd.
8. Brenner, D. J., Steigerwalt, P. E., Epple, P. et al. (1998). *J. Clin. Microbiol.*, **26**: 1695–1703.
9. Luck, P. C., Helbig, J. H., Ehret, W. and Ott, M. (1995). *Int. J. Med. Microbiol. Virol. Parasitol. Infect. Dis.*, **282**: 35–39.
10. Brenner et al. *Legionella pneumophila* Serogroup Lansing 3 Isolated from a Patient with Fatal Pneumonia, and Descriptions of *L. pneumonia* subsp. *pneumophila* subsp. nov., *L. pneumonia* subsp. *fraseri* subsp. nov., and *L. pneumonia* subsp. *pascullei* subsp. nov. *Journal of Clinical Microbiology*; 1988; 26: 1695-1703.
11. Helbig, J. H et al. Pan-European Study on Culture Proven Legionnaires' Disease: Distribution of *Legionella pneumophila* Serogroups and Monoclonal Subgroups. *European Journal of Clinical Microbiology and Infectious Disease*; 2002; 21: 710-716.
12. Lück, C et al. Isolation of a *Legionella pneumophila* Strain Serologically Distinguishable from all known Serogroups. *Zentralblatt fuer Bakteriologie*. 1995; 282: 35-39.
13. Fry, N. K and Harrison, T. G. An evaluation of intergenic rRNA gene sequence length polymorphism analysis for the identification of Legionella species. *Molecular Identification and Epidemiology*; 1998; 47: 667-678.
14. Lück, P. C et al. DNA Polymorphisms in Strains of *Legionella pneumophila* Serogroups 3 and 4 Detected by Macrorestriction Analysis and Their Use for Epidemiological Investigation of Nosocomial Legionellosis. *Applied and Environmental Microbiology*; 1995; 61: 2000-2003.
15. Helbig, J. H. 2003. Personal Communication.

17. SYMBOL LEGEND

	Catalogue Number
	In Vitro Diagnostic Medical Device
	Consult Instructions for Use (IFU)
	Temperature Limitations (Storage temp.)
	Contains sufficient for <N> tests
	Do not use if package is damaged
	Do not re-use
	Batch Code (Lot Number)
	Use By (Expiration Date)
	Importer
	Unique Device Identifier
	Authorized representative in the European Community
	UK Conformity Assessed
	European Conformity Assessment
	Manufacturer

 Oxid Limited, Wade Road, Hampshire
Basingstoke, RG24 8PW, UK
+44 (0) 1256 841144 thermofisher.com

For technical assistance www.thermofisher.com

Version	Date of modifications introduced
X5782F	May 2025 Updated to conform to IVDR

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16. REFERENCE

- Sedgwick, A. K. and Tilton, R. C. (1983). *J. Clin. Microbiol.*, **17**: 365–368.
- Ciesielski, C. A., Blaser, M. J. and Wang, W. L. (1986). *Infect. Immun.*, **51**: 397–404.
- Dennis, P. J. L. (1988). Isolation of Legionellae from Environmental Specimens str. 31–34. In Harrison, T. G. and Taylor, A. G. (eds). A Laboratory Manual for Legionellae. John Wiley & Sons Ltd., Chichester, UK.
- Dournon, E. (1988). Isolation of Legionellae from Clinical Specimen str. 13–30. In Harrison, T. G. and Taylor, A. G. (eds). A Laboratory Manual for Legionella. John Wiley & Sons Ltd., Chichester, UK.
- Cowan and Steel's Manual for the Identification of Medical Bacteria 3rd ed. Barrow, G. I. and Feltham, R. K. A. (eds) (1993). str. 161–163. University Press, London.
- Harrison, T. G. and Taylor, A. G. (1988). Identification of Legionellae by Serological Methods. In Harrison, T. G. and Taylor, A. G. (eds). A Laboratory Manual for *Legionella*. John Wiley & Sons Ltd., Chichester, UK.
- Data on file, Oxoid Ltd.
- Brenner, D. J., Steigerwalt, P. E., Epple, P. et al. (1998). *J. Clin. Microbiol.*, **26**: 1695–1703.
- Luck, P. C., Helbig, J. H., Ehret, W. and Ott, M. (1995). *Int. J. Med. Microbiol. Virol. Parasitol. Infect. Dis.*, **282**: 35–39.
- Brenner et al. *Legionella pneumophila* Serogroup Lansing 3 Isolated from a Patient with Fatal Pneumonia, and Descriptions of *L. pneumonia* subsp. *pneumophila* subsp. nov., *L. pneumophila* subsp. *fraseri* subsp. nov., and *L. pneumophila* subsp. *pascullei* subsp. nov. *Journal of Clinical Microbiology*; 1988; 26: 1695–1703.
- Helbig, J. H et al. Pan-European Study on Culture Proven Legionnaires' Disease: Distribution of *Legionella pneumophila* Serogroups and Monoclonal Subgroups. *European Journal of Clinical Microbiology and Infectious Disease*; 2002; 21: 710–716.
- Lück, C et al. Isolation of a *Legionella pneumophila* Strain Serologically Distinguishable from all known Serogroups. *Zentralblatt fuer Bakteriologie*. 1995; 282: 35–39.
- Fry, N. K and Harrison, T. G. An evaluation of intergenic rRNA gene sequence length polymorphism analysis for the identification of Legionella species. *Molecular Identification and Epidemiology*; 1998; 47: 667–678.
- Lück, P. C et al. DNA Polymorphisms in Strains of *Legionella pneumophila* Serogroups 3 and 4 Detected by Macrorestriction Analysis and Their Use for Epidemiological Investigation of Nosocomial Legionellosis. *Applied and Environmental Microbiology*; 1995; 61: 2000–2003.
- Helbig, J. H. 2003. Osobní sdělení.

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16. VIITED

- Sedgwick, A. K. ja Tilton, R. C. (1983). *J. Clin. Microbiol.*, **17**: 365–368.
- Ciesielski, C. A., Blaser, M. J. ja Wang, W. L. (1986). *Infect. Immun.*, **51**: 397–404.
- Dennis, P. J. L. (1988). Isolation of Legionellae from Environmental Specimens, lk 31–34. In Harrison, T. G. ja Taylor, A. G. (toim). A Laboratory Manual for Legionellae. John Wiley & Sons Ltd., Chichester, UK.
- Dournon, E. (1988). Isolation of Legionellae from Clinical Specimen, lk 13–30. In Harrison, T. G. ja Taylor, A. G. (toim). A Laboratory Manual for Legionella. John Wiley & Sons Ltd., Chichester, UK.
- Cowan and Steel's Manual for the Identification of Medical Bacteria, 3rd trükk, Barrow, G. I. ja Feltham, R. K. A. (toim) (1993). Lk 161–163. University Press, London.
- Harrison, T. G. ja Taylor, A. G. (1988). Identification of Legionellae by Serological Methods. In Harrison, T. G. ja Taylor, A. G. (toim). A Laboratory Manual for *Legionella*. John Wiley & Sons Ltd., Chichester, UK.
- Uuringuandmed, Oxoid Ltd.
- Brenner, D. J., Steigerwalt, P. E., Eppler, P. et al. (1998). *J. Clin. Microbiol.*, **26**: 1695–1703.
- Luck, P. C., Helbig, J. H., Ehret, W. ja Ott, M. (1995). *Int. J. Med. Microbiol. Virol. Parasitol. Infect. Dis.*, **282**: 35–39.
- Brenner et al. *Legionella pneumophila* Serogroup Lansing 3 Isolated from a Patient with Fatal Pneumonia, and Descriptions of *L. pneumonia* subsp. *pneumophila* subsp. nov., *L. pneumonia* subsp. *fraseri* subsp. nov., and *L. pneumonia* subsp. *pascullei* subsp. nov. *Journal of Clinical Microbiology*; 1988; 26: 1695–1703.
- Helbig, J. H et al. Pan-European Study on Culture Proven Legionnaires' Disease: Distribution of *Legionella pneumophila* Serogroups and Monoclonal Subgroups. European Journal of Clinical Microbiology and Infectious Disease; 2002; 21: 710–716.
- Lück, C et al. Isolation of a *Legionella pneumophila* Strain Serologically Distinguishable from all known Serogroups. Zentralblatt fuer Bakteriologie. 1995; 282: 35–39.
- Fry, N. K and Harrison, T. G. An evaluation of intergenic rRNA gene sequence length polymorphism analysis for the identification of Legionella species. Molecular Identification and Epidemiology; 1998; 47: 667–678.
- Lück, P. C et al. DNA Polymorphisms in Strains of *Legionella pneumophila* Serogroups 3 and 4 Detected by Macrorestriction Analysis and Their Use for Epidemiological Investigation of Nosocomial Legionellosis. *Applied and Environmental Microbiology*; 1995; 61: 2000–2003.
- Helbig, J. H. 2003. Islikultus.

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4. Dournon, E. (1988). Isolation of Legionellae from Clinical Specimen p. 13–30. In Harrison, T. G. and Taylor, A. G. (eds). A Laboratory Manual for Legionella. John Wiley & Sons Ltd., Chichester, UK.
5. Cowan and Steel's Manual for the Identification of Medical Bacteria 3rd ed. Barrow, G. I. and Feltham, R. K. A. (eds) (1993). p. 161–163. University Press, London.
6. Harrison, T. G. and Taylor, A. G. (1988). Identification of Legionellae by Serological Methods. In Harrison, T. G. and Taylor, A. G. (eds). A Laboratory Manual for *Legionella*. John Wiley & Sons Ltd., Chichester, UK.
7. Dostupni podaci, Oxoid Ltd.
8. Brenner, D. J., Steigerwalt, P. E., Eppe, P. et al. (1998). *J. Clin. Microbiol.*, **26**: 1695–1703.
9. Luck, P. C., Helbig, J. H., Ehret, W. and Ott, M. (1995). *Int. J. Med. Microbiol. Virol. Parasitol. Infect. Dis.*, **282**: 35–39.
10. Brenner et al. *Legionella pneumophila* Serogroup Lansing 3 Isolated from a Patient with Fatal Pneumonia, and Descriptions of *L. pneumonia* subsp. *pneumophila* subsp. nov., *L. pneumophila* subsp. *fraseri* subsp. nov., and *L. pneumophila* subsp. *pascuelli* subsp. nov. *Journal of Clinical Microbiology*; 1988; 26: 1695-1703.
11. Helbig, J. H et al. Pan-European Study on Culture Proven Legionnaires' Disease: Distribution of *Legionella pneumophila* Serogroups and Monoclonal Subgroups. *European Journal of Clinical Microbiology and Infectious Disease*; 2002; 21: 710-716.
12. Lück, C et al. Isolation of a *Legionella pneumophila* Strain Serologically Distinguishable from all known Serogroups. *Zentralblatt fuer Bakteriologie*. 1995; 282: 35-39.
13. Fry, N. K and Harrison, T. G. An evaluation of intergenic rRNA gene sequence length polymorphism analysis for the identification of Legionella species. *Molecular Identification and Epidemiology*; 1998; 47: 667-678.
14. Lück, P. C et al. DNA Polymorphisms in Strains of *Legionella pneumophila* Serogroups 3 and 4 Detected by Macrorestriction Analysis and Their Use for Epidemiological Investigation of Nosocomial Legionellosis. *Applied and Environmental Microbiology*; 1995; 61: 2000-2003.
15. Helbig, J. H. 2003. Personal Communication.

17. TUMAČ SIMBOLA

	Kataloški broj
	In vitro dijagnostički medicinski proizvod
	Pogledajte upute za uporabu (IFU)
	Ograničenja temperature (temp. skladištenja)
	Sadržava dovoljnu količinu za <N> testova
	Nemojte upotrebljavati ako je pakiranje oštećeno
	Nemojte ponovno upotrebljavati
	Kôd serije (broj partije)
	Upotrijebiti do (datum isteka roka trajanja)
	Uvoznik
	Jedinstvena identifikacija proizvoda
	Ovlašteni zastupnik u Europskoj zajednici
	Ocijenjena sukladnost za UK
	Europska ocjena sukladnosti
	Proizvođač

Oxoid Limited, Wade Road, Hampshire Basingstoke, RG24 8PW, UK

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16. LITERATŪRA

- Sedgwick, A. K. and Tilton, R. C. (1983). *J. Clin. Microbiol.*, **17**: 365-368.
- Ciesielski, C. A., Blaser, M. J. and Wang, W. L. (1986). *Infect. Immun.*, **51**: 397-404.
- Dennis, P. J. L. (1988). Isolation of *Legionellae* from Environmental Specimens p. 31–34. In Harrison, T. G. and Taylor, A. G. (eds). A Laboratory Manual for *Legionellae*. John Wiley & Sons Ltd., Chichester, UK.
- Dournon, E. (1988). Isolation of *Legionellae* from Clinical Specimen p. 13–30. In Harrison, T. G. and Taylor, A. G. (eds). A Laboratory Manual for *Legionella*. John Wiley & Sons Ltd., Chichester, UK.
- Cowan and Steel’s Manual for the Identification of Medical Bacteria 3rd ed. Barrow, G. I. and Feltham, R. K. A. (eds) (1993). p. 161–163. University Press, London.
- Harrison, T. G. and Taylor, A. G. (1988). Identification of *Legionellae* by Serological Methods. In Harrison, T. G. and Taylor, A. G. (eds). A Laboratory Manual for *Legionella*. John Wiley & Sons Ltd., Chichester, UK.
- Faile pateikti duomenys, Oxoid Ltd.
- Brenner, D. J., Steigerwalt, P. E., Epple, P. et al. (1998). *J. Clin. Microbiol.*, **26**: 1695-1703.
- Luck, P. C., Helbig, J. H., Ehret, W. and Ott, M. (1995). *Int. J. Med. Microbiol. Virol. Parasitol. Infect. Dis.*, **282**: 35-39.
- Brenner et al. *Legionella pneumophila* Serogroup Lansing 3 Isolated from a Patient with Fatal Pneumonia, and Descriptions of *L. pneumonia* subsp. *pneumophila* subsp. nov., *L. pneumonia* subsp. *fraseri* subsp. nov., and *L. pneumonia* subsp. *pascuelli* subsp. nov. *Journal of Clinical Microbiology*; 1988; 26: 1695-1703.
- Helbig, J. H et al. Pan-European Study on Culture Proven Legionnaires’ Disease: Distribution of *Legionella pneumophila* Serogroups and Monoclonal Subgroups. *European Journal of Clinical Microbiology and Infectious Disease*; 2002; 21: 710-716.
- Lück, C et al. Isolation of a *Legionella pneumophila* Strain Serologically Distinguishable from all known Serogroups. *Zentralblatt fuer Bakteriologie*. 1995; 282: 35-39.
- Fry, N. K and Harrison, T. G. An evaluation of intergenic rRNA gene sequence length polymorphism analysis for the identification of *Legionella* species. *Molecular Identification and Epidemiology*; 1998; 47: 667-678.
- Lück, P. Cet al. DNA Polymorphisms in Strains of *Legionella pneumophila* Serogroups 3 and 4 Detected by Macrorestriction Analysis and Their Use for Epidemiological Investigation of Nosocomial Legionellosis. *Applied and Environmental Microbiology*; 1995; 61: 2000-2003.
- Helbig, J. H. 2003. Personal Communication.

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 Oxoid Limited, Wade Road, Hampshire Basingstoke, RG24 8PW, JK
+44 (0) 1256 841144 thermofisher.com

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Spausdinta JK

Key Code TSMX5782F
www.thermofisher.com

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5. Cowan and Steel’s Manual for the Identification of Medical Bacteria 3rd ed. Barrow, G. I. and Feltham, R. K. A. (eds) (1993). p. 161–163. University Press, London.

6. Harrison, T. G. and Taylor, A. G. (1988). Identification of Legionellae by Serological Methods. In Harrison, T. G. and Taylor, A. G. (eds). A Laboratory Manual for *Legionella*. John Wiley & Sons Ltd., Chichester, UK.

7. Data on file, Oxoid Ltd.

8. Brenner, D. J., Steigerwalt, P. E., Epple, P. et al. (1998). *J. Clin. Microbiol.*, **26**: 1695–1703.

9. Luck, P. C., Helbig, J. H., Ehret, W. and Ott, M. (1995). *Int. J. Med. Microbiol. Virol. Parasitol. Infect. Dis.*, **282**: 35–39.

10. Brenner et al. Legionella pneumophila Serogroup Lansing 3 Isolated from a Patient with Fatal Pneumonia, and Descriptions of L. pneumonia subsp. pneumonia subsp. nov., L. pneumophila subsp. fraseri subsp. nov., and L. pneumophila subsp. pasculli subsp. nov. Journal of Clinical Microbiology; 1988; 26: 1695-1703.

11. Helbig, J. H et al. Pan-European Study on Culture Proven Legionnaires’ Disease: Distribution of Legionella pneumophila Serogroups and Monoclonal Subgroups. European Journal of Clinical Microbiology and Infectious Disease; 2002; 21: 710-716.

12. Lück, C et al. Isolation of a Legionella pneumophila Strain Serologically Distinguishable from all known Serogroups. Zentralblatt fuer Bakteriologie. 1995; 282: 35-39.

13. Fry, N. K and Harrison, T. G. An evaluation of intergenic rRNA gene sequence length polymorphism analysis for the identification of Legionella species. Molecular Identification and Epidemiology; 1998; 47: 667-678.

14. Lück, P. C et al. DNA Polymorphisms in Strains of Legionella pneumophila Serogroups 3 and 4 Detected by Macrorestriction Analysis and Their Use for Epidemiological Investigation

of Nosocomial Legionellosis. Applied and Environmental Microbiology; 1995; 61: 2000-2003.

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Oxoid Limited, Wade Road, Hampshire
Basingstoke, RG24 8PW, UK

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16. REFERENCIAS

- Sedgwick, A. K. and Tilton, R. C. (1983). *J. Clin. Microbiol.*, **17**: 365–368.
- Ciesielski, C. A., Blaser, M. J. and Wang, W. L. (1986). *Infect. Immun.*, **51**: 397–404.
- Dennis, P. J. L. (1988). Isolation of Legionellae from Environmental Specimens p. 31–34. In Harrison, T. G. and Taylor, A. G. (eds). A Laboratory Manual for Legionellae. John Wiley & Sons Ltd., Chichester, UK.
- Dournon, E. (1988). Isolation of Legionellae from Clinical Specimen p. 13–30. In Harrison, T. G. and Taylor, A. G. (eds). A Laboratory Manual for Legionella. John Wiley & Sons Ltd., Chichester, UK.
- Cowan and Steel’s Manual for the Identification of Medical Bacteria 3rd ed. Barrow, G. I. and Feltham, R. K. A. (eds) (1993). p. 161–163. University Press, London.
- Harrison, T. G. and Taylor, A. G. (1988). Identification of Legionellae by Serological Methods. In Harrison, T. G. and Taylor, A. G. (eds). A Laboratory Manual for *Legionella*. John Wiley & Sons Ltd., Chichester, UK.
- Data on file, Oxoid Ltd.
- Brenner, D. J., Steigerwalt, P. E., Epple, P. et al. (1998). *J. Clin. Microbiol.*, **26**: 1695–1703.
- Luck, P. C., Helbig, J. H., Ehret, W. and Ott, M. (1995). *Int. J. Med. Microbiol. Virol. Parasitol. Infect. Dis.*, **282**: 35–39.
- Brenner et al. Legionella pneumophila Serogroup Lansing 3 Isolated from a Patient with Fatal Pneumonia, and Descriptions of L. pneumonia subsp. pneumophila subsp. nov., L. pneumophila subsp. fraseri subsp. nov., and L. pneumophila subsp. pascullei subsp. nov. Journal of Clinical Microbiology; 1988; 26: 1695-1703.
- Helbig, J. H et al. Pan-European Study on Culture Proven Legionnaires’ Disease: Distribution of Legionella pneumophila Serogroups and Monoclonal Subgroups. European Journal of Clinical Microbiology and Infectious Disease; 2002; 21: 710-716.
- Lück, C et al. Isolation of a Legionella pneumophila Strain Serologically Distinguishable from all known Serogroups. Zentralblatt fuer Bakteriologie. 1995; 282: 35-39.
- Fry, N. K and Harrison, T. G. An evaluation of intergenic rRNA gene sequence length polymorphism analysis for the identification of Legionella species. Molecular Identification and Epidemiology; 1998; 47: 667-678.
- Lück, P. C et al. DNA Polymorphisms in Strains of Legionella pneumophila Serogroups 3 and 4 Detected by Macrorestriction Analysis and Their Use for Epidemiological Investigation of Nosocomial Legionellosis. Applied and Environmental Microbiology; 1995; 61:

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 Oxoid Limited, Wade Road, Hampshire Basingstoke, RG24 8PW, UK
+44 (0) 1256 841144 thermofisher.com

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16. RÉFÉRENCES

- 1. Sedgwick, A. K. and Tilton, R. C. (1983). J. Clin. Microbiol., 17: 365–368.
- 2. Ciesielski, C. A., Blaser, M. J. and Wang, W. L. (1986). Infect. Immun., 51: 397–404.
- 3. Dennis, P. J. L. (1988). Isolation of Legionellae from Environmental Specimens p. 31–34. In Harrison, T. G. and Taylor, A. G. (eds). A Laboratory Manual for Legionellae. John Wiley & Sons Ltd., Chichester, UK.
- 4. Dournon, E. (1988). Isolation of Legionellae from Clinical Specimen p. 13–30. In Harrison, T. G. and Taylor, A. G. (eds). A Laboratory Manual for Legionella. John Wiley & Sons Ltd., Chichester, UK.
- 5. Cowan and Steel’s Manual for the Identification of Medical Bacteria 3rd ed. Barrow, G. I. and Feltham, R. K. A. (eds) (1993). p. 161–163. University Press, London.
- 6. Harrison, T. G. and Taylor, A. G. (1988). Identification of Legionellae by Serological Methods. In Harrison, T. G. and Taylor, A. G. (eds). A Laboratory Manual for Legionella. John Wiley & Sons Ltd., Chichester, UK.
- 7. Data on file, Oxoid Ltd.
- 8. Brenner, D. J., Steigerwalt, P. E., Epple, P. et al. (1998). J. Clin. Microbiol., 26: 1695–1703.
- 9. Luck, P. C., Helbig, J. H., Ehret, W. and Ott, M. (1995). Int. J. Med. Microbiol. Virol. Parasitol. Infect. Dis., 282: 35–39.
- 10. Brenner et al. Legionella pneumophila Serogroup Lansing 3 Isolated from a Patient with Fatal Pneumonia, and Descriptions of L. pneumonia subsp. pneumophila subsp. nov., L. pneumophila subsp. fraseri subsp. nov., and L. pneumophila subsp. pascuelli subsp. nov. Journal of Clinical Microbiology; 1988; 26: 1695-1703.
- 11. Helbig, J. H et al. Pan-European Study on Culture Proven Legionnaires’ Disease: Distribution of Legionella pneumophila Serogroups and Monoclonal Subgroups. European Journal of Clinical Microbiology and Infectious Disease; 2002; 21: 710-716.
- 12. Lück, C et al. Isolation of a Legionella pneumophila Strain Serologically Distinguishable from all known Serogroups. Zentralblatt fuer Bakteriologie. 1995; 282: 35-39.
- 13. Fry, N. K and Harrison, T. G. An evaluation of intergenic rRNA gene sequence length polymorphism analysis for the identification of Legionella species. Molecular Identification and Epidemiology; 1998; 47: 667-678.
- 14. Lück, P. C et al. DNA Polymorphisms in Strains of Legionella pneumophila Serogroups 3 and 4 Detected by Macrorestriction Analysis and Their Use for Epidemiological Investigation of Nosocomial

Legionellosis. Applied and Environmental Microbiology; 1995; 61: 2000-2003.

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Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, UK

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16. BIBLIOGRAFIA:

- 1. Sedgwick, A. K. and Tilton, R. C. (1983). *J. Clin. Microbiol.*, **17**: 365–368.
- 2. Ciesielski, C. A., Blaser, M. J. and Wang, W. L. (1986). *Infect. Immun.*, **51**: 397–404.
- 3. Dennis, P. J. L. (1988). Isolation of Legionellae from Environmental Specimens p. 31–34. In Harrison, T. G. and Taylor, A. G. (eds). A Laboratory Manual for Legionellae. John Wiley & Sons Ltd., Chichester, UK.
- 4. Dournon, E. (1988). Isolation of Legionellae from Clinical Specimen p. 13–30. In Harrison, T. G. and Taylor, A. G. (eds). A Laboratory Manual for Legionella. John Wiley & Sons Ltd., Chichester, UK.
- 5. Cowan and Steel's Manual for the Identification of Medical Bacteria 3rd ed. Barrow, G. I. and Feltham, R. K. A. (eds) (1993). p. 161–163. University Press, London.
- 6. Harrison, T. G. and Taylor, A. G. (1988). Identification of Legionellae by Serological Methods. In Harrison, T. G. and Taylor, A. G. (eds). A Laboratory Manual for *Legionella*. John Wiley & Sons Ltd., Chichester, UK.
- 7. Data on file, Oxoid Ltd.
- 8. Brenner, D. J., Steigerwalt, P. E., Epple, P. et al. (1998). *J. Clin. Microbiol.*, **26**: 1695–1703.
- 9. Luck, P. C., Helbig, J. H., Ehret, W. and Ott, M. (1995). *Int. J. Med. Microbiol. Virol. Parasitol. Infect. Dis.*, **282**: 35–39.
- 10. Brenner et al. *Legionella pneumophila* Serogroup Lansing 3 Isolated from a Patient with Fatal Pneumonia, and Descriptions of L. pneumonia subsp. pneumophila subsp. nov., L. pneumophila subsp. fraseri subsp. nov., and L. pneumophila subsp. pascullei subsp. nov. *Journal of Clinical Microbiology*; 1988; 26: 1695-1703.
- 11. Helbig, J. H et al. Pan-European Study on Culture Proven Legionnaires' Disease: Distribution of *Legionella pneumophila* Serogroups and Monoclonal Subgroups. *European Journal of Clinical Microbiology and Infectious Disease*; 2002; 21: 710-716.
- 12. Lück, C et al. Isolation of a *Legionella pneumophila* Strain Serologically Distinguishable from all known Serogroups. *Zentralblatt fuer Bakteriologie*. 1995; 282: 35-39.
- 13. Fry, N. K and Harrison, T. G. An evaluation of intergenic rRNA gene sequence length polymorphism analysis for the identification of *Legionella* species. *Molecular Identification and Epidemiology*; 1998; 47: 667-678.
- 14. Lück, P. C et al. DNA Polymorphisms in Strains of *Legionella pneumophila* Serogroups 3 and 4 Detected by Macrorestriction Analysis and Their Use for Epidemiological Investigation of Nosocomial Legionellosis. *Applied and Environmental Microbiology*; 1995; 61: 2000-2003.
- 15. Helbig, J. H. 2003. Personal Communication.

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Oxoid Ltd Wade Road,
Basingstoke, Hampshire,
RG24 8PW, Regno Unito

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2

16. REFERENCIAS:

- Sedgwick, A. K. and Tilton, R. C. (1983). *J. Clin. Microbiol.*, **17**: 365-368.
- Ciesielski, C. A., Blaser, M. J. and Wang, W. L. (1986). *Infect. Immun.*, **51**: 397-404.
- Dennis, P. J. L. (1988). Isolation of Legionellae from Environmental Specimens p. 31-34. In Harrison, T. G. and Taylor, A. G. (eds). A Laboratory Manual for Legionellae. John Wiley & Sons Ltd., Chichester, UK.
- Dournon, E. (1988). Isolation of Legionellae from Clinical Specimen p. 13-30. In Harrison, T. G. and Taylor, A. G. (eds). A Laboratory Manual for Legionella. John Wiley & Sons Ltd., Chichester, UK.
- Cowan and Steel's Manual for the Identification of Medical Bacteria 3rd ed. Barrow, G. I. and Feltham, R. K. A. (eds) (1993). p. 161-163. University Press, London.
- Harrison, T. G. and Taylor, A. G. (1988). Identification of Legionellae by Serological Methods. In Harrison, T. G. and Taylor, A. G. (eds). A Laboratory Manual for *Legionella*. John Wiley & Sons Ltd., Chichester, UK.
- Data on file, Oxoid Ltd.
- Brenner, D. J., Steigerwalt, P. E., Epple, P. et al. (1998). *J. Clin. Microbiol.*, **26**: 1695-1703.
- Luck, P. C., Helbig, J. H., Ehret, W. and Ott, M. (1995). *Int. J. Med. Microbiol. Virol. Parasitol. Infect. Dis.*, **282**: 35-39.
- Brenner et al. Legionella pneumophila Serogroup Lansing 3 Isolated from a Patient with Fatal Pneumonia, and Descriptions of *L. pneumonia* subsp. *pneumophila* subsp. nov., *L. pneumophila* subsp. *fraseri* subsp. nov., and *L. pneumophila* subsp. *pascullei* subsp. nov. *Journal of Clinical Microbiology*; 1988; 26: 1695-1703.
- Helbig, J. H et al. Pan-European Study on Culture Proven Legionnaires' Disease: Distribution of Legionella pneumophila Serogroups and Monoclonal Subgroups. *European Journal of Clinical Microbiology and Infectious Disease*; 2002; 21: 710-716.
- Lück, C et al. Isolation of a Legionella pneumophila Strain Serologically Distinguishable from all known Serogroups. *Zentralblatt fuer Bakteriologie*. 1995; 282: 35-39.
- Fry, N. K and Harrison, T. G. An evaluation of intergenic rRNA gene sequence length polymorphism analysis for the identification of Legionella species. *Molecular Identification and Epidemiology*; 1998; 47: 667-678.
- Lück, P. C et al. DNA Polymorphisms in Strains of Legionella pneumophila Serogroups 3 and 4 Detected by Macrorestriction Analysis and Their Use for Epidemiological Investigation of Nosocomial Legionellosis. *Applied and Environmental Microbiology*; 1995; 61:

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Oxoid Ltd Wade Road,
Basingstoke, Hampshire,
RG24 8PW, Reino Unido

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16. REFERENSER:

1. Sedgwick, A. K. and Tilton, R. C. (1983). *J. Clin. Microbiol.*, **17**: 365–368.
2. Ciesielski, C. A., Blaser, M. J. and Wang, W. L. (1986). *Infect. Immun.*, **51**: 397–404.
3. Dennis, P. J. L. (1988). Isolation of *Legionellae* from Environmental Specimens p. 31–34. In Harrison, T. G. and Taylor, A. G. (eds). A Laboratory Manual for *Legionellae*. John Wiley & Sons Ltd., Chichester, UK.
4. Dournon, E. (1988). Isolation of *Legionellae* from Clinical Specimen p. 13–30. In Harrison, T. G. and Taylor, A. G. (eds). A Laboratory Manual for *Legionella*. John Wiley & Sons Ltd., Chichester, UK.
5. Cowan and Steel's Manual for the Identification of Medical Bacteria 3rd ed. Barrow, G. I. and Feltham, R. K. A. (eds) (1993). p. 161–163. University Press, London.
6. Harrison, T. G. and Taylor, A. G. (1988). Identification of *Legionellae* by Serological Methods. In Harrison, T. G. and Taylor, A. G. (eds). A Laboratory Manual for *Legionella*. John Wiley & Sons Ltd., Chichester, UK.
7. Data on file, Oxoid Ltd.
8. Brenner, D. J., Steigerwalt, P. E., Epple, P. et al. (1998). *J. Clin. Microbiol.*, **26**: 1695–1703.
9. Luck, P. C., Helbig, J. H., Ehret, W. and Ott, M. (1995). *Int. J. Med. Microbiol. Virol. Parasitol. Infect. Dis.*, **282**: 35–39.
10. Brenner et al. *Legionella pneumophila* Serogroup Lansing 3 Isolated from a Patient with Fatal Pneumonia, and Descriptions of *L. pneumonia* subsp. *pneumophila* subsp. nov., *L. pneumophila* subsp. *fraseri* subsp. nov., and *L. pneumophila* subsp. *pascuelli* subsp. nov. *Journal of Clinical Microbiology*; 1988; 26: 1695-1703.
11. Helbig, J. H et al. Pan-European Study on Culture Proven Legionnaires' Disease: Distribution of *Legionella pneumophila* Serogroups and Monoclonal Subgroups. *European Journal of Clinical Microbiology and Infectious Disease*; 2002; 21: 710-716.
12. Lück, C et al. Isolation of a *Legionella pneumophila* Strain Serologically Distinguishable from all known Serogroups. *Zentralblatt fuer Bakteriologie*. 1995; 282: 35-39.
13. Fry, N. K and Harrison, T. G. An evaluation of intergenic rRNA gene sequence length polymorphism analysis for the identification of *Legionella* species. *Molecular Identification and Epidemiology*; 1998; 47: 667-678.
14. Lück, P. C et al. DNA Polymorphisms in Strains of *Legionella pneumophila* Serogroups 3 and 4 Detected by Macrorestriction Analysis and Their Use for Epidemiological Investigation of Nosocomial Legionellosis. *Applied and Environmental Microbiology*; 1995; 61: 2000-2003.
15. Helbig, J. H. 2003. Personal Communication.

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