www.thermofisher.com

# Oxoid Salmonella Test Kit

REF DR1108A..... 100 EN

## 1. INTRODUCTION

The genus *Salmonella* consists of over 2,400 serotypes, which may cause acute infectious gastroenteritis, frequently associated with food poisoning. Since many *Salmonella* species infect domestic animals, either clinically or subclinically, cases of *Salmonella* food poisoning usually originate from animal sources. In contrast, *S. typhi* and *S. paratyphi* only infect humans and cause enteric fever. Salmonellae are often present in a sample in small numbers, together with other Enterobacteriaceae. In order to promote the growth of *Salmonella*, enrichment and/or selective media are used, such as Rappaport Medium, broths for inducing motility (e.g. M broth), and those containing selenite, (RV- CM0669, RVS-CM0866) deoxycholate citrate (DCA – CM0035) or bismuth sulphite.<sup>1,2</sup>

IVD For in vitro diagnostic use only

## 2. INTENDED USE

The Oxoid™ Salmonella Test Kit is a qualitative latex slide agglutination test for the confirmatory identification of *Salmonella* species isolates by the detection of *Salmonella* antigens from isolates grown in selective and/or enrichment broths or solid media. The test can be used in a diagnostic workflow to aid clinicians in deciding treatment options for patients suspected of having *Salmonella* infections. The device is not automated, is for professional use only and is not a companion diagnostic or intended for self-testing. The test can be used for screening samples from the food industry for the presence of *Salmonella* species.

## 3. PRINCIPLE OF THE TEST

The latex agglutination test is a traditional, long established method used for simple presumptive identification of microorganisms. Polyvalent antisera are prepared against a wide range of *Salmonella* flagellar antigens and used to coat latex particles. When mixed with a suspension of *Salmonella* containing these antigens, the latex particles rapidly agglutinate to form visible clumps. The Oxoid Salmonella Kit detects the majority of common *Salmonella* species including *S. typhimurium* and *S. enteritidis*. In order to minimise cross reactions with other Enterobacteriaceae, antibodies to the principal somatic antigens are removed during preparation of the antisera. Consequently, the Oxoid Salmonella Kit reacts predominantly with flagellar antigens although the reagent will react with the non-motile species *S. pullorum* and *S. gallinarum*.

## 4. MATERIALS SUPPLIED

Salmonella Latex Reagent (1101): Latex particles coated with rabbit antiserum against *Salmonella* antigens, preserved with 0.09% merthiolate. 5.0ml

0.85% isotonic saline (1098), preserved with 0.09% sodium azide. 5ml

Positive Control Suspension (IN1097): A suspension of non-viable *Salmonella* preserved with 1% formalin (0.5ml)

Disposable reaction cards

## Materials Required but Not Supplied

Disposable mixing sticks

Pasteur pipettes

## 5. STORAGE

The Oxoid Salmonella Test Kit should be stored at 2-8°C. Do not freeze. The kit should not be used after the expiry date printed on the outside of the carton.

Return kit to normal storage conditions after use.

## 6. PRECAUTIONS

- Kit reagents may settle during storage and should be mixed by gentle inversion prior to use.
- Appropriate precautions should be taken when working with potential pathogens.
- Do not mouth pipette.
- Discarded reagents and materials must be autoclaved before disposal.
- Do not cross contaminate solutions or samples.
- Do not use after expiry date printed on outside of carton.
- Ensure the reaction card is clean and dry before use.
- The isotonic saline (1098) is preserved with sodium azide which can form potentially explosive compounds on contact with copper and lead plumbing. On disposal of saline, flush with copious quantities of water to prevent azide build up.
- Please refer to the Safety Data Sheet (SDS), available on the ThermoFisher website, and product labelling for information on potentially hazardous components.
- Disposable reaction cards are for single use only and must not be re-used.
- In the event of a malfunction do not use the device.
- Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

## 7. METHOD OF USE

Allow the Oxoid Salmonella Test Kit reagents to equilibrate to room temperature before use.

## Specimens

Selective and/or enrichment broths or solid media can be tested with the Oxoid Salmonella Test Kit. The recommended procedure for optimum performance of the Kit is shown below:

Food samples

- Homogenize the food sample in 225ml Resuscitation broth (e.g BPW (ISO) CM1049).
- Incubate the broth aerobically for 18 hours at 35°C.
- Transfer 1ml of incubated broth in either:
  - 10ml of Selective broth (CM1048 + SR0181, or CM0399) or
  - 100ml (or 0,1ml of incubated broth in 10 ml) of Rappaport Broth (CM0669, or CM0399)
- Incubate the broth aerobically for 18 hours at 35°C.
- Subculture 1ml of the broth into 10ml/M broth

- Further incubate aerobically for 6 hours at 35°C.

- Use with the Oxoid Salmonella kit as directed in Section 8

## Controls

The following controls should be performed weekly, and on receipt of a new product batch..

- Reagent Control

Add 1 drop of Oxoid Salmonella Latex Reagent (1101) to 1 drop of saline (1098) in the same circle on a slide, mix for 2 minutes observing for agglutination. No agglutination should occur. If this control shows agglutination either the Latex Reagent (1101) or the saline (1098) is probably contaminated and should be discarded.

## 8. PROCEDURE

### A. Method for Testing Broths

- Mix the broth culture by gently inverting, and remove 1 drop of culture using a Pasteur pipette.
- Place a single drop of the broth culture within one circle on the reaction card.
- Mix the Oxoid Salmonella Latex Reagent (1101) by gently inverting and add 1 drop to the broth on the card. DO NOT ALLOW THE DROPPER TO TOUCH THE DROP OF CULTURE AS THIS MAY IMPACT THE PRODUCTS PERFORMANCE. Mix thoroughly using a clean mixing stick or inoculating loop. Gently rock the reaction card 2 or 3 times. Extensive rocking of the slide is not necessary.
- Examine for agglutination within a maximum of 2 minutes.
- After reading dispose of the reaction cards in disinfectant.

**Note:** When testing food powders cultured in selective enrichment broths, ensure that the food product alone in enrichment broth does not cause non-specific agglutination. Place one drop of food product in enrichment broth on a reaction card. Observe for any agglutination or clumping which would indicate autoagglutination. Such products cannot be tested with the Oxoid Salmonella Test Kit.

### B. Method for Identification from Selective Solid Media

- Place a reaction card on the work bench.
- Add 1 drop of saline (1098) within one circle on the reaction card.
- Using a mixing stick or inoculating loop emulsify 2 - 3mm of suspect colonies in the drop of saline to produce a thick smooth suspension. Suspensions should be made from colonies with morphologies resembling *Salmonella* species.
- Observe the suspension for any agglutination or clumping which would indicate autoagglutination. If the suspension remains smooth, proceed to Step 5. (See Limitations of Use, point 1.)
- Mix the Oxoid Salmonella Latex Reagent (1101) gently by inverting and add a drop to the saline suspension. DO NOT ALLOW THE DROPPER TO TOUCH THE ORGANISM SUSPENSION. Mix the Latex Reagent and organism suspension with a clean mixing stick for 30 seconds and gently rock the reaction card two or three times. Excessive rocking of the slide is not necessary. Examine for agglutination within a maximum of two minutes.
- After reading, discard the used reaction cards into a suitable disinfectant and follow local disposal procedures.

## 9. INTERPRETATION

Agglutination within two minutes is a positive result and indicates the presence of *Salmonella* in the sample. Care must be taken when testing broth cultures that adherence of latex to particulate matter is not interpreted as agglutination. Absence of agglutination indicates that *Salmonella* is not present in the test culture.

## 10. LIMITATIONS OF USE

- Rough strains of *Salmonella* are known to cause non-specific autoagglutination in saline alone and therefore cannot be tested with the Oxoid Salmonella Test Kit.
- Some non-motile strains may not be detected by the Oxoid Salmonella Test Kit.
- Some oxidase-positive organisms may give false positive reactions.
- Old stock cultures of Enterobacteriaceae on Nutrient Agar slopes may cause non-specific agglutination, whereas old stocks of *Salmonella* may give false negative results. Fresh subcultures should be prepared for testing.
- Identification with the Oxoid Salmonella Test Kit is presumptive and all positive results should be confirmed by subculture of broths, further confirmatory identification tests and serotyping of pure cultures.

## 11. POSITIVE CONTROL SUSPENSION

The positive control suspension should be performed weekly, and on receipt of a new product batch to check the correct working of the latex reagents. Use a single drop of the non-viable positive control suspension as the test sample. Gently invert the positive control suspension before use. The test latex must give a positive reaction within 30 seconds.

## 12. PERFORMANCE CHARACTERISTICS

### Centre 1

561 *Salmonella* isolates were thoroughly characterised by serotyping. All these organisms, as well as 359 isolates representing 13 different non-*Salmonella* bacteria, were cultured prior to testing with the Oxoid Salmonella Test Kit at an independent laboratory.

Results are as shown below:

|              | Oxoid Salmonella Test Kit |     | Total |
|--------------|---------------------------|-----|-------|
|              | +                         | -   |       |
| Confirmed +  | 561                       | 0   | 561   |
| Salmonella - | 10                        | 349 | 359   |
| Total        | 571                       | 349 | 920   |

The performance characteristics of The Oxoid Salmonella Test Kit were; Sensitivity 100.0% and Specificity 97.2%.

The predictive value of a positive result was 98.2% and the predictive value of a negative result was 100.0%.

### Centre 2

64 food samples suspected of *Salmonella* contamination were tested for *Salmonella* at a Food Hygiene laboratory. Homogenised food was incubated for 18 hours in Resuscitation Broth, followed by 18 hours in selective enrichment broth (Cystine mannitol selenite or Rappaport Broth). Broths were subcultured onto selective agar plates for 24-48 hours. Colonies resembling *Salmonella* were tested with the Oxoid Salmonella Test Kit, or by classical biochemical tests and serotyping.

Results are shown below:

| Oxoid Salmonella Test Kit                                |      |   |
|----------------------------------------------------------|------|---|
|                                                          | +    | - |
| Salmonella confirmed by biochemical tests and serotyping | 53   | 0 |
| % Sensitivity                                            | 100% |   |
| % Specificity                                            | 100% |   |

The predictive value of both negative and positive results was 100%

13. REFERENCES

1. Data on file at Oxoid.
2. Personal Communication from Dr A. C. Baird-Parker.
3. Poster presentation: "Validation of a Latex Agglutination Method for Confirmation of *Salmonella* colonies from Six Selective Media", I.Fairlamb, T. Organ, Oxoid Ltd., Wade Road, Basingstoke, Hampshire RG24 8PW.
4. Poster presentation: "Evaluation of the Oxoid Salmonella Latex Test as a Screening Test for *Salmonella*", A. M. Paccagnella, Provincial Laboratory, BC Centre for Disease Control, Vancouver, Canada.

14. SYMBOL LEGEND

|                                                                                    |                                                     |
|------------------------------------------------------------------------------------|-----------------------------------------------------|
|    | Catalogue Number                                    |
|    | In Vitro Diagnostic Medical Device                  |
|    | Consult Instructions for Use (IFU)                  |
|    | Temperature Limitations (Storage temp.)             |
|    | Contains sufficient for <N> tests                   |
|    | Do not use if package is damaged                    |
|    | Do not re-use                                       |
|    | Batch Code (Lot Number)                             |
|    | Use By (Expiration Date)                            |
|  | Importer                                            |
|  | Unique Device Identifier                            |
|  | Authorized representative in the European Community |
|  | UK Conformity Assessed                              |
|  | European Conformity Assessment                      |
|  | Made in Australia                                   |
|  | Manufacturer                                        |

Oxoid Limited, Wade Road, Basingstoke,  
Hampshire, RG24 8PW, UK  
www.thermofisher.com

For all enquiries please contact your local distributor.

| Version | Date of modifications introduced                   |
|---------|----------------------------------------------------|
| X7990C  | November 2023<br>Updated to meet IVDR requirements |

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