

Technical Data

Cap VS KTH 45 OV grey DIA PE gasket

15.12.2016 TH

Type	Cap VS KTH 45 OV grey DIA PE gasket
Raw material	HD-PE
Color	Grey
Stamping	Description, recycling symbol
Total weight	18,3 g
Dimensions:	
Closure depth	20,20 mm ± 0,30 mm
Thread inner diameter	41,30 mm + 0,50 mm
Thread outer diameter	44,80 mm + 0,50 mm
Gasket	
Material	PE
Density	300 kg/m ³
Dimensions	Ø 41,00 mm x ø 27,80 x 3 mm
Storage:	protect the caps against longer sun light irradiation.

This information is supplied without liability. Subject to modification.