

## Deltalab Sterile Micro Tip Pasteur Pipettes (5.5ml)

**Product Code:** 210024

**Brand:** Deltalab

Made of transparent and inert polyethylene, Deltalab's Sterile Micro Tip Pasteur Pipettes (5.5ml) can be used for fragile, small drop work.

Each pipette is sterilised by ethylene oxide.



### Specifications:

|                               |              |
|-------------------------------|--------------|
| <b>Volumetrical Tolerance</b> | ± 10%        |
| <b>Material</b>               | Polyethylene |
| <b>Sterility</b>              | Sterile      |
| <b>Length (mm)</b>            | 153          |
| <b>Stem (Ø mm)</b>            | 3            |
| <b>Total Capacity (ml)</b>    | 5.5          |
| <b>Bulb Draw (ml)</b>         | 3            |
| <b>Quantity</b>               | 500          |

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)