


CATALOG NO.	NOMINAL CAPACITY (GAL)	MATERIALS	
		CARBOY	CLOSURE
002214	5	HDPE	HDPE

SIZE CODE.....	-0050
SIZE (NOM VOL L/GAL):	20/5
NOM BRIM CAP (L).....	20.2
NECK FINISH (mm).....	70
THREAD (T.P.I.).....	4
NOM. WEIGHT W/ CLOSURE(gms).....	2300
B (OVERALL HEIGHT)...	15.6 ±.4
D (WIDTH).....	12.56 ±.25
E (DEPTH).....	9.08 ±.20
W.T.(NOM WALL THK)...	0.10
T (DIA).....	2.75 ±.03
E (DIA).....	2.60 ±.03
I (DIA).....	2.25 ±.05
S (START OF THREAD)..	0.04 ±.03
L (DIM).....	0.70 ±.03
H (DIM).....	1.93 ±.04
P (PITCH).....	.250
F (FACE).....	0.06
ANGLE A.....	10°
ANGLE B.....	45°
G DIA (CLOSURE O.D.):	3.12 ±.03
J (CLOSURE HEIGHT)...	1.00 ±.02



SPECS ON BACK

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	MCC	APP'D	PVC	PSS 00-2214-0050-S	REV. 1291
CK'D	PVC	12/5/91			
DESCRIPTION:					