


CATALOG NO.	NOMINAL CAPACITY (oz)	MATERIALS			CLOSURE COLOR
		BOTTLE	CLOSURE	DRAW TUBE	
002422-2500	16	LOPE	PP	PPCO	YELLOW
002422-3500	16	LOPE	PP	PPCO	BLUE
002422-4500	16	LOPE	PP	PPCO	GREEN
002422-5500	16	LOPE	PP	PPCO	RED

SIZE (NOM VOL mL/oz)	500/16	
NOM. BRIM CAP. (mL):	570	
NECK FINISH (mm):	38-415	
THREAD (T.P.I.):	6	
NOM. WEIGHT W/CLOSURE (gms):	116	
	INCH	MM
A (NECK I.D.):	1.02 ±0.03	25.9 ±0.8
C (HT. W/O CLOSURE):	7.00 ±0.09	177.8 ±2.3
D (DIA):	2.87 ±0.08	72.9 ±2.0
H (NOM OVERALL HT):	8.8	224
W.T.(NOM. WALL THK):	0.040	1.0

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	NS	APP'D	PVC	PSS 00-2422-XXXX-1	REV. -
CK'D	PVC		5/13/93		
DESCRIPTION: (NARROW MOUTH) COLOR-CODED WASH BOTTLES					