


| ITEM         | MATERIALS    |
|--------------|--------------|
| VACUUM PLATE | PC ( WHITE ) |
| GASKETS      | NEOPRENE     |
| ADAPTER      | HDPE         |

|                        |           |       |
|------------------------|-----------|-------|
| CAT. NO. ....          | 5306-0130 |       |
| NOM. WEIGHT ( gms )... | 900       |       |
|                        | INCHES    | MM    |
| A ( DIA ).....         | 13.5      | 343   |
| B ( DIA ).....         | 12.80     | 325.1 |
| C ( DIA ).....         | 11.52     | 292.6 |
| D ( DIA ).....         | 8.32      | 211.3 |
| E ( DIA ).....         | 6.31      | 160.3 |
| F .....                | 0.10      | 2.5   |
| G .....                | 1.07      | 27.2  |
| H .....                | 14.2      | 361   |



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

**ThermoFisher**  
**SCIENTIFIC**

Thermo Fisher Scientific  
75 Panorama Creek Drive  
Rochester, NY 14625-2385, USA

|             |              |                  |           |
|-------------|--------------|------------------|-----------|
| MADE<br>NS  | APP'D<br>PVC | PSS-00-5306-0130 | REV.<br>— |
| CK'D<br>PVC | 2/28/94      |                  |           |

DESCRIPTION:  
VACUUM PLATE

NOTES:

- ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
- THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.