


CATALOG NO.....:	6120-0010	
NOM. WT.(gn)....:	3	
	INCHES	MM
A (OVERALL)....:	2.20 ±.03	55.9 ±0.8
B DIA.....:	0.57 ±.02	14.5 ±0.5

COMPONENTS

ITEM	MATERIAL
FEMALE HALF	HDPE
MALE HALF	HDPE
VALVE DISC	SILICONE RUBBER

NOTE:

- FITS 1/4" THRU 5/16" I.D. TUBING
- NOT FOR USE WITH LIQUIDS



SPECS ON BACK

NOT TO SCALE

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	MC	APP'D	MP	PSS 00-6120-0010	REV. 1002
CK'D	MP		10/28/02		

DESCRIPTION:

CHECK VALVE AND POSITIVE CONNECTOR