


CATALOG NUMBER.....	7211-0500	
MATERIAL.....	PP	
	INCH	MM
NOM. HOSE BARB SIZE:	1/2	13
A HOSE BARB I.D.....	0.31	7.8
B DIA.....	0.53	13.5
C (BARB LENGTH).....	0.93	23.6
D (OVERALL LENGTH)..:	2.02	51.3
E DIA.....	1.98	50.3

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	MC ²	APP'D PVC	PSS 00-7211-0500	REV. 0502
CK'D	PVC	5/1/02		
DESCRIPTION: 1" LADISH TRI-CLAMP x HOSE BARB SANITARY ADAPTER				