


CATALOG NO.	NOMINAL CAPACITY, ( oz )	MATERIALS	
		BOTTLE	CLOSURE
312002	6	HDPE	PP

SIZE CODE.....:	-0006	
SIZE (NOM VOL mL/oz):	175/6	
NOM. BRIM CAP. (mL):	200	
NECK FINISH (mm)....:	24-415	
THREAD (T.P.I.).....:	8	
NOM. WEIGHT W/CLOSURE (gms).....:	27	
	INCH	MM
A (NECK I.D.).....:	0.69 ±0.02	17.5 ±0.5
B (HT. W/CLOSURE)....:	4.90 ±0.08	124.5 ±2.0
C (HT. W/O CLOSURE):	4.81 ±0.08	122.2 ±2.0
D (DIA).....:	2.13 ±0.06	54.1 ±1.5
W.T.(NOM. WALL THK):	0.040	1.02

NOTES:

- ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
- THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

**ThermoFisher**  
**SCIENTIFIC**

Thermo Fisher Scientific  
75 Panorama Creek Drive  
Rochester, NY 14625-2385, USA

MADE MCC		APP'D RAJ	PSS 31-2002-0006-S	REV. 0995
CK'D DM	9/28/95	9/28/95		
DESCRIPTION: ROUND NARROW MOUTH LAB QUALITY BOTTLE				