

# 7900HT System Spectral Calibration Kit, Fast 96-well

Catalog Number 4351653

Pub. No. 4351651 Rev. B

**WARNING!** Read the Safety Data Sheets (SDSs) and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves. Safety Data Sheets (SDSs) are available from [thermofisher.com/support](http://thermofisher.com/support).

## Kit contents and storage

Contents	Quantity	Storage
Background Plate sealed with an optical cover	1	Store at -25°C to -15°C
Pure Dye Plates sealed with optical covers	2	

## Pure Dye Plate layout

	FAM <sup>TM</sup> STD			JOE <sup>TM</sup> STD			NED <sup>TM</sup> STD			ROX <sup>TM</sup> STD		
	1	2	3	4	5	6	7	8	9	10	11	12
A												
B												
C												
D												
E												
F												
G												
H												

Pure Dye Plate 1

	SYBR <sup>TM</sup> GREEN STD			TAMRA <sup>TM</sup> STD			TET <sup>TM</sup> STD			VIC <sup>TM</sup> STD		
	1	2	3	4	5	6	7	8	9	10	11	12
A												
B												
C												
D												
E												
F												
G												
H												

Pure Dye Plate 2

## Related documentation

For detailed information on instrument setup and the calibration process, refer to the *Applied Biosystems™ 7900HT Fast Real-Time PCR System User Bulletin* (Pub. no. 4352533).

## Limited product warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale found on Life Technologies' website at [www.thermofisher.com/us/en/home/global/terms-and-conditions.html](http://www.thermofisher.com/us/en/home/global/terms-and-conditions.html). If you have any questions, please contact Life Technologies at [www.thermofisher.com/support](http://www.thermofisher.com/support).

The information in this guide is subject to change without notice.

**DISCLAIMER**

TO THE EXTENT ALLOWED BY LAW, LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE, OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.

Important Licensing Information: These products may be covered by one or more Limited Use Label Licenses. By use of these products, you accept the terms and conditions of all applicable Limited Use Label Licenses.

Corporate entity: Life Technologies | Carlsbad, CA 92008 USA | Toll Free in USA 1.800.955.6288

©2015 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.

For support visit [thermofisher.com/support](http://thermofisher.com/support) or email [techsupport@lifetech.com](mailto:techsupport@lifetech.com)  
thermofisher.com

**ThermoFisher**  
S C I E N T I F I C