
InstantOne ELISA Enhancer Solution

Catalog Number: IOES

RUO: For Research Use Only. Not for use in diagnostic procedures.

Product Information

Contents: InstantOne ELISA Enhancer Solution



Catalog Number: IOES

Handling Conditions: Store at 2-8C



Temperature Limitation: Store at 2-8C



Batch Code: Refer to vial



Use By: Refer to vial

Description

The InstantOne ELISA Enhancer Solution is a supplement added to the 5X Cell Lysis Buffer (Cat. No. IOCLB1) to yield the Cell Lysis Mix solution that is applied to cultured cells or tissues prior to running the assays. The 1X Cell Lysis Mix is prepared by simply diluting Cell Lysis Mix (5X) (a mixture of the Cell Lysis Buffer (5X) and the Enhancer Solution) 5-fold with deionized water. This buffer is used to lyse cells after the removal of culture medium, and is typically used to lyse adherent cells or non-adherent cells that have been harvested by centrifugation. Lysis Mix should be used as the diluent for any dilution of cellular lysates that are required.

Not for further distribution without written consent.

Copyright © 2016 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.

Tel: 888.999.1371 or 858.642.2058 • Fax: 858.642.2046 • thermofisher.com/ebioscience • info@ebioscience.com