

# Recombinant Human Noggin

Catalog Number PHC1506 (20 µg)

Pub. No. MAN0003548 Rev. A.0








## Product specifications

<b>Lot number</b>	See product label.
<b>Molecular weight</b>	Non-glycosylated, non-disulfide-linked homodimer consisting of two 206 amino acid polypeptide chains, with a total molecular weight of 46.2 kDa (each subunit having a molecular weight of 23.1 kDa).
<b>Purity</b>	>95% as determined by SDS PAGE analysis.
<b>Biological activity</b>	ED <sub>50</sub> <3 ng/mL, determined by measuring the dose dependent inhibition of alkaline phosphatase production induced by 5 ng/mL of BMP-4 in ATDC-5 chondrogenic cells. Determine the optimal concentration for each specific application using an initial dose response assay.
<b>Formulation</b>	Lyophilized, carrier free.
<b>Sterility</b>	Filtered before lyophilization through a 0.22 micron sterile filter.
<b>Endotoxin</b>	<0.1 ng/µg
<b>Production</b>	Produced in <i>E. coli</i> and purified via sequential chromatography.
<b>Reconstitution recommendation</b>	Centrifuge the vial briefly, before opening to bring the contents to the bottom. Reconstitute the lyophilized protein in 10 mM acetic acid to a concentration of 0.1–1.0 mg/mL. Apportion the reconstituted protein into working aliquots and store at ≤ –20°C. Make any further dilutions of the reconstituted protein in appropriate buffered solution containing a carrier protein, such as 0.1% HSA or BSA.
<b>Suggested working dilutions</b>	The optimal concentration should be determined for each specific application.
<b>Storage</b>	The lyophilized protein is stable at room temperature for 3 weeks. Store at –20°C, preferably desiccated. Upon reconstitution, apportion into working aliquots and store at 4°C for 2–7 days, or at –20°C. Avoid repeated freeze-thaw cycles.
<b>Expiration date</b>	Expires one year from date of receipt when stored as instructed.
<b>References</b>	Que, J, et al. (2006) Morphogenesis of the trachea and esophagus: current players and new roles for noggin and Bmps. <i>Differentiation</i> 74(7):422-437. Yanagita, M. (2005) BMP antagonists: their roles in development and involvement in pathophysiology. <i>Cytokine Growth Factor Rev.</i> 16(3):309-317. Chen, D, et al. (2004) Bone morphogenetic proteins. <i>Growth Factors</i> 22(4):233-241. Sebald, W, et al. (2004) Molecular recognition in bone morphogenetic protein (BMP)/receptor interaction. <i>Biol. Chem.</i> 385(8):697-710.

## Limited product warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale found on Life Technologies' website at [www.thermofisher.com/us/en/home/global/terms-and-conditions.html](http://www.thermofisher.com/us/en/home/global/terms-and-conditions.html). If you have any questions, please contact Life Technologies at [www.thermofisher.com/support](http://www.thermofisher.com/support).

## Explanation of Symbols

Symbol	Description	Symbol	Description	Symbol	Description
	Manufacturer		Catalog number		Batch code
	Use by		Temperature limitation		
	Consult instructions for use		Caution, consult accompanying documents		



**Manufacturer:** Life Technologies Corporation | 5781 Van Allen Way | Carlsbad, CA 92008

The information in this guide is subject to change without notice.

**DISCLAIMER:** TO THE EXTENT ALLOWED BY LAW, THERMO FISHER SCIENTIFIC INC. AND/OR ITS AFFILIATE(S) WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE, OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.

**Important Licensing Information:** This product may be covered by one or more Limited Use Label Licenses. By use of this product, you accept the terms and conditions of all applicable Limited Use Label Licenses.

©2018 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.