

# Recombinant Human Platelet Derived Growth Factor-BB (PDGF-BB)

Catalog Number PHG0044 (5 µg), PHG0045 (10 µg), PHG0046 (50 µg), PHG0041 (100 µg), PHG0043 (1 mg)

Pub. No. MAN0003588 Rev. B.0








## Product specifications

<b>Lot number</b>	See product label.
<b>Molecular weight</b>	25.4 kDa, homodimer
<b>Purity</b>	>95% as determined by SDS PAGE analysis.
<b>Amino acid sequence</b>	SLGSLTIAEP AMIAECKTRT EVFEISRRLI DRTNANFLWW PPCVEVQRCS GCCNNRNVQC RPTQVQLRPV QVRKIEIVRK KPIFKKATVT LEDHLACKCE TVAAARPVT
<b>Biological activity</b>	ED <sub>50</sub> <5.00 ng/mL, determined by the dose dependent proliferation of BALB/c 3T3 cells. Determine the optimal concentration for each specific application using an initial dose response assay.
<b>Formulation</b>	Lyophilized, carrier free.
<b>Sterility</b>	The protein is eluted in acetonitrile and then lyophilized under aseptic conditions.
<b>Endotoxin</b>	<0.1 ng/µg
<b>Production</b>	Produced in <i>E. coli</i> and purified via sequential chromatography.
<b>Reconstitution recommendation</b>	Centrifuge the vial briefly, before opening to bring the contents to the bottom. Reconstitute the lyophilized protein in 100 mM acetic acid containing 0.1% BSA to 0.1–1.0 mg/mL to regain full activity. Apportion the reconstituted protein into working aliquots and store at ≤ –20°C. Make any further dilutions of the reconstituted protein in low endotoxin medium or buffered solution with FBS or tissue culture grade BSA.
<b>Suggested working dilutions</b>	The optimal concentration should be determined for each specific application.
<b>Storage</b>	Store the lyophilized protein at 2–8°C, preferably desiccated. Upon reconstitution, apportion into working aliquots and store at ≤ –20°C. Avoid repeated freeze-thaw cycles.
<b>Expiration date</b>	Expires one year from date of receipt when stored as instructed.
<b>References</b>	<p>Jin, P, et al. (1991) Recombinant platelet-derived growth factor-BB stimulates growth and inhibits differentiation of rat L6 myoblasts. <i>J. Biol. Chem.</i> 266(2):1245-1249.</p> <p>Dong, F, et al. (1998) Cyclin D3-associated kinase activity is regulated by p27kip1 in BALB/c 3T3 cells. <i>Mol. Biol. Cell.</i> 9(8):2081-2092.</p> <p>Dong, F, et al. (1998) The role of cyclin D3-dependent kinase in the phosphorylation of p130 in mouse BALB/c 3T3 fibroblasts. <i>J. Biol. Chem.</i> 273(11):6190-6195.</p> <p>Hauser, PJ, et al. (1997) p107 and p130 associated cyclin A has altered substrate specificity. <i>J. Biol. Chem.</i> 272(36):22954-22959.</p> <p>Zhang, X, Wharton, W, Donovan, M, Coppola, D, Croxton, R, Cress, WD, and Pledger, WJ. (2000) Density dependent growth inhibition of fibroblasts ectopically expressing p27kip1. <i>Mol. Biol. Cell</i> 11(6):2117-2130.</p> <p>Wharton, W, et al. (2000) Inhibition of mitogenesis in BALB/c-3T3 cells in trichostatin A. <i>J. Biol. Chem.</i> 275(43):33981-33987.</p> <p>Koh, JS, et al. (2000) Cytokine dysregulation induced by apoptotic cells is a shared characteristic of murine lupus. <i>J. Immunol.</i> 165:4190-4201.</p>

## Limited product warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale at [www.thermofisher.com/us/en/home/global/terms-and-conditions.html](http://www.thermofisher.com/us/en/home/global/terms-and-conditions.html). If you have any questions, please contact Life Technologies at [www.thermofisher.com/support](http://www.thermofisher.com/support).

## Explanation of Symbols

Symbol	Description	Symbol	Description	Symbol	Description
	Manufacturer		Catalog number		Batch code
	Use by		Temperature limitation		
	Consult instructions for use		Caution, consult accompanying documents		



Life Technologies Corporation | 5781 Van Allen Way | Carlsbad, CA 92008

For descriptions of symbols on product labels or product documents, go to [thermofisher.com/symbols-definition](http://thermofisher.com/symbols-definition).

The information in this guide is subject to change without notice.

**DISCLAIMER:** TO THE EXTENT ALLOWED BY LAW, THERMO FISHER SCIENTIFIC INC. AND/OR ITS AFFILIATE(S) WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE, OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.

**Important Licensing Information:** This product may be covered by one or more Limited Use Label Licenses. By use of this product, you accept the terms and conditions of all applicable Limited Use Label Licenses.

©2020 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.