

ThermoFisher S C I E N T I F I C	Document Owner Department: QC	BT-SPEC-0641
		Page 1 of 2
OXOID QUALITY ASSURANCE PRODUCT SPECIFICATION		
SPECTINOMYCIN SH10 CT0046B		

Product Details

Product Name	SPECTINOMYCIN
Disc Code	SH10
Disc Concentration	10µg
CT Code	CT0046B
Shelf Life (years)	3
CLSI Recommended Concentration	No
EUCAST Recommended Concentration	No

Microbiological Performance Criteria


Organism	ATCC/NCTC Number	Range (mm)
<i>Escherichia coli</i>	ATCC® 25922	11 – 14
<i>Klebsiella pneumoniae</i>	ATCC® 29665	11 – 17
<i>Haemophilus influenzae</i>	ATCC® 49247	9 – 15
<i>Streptococcus pneumoniae</i>	ATCC® 49619	9 – 14
<i>Staphylococcus aureus</i>	ATCC® 25923	0 – 0
<i>Pseudomonas aeruginosa</i>	ATCC® 27853	0 – 0

References

Tested in accordance with the current CLSI methodology, the control ranges are internally derived based on historical data.

Other additional challenging strains may be employed.

The information given above is believed to be correct, however both the information and the product are offered without warranty for any specific use, other than that specified in the current Oxoid manual.

	Document Owner Department: QC	BT-SPEC-0641
		Page 2 of 2
OXOID QUALITY ASSURANCE PRODUCT SPECIFICATION		
SPECTINOMYCIN SH10 CT0046B		

Revision History

Section / Step	Description of Change	Reason for Change	Reference
Entire document/ References	Update to current format / Update of current performance statement	Change control	BT-CC-2887