	Document Owner Department: QC	BT-SPEC-0708
		Page 1 of 2
OXOID QUALITY ASSURANCE PRODUCT SPECIFICATION		
CEFTIZOXIME ZOX30 CT0477B		

Product Details

Product Name	CEFTIZOXIME
Disc Code	ZOX30
Disc Concentration	30µg
CT Code	CT0477B
Shelf Life (years)	3
CLSI Recommended Concentration	Yes
EUCAST Recommended Concentration	No

Microbiological Performance Criteria

CLSI		
Organism	ATCC/NCTC Number	Range (mm)
<i>Escherichia coli</i>	ATCC®25922	30 – 36
<i>Staphylococcus aureus</i>	ATCC®25923	27 – 35
<i>Pseudomonas aeruginosa</i>	ATCC®27853	12 – 17
<i>Haemophilus influenzae</i>	ATCC®49247	29 – 39
<i>Neisseria gonorrhoeae</i>	ATCC®49226	42 – 51
<i>Streptococcus pneumoniae</i>	ATCC®49619	28 – 34


References

For recommended concentrations, the zone limits are in accordance with the current performance standards for antimicrobial disc susceptibility tests;

- CLSI: Performance standards for antimicrobial susceptibility testing

Other additional challenging strains may be employed.

The information given above is believed to be correct, however both the information and the product are offered without warranty for any specific use, other than that specified in the current Oxoid manual.

	Document Owner Department: QC	BT-SPEC-0708
		Page 2 of 2
OXOID QUALITY ASSURANCE PRODUCT SPECIFICATION		
CEFTIZOXIME ZOX30 CT0477B		

Revision History

Section / Step	Description of Change	Reason for Change	Reference
Entire document/ / References	Update to current format and minor typographical errors / Update of current performance standards	Change control	BT-CC-2856