

ThermoFisher S C I E N T I F I C	Document Owner Department: QC	BT-SPEC-0749
		Page 1 of 2
OXOID QUALITY ASSURANCE PRODUCT SPECIFICATION		
LINEZOLID LZD30 CT1650B		

Product Details

Product Name	LINEZOLID
Disc Code	LZD30
Disc Concentration	30µg
CT Code	CT1650B
Shelf Life (years)	2
CLSI Recommended Concentration	Yes
EUCAST Recommended Concentration	No

Microbiological Performance Criteria

CLSI:		
Organism	ATCC/NCTC Number	Range (mm)
<i>Staphylococcus aureus</i>	ATCC®25923	25 – 32
<i>Streptococcus pneumoniae</i>	ATCC®49619	25 – 34


References

For recommended concentrations, the zone limits are in accordance with the current performance standards for antimicrobial disc susceptibility tests;

- CLSI: M100 Performance standards for antimicrobial susceptibility testing

Other additional challenging strains may be employed.

The information given above is believed to be correct, however both the information and the product are offered without warranty for any specific use, other than that specified in the current Oxoid manual.

	Document Owner Department: QC	BT-SPEC-0749
		Page 2 of 2
OXOID QUALITY ASSURANCE PRODUCT SPECIFICATION		
LINEZOLID LZD30 CT1650B		

Revision History

Section / Step	Description of Change	Reason for Change	Reference
Entire document/ References	Update to current format/ Update to current performance standards	Change control	BT-CC-2820