


Certificate of Quality Thermo Scientific ClipTip™

Product Code: 94410150
Product Name: ClipTip 384 125, Non-filtered, Non-sterile, Rack
Lot Number: 1360266
Expiry date: 2027-06

Thermo Fisher Scientific guarantees that Thermo Scientific ClipTip™ Pipette tips meet the industry's highest standards for molding precision and clarity. Our operation system is certified and conforms to ISO 9001 and ISO 14001. ClipTip pipette tips are manufactured under controlled environment and are certified by validation of the manufacturing process to be free of RNase, DNase, ATP and endotoxin contamination.

DNA	< 30 pg
RNase	< 10 ⁻⁹ Kunitz units/μl
DNase	< 10 ⁻⁷ Kunitz units/μl
Endotoxins	< 0.06 EU/ml
ATP	< 10 ⁻¹³ mg/μl

In Monterrey, 2022



Ana Edith Martinez
Quality Manager

Made in Mexico



Thermo Fisher Scientific Oy
Ratastie 2
01620 Vantaa
Finland